

## Adhesion (Crystal Violet assay)

**Warning: Crystal Violet is a suspected carcinogen. Crystal Violet will stain your skin or clothes by contact and cannot easily be removed. Use gloves at all times when handling Crystal Violet.**

With this protocol you can compare adhesion ability/biofilm forming ability of isolated clones/mutants. The adhesion assay can be used as the first screening of a series of strains before investigating the most interesting ones in more details in our flow-chamber system.

### Day 1:

Inoculate strains you want to analyze in 10 ml test tubes. Use dilute media to avoid excessive growth.

### Day 2 (or when the plats have suitable growth):

Inoculate microtiter plate with the different strains of interest. Use 150 µl medium for each well. Use suitable microtiter plates. U-bottom or flat-bottom are suited. 0.1% Crystal Violet solution [0.1 g Crystal Violet in 100 ml water].

1. Use 2 microtiter plates for each set of clones
2. Place the inoculated microtiter plate in a thick plastic bag and incubate over night at 37°C or 30°C (depending on strain etc., also, some clinical isolates needs two/three days' incubation).

### Day 3:

3. After incubation empty all the wells by simply throwing out the liquid in a clinical waste bag (without using a pipette).
4. Transfer pre-warmed (37°C or 30°C) medium or 0.9% NaCl to a sterile Petri dish and use a multi-pipette to gently add 200 µl (avoid overloading) medium/NaCl to each well, and discard it by throwing. Repeat this washing step twice.
5. Add 200 µl (avoid overloading) Crystal Violet solution to each well (**use gloves**) and let stand on the table for 15-20 min.
6. Discard the Crystal Violet solution by throwing (**use gloves, beware of splashes**).
7. Wash three times with water (3x200 µl (avoid overloading)). Again discard by throwing.
8. Add 200 µl (avoid overloading) 96% ethanol; pipette up and down thoroughly in order to dissolve the Crystal Violet (use multi-pipette).
9. Read the plate in an Elisa reader using an excitation of 585 nm and make a column diagram data using spreadsheet software.

### Reference:

Reisner A, Krogfelt KA, Klein BM et al (2006) In vitro biofilm formation of commensal and pathogenic *Escherichia coli* strains: impact of environmental and genetic factors. *J Bacteriol* 188: 3572-3581