

Exercise 1: Molecular Diagnostics

Teachers:

Søren Damkiær

Pedersen

Claus Sternberg



Exercise 1

- Purpose
 - To test several different diagnostic methods:
 - FISH
 - 16S rRNA RFLP
 - Sequencing of 16s rRNA
 - Biochemical tests (BIOLOG and antibiotic sensitivity)
- Each team will identify 2 clinically isolates + one isolate from your body or the environment

Introduction

- Why identify bacteria:
 - Targeted treatment of infected patients
 - Test and control of food products
- How bacteria are identified:
 - Phenotypic
 - Genotypic

Phenotypic identification

- Phenotypic identification:
 - Microscopy
 - Biochemical assays
 - testing for specific metabolisms or other traits
 - Staining
 - E. g. Gram staining


Genotypic identification

- Methods based on bacterial DNA/RNA
 - Fluorescent in situ hybridisation (FISH)
 - 16S rRNA restriction fragment length polymorphism (RFLP)
 - rRNA sequencing
- Advantages:
 - Fast, do not require cultivation, precise

FISH

- Based on 16S rRNA
- All bacteria have 16S rRNA
 - Conserved regions
 - Variable regions

16S - rRNA *Escherichia coli*




FISH

- Hybridisation of fluorescent probes complimentary to 16s rRNA

● GATTCCTAGATTAGGTCGT

- Different levels:


– Domain	Bacteria
– Phylum	Proteobacteria
– Class	Gamma Proteobacteria
– Order	Pseudomonadales
– Family	Pseudomonadaceae
– Genus	<i>Pseudomonas</i>
– Species	<i>aeruginosa</i>
– Subspecies	
– Strain	PAO1



FISH

- Different levels of probes:


- Domain
- Phylum
- Class
- Order
- Family
- Genus
- Species
- Subspecies
- Strain



Probe	Sequence 5' → 3'	Label
β-proteobacteria	GCC TTC CCA CTT CGT TT	CY5 (Blue)
γ-proteobacteria	GCC TTC CCA CAT CGT TT	CY3 (Red)
Bacteria	GCT GCC TCC CGT AGG AGT	FITC (green)


FISH

- Cells are fixed and made permeable
- Probes are hybridised
- The result is viewed using epifluorescence microscopy



16S rRNA RFLP

- Restriction Fragment Length Polymorphism
 - The gene encoding 16s rRNA is amplified by PCR
 - The PCR fragment is cut with restriction enzymes
 - Cut DNA is visualised using gel electrophoresis
 - Fragment pattern = fingerprint of the organism



rDNA sequencing

- Using 16S rDNA sequencing, the strain or its relation can be determined specifically
- The sequencing usually takes place externally
- The sequence is compared with sequence databases on the internet
- Search tool: BLAST
- use: <http://www.ncbi.nlm.nih.gov/>

NCBI - Blast

Basic BLAST

Choose a BLAST program to run.



nucleotide blast

Search a **nucleotide** database using a **nucleotide** query

Algorithms: blastn, megablast, discontinuous megablast

protein blast

Search **protein** database using a **protein** query

Algorithms: blastp, psi-blast, phi-blast

blastx

Search **protein** database using a **translated nucleotide** query

tblastn

Search **translated nucleotide** database using a **protein** query

tblastx

Search **translated nucleotide** database using a **translated nucleotide** query

Enter Query Sequence

Enter accession number, gi, or FASTA sequence [Clear](#)

[?](#) Query subrange

From

To

Or, [?](#) load file

Browse...

[?](#) Job title

Enter a descriptive title for your BLAST search

Choose Search Set

Database

☐ Human genomic + transcript ☐ Mouse genomic + transcript ☒ Others (nr etc.):

Nucleotide collection (nr/nt)

Organism

Enter organism name or id—completions will be suggested

[Optional](#)

Enter organism common name, binomial, or tax id. Only 20 top taxa will be shown.

[?](#) Entrez Query

[Optional](#)

Enter an Entrez query to limit search

Program Selection

BLAST

Search database nr using **Megablast** (Optimize for highly similar sequences)

Blast


Sequences producing significant alignments:		Score (bits)	E Value
gi 41562 emb X05173.1 ECGLN	E.coli glnALG operon with glnA ...	952	0.0
gi 304961 gb L19201.1 ECOUW87	E. coli chromosomal region fr...	952	0.0
gi 1790295 gb AE000462.1 AE000462	Escherichia coli K12 MG16...	952	0.0
gi 13364198 dbj AP002567.1 	Escherichia coli O157:H7 DNA, c...	928	0.0
gi 12518741 gb AE005617.1 AE005617	Escherichia coli O157:H7...	928	0.0
gi 146156 gb M13746.1 ECOGLNAB	E.coli glnA gene encoding gl...	920	0.0
gi 24054430 gb AE015401.1 	Shigella flexneri 2a str. 301 se...	912	0.0
gi 30042943 gb AE016990.1 	Shigella flexneri 2a str. 2457T ...	912	0.0
gi 26111017 gb AE016770.1 	Escherichia coli CFT073 section ...	864	0.0
gi 16422561 gb AE008887.1 	Salmonella typhimurium LT2, sect...	500	e-138
gi 154089 gb M14536.1 STYGLNA	S.typhimurium glnA gene encod...	500	e-138
gi 29139451 gb AE016846.1 	Salmonella enterica subsp. enter...	492	e-136
gi 16504729 emb AL627280.1 	Salmonella enterica serovar Typ...	492	e-136
gi 3808289 gb AF072440.1 AF072440	Enterobacter gergoviae GT...	468	e-129
gi 23630224 gb AY144623.1 	Pantoea agglomerans glutamine sy...	389	e-105
gi 21960789 gb AE013983.1 	Yersinia pestis KIM section 383 ...	309	3e-81
gi 15978115 emb AJ414141.1 	Yersinia pestis strain CO92 com...	309	3e-81
gi 45434720 gb AE017127.1 	Yersinia pestis biovar Mediaevai...	309	3e-81
gi 146158 gb J01618.1 ECOGLNACR	Escherichia coli strain K-1...	256	4e-65
gi 33087906 gb AY333022.1 	Yersinia bercovieri strain WS 39...	230	2e-57
gi 33087902 gb AY333020.1 	Yersinia bercovieri strain WS 19...	226	3e-56
gi 33087898 gb AY333018.1 	Yersinia bercovieri strain WA 17...	226	3e-56
gi 33087896 gb AY333017.1 	Yersinia bercovieri strain ATCC ...	226	3e-56
gi 33087952 gb AY333045.1 	Yersinia kristensenii strain WA ...	202	5e-49
gi 33087916 gb AY333027.1 	Yersinia enterocolitica strain 1...	202	5e-49

Biochemical Tests

- BIOLOG test
 - Based on the ability of the bacteria to utilize different carbon sources
- Antibiotic sensitivity
 - Important for treatment and prognosis
 - Otherwise "identical" strains can differ in AB pattern

Biolog GEN III system

96 wells microplate assay containing different C-sources and antibiotics



Growth pattern
identify the bacteria
down to species level

GEN III MicroPlate™


A1 Negative Control	A2 Dextrin	A3 D-Maltose	A4 D-Trehalose	A5 D-Cellobiose	A6 Gentibiose	A7 Sucrose	A8 D-Turanose	A9 Stachyose	A10 Positive Control	A11 pH 6	A12 pH 5
B1 D-Raffinose	B2 α-D-Lactose	B3 D-Melibiose	B4 β-Methyl-D-Glucoside	B5 D-Salicin	B6 N-Acetyl-D-Glucosamine	B7 N-Acetyl-β-D-Mannosamine	B8 N-Acetyl-D-Galactosamine	B9 N-Acetyl Neuraminic Acid	B10 1% NaCl	B11 4% NaCl	B12 8% NaCl
C1 α-D-Glucose	C2 D-Mannose	C3 D-Fructose	C4 D-Galactose	C5 3-Methyl Glucose	C6 D-Fucose	C7 L-Fucose	C8 L-Rhamnose	C9 Inosine	C10 1% Sodium Lactate	C11 Fusidic Acid	C12 D-Serine
D1 D-Sorbitol	D2 D-Mannitol	D3 D-Arabitol	D4 myo-Inositol	D5 Glycerol	D6 D-Glucose-6-PO4	D7 D-Fructose-6-PO4	D8 D-Aspartic Acid	D9 D-Serine	D10 Troleandomycin	D11 Rifamycin SV	D12 Minocycline
E1 Gelatin	E2 Glycyl-L-Proline	E3 L-Alanine	E4 L-Arginine	E5 L-Aspartic Acid	E6 L-Glutamic Acid	E7 L-Histidine	E8 L-Pyrogutamic Acid	E9 L-Serine	E10 Lincomycin	E11 Guanidine HCl	E12 Niaproof 4
F1 Pectin	F2 D-Galacturonic Acid	F3 L-Galactonic Acid Lactone	F4 D-Gluconic Acid	F5 D-Glucuronic Acid	F6 Glucuronamide	F7 Mucic Acid	F8 Quinic Acid	F9 D-Saccharic Acid	F10 Vancomycin	F11 Tetrazolium Violet	F12 Tetrazolium Blue
G1 p-Hydroxy-Phenylacetic Acid	G2 Methyl Pyruvate	G3 D-Lactic Acid Methyl Ester	G4 L-Lactic Acid	G5 Citric Acid	G6 α-Keto-Glutaric Acid	G7 D-Malic Acid	G8 L-Malic Acid	G9 Bromo-Succinic Acid	G10 Nalidixic Acid	G11 Lithium Chloride	G12 Potassium Tellurite
H1 Tween 40	H2 n-Amino-Butyric Acid	H3 α-Hydroxy-Butyric Acid	H4 β-Hydroxy-D-L-Butyric Acid	H5 α-Keto-Butyric Acid	H6 Acetoacetic Acid	H7 Propionic Acid	H8 Acetic Acid	H9 Formic Acid	H10 Aztreonam	H11 Sodium Butyrate	H12 Sodium Bromate

Antibiotic Testing

There are several ways of determining the sensibility to different antibiotics.

They will be addressed in Ex 2.

In Ex 1 we use "Rosco- tabs"



Exercise 1

- Remember:
 - Plan the experiments
 - Use GLOVES
 - Work sterile
 - All cell-containing liquids, supernatants etc. must be disposed in the culture disposal flasks in the fume hood.
- Fixation and hybridisation of cells: **MUST** be in the **FUME HOOD**

Exercise 1. Flow chart

Day	1 A FISH	1 B PCR	1 C RFLP	1 D Antibiotic test	1 E BIOLOG ID
Monday 10/1	Inoculate liquid media (also used in 1B)			Inoculate plate media (BA). Sample own isolate.	
Tuesday 11/1	Fix cells (do this first today). FISH Procedure.	Prepare Chromosomal DNA			
Wednesday 12/1	Microscopy of FISH samples. Can be continued Thursday	16S PCR amplification. (also used in 1C)		AB test plate preparation.	
Thursday 13/1			Restriction enzyme cut PCR products	AB and Blue plate inoculation	Inoculate BIOLOG plates
Friday 14/1			Run agarose gel Take Gel photo	Read Plates	Read BIOLOG plates

Preparation of Chromosomal DNA

1. Harvest 500 μ l culture : Spin down cells (5 min at 7000 g) and decant.
2. Wash with 1 ml TNE : Resuspend in 1 ml TNE; vortex; spin down cells (5 min at 7000 g); decant.
3. Resuspend pellet in 270 μ l TNEX.
4. Add 30 μ l of a freshly prepared lysozyme solution (5 mg/ml in H₂O).
5. Add 7.5 μ l Proteinase K solution (20 mg/ml in H₂O).
6. Incubate for 90 min at 37°C and then at 65°C for another 90 min.
7. Add 15 μ l 5 M NaCl and mix gently by inverting the tube.
8. Add gently 1 ml 96% EtOH. After 2-3 min mix gently and incubate 15 min at -20°C.
9. Spin at 15000 g for 10 min at 4°C. —————> DNA visible
10. Wash pellet with 1 ml ice cold 96% EtOH. Be careful! Don't loose the DNA!
11. Spin at 15000 g for 10 min at 4°C.
12. Remove remaining EtOH by pipetting followed by incubation at 37°C for 3-5 min with the lid open (do not over dry pellet).
13. Add 500 μ l dH₂O and mix gently by inverting the tube. Store at 20°C o.n. to dissolve the DNA; then at -20°C.


Don't vortex chromosomal DNA or it will shear.

FISH day 2

- 2 Fume hoods + 1 incubator chamber 37 degree
- Fixation of bacteria:
 - Fixative waste in the container
 - Use blue (nitril) gloves
 - Recommended to continue right away with the hybridization step

FISH day 2

- Hybridization step:
 - Probes are ready for use
 - 3 probes separately, so three wells per strain!
 - One 37 degree incubator for the 3h hybridization



- Slides are collected when finished (remember name)

FISH day 3

- Visualization of hybridized cells using epifluorescence microscopy (1. floor, room 125)
- Time schedule starting from 9:30 and running for the rest of the day
- New team roughly every 30-45 minutes

Team 1: 9:30

Team 2: 10:15

Team 3: 11:00

Team 4: 11:45

Team 5: 13:00

Team 6: 13:45

Team 7: 14:30

Team 8: 15:00

Team 9: 15:45

Team 10: 16:30