

BiOLOG

MicroStation[™] System
MicroLog
Version 4.2

User Guide



© 2007 Biolog, Inc.
All rights reserved

Trademarks

OmniLog is a registered trademark, MicroStation, MicroLog, MicroPlate, Streakerz, LongSwabs, RetroSpect are trademarks of Biolog, Inc.

Microsoft Windows, Windows XP, Excel and Access are either registered trademarks of Microsoft Corporation in the United States and/or other countries.

Pentium is a registered trademark of Intel Corporation.

Macintosh is a registered trademark of Apple, Inc.

InstallShield is a registered trademark of Macrovision Corporation.

BDFalcon is a trademark of Becton Dickenson Co.

Phytigel is a trademark of Sigma-Aldrich Biotechnology LP.

Technical Service

For technical and sales assistance, contact Biolog at:

| | |
|------------------|---|
| Address: | 21124 Cabot Blvd Hayward, CA 94545-1130 U.S.A. |
| Tel: | 510-785-2564 (M-F, 7:30 a.m. to 5:00 p.m., PST) |
| Fax: | 510-782-4639 |
| email: | tech@biolog.com |
| web page: | http://www.biolog.com |

MicroStation System, MicroLog Version 4.2 User Guide

Part number: 35510

Revisions: 2003; 2004; November 2007

Not for human *in vitro* diagnostic use.

LICENSE AGREEMENT

Notice: This is a legal agreement between you, an individual or entity ("Licensee"), and Biolog, Inc. ("Licensor"). By purchasing or using this software, whether it is on a diskette or preinstalled on a computer, you agree to be bound by this agreement. If you do not accept this agreement, please immediately return the unopened software case containing the complete contents, as well as any computer on which such contents were preinstalled, to Licensor or other place of purchase.

License: Licensor grants Licensee the right to use the enclosed copy of Licensor's MicroLog software (the "Software") at the single microbiology laboratory (at one physical location) for which this license is granted (the "Site"). The Software shall include all related materials and documentation provided by Licensor. The Software shall be used only in the conduct of Licensee's own business and Licensee shall not permit any third party to use the Software. Licensee shall make available implementation computer equipment and software configurations approved by Licensor as adequate for implementation of the Software.

Ownership: The Software and all modifications thereto and all copies thereof are proprietary to Licensor and title thereto remains in Licensor. All applicable rights to patents, copyrights, know-how, trademarks, and trade secrets in the Software and any modifications made thereto are and shall remain in Licensor. Licensee shall not sell, transfer, publish, disclose, display or otherwise make available the Software or copies thereof to others. Licensee agrees to secure and protect each module, software product, and related documentation in a manner consistent with the maintenance of Licensor's rights therein and to take appropriate actions by instruction or agreement with its employees or consultants who are permitted access to the Software to satisfy its obligations hereunder. Any copies made by Licensee of the Software and other programs developed hereunder, including translations, compilations, partial copies with modifications and updated works, are the property of Licensor. Violation of any provision of this paragraph shall be the basis for immediate termination of this Agreement. No copying or transmission of the Software is permitted.

LIMITED WARRANTY

For one year after the date of the invoice from Seller, Seller will either, at its sole discretion, repair or replace any product covered by such invoice purchased by you from Seller or Seller's authorized distributor which, in Seller's judgment, is defective in materials or workmanship, or refund the amount paid by you for such product. If the product is dated as to shelf life, this limited warranty shall in no event extend beyond such shelf life. In the case of third-party computer operating system software, this limited warranty shall extend only for thirty (30) days after the date of invoice and shall entitle you only to an exchange, with no right to a refund. This limited warranty does not cover: (a) computer applications software (which is governed by the Biolog, Inc. Software License Agreement); (b) consumable or disposable items not dated by Seller as to shelf life, including, without limitation, lamp assemblies, bulbs, fuses, and batteries; or (c) damage caused by batteries, probes or electrodes. This limited warranty does not apply to any product that has been misused, neglected, modified, or repaired by anyone other than an authorized service facility, or to any product requiring refrigeration or other special handling that is not properly refrigerated or so handled. All warranty registration cards included with the products must be completed and returned by you within ten (10) days after the date of purchase. Any product which is or becomes defective during said one year period shall be returned by you, transportation, custom duties and other fees prepaid, in the original carton to Seller's principal facility in Hayward, California, or such other place which may be designated by Seller. You must obtain prior approval from Seller before returning any product to Seller. Transportation charges, custom duties and other fees for return to you shall be paid by you, and you shall bear all risk of loss or damage during transportation to and from Seller. Seller shall reimburse you transportation charges, custom duties and other charges incurred by you in returning product pursuant to this paragraph only if: (a) you have complied with all of these Terms and Conditions; and (b) the returned product is defective in materials or workmanship and is covered by the limited warranty provided in this paragraph. THERE ARE NO WARRANTIES WHICH EXTEND BEYOND THIS DESCRIPTION ON THE FACE HEREOF, ANY COURSE OF DEALING, CUSTOM OR USAGE OF TRADE OR COURSE OF PERFORMANCE NOTWITHSTANDING. THE FOREGOING IS SPECIFICALLY IN LIEU OF ALL OTHER WARRANTIES EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR USE FOR A PARTICULAR PURPOSE AND ANY WARRANTIES THAT THE PRODUCTS DO NOT INFRINGE THE PATENT OR OTHER INTELLECTUAL PROPERTY RIGHTS OF THIRD PARTIES. NO AGENT OR REPRESENTATIVE OF SELLER HAS ANY AUTHORITY TO BIND SELLER TO ANY AFFIRMATION, REPRESENTATION OR WARRANT CONCERNING PRODUCTS MADE BY SELLER. IN NO EVENT SHALL SELLER BE LIABLE FOR LOSS OF PROFIT, LOSS OF DATA, DAMAGE TO OTHER EQUIPMENT USED IN CONJUNCTION WITH THE PRODUCTS WHETHER OR NOT USED PROPERLY, OR ANY OTHER INCIDENTAL, CONSEQUENTIAL, OR SPECIAL DAMAGES. YOU AGREE THAT SELLER'S LIABILITY FOR DAMAGES, IF ANY, SHALL NOT EXCEED THE PAYMENTS MADE BY YOU FOR ANY DEFECTIVE PRODUCT.

Schedule of Current Pages/Sections

The current MicroStation System, MicroLog Version 4.2 User Guide is composed of the following pages or sections, dated as shown:

| Page/Section | Date |
|------------------------------------|-------------|
| Cover Page | 15 Apr 07 |
| Technical Service | 1 Nov 07 |
| License Agreement | 1 Nov 07 |
| Schedule of Current Pages/Sections | 1 Nov 07 |
| Table of Contents | 15 Apr 07 |
| Section 1 | 15 Apr 07 |
| Section 2 | 15 Apr 07 |
| Section 3 | 1 Nov 07 |
| Section 4 | 15 Apr 07 |
| Section 5 | 15 Apr 07 |
| Section 6 | 15 Apr 07 |
| Section 7 | 15 Apr 07 |
| Section 8 | 15 Apr 07 |
| Section 9 | 15 Apr 07 |
| Section 10 | 15 Apr 07 |
| Section 11 | 15 Apr 07 |
| Section 12 | 15 Apr 07 |
| Appendix Cover Page | 15 Apr 07 |
| Appendix 1 | 15 Apr 07 |
| Appendix 2 | 15 Apr 07 |
| Appendix 3 | 15 Apr 07 |
| Appendix 4 | 15 Apr 07 |
| Appendix 5 | 15 Apr 07 |
| Appendix 6 | 1 Nov 07 |

Table of Contents

| Section Title | Sec. Page |
|--|----------------|
| 1. Welcome..... | 1.1 |
| System Requirements..... | 1.1 |
| MicroStation/MicroLog Software Configurations..... | 1.2 |
| MicroStation/MicroLog User Guide Access..... | 1.3 |
| 21 CFR Part 11 Compliance..... | 1.3 |
| 2. Installation and Registration | 2.1 |
| Installing the MicroStation/MicroLog Software..... | 2.1 |
| First Log-In and Setting up the Administrator..... | 2.2 |
| Registering your Software..... | 2.4 |
| Installing Biolog Databases..... | 2.5 |
| 3. MicroStation System Overview..... | 3.1 |
| How It Works..... | 3.1 |
| The Identification Process..... | 3.2 |
| Easy-to-Use Software..... | 3.3 |
| Logging In and Out..... | 3.4 |
| The Math Behind the Software..... | 3.5 |
| Identifying Microbes and Managing Data..... | 3.6 |
| 4. Preparing Samples..... | 4.1 |
| Isolating a Pure Culture..... | 4.1 |
| Gram Staining..... | 4.2 |
| Wet Prep..... | 4.3 |
| Characterizing Aerobic Bacteria..... | 4.3 |
| Characterizing Anaerobic Bacteria..... | 4.5 |
| Characterizing Filamentous Fungi & Yeast..... | 4.5 |
| Culturing Your Microbe..... | 4.6 |
| Preparing Inocula..... | 4.7 |
| Special Procedures for Spore-Forming Gram-Positive Rods..... | 4.9 |
| Special Procedures for Filamentous Fungi..... | 4.13 |
| Inoculating MicroPlates..... | 4.14 |
| Incubating MicroPlates..... | 4.15 |
| Special Procedures for Incubating Anaerobic MicroPlates..... | 4.16 |
| Sample Preparation Process for Filamentous Fungi..... | 4.17 |
| Sample Preparation Process for GN,GP,AN, YT | 4.17 |

| | |
|---|----------------|
| 5. Reading MicroPlates..... | 5.1 |
| Logging In..... | 5.1 |
| Choosing Manual, Reader, or File Mode..... | 5.1 |
| Choosing Printer Method..... | 5.3 |
| Reading Plates Manually..... | 5.3 |
| Reading Plates Using the Plate Reader..... | 5.7 |
| Reading from a Saved File..... | 5.10 |
| Setting Up a Worksheet..... | 5.11 |
| Setting Up to Save Files..... | 5.14 |
| Using Worksheets to Read Multiple MicroPlates..... | 5.16 |
| Adjusting Thresholds Manually..... | 5.19 |
| Choosing a Database to Search..... | 5.20 |
| Exiting..... | 5.22 |
| 6. Interpreting Results..... | 6.1 |
| Reading Printouts..... | 6.1 |
| Assessing the Accuracy of Your ID..... | 6.1 |
| Pinning Down an Uncertain ID..... | 6.2 |
| Species Comparison..... | 6.5 |
| Searching and Assessing Species Information..... | 6.6 |
| Using the View Database Button | 6.6 |
| Using Macroscopic and Microscopic Pictures (FF)..... | 6.12 |
| Using the “Other” Row..... | 6.13 |
| Assessing Raw OD Values..... | 6.14 |
| Assessing Data Values..... | 6.15 |
| 7. Data Management and Compiling Databases | 7.1 |
| Backup Data File..... | 7.2 |
| Combining Data Files..... | 7.3 |
| Viewing and Editing Files..... | 7.4 |
| Compiling Data Files..... | 7.7 |
| Exporting ASCII Data..... | 7.16 |
| Exporting to Microsoft Access..... | 7.21 |
| Cluster Analysis..... | 7.23 |
| 8. Administration and Security..... | 8.1 |
| What is Restricted Access Mode?..... | 8.1 |
| Administrator Options | 8.2 |
| Creating a User List | 8.4 |
| Viewing the Log-In Log | 8.5 |
| Interpreting the Log-In Log..... | 8.6 |
| Changing or Lost Passwords | 8.6 |
| Unauthorized Log-In Attempts..... | 8.8 |
| 9. Technical Notes..... | 9.1 |
| Materials List | 9.1 |
| Media Preparation | 9.2 |
| Specialized MicroPlates | 9.6 |

| | |
|--|-------------|
| 10. Frequently Asked Questions..... | 10.1 |
| 11. Troubleshooting..... | 11.1 |
| Gram Stain Identification..... | 11.1 |
| Culturing Microorganisms..... | 11.3 |
| Preparing Inocula..... | 11.4 |
| Inoculating MicroPlates..... | 11.5 |
| Incubating MicroPlates..... | 11.5 |
| MicroStation Reader..... | 11.6 |
| 12. Glossary..... | 12.1 |
| 13 Appendices..... | 13.1 |
| Appendix 1: Data Statistics for Various Plate Types..... | 13.1 |
| Appendix 2: Release 4.2 Setup Flowchart..... | 13.2 |
| Appendix 3: FF Setup Flowchart..... | 13.3 |
| Appendix 4: Database Species List and Their Characteristics..... | 13.5 |
| Appendix 5: Program Printouts..... | 13.40 |
| Appendix 6: Dangerous Pathogen Database..... | 13.52 |

1. Welcome

In this section:

- ➔ System Requirements
- ➔ Installing MicroStation/MicroLog Software and Biolog Databases
- ➔ MicroStation/MicroLog Software Configurations

Welcome to Biolog's MicroStation™ System. The MicroStation System includes everything you need to use all the features of MicroStation software. If you ordered a MicroStation System with a computer, it should include the following:

- Multimedia-compatible computer with color monitor
- Windows® operating system installed
- MicroStation Reader
- Turbidimeter
- MicroStation/MicroLog software on CD ROM
- Multichannel repeating pipettor

Note: MicroLog 1 and MicroLog 2 systems do not include the MicroStation Reader or pipettor.

System Requirements

If you are installing MicroStation/MicroLog software on a computer system you already have, you will need the following minimum requirements:

- Pentium® PC with 512 Cache, 700 mHz or comparable
- 256 MB RAM
- 80 GB Hard Disk Space
- CD ROM drive
- Microsoft two-button mouse
- Monitor: 800x600 pixel resolution; 6 bit, 256 color
- Windows XP®, Service Pack 2

Note: The software is not designed to run on a Macintosh®.

NOTE:

Throughout this guide, references are made to other sections in the following manner:

- Section 5 page 3, referring to the third page of the fifth section.

Or,

- Page 3.5 referring to the fifth page of the third section.

Connecting the MicroStation System

1. Unpack all materials. Your MicroStation arrives in many boxes. Check that all items have arrived against the packing slip.
2. Connect the computer, monitor, keyboard, and mouse.
3. Connect the reader, per enclosed instructions.
4. Set up the turbidimeter, per enclosed instructions.
5. Plug in the pipetter to fully charge the battery.

MicroStation/MicroLog Software Configurations

In order to make the power of the MicroStation System available to a variety of microbiology laboratories with different needs and budgets, Biolog sells three configurations of MicroStation/MicroLog software.

This user guide contains instructions and information for all basic and advanced features of MicroStation/MicroLog software. Depending on the software configuration you purchased, some sections of the manual may not apply to your lab.

Table 1.1 shows the differences between configurations.

| Software Configuration | Specifications |
|-----------------------------|---|
| MicroLog 1 | <ul style="list-style-type: none"> ➔ Identification only ➔ Visual reading only (no MicroStation Reader) |
| MicroLog 2 | <ul style="list-style-type: none"> ➔ All advanced features, Restricted or Unrestricted Mode ➔ Visual and file reading only (no MicroStation Reader) |
| MicroStation/ MicroLog 3 | <ul style="list-style-type: none"> ➔ All advanced features, Restricted or Unrestricted Mode ➔ MicroStation Reader (visual and file reading also possible) |

TABLE 1.1. MICROSTATION/MICROLOG SOFTWARE CONFIGURATIONS

If you bought MicroLog1

- ➔ Skip material dealing with the MicroStation Reader, Worksheets, Adjusting Thresholds, and Data Management, Compiling Databases and Restricted Access.

If you bought MicroLog2

- ➔ Skip material dealing with MicroStation Reader.

If you bought MicroStation/ MicroLog 3

- ➔ Read the entire User Guide.

MicroStation/MicroLog User Guide Access

In order to access the MicroStation/ MicroLog 4.2 User Guide, Acrobat Reader must be installed in your computer. If you do not have Acrobat Reader already installed in your computer, then a set-up program called “Acrobat Reader 7.0.exe exists in the Manual sub-directory under ...**Program Files\Biolog\MLN_XX_XX (e.g. ML3 version 4.20.05 is 42_05)**. Run the set-up application and follow its instructions.

21 CFR Part 11 Compliance

The MicroStation/MicroLog system is designed to exist within a 21 CFR Part 11 environment. We provide the basic functionality to support our customer’s compliance efforts.

The objective of the original regulation is to ensure integrity of electronic records by:

1. Limiting Access
2. Ensuring that original record files can not be modified.
3. Providing documentation of what was changed, by whom and when.

Within “Restricted Access Mode” the software will enforce all of the security Microbial identification tasks accomplished in accordance with the guidelines of 21 CFR Part 11, ensuring data integrity and security controls. Database management (audit control) will need user procedural controls to manually freeze files after edits are made

The software is both password protected and encrypted.

The MicroStation/MicroLog Software combined with the RetroSpect™ Trending and Tracking Software meets the needs of facilities that must adhere to strict regulatory guidelines.


2. Installation and Registration

In this section:

- System Requirements
- Installing MicroStation/MicroLog Software and Biolog Databases
- MicroStation/MicroLog Software Configurations

Installing the MicroStation/MicroLog Software

1. Insert the software CD into the CD ROM drive.
2. The software setup Installer will initialize. This will take a few moments to prepare.
3. The software – **InstallShield® Wizard** will appear. Click **Next**.



4. The **License Agreement** will appear. Review the terms and select the **I accept the terms in the license agreement** radio button. Click **Next**.
5. The next screen that appears gives you the option of where to install the software root directory folders:

Destination Folder - The root directory (for the software) may be installed anywhere on a stand alone PC or on a network drive. If you do install the directory on a network drive, please remember that this program is currently not designed for multiple users at one time.

- The default/suggested installation location is ...\\Program Files\\Biolog\\MLN_XX_XX (e.g. **ML3 version 4.20.05 is 42_05**). Click **Next** to accept this default location.
- To choose your own location, click **Change**. Select the **Look In** pull down menu to browse for a location. Select the location and click **OK**.


Important!

- The person installing the software must be logged in through the Windows operating system as the **Administrator**.

6. Now the **InstallShield Wizard** is ready to begin installation. Choose **Back** to make any changes, **Cancel** to exit the InstallShield Wizard and not install the program, or **Install** to proceed with installation.
7. The next screen that appears asks you to enter the serial number of your program (found on the spine of the software jewel case). Enter the **Serial Number** and click **OK**.

Note:

The Readme Information is not printable at the InstallShield Wizard window. The file is located in the MicroStation\MicroLog Software root directory and is named *readme.rtf*. Go to this location to open and print the Readme information.



8. The **Installing Software** screen will appear while the installation is taking place. Click **Finish** when the InstallShield Wizard Completed screen appears.

The following shortcut icons will now be installed on your desktop:

- **MLN_XX_XX.exe** is for the program
- **MicroLog ID User Guide**

Please remember that the program will be accessible for up to 30 days. After 30 days the software must be registered to use!

First Log-In and Setting up an Administrator

The program is ready to be opened for the first time. The program operates in the Restricted Access Mode.


Restricted Access Mode requires that only the **Program Administrator** oversees and controls who has access to the program. The person who is designated as the Original Program Administrator should perform the first Log-in. This individual then manages who has access to the application and what tasks they can perform. If desired, more than one user can be assigned Administrator Privileges on the User List.

The Administrator will:

- Assign User names and passwords for those who will use the system.
- Assign access privileges to each user.

The following steps must be implemented by the user who will act as the program Administrator.

1. Click on the '**MicroStation/MicroLog**' shortcut icon on your desktop.
2. The **Welcome** screen will appear. Click the **Administrator Login** key in the lower right hand corner.



3. The **Administrator Dialog** box will open.
 - Enter a Username that is at least 1 character in length.
 - Enter a Password that is at least 6 characters in length and contains at least one number. *The password is case sensitive.*
4. Click **OK**.
5. When the message "Click on the Log In button at the top right to log in to the program" appears, click **OK**.
6. Click on the **Log In** box located in the upper right hand corner of the **Welcome** window. A **Password Dialog** box will appear.

Note:

Not all users will have full administrative privileges. To learn more about assigning user names and passwords, as well as setting up different levels of user access, please refer to *Section 9*.



7. Enter the Administrator username and password you set in **Step 3**. Click **OK**.
8. You will now be logged in with full user privileges.

Registering your Software

After initial installation, the **Welcome** tab will show “**Temporary Registration Days Left: 3**” in red. The software will count down how many days you have left to register. You must click the **Registration button** to start the registration process.

Registration Process

1. Generate a **User Key** and send to Biolog
2. Load the **Registration Key** from Biolog

Follow the steps outlined below to generate a User Key:

1. At the **Welcome** tab, click on the **Registration** button.
2. The **Registration Form** window appears.

Note:

There is only 1 registration button access per session. Log-Out and Log-In for additional access.

• Each computer requires a separate registration key.


• Your registration key should arrive within 48 business hours, providing it is received Monday to Thursday during regular Biolog business hours (M-F 8:30 A.M. – 5:00 P.M. PST).

3. Fill out every line of the **Registration Form**. (* Required Fields)
4. Click the **Save User Key** button.
5. A **Save As** window appears. Type a file name for the User Key File in the **File Name** field. The **Save as Type** field should show Text Files as the default file type.
6. Click **Save**. The User Key will be saved as a text file.
7. Click **Close**.

8. E-mail the User Key File to tech@biolog.com.

Follow the steps outlined below to load the Registration Key:

1. Biolog will e-mail you a Registration Key. Save the attached registration key file on your hard drive (...\\MLN_XX_XX\\Certificates).
2. Once you have saved your Registration key, open the MicroStation/MicroLog Software program, and click on the **Registration** button to open the **Registration Form**.
3. Click on the **Registration Key** tab that is on the registration form. Then click the **Load Registration Key** button to load the registration key file and process the registration.



4. The Temporary Registration box on the Welcome screen should now read **“Registered”**, and turn from red to gray.
5. The Registration Button is no longer present.

Installing Biolog Databases

Installing a MicroLog database

Install databases in any order. There are six available databases each on a separate CD.

- Gram Negative (GN)
- Gram Positive (GP)
- Anaerobe (AN)
- Yeast (YT)
- Filamentous Fungi (FF)
- Dangerous Pathogens (DP)

Note: Yeast or Filamentous Fungi databases requires use of a MicroStation Reader.

1. Check to make sure you have the MicroLog database CD ROM.
2. Put the database CD ROM into the CD ROM drive.
3. Click **Start**.
4. Select **Run**. Click on Browse. Click on the CD ROM drive.
5. The field should read **D:\SetUp(database).exe**. If drive D is not your CD ROM drive, then type in the correct drive letter.
6. Click **OK**.
7. Follow the instructions on the pop-up dialogs. Make note of the drive and directory path into which the database is installed.
8. To install the database into the default directory (...**Program Files\Biolog\MLN_XX_XX**), click **Finish**. If you chose a different directory for the ML software, then type in the name of that same directory, before clicking **Finish**.

Running MicroLog software

Double-click the MicroLog desktop icon.

3. MicroStation System/MicroLog Overview

In this section:

- ➔How It Works
- ➔The Identification Process
- ➔Easy-to-Use Software
- ➔The Math Behind the Software
- ➔Identifying Microbes and Managing Data

The MicroStation System/MicroLog is an easy-to-use yet advanced tool for identifying and characterizing microorganisms. Our combined databases include over 1,900 species of aerobic bacteria, anaerobic bacteria, fungi and yeasts. They contain nearly all of the significant species encountered in diverse practices of microbiology, including pharmaceutical, biotechnology, cosmetic, and medical device companies; veterinary medicine; agriculture and environmental science; food processing, spoilage, and safety; reference laboratories; industrial microbiology; and research and education.

MicroStation/MicroLog continues to expand and evolve to keep pace with progress in the field of microbiology. Every month researchers discover new species of microorganisms and recognize their importance. Biolog's patented microbial identification technology with 95 carbon source utilization tests in a microtiter plate format (MicroPlate™) can recognize over 4×10^{28} possible metabolic patterns. This provides room for future growth of the system, so that the technology will remain state-of-the-art.

How It Works

Biolog's innovative, patented technology uses each microbe's ability to use particular carbon sources to produce a unique pattern or "fingerprint" for that microbe. As a microorganism begins to use the carbon sources in certain wells of the MicroPlate, it respire. For bacteria, this respiration process reduces a tetrazolium redox dye and those wells change color to purple. The end result is a pattern of colored wells on the MicroPlate that is characteristic of that microorganism. For fungi, respiration and assimilation are detected. The color change in the wells is reddish orange due to the respiration of fungi, which reduces the dye. Assimilation or growth is detected by the turbidity of the well.

A bacterial pattern is readable either visually or by a fiber optic reading instrument – the MicroStation Reader. This reader is required to read a yeast or fungal pattern. The fingerprint data is fed into the software, which searches its extensive database and makes a identification in seconds. By developing a simple tool that allows 95 simultaneous carbon source utilization tests, Biolog has accomplished its goal of producing an efficient, easy-to-use, powerful, and reliable microbe identification system.

The Identification Process

Microbial identification involves five basic steps. These steps apply to all identifications. A small number of species have peculiarities that may require an extra step or special handling techniques.

The Microbial Identification Process

Step 1

Isolate a pure culture on Biolog media



Step 2

Do a Gram stain (or wet prep for FF)
and determine testing protocol



Step 3

Prepare inoculum at specified cell density



Step 4

Inoculate and incubate MicroPlate



Step 5


Read MicroPlate and determine ID

Step 1: Isolate a pure culture on Biolog media

Isolating a pure culture is not always easy. For example: Bacteria often have sticky surfaces and cells sometimes stick together in clumps. As a first step to accurate microbe identification, streak agar plates using correct techniques to generate well isolated colonies. Always use Biolog-recommended culture media and growth conditions. See Section 4 for full culturing and incubation instructions.

Step 2: Do a Gram stain and determine testing protocol

For bacteria, proper Gram stain technique and interpretation are the important second step in the ID process. See pages 4.2 and 10.1-10.2. For FF identification, use the wet prep test as necessary to differentiate yeasts from filamentous fungi.



*Follow
directions
closely to
obtain accurate
results.*

Step 3: Prepare inoculum at specified cell density

Microbiologists are sometimes trained to prepare cell suspensions by judging cell density by eye. This practice will not yield accurate and reproducible results. Cell density determines oxygen concentration – a key parameter to control when testing microorganisms in MicroPlates. In addition, Biolog has carefully optimized the required inoculating fluids. Don't deviate from Biolog's inocula preparation directions See pages 4.7 - 4.10.

Step 4: Inoculate and incubate MicroPlate

Pipet the specified amount of cell suspension into the MicroPlate, put the lid on, and incubate under the same conditions of temperature and atmosphere used to culture the microorganism. Biolog MicroPlates do not need oil overlays or color-developing chemicals.

Step 5: Read MicroPlate and determine ID

After an appropriate incubation time, read MicroPlates either by eye or using the MicroStation Reader. In either case, the pattern formed in the wells is entered into the software, which searches the database and provides an identification in seconds.

Easy-to-Use Software

MicroStation/MicroLog provides an easy-to-follow visual software interface to lead you through the identification process. The software can

be run in a restricted access mode that provides lists of user names, passwords, and a hierarchy of log-in privileges. This ensures the integrity of electronic records by limiting access, providing documentation of changes made, creating audit trails, allowing you to freeze data, and ensuring that original file records have not been modified or deleted.

You'll spend most of your time using the Data window, with its areas for entering sample identifiers (strain name, strain number, etc.), viewing and entering MicroPlate reactions (positive, negative, and borderline), and viewing microbial identifications.

Once you're familiar with the system, you'll be able to prepare your samples using proper techniques, read MicroPlates either singly or in batches, prepare worksheets to save time, build your own database, and analyze data using the advanced functions built into the software.

The screenshot shows the Biolog MicroLog3 4.20 Data window. It is divided into several sections:

- Sample identifiers:** A section on the left with fields for Current Time (May 04 2001 21:17), Plate Number (2), Sample Number, Plate Type (FF), Strain Type (Penicillium), Strain Name (PEN.CITE BGC), Strain Number (CBS 454.93), Incubation Time (48 HR), and Other.
- MicroPlate reactions:** A grid on the right showing reactions for 12 wells (1-12) across 8 rows (A-H). The reactions are indicated by symbols like circles and dots.
- List of ranked choices:** A table at the bottom showing species identification results. The table has columns: ML4, NAME, PROB, SIM, DIST, and TYPE. The first row is highlighted, showing 'Penicillium citreonigrum Dierckx BGC' with a probability of 99%.
- ID box with first-ranked choice:** A box on the right pointing to the first row of the table, indicating the first-ranked choice.

Buttons at the top include Welcome, Set Up, Data, Data File, and Exit. Buttons at the bottom include Select Read, Read Next, Print, Save, View Database, and Cluster Menu.

DATA WINDOW, WITH RESULTS

Logging In and Out

First, follow these steps to Log In and Out of the program, as well as switch users. Remember, every person who wants to use the program must always enter a User ID and password to access the program. This provides important security to the program by putting controls on and keeping a record of who accesses the program.

Logging In

1. Click on the 'MicroStation/MicroLog' shortcut icon on your desktop
2. The **Welcome** screen will appear, along with a **Password Dialog** box.

3. Enter your Username and Password. Click **OK**.
4. You will now be Logged In, and the Welcome interface will appear.

Timed Log Out

The software program has a default setting to automatically **Log Out** if the program has not been used for 15 minutes. The **Login Time Out** screen will appear. Simply select the **Click Here** button to stop the Log Out process and continue using the program. Otherwise the program will log off, return to the **Welcome** window, and display a new **Log In** dialog box.

Logging Out and Switching Users

There are 3 ways to end your use of the program:

1. Select the **Exit Tab** and then click the **Exit** button.
2. Click on the **Change Users** button in the center of the **Welcome** window. The **Password Dialog** box will be displayed for the new user to enter their username and password.
3. Click on the windows page “**X**” button.

The Math Behind the Software

MicroStation/MicroLog software uses extensive computer algorithms to take the information from the observed pattern and compare it to the database.


In simple terms, the software rapidly compares the positive-negative-borderline color pattern in the MicroPlate to the species pattern in the appropriate database. The patterns that most closely match your microbe’s pattern are shown on the screen in ranked order. Before making a decision on the result, the software considers the possibility that even the first-ranked choice may not be a good match. It looks to see whether the first choice match is really “close enough” to be acceptable. If not, a “No ID” designation will result.

For GN, GP, and AN databases, MicroStation/MicroLog uses a newly developed pattern matching method called Progressive ID (PID). This method more accurately identifies species patterns by considering the progressive sequence in which purple wells are formed. Typically, microorganisms will use their favorite carbon sources most rapidly and completely, resulting in dark purple wells that form quickly. Less-preferred carbon sources are consumed slowly or incompletely, resulting in slower-forming or lighter purple color. The extra information considered by the PID matching method brings a higher level of consistency and accuracy, representing another innovation in Biolog’s technology.

For YT and FF databases, MicroStation/MicroLog uses the developed pattern matching method called Endpoint ID (EID). This method more accurately identifies species patterns by considering the daily endpoint determinations. For the YT MicroPlate, turbidity may be caused by either growth or color formation. For the FF MicroPlate, both growth (measured by turbidity) and color formation are measured independently for each well.

Identifying Microbes and Managing Data


The software is structured to move you easily through identifying microbes and managing data. Navigation tools include a menu bar at the top of the window, menu tabs to move from window to window, drop-down lists to choose from pre-set choices, selection bars to reach detailed setup windows, and fields to type in specific data.



The software will guide you in answering the following questions:

- Do I plan to read MicroPlates visually?
- Will I use the reader or will I read back from a file?
- Do I want an automatic printout of results?
- Do I want to select a specific printer?

If you plan to enter MicroPlate reactions manually you'll enter sample identifiers using the mouse, keyboard, or on-screen methods to indicate positive, negative, or borderline MicroPlate reactions.



ENTERING SAMPLE IDENTIFIERS AND MICROPLATE REACTIONS

In addition, the software allows you to:

- Use worksheets
- Save data
- Edit files and build databases
- Adjust thresholds
- View a database
- Perform cluster analysis

Using worksheets

Worksheets allow you to enter all of your sample identifiers ahead of time. This speeds the reading process and enables you to read high MicroPlate volumes more efficiently.

| | Plate # | Sample # | Plate | Strain Type | Strain Name | Strain # | Other |
|------|---------|----------|-------|-------------|-------------|----------|-------|
| => 1 | - | - | - | - | - | - | - |
| 2 | - | - | - | - | - | - | - |
| 3 | - | - | - | - | - | - | - |
| 4 | - | - | - | - | - | - | - |
| 5 | - | - | - | - | - | - | - |
| 6 | - | - | - | - | - | - | - |
| 7 | - | - | - | - | - | - | - |
| 8 | - | - | - | - | - | - | - |
| 9 | - | - | - | - | - | - | - |
| 10 | - | - | - | - | - | - | - |

WORKSHEET (WITH SAMPLE ENTRIES) FOR BATCHING MICROPLATES

Saving data

You can save your results into files and manage them by using the Output options on the **Set Up** window.

Input / Output Database

Input

Data Input Mode: Manual

Use Worksheet: No

Output

Automatic Printout: No

Print As Text: No

Printer Setup

Write Data To Data File: Yes

Output File Name:


Auto Write to Data File: No

File output area

FIELDS FOR NAMING AND SAVING FILES

Editing files and building databases


You can also select which database you want to use – Biolog's or your own. Advanced Data Management functions allow you to view data by strain or plate type, edit data files, build your database, and perform a number of other database management techniques.



DATABASE MANAGEMENT WINDOW (FOR ADVANCED DATABASE FUNCTIONS)

Adjusting thresholds

Another software advanced feature for bacteria is the capability of adjusting thresholds by using histograms. This function may be appropriate if you feel that the automatic reading of your MicroPlate does not match your visual interpretation.




HISTOGRAM AND MICROPLATE DATA (FOR MANUAL ADJUSTMENT OF THRESHOLDS)

Viewing databases

The **View Database** function allows you to access end-point and/or progressive data of specific organisms in each of Biolog's databases.


The function also allows you to view the Macroscopic and Microscopic image libraries included in the Filamentous Fungi (FF) Database.





EXAMPLE OF MICRO PICTURES FROM THE FF DATABASE IMAGE LIBRARY

Performing cluster analysis

Cluster analysis is an advanced feature that offers visual ways to simplify interpretation of data sets. You'll be able to view dendrograms, two-dimensional, and three-dimensional cluster diagrams.





4. Preparing Samples

In this section:

- Isolating a Pure Culture
- Gram Staining
 - Wet Prep
- Characterizing Aerobic Bacteria
- Characterizing Anaerobic Bacteria
- Characterizing Filamentous Fungi and Yeast
- Culturing your Microbe
- Preparing Inocula
- Special Procedures for Spore-Forming Gram-Positive Rods
- Special Procedures for Filamentous Fungi
 - Inoculating MicroPlates
 - Incubating MicroPlates
- Special Procedures for Incubating Anaerobic MicroPlates
- Sample Preparation Process

As with any system, the precision and accuracy of MicroStation/MicroLog results require proper sample preparation. The most common identification problems result from improper lab technique and using non-recommended media. Using good sterile technique and the correct media greatly increase the likelihood of problem-free microbe identification.

If you're thinking of using non-recommended media – *don't*. The extensive MicroLog databases were specifically developed using Biolog's carefully selected media. Using non-specified media causes changes in the physiology of microbes. Luxuriant growth and full metabolic activity are required for the system to work well.

The flowchart on page 2.2 provides an outline of the identification process. The tables on page 4.16-17 give a comprehensive overview of the sample preparation process, as does the flowchart in the Appendices. See pages 9.2-9.6 for detailed media preparation instructions. Refer to Appendix 6 when using the DP (Dangerous Pathogen) database.

For additional information refer to Instructions for Use specific to the plate type (supplied with each purchase of MicroPlates).

Isolating a Pure Culture

Culture your sample on Biolog's general-purpose culture medium. Nearly all the bacterial species in the MicroLog databases will grow on Biolog Universal Growth media (BUG or BUA), generally with the addition of either 5% sheep blood or maltose. Prepare media according to package insert. Yeast are cultured on BUY (YT), while filamentous fungi (FF) are cultured on 2% Malt Extract Agar. .

Each shipment of Biolog MicroPlates comes with instructions. Refer to these for detailed directions about how to grow a healthy culture. In general:

- Make certain you have a pure culture (a single strain).
- Incubate at optimal temperature. Most bacteria should be grown at 35-37° C or 30° C depending on genera. (A few gram positive, anaerobic organisms, yeasts and filamentous fungi are grown at 26° C.)

Biolog special agar:
BUG = Biolog Universal Growth
BUA = Biolog Universal Anaerobe
BUY = Biolog Universal Yeast
2% ME = 2% Malt Extract Agar

Preparing an inoculum directly from a mixed-growth plate will cause identification problems. Use the colony magnifier lamp to closely examine colonies.

- Do not allow cultures to grow for too long. Maximum growth is 24 hours for most bacteria, 72 hours for yeasts, and 10 days for most fungi. Some exceptionally slow-growing or fastidious bacteria may require 48 hours of growth or the use of multiple growth plates.

Checking to see if your culture is pure

Colonies on the plate that seem to be isolated, may in fact be the result of mixed growth. This is especially true with *Staphylococcus* species. Careful visual examination is essential to ensuring that a culture is pure. It may be useful to use a colony magnifier lamp to assist in the examination. If a colony shows any hint of pleomorphism, it is probably not a pure culture and requires additional re-streaking and isolation.

Carefully examine areas of confluent growth. If the lawn is not uniform in texture and color, this may indicate that the culture is not pure. Once again, restreak for isolation.

Note: Unless your organism is a very slow grower, we don't recommend using lawn growth.

The opposite problem can also occur; sometimes a culture may be pure, but give the appearance of heterogeneity. This is due to a rather common phenomenon whereby microorganisms produce more than one colony type. To be certain of its identity, purify and test each colony type individually. Be aware that different colony types of the same strain may ID but give different phenotype profiles.

Gram Staining

A Gram stain is essential in order to select the proper media and MicroPlates for testing bacteria. Once your initial culture has been incubated sufficiently, perform a standard Gram stain. Determine the following:

- Is the microbe a bacterium or a yeast?
- Is the microbe gram positive or negative?
- Are the cells cocci or rods?
- Do the cells form spores?

Based on the source of the sample, the initial growth conditions, and the Gram stain, you should know what basic type of microbe you have. See pages 10.1-10.2 for specific Gram stain information. Refer to a basic microbiology textbook if you are unsure of these procedures.

Wet Prep

A wet prep is performed to determine if the microbe is a yeast or filamentous fungi. This is an optional test, required only if the form is uncertain. Refer to a basic microbiology textbook if you are unsure of this procedure.

Characterizing Aerobic Bacteria

Gram-negative and gram-positive aerobes require further investigation to determine proper handling and identification.

Attention!

Please refer to Appendix 6 when using the DP (Dangerous Pathogen) Database.

Characterizing gram-negative microbes

If your microbe is gram negative, several additional tests will help determine whether your microbe is non-enteric (GN-NENT), enteric (GN-ENT) or fastidious (GN-FAS) as well as determine the proper setup protocol. This has important implications in regards to choice of incubation temperature, inoculum density, and inoculating fluid.

Performing an oxidase test

- Do an oxidase test on all gram-negative organisms. Refer to a basic microbiology textbook if you are unsure of this procedure.
- If the oxidase test is positive, you have a non-enteric (GN-NENT) microbe.
- If the oxidase test is negative you most likely have an enteric (GN-ENT) microbe.
- If you are uncertain about whether the microbe is enteric or non-enteric, set up a TSI slant.

Setting up a TSI slant

- Prepare a TSI slant on oxidase-negative microbes. Refer to a basic microbiology textbook if you are unsure of this procedure.
- If the TSI slant shows a K/K (alkaline/alkaline) or K/A^w (alkaline/weak acid) reaction, your microbe is non-enteric (GN-NENT).
- If the TSI slant shows an A/A (acid/acid) or K/A (alkaline/acid) reaction, your microbe is enteric (GN-ENT). This is true with a few exceptions (such as *Pasteurella species*=A/A, or *Acinetobacter A/late A*), which can be weak oxidase positive.

Enteric: gram-negative bacteria belonging to Enterobacteriaceae group

Non-enteric: gram-negative bacteria not belonging to Enterobacteriaceae group

Recognizing fastidious gram negatives

- Fastidious gram-negative bacteria are primarily isolated from the respiratory tracts of mammals.
- These bacteria grow poorly or not at all on BUG + B medium. Some *Moraxella species* do grow well on BUG + B, but should be set up as fastidious gram negatives. They generally stain as gram-negative diplococci.
- They grow much better on chocolate agar at 35-37° C in an atmosphere of 6.5% CO₂.

Characterizing gram-positive microbes

If your microbe is gram positive, pay special attention to the cell morphology in the Gram stain. Note the following:

- Determine whether the microbe is a coccus or a rod.
- If you are unsure if your organism is a coccus or a rod (e.g., coccobacillus), set up a MicroPlate and select the **GP-COCCUS-ROD** strain type option to determine the correct identification.
- If the microbe is a rod (GP-ROD), determine whether it is spore-forming bacillus (GP-ROD SB, consisting of *Bacillus* and species formerly called *Bacillus*). Spore-formers require special handling and testing. See pages 4.10- 4.12.
- It can be difficult to recognize *Bacillus* species from a Gram stain, but most can be recognized by colony morphology. They are often fast growers, forming colonies that have unusual characteristic textures (e.g., slimy, crusty, dry, embedded, or forming skin-like pellicles). Prepare a wet mount to look for spores, retractile ovals or spheres.

Performing a catalase test

It is advisable to perform a catalase test (3% hydrogen peroxide) on GP-COCCUS and GP-ROD isolates. The catalase reaction can help you verify or narrow down the species identification. Refer to a basic microbiology textbook if you are unsure of this procedure.

If the organism is catalase-negative, this generally indicates a slow grower (CO₂ may be required). Use a swab or loop to subculture the organism onto the agar to obtain a lawn of growth. You may need more than one agar plate.

Do not perform a catalase test from a blood plate, it may give false positive results.

Characterizing Anaerobic Bacteria

Gram-negative anaerobic rods require further investigation to determine proper handling and identification.

Characterizing rapidly-growing, gram-negative microbes

If your microbe is rod-shaped, anaerobic, gram negative, and grows rapidly, a kanamycin disk test will help differentiate *Fusobacterium* from *Bacteroides/Prevotella*. This has important implications in regards to choosing the appropriate inoculation protocol.

Performing a kanamycin disk test

- Streak rapidly-growing, gram-negative rods on BUA plates.
- Add a 1 mg kanamycin disk to the first quadrant (immediately after streaking, before incubating the organism).
- Incubate the plate overnight.
- If the bacterium is sensitive (a 10 mm zone of clearing will surround the kanamycin disk), you have a presumptive *Fusobacterium*.
- If the bacterium is resistant (no or <10 mm zone of clearing around the kanamycin disk), you have a presumptive *Bacteroides* or *Prevotella*.
- This distinction is important when inoculating AN MicroPlates. See pages 4.15-4.16.

Characterizing Filamentous Fungi and Yeast

The Filamentous Fungi database is unlike any of Biolog's other databases. It consists of eight separate databases that are accessed by selecting the **Strain Type**. It is important to correctly determine the database you wish to search by entering the appropriate **Strain Type** on the worksheet or in the **Data** window for single FF MicroPlate reads.

- The Air database includes the most commonly isolated filamentous fungi and yeasts from the air.
- The Food database includes the most commonly isolated filamentous fungi and yeast from food.
- The Yeast database contains a limited number of yeasts.*
- The Aspergillus, Colletotrichum, Fusarium, Penicillium, and Trichoderma databases contain anamorphs and teleomorphs related to the specific genera.*

**Note: These databases are intended for use by researchers working on a specific genera, not for routine food and air isolates.*

Culturing Your Microbe

From a single isolated colony subculture your microbe onto the proper agar medium. Most of the species in Biolog's database are relatively fast growers. However, if your microbe grows slowly, you may need to streak (subculture) more than one agar plate.

Incubate according to instructions in package insert. In general:

- Incubate most organisms for 16-24 hours. Very slow-growing species and yeasts may require 48 hours. Slow-growing anaerobic bacteria may require even longer incubation. Incubate filamentous fungi for 5-10 days (sporulation of isolate). Slower growing filamentous fungi may require longer incubation for sporulation.
- Don't over-incubate your culture. Excess incubation can cause microbes to enter a stationary phase, during which they lose viability and metabolic activity.
- For best results, subculture anaerobic bacteria twice on BUA. Prepare the suspension from the second plate only.
- For filamentous fungi, subculture less hyphae for fast growing fungi and more hyphae for slow growing fungi.

Table 4-1 and Appendix 2 will help you select the correct culture medium.

Media shorthand:

B = blood
M = maltose
T = thioglycolate
ME = malt extract

| Organism Type | Organism Abbreviations | Culture Media |
|--|------------------------|---------------|
| Gram Negative Non-Enteric** | GN-NENT | BUG + B |
| Gram Negative Enteric | GN-ENT | BUG + B |
| Gram Negative Fastidious | GN-FAS | Chocolate |
| Gram Positive Cocci | GP-COCCUS | BUG + B |
| Gram Positive Rods (non-spore forming) | GP-ROD | |
| Gram Positive Rods (spore forming bacillus) | GP-ROD SB | BUG + M + T |
| Anaerobes | AN | BUA + B |
| Yeasts | YT | BUY |
| Filamentous Fungi (and select yeast species) | FF | 2% ME |

TABLE 4-1: SELECTING THE CORRECT CULTURE MEDIUM

**Many agriculture organisms should be grown on BUG (w/o Blood). For listing, see Appendix 4.

Preparing Inocula

STOP!

Is your suspension density accurate? Did you calibrate using Biolog's Turbidity Standards?

Once your microbe is isolated and cultured, prepare a liquid inoculum. Refer to instructions enclosed with Biolog MicroPlates for detailed directions. In general:

- The inoculum **MUST** be within the range specified by the turbidity standards accompanying your database ($\pm 2\%$ T).
- The inoculum **MUST** be uniformly suspended. If a bacterial organism forms clumps, you will need to use special techniques to achieve a homogenous suspension. (see "Dry Tube Method" page 4.11)

Table 4-2 will help you decide which suspension medium to use:

| Organism Type | Inoculating Fluid | Turbidity Standards | Inoculum Density |
|--|-------------------|--------------------------|------------------|
| Gram Negative Non-Enteric | GN/GP-IF | GN-NENT | 52% T |
| Gram Negative Enteric | GN/GP-IF + T | GN-ENT | 61% T |
| Gram Negative Fastidious | GN/GP-IF + T | GP-COC & GP-ROD & GN-FAS | 20% T |
| Gram Positive Cocci | GN/GP-IF + T | GP-COC & GP-ROD & GN-FAS | 20% T |
| Gram Positive Rods (non-spore forming) | | | |
| Gram Positive Rods (spore forming bacillus) | GN/GP-IF | GP-ROD SB | 28% T |
| Anaerobes | AN-IF | AN | 65% T |
| Yeasts | Sterile water | YT | 47% T |
| Filamentous Fungi (and select yeast species) | FF-IF | FF | 75% T |

TABLE 4-2: SELECTING THE CORRECT INOCULATING FLUID, TURBIDITY STANDARDS, AND CELL DENSITY

Adding thioglycolate to inoculating fluid

Some microbes (see Table 4-2) require the addition of thioglycolate (T) to the inoculum. Thioglycolate acts as an anticapsule agent; it decreases production of bacterial capsules so that strains give more consistent patterns. Biolog provides you with droppers of thioglycolate. To add thioglycolate:

*Don't forget
to add thioglycolate
when working with:*

**GN Enteric
GN Fastidious
GP Cocci & Rods**

*Remember,
your microbes
are alive.
Treat them
with care. Use
fresh cultures.*

Caution!
Use special care!

**Do not
excessively
vortex FF-IF
suspensions.
The FF-IF is
viscous.**

**It will trap air
bubbles that will
interfere with
absorbance
readings.**

1. Hold reagent dropper upright and point tip away from you. Using a dissecting hemostat fully crush ampule close to its center one time only. Tap the bottom on benchtop a few times.

Caution: *Do not manipulate dropper any further, as the plastic may puncture, causing injury.*

2. Invert the dropper for convenient drop-by-drop dispensing of reagent.
3. Dispense precisely 3 drops of concentrated thioglycolate into your inoculating fluid. Do not use more than 3 drops per 18-20 ml. This gives you a 5 mM final concentration.
4. Each dropper contains about 15 drops. You can continue to use an open dropper for the rest of the day, and then discard the remainder.

Adding sodium salicylate to inoculating fluid

If gram-positive bacteria are false positive even after you've added thioglycolate to the inoculating fluid, try adding sodium salicylate:

1. Add thioglycolate as described above.
2. Add 1 ml of sterile 100 mM sodium salicylate to 18-20 ml of GN/GP-IF.
3. Proceed with sample preparation.

Preparing inocula

Once you've added thioglycolate (if necessary), prepare the inocula:

1. Establish the appropriate turbidity range on your turbidimeter by adding and subtracting 2% T to the percent transmittance measured with the appropriate turbidity standard.
2. Dip a sterile swab into the inoculating fluid to moisten it.
3. Lift cells from the agar by rolling the swab over the colonies, rather than sliding across them. Be sure not to pick up any agar.
4. Twirl the swab against the inside surface of the tube (above the fluid line) to gently dissolve colonies.

5. Dip the swab into the fluid and stir with a up-down motion to the bottom of the tube to create a uniform suspension. You can also use a sterile transfer pipette to mix without creating an aerosol. Recap the tube and invert if using Biolog's GN/GP-IF.
6. Adjust the inoculum density so it is within the specified turbidity range. You can change the density by adding more cells (to increase density) or more inoculating fluid (to lower density). Be sure to vortex or mix gently before taking the final reading.
7. Examine your suspension and make certain that it is homogeneously suspended and free of clumps. It is essential that the cell suspension does not contain clumps. If there are only a few clumps, allow them to settle to the bottom, pouring off the supernatant. *If the bacterial suspension is not homogenous, use the special dry tube procedure described in Step 4 on page 4.11 and in Section 10.*
8. For anaerobic bacteria, prepare the suspension in batches of no more than six at a time. The time from beginning the first suspension to finishing the last suspension should not exceed 5 minutes.

Special Procedures for Spore-Forming Gram-Positive Rods

This section describes the procedure and details the techniques found to provide consistent identification of spore-forming gram-positive rods (*Bacillus* species). *Bacillus* species require a special agar media and streaking technique to keep cells in a vegetative growth phase and minimize spore formation. Most species are strong producers of capsular polysaccharide. This polysaccharide material interferes with cell testing and causes cells to clump, making it difficult to prepare cell suspensions. Thioglycolate is used in this special procedure to reduce polysaccharide production. Preparing a uniform cell suspension free from clumps is important for accurate identification.


Always use the following special swabbing and streaking technique with *Bacillus* species:

1. Isolate a pure culture of the Spore-Forming Gram-Positive Rod.
2. Subculture the spore-forming gram-positive rod on BUG plus Maltose plus Thioglycolate (BUG + M + T). This will stimulate cells into vegetative growth phase and decrease cell clumping. Prepare the agar media according to the following procedure.
 - a. Add 8 drops of thioglycolate (Biolog P/N 73011) from the dropper ampoule into 3 mL of sterile water. Mix thoroughly.

- b. Use a sterile swab (Biolog P/N 3023 or 3021) wet with the thioglycolate solution to spread a thin film across the entire surface of the BUG + M agar plate (Biolog P/N 71103).
 - c. Allow the thioglycolate to dry on the agar before streaking the spore-forming gram-positive rod strain (takes about 5 minutes).
 3. With a Biolog sterile Streakerz™ stick (Biolog P/N 3025 or 3026), touch an isolated colony and transfer to the BUG + M + T agar media. Make a plus sign (+) on the center of the agar media. See Figure 1.
 - a. If the organism is a fast grower and spreads, make a single very thin plus sign (+). This will avoid overgrowth and sporulation.
 - b. If the organism is a very slow grower and does not spread then make up to three (3) streaks in the shape of the letter N in each direction of the plus. Leave a small distance between the three streaks to promote vegetative growth on all three streaks. Prepare multiple agar plates (2 to 3) to insure that enough cells are available to prepare the inoculation suspension.


Figure 1: Plus Format to Streak

A



PLUS "+" STREAKING METHOD

B



CELL GROWTH AFTER INCUBATION

- c. The goal of this streaking technique is to limit cell growth to one (1) to three (3) lines so that the cells along the edges have a good supply of nutrients. This keeps the cells in an active state and decreases sporulation.
 4. Incubate at 30°C for 16-20 hours; do not exceed 16 hours for fast-growing strains.

5. Prepare a cell suspension with the following “**Dry Tube Method**” for mucoid or dry colony types.
 - a. Use a sterile, dry glass test tube.
 - b. Using a Streakerz stick, take cells starting at the ends and edges of the “Plus” sign as shown by the arrows in the Figure 1. The cells located outside the boxed area illustrated in Figure 1 are the most metabolically active cells. Avoid picking cells from the boxed area that are more likely to have sporulated. Be careful not to scrape up agar from the plate while harvesting the cells.
 - c. Some mucoid strains as well as strains that become embedded in the agar may not adhere to the Streakerz stick making it difficult to transfer cell mass to the dry test tube. In this case, use the edge of a sterile microscope slide to scrape across the inoculation lines on the agar plate and recover cell mass. Be careful not to dig in to the agar. Transfer the cell mass from the slide to the inner wall of a dry test tube with the aid of a Streakerz stick.
 - d. Crush and break up the cell mass against the side of the tube. Use an up and down motion in combination with a circular motion to crush and spread the cells against the inner wall of the dry glass tube. Continue until a sufficient mass of cells adheres to the wall of the glass tube forming a uniform thin film on the inside surface free of clumps.
 - e. Transfer approximately 5 mL from a tube of GN/GP-IF (Biolog P/N 72101) into the glass tube containing the cells. Use an up and down rubbing motion to slide the bacterial film from the tube wall downward, into the GN/GP-IF. Continue until all cells are suspended in the inoculation fluid.
 - f. Examine the cell suspension. Hold the tube up to the light and see if the cells are completely suspended or if clumps and strings remain. If clumps or strings remain, follow one these methods to remove them:

Method 1: If the suspension is a mixture of single cells and clumps, let the mixture stand at room temperature until the clumps and debris have settled to the bottom. Use the homogenous cell suspension layer free of debris. If the suspension contains only strings and clumps (observed with highly mucoid strains like *Bacillus amyloliquefaciens* and *licheniformis*), place the suspension test tube in a warm ($36 \pm 1^\circ\text{C}$) incubator or water bath for up to 20 minutes. The heat

aids in liquefying the cell mucous. Swirl the contents of the tube occasionally during the incubation. Look for the disappearance of the distinct mucoid strings as evidence that the sample is ready.

Settling of the debris is critical to eliminating clumps and obtaining a homogeneous cell suspension for accurate identification. This step may take up to 20 minutes.

OR

Method 2: Filter the cell mixture through a 70 µm cell strainer (BD Falcon™ 352350) to pull out the cell debris from the homogeneous cell suspension.

- g. Prepare a homogeneous cell suspension. Place a new tube of GN/GP-IF into the Biolog Turbidimeter. Adjust the zero of the turbidimeter for the tube. Carefully remove the homogenous cell suspension using a transfer pipette (Biolog P/N 3019) making sure that the cell debris in the bottom of the tube is not taken. Add the cell suspension drop wise, with mixing to the tube of GN/GP-IF in the turbidimeter. Use a swab to stir and mix the cell suspension into the GN/GP IF. Continue to add cells to achieve a cell turbidity of $28\% \pm 2\%$. Uniformly suspended cells typically give stable turbidimeter readings whereas remaining clumps and strings cause the turbidimeter to fluctuate as they pass across the light path. Remaining clumps can be removed by letting the tube stand until the clumps settle. Remove the uniformly suspended cells to inoculate the MicroPlate.
 - h. If more cells are needed to get the required cell density, repeat the dry tube method preparation using a new empty, sterile dry glass tube. Do not use the same tube that is wet with GN/GP-IF.
6. Inoculate a GP2 MicroPlate™ following the procedure in the Instructions for Use.
7. Incubate at 30°C for 4 – 6 hours and/or 16 – 24 hours to allow the pattern to form.

General hints for handling *Bacillus*

- ☒ Always grow *Bacillus* species on BUG + M + T for identification.
- ☒ Keep the incubation time on BUG + M + T as short as possible to reduce sporulation (i.e., incubate for 16 rather than 24 hours if the organism grows fast).
- ☒ For slow growing strains you may need to use multiple (2 to 3) BUG + M + T agar plates to obtain enough cell mass.
- ☒ Use the dry tube method described above to properly prepare the inoculum. Allow the suspension to stand until all clumps settle to the bottom of the tube. Be patient.
- ☒ *Bacillus cereus* and *Bacillus thuringiensis* are considered indistinguishable by biochemical assays. These two species are listed in the Biolog database as *Bacillus cereus/thuringiensis*. *Bacillus anthracis* is also very closely related to these taxa.

Special Procedures for Filamentous Fungi

Because of their potentially pathogenic nature, filamentous fungi need special care and handling.

Laboratory Safety¹

Sporulating cultures require extra handling precautions to prevent spores or conidia from escaping into the air. Manipulating sporulating cultures can sometimes present a risk of infection or allergy, contaminate laboratory air, and contaminate bacterial cultures. We recommend that filamentous fungi be handled in a biological safety cabinet. This practice also diminishes the risk of infection by dimorphic fungal pathogens.

Follow the guidelines in this chapter for all specimen preparation, characterization, inoculum preparation, and MicroPlate inoculation.

General hints for handling filamentous fungi

1. Incubate subcultured filamentous fungi in a closed container. Place the agar plates in a plastic container or in sealed plastic bags on trays.

¹ Several organizations have published guidelines for laboratory safety in microbiology – the Centers for Disease Control and National Institutes of Health (1988), the Medical Research Council of Canada (1990), and the World Health Organization (1993).

2. Do **NOT** add a source of moisture to this container. Mycelial growth will overwhelm the agar plate and cause crossover contamination of other cultures.
3. Incubate the 2% ME agar at 26° C for 5-10 days to induce conidia formation (incubate longer if required). Incubate yeasts for 1-2 days.

Inoculating MicroPlates

Inoculate the suspension into one of the following Biolog MicroPlates:

- GN2: For all gram-negative aerobic bacteria
- GP2: For all gram-positive aerobic bacteria
- AN: For all anaerobic bacteria
- YT: For all yeasts visceral
- FF: For all filamentous fungi and select yeast strains

Pipet the inoculum into a Biolog MicroPlate within 20 minutes. Waiting any longer may cause inaccurate identification. Anaerobic bacteria are especially sensitive to delays. If you are running a batch of MicroPlates, set them up (from preparing the inoculum to pipetting into the MicroPlate) so you will not exceed the 20 minute deadline.

Caution!
Pipet the
inoculum into a
MicroPlate
within
20 minutes!
Use sterile
technique.

Inoculating protocol

1. Label the MicroPlate with the organism name/number and plate type (e.g., GN, GP). Label the side of the MicroPlate itself, not the lid.
2. Using aseptic technique, pour the cell suspension into a multichannel pipette reservoir.
3. Firmly attach eight sterile tips to the pipettor.
4. Fill the tips with the suspension. Check to see that all tips fill equally.
5. Prime the tips by dispensing once back into the reservoir. The electronic pipettor will do this automatically.
6. Fill all MicroPlate wells by placing the tips at an angle inside the top of the wells. Take care not to splash from one well to another. Avoid contamination. Avoid touching the bottom of the wells, which could transfer carbon sources.
 - Gram negatives and gram positives → 150 µl per well
 - Anaerobes, yeasts and filamentous fungi → 100 µl per well

Caution!
Accurately
pipet
recommended
volumes.

7. If the fluid level in the tips gets low, refill and continue dispensing until all wells are full.
8. Cover the MicroPlate with its lid.

Incubating MicroPlates

As soon as you dispense the suspension, incubate the MicroPlate using the appropriate temperature, atmosphere, and time conditions. Table 4-3 (page 4-15) and Appendix 2 will help you select the correct conditions.

- For aerobes: Use a moist chamber during incubation to prevent evaporation. A plastic container with wet paper towels on the bottom should be sufficient.
- For filamentous fungi: MicroPlates must be incubated in a closed container or in sealed plastic bags on trays. Do not add a source of moisture. Mycelial growth will overwhelm the MicroPlate and cause reaction crossover within the MicroPlate.

MicroPlates of aerobic bacteria can be read after 4-6 hours of incubation. If reactions are insufficient, read again at 16-24 hours.

| Organism Type | Temperature | Atmosphere | Incubation Time |
|---|--|---|------------------------|
| Gram Negative Non-Enteric | 30° C | Air | 4-6 and 16-24 hours |
| Gram Negative Enteric | 35-37° C | Air | 4-6 and 16-24 hours |
| Gram Negative Fastidious | 35-37° C | 6.5% CO ₂ | 4-6 and 16-24 hours |
| Gram Positive Cocci Gram Positive Rods (non-spore forming) | 35-37° C or 30° C or 26° C (if required) | Air or 6.5% CO ₂ (if required) | 4-6 and 16-24 hours |
| Gram Positive Rods (spore forming bacillus) | 30° C or 55° C (if required) | Air | 4-6 and 16-24 hours |
| Anaerobes | 35° C or 30° C or 26° C (if required) | Anaerobic-H ₂ free | 20-24 hours |
| Yeasts | 26° C | Air | 24, 48, or 72 hours |
| Filamentous Fungi (and select yeast species) | 26° C | Air | 24, 48, 72 or 96 hours |

TABLE 4-3: SELECTING CORRECT INCUBATING CONDITIONS

Special Procedures for Incubating Anaerobic MicroPlates

1. In most cases, leave AN MicroPlates in their foil pouch until immediately before inoculation. The only exception is if you presume your organism is either *Bacteroides* or *Prevotella*. If your gram-negative rod is resistant to kanamycin (see page 4.5) and has grown well enough overnight to prepare a suspension, open the AN MicroPlate foil pouch and expose it to air for 20 minutes before inoculation.
2. For all species, keep the air exposure time consistent. Do **NOT** exceed 5 minutes between the time you begin and end inoculating AN MicroPlates for a single run. You should be able to inoculate six MicroPlates in 5 minutes. If the run requires more than six MicroPlates, use two people or divide the number of suspensions into smaller batches.
3. Wait 10 minutes after inoculating the AN MicroPlate before placing it in an anaerobic jar. This is necessary to slightly oxidize the buffer.
 - Rectangular jars can hold six or 24 MicroPlates, are easy to work with, and are inexpensive
 - 3.5 liter jars hold up to nine MicroPlates
 - 10.5 liter jars hold up to 28 MicroPlates and can stand upright
4. Use an anaerobic, atmosphere-generating system that does not produce hydrogen. Organisms with strong hydrogenases will reduce the tetrazolium in all wells if hydrogen is present. Add an anaerobic indicator strip to monitor the atmosphere. Call Biolog Technical Service for suggestions on where to find anaerobic supplies.
5. Incubate most anaerobes at 35° C for 20-24 hours.
6. After incubation, check the color of the MicroPlates before you open the anaerobic jar. MicroPlates should be clear with purple wells. Once you expose them to air, negative wells will slowly take on a faint green-blue color. If this green-blue color is present before you open the jar, the atmosphere was probably not fully anaerobic.

Caution!

If you're working with anaerobes:
→ Use anaerobically prepared media.
→ Inoculate MicroPlates in air.
→ After inoculation, wait 10 minutes before moving the MicroPlates into a hydrogen-free anaerobic jar.

Sample Preparation Process

For Filamentous Fungi

| | | |
|---|-------------------|----------------|
| Initial Culture Medium | 2% ME | 2% ME |
| Wet Prep Results | Hyphal elements | Yeast cells |
| Microbe Type | Filamentous Fungi | Yeast |
| Culture Medium | 2% ME | 2% ME |
| Temperature | 26° C | 26° C |
| Atmosphere | Air | Air |
| Culture Incubation Time | 5 - 10 days | 2 days |
| Inoculating Fluid | FF-IF | FF-IF |
| Inoculum Turbidity/ Turbidity Standard | 75% T FF | 75% T FF |
| MicroPlate Type/μl per well | FF 100 | FF 100 |
| Incubation time (hours) | 24, 48, 72, 96 | 24, 48, 72, 96 |

ABBREVIATIONS

2% ME = 2% Malt Extract Agar
 FF = Filamentous Fungi
 %T = percent transmittance
 IF = inoculating fluid

Sample Preparation Process For Gram Negatives, Gram Positives, Anaerobes and Yeasts

| | | | | | | | |
|---|--|-------------------------------------|--|--|---|--|---|
| Initial culture medium | BUG + B* | | | | | BUA + B | BUY |
| | Gram stain and observe cell morphology | | | | | | |
| Gram stain results | <div> <div>↓</div> <div>Gram Negatives</div> <div>↓</div> </div> | | | <div> <div>↓</div> <div>Gram Positives</div> <div>↓</div> </div> | | <div> <div>↓</div> <div>Anaerobes</div> </div> | <div> <div>↓</div> <div>Yeasts</div> </div> |
| Characterization test | Oxidase positive requires 30° C | Oxidase negative and TSI=A/A or K/A | Requires CO ₂ or chocolate agar or < 1 mm colonies on BUG + B | | | Kanamycin test for rapidly-growing, gram-negative rods | |
| Characterization test | Oxidase negative and TSI=K/K or K/A ^w | | | | | | |
| Microbe type | GN-NENT | GN-ENT | GN-FAS | GP-COCCUS-ROD, GP-COCCUS, GP-ROD | GP-ROD SB (spore-forming bacillus) | AN | YT |
| Culture medium | BUG + B or TSA + B | BUG + B or TSA + B | Chocolate | BUG + B | BUG + M + T (8drops in 3ml H ₂ O) “+”streaking | BUA + B | BUY |
| Temperature | Typically 30° C | Typically 35-37° C | Typically 35-37° C | Typically 35-37° C | Typically 30° C | Typically 35-37° C | 26° C |
| Atmosphere | Air | Air | 6.5% CO ₂ | Air or 6.5% CO ₂ (if required) | Air | ANA-H ₂ free | Air |
| Inoculating fluid | GN/GP-IF (if A1 well is pos, add thioglycolate) | GN/GP-IF + T (3 drops) | GN/GP-IF + T (3 drops) | GN/GP-IF + T (3 drops) (if A1 well is pos, add 1 ml salicylate) | GN/GP-IF (NO THIO) | AN-IF | Water |
| Inoculum turbidity/ Turbidity standard | 52% T GN-NENT | 61% T GN-ENT | 20% T GP-COC & GP-ROD & GN-FAS | 20% T GP-COC & GP-ROD & GN-FAS | 28% T GP-ROD SB | 65% T AN | 47% T YT |
| MicroPlate type/ µl per well | GN2 150 | GN2 150 | GN2 150 | GP2 150 | GP2 150 | AN 100 | YT 100 |
| Incubation time (hours) | 4-6, 16-24 | 4-6, 16-24 | 4-6, 16-24 | 4-6, 16-24 | 4-6, 16-24 | 20-24 | 24, 48, 72 |

*Note: Agricultural bacteria may be grown on BUG without blood.and at 30° C for enterics.

ABBREVIATIONS

TSI = Triple Sugar Iron Slant; A = Acid, A^w = Weak Acid, K = Alkaline; BUG + B = Biolog Universal Growth Medium + blood; BUG + M = Biolog Universal Growth Medium + maltose; T = thioglycolate; %T = percent transmittance; GN-NENT = gram-negative non-enteric; GN-ENT = gram-negative enteric; GN-FAS = gram-negative fastidious; GP-COCCUS & GP-COC = gram-positive cocci; GP-ROD = gram-positive rod; GP-ROD SB = gram-positive rod, spore-forming bacillus; AN = anaerobe; YT = yeast; IF = inoculating fluid

5. Reading MicroPlates

In this section:

- Logging-In
- Choosing Manual, Reader or File Mode
- Reading Plates Manually or using the Plate Reader
- Reading from a Saved File
- Setting Up a Worksheet
- Setting Up to Save Files
- Saving Files
- Using a Worksheet to Read Multiple MicroPlates
- Adjusting Thresholds Manually
- Choosing a Database to Search
- Exiting

Depending on which MicroStation/ MicroLog software version you are using, the process of reading MicroPlates can be done in a simplified way or using the detailed method of creating worksheets.

- MicroLog 1 software offers only the simplified reading process (entering MicroPlate reactions manually’ “by eye”, printing the results, and interpreting the results.
- MicroLog 2 software uses the same manual reading process as the MicroLog 1 but also allows you to save results, as well as create a user built database.
- MicroStation/ MicroLog 3 software allows you to create worksheets to read multiple MicroPlates, save files, work with databases, and use advanced functions (such as cluster analysis and threshold adjustment)

Logging-In

Unless the software has been set up to run in Unrestricted mode, the MicroStation system/ MicroLog software will always require user log-in before beginning the process of entering information and reading MicroPlates.

Follow the instructions for logging-in on page 3.2.


Choosing Manual, Reader or File Mode

You can run the identification process in three modes:

- **Manual** mode allows you to enter MicroPlate reactions by hand SEE Section 5, Page 3
- **Reader** mode is used in conjunction with the MicroStation reader, where MicroPlate reactions are read and entered automatically. SEE Section 5, Page 7
- **File** mode lets you recall and view MicroPlate data already stored in a file. SEE Section 5, Page 10

Note: While running in Restricted mode, you have access to these features only if you have “Log-In” and “Set Up” privileges.


To select whether you want to work in manual, reader, or file mode, click the **Set Up** tab on the **Welcome** window.



WELCOME WINDOW

Selecting Data Input (Read) mode

1. Use the **Data Input Mode** drop-down list to select **Manual**, **Reader**, or **File**.



SET UP WINDOW (WITH MANUAL MODE SELECTED)

2. Select the correct MicroStation reader. **Molecular Devices** or **Bio-Tek**.
3. Select **Yes** or **No** in the **Use Worksheet** drop-down list. For further information on using a worksheet, see Section 5, page 12-16.

4. Select **Yes** or **No** in the **Write Data to Data File** drop-down list. For further information on writing to a data file, see Section 5 pages 15-16. **Note:** If the program is operating in **Restricted Access Mode**, this option defaults to **Yes**, and only users with **SET UP** privileges have the ability to change this option.

Choosing printer method

1. Click the **Automatic Printout** drop-down list on the **Set Up** window to select **Yes** (if you want to automatically print out your result) or **No** (if you don't want automatic printout).
2. Click the **Print in Text Mode** drop-down list to select **Yes** (if you want speedier printing) or **No** (if you want higher resolution printing).
3. Click the **Printer Setup** bar if you need to select a printer. This bar is not necessary for text printing and will disappear if you select text mode.

Reading Plates Manually

1. Use the Data Input Mode drop down menu to select **Manual** as the Data Input Mode.
2. Click the **Data** tab on the Menu Bar. The **Data** window appears.

Caution!
Entering reactions manually is NOT recommended for YT and FF MicroPlates.

The screenshot shows the MicroLog 4.20.05 interface. The 'Data' tab is selected. The 'Plate Info' section contains the following fields:

- Current Time: Apr 15 2007 21:42
- Plate Number: 1
- Sample Number: (empty)
- Plate Type: GN2
- Strain Type: NOT SELECTED
- Strain Name: (empty)
- Strain Number: (empty)
- Incubation Time: NOT SELECTED
- Other: (empty)

The 'Pos/Neg' grid shows a 12x12 layout of wells. The grid is currently empty, with all wells showing a blank circle.

Below the grid are buttons for Re-Read, Read Next, Print, Save, View Database, and Cluster Menu.

At the bottom, a table lists data for wells 1 through 12, with columns for NAME, PROB, SIM, DIST, and TYPE.

| ML4 | NAME | PROB | SIM | DIST | TYPE |
|-------|------|------|-----|------|------|
| 1 | .. | .. | .. | .. | .. |
| 2 | .. | .. | .. | .. | .. |
| 3 | .. | .. | .. | .. | .. |
| Other | .. | .. | .. | .. | .. |

DATA WINDOW (IN MANUAL MODE)

When you're entering sample identifiers, make sure you select the correct strain type. The default entry in that field is "NOT SELECTED." If default strain type is selected, no database will be searched

Entering sample identifiers

1. **Current Time:** Your computer's clock automatically tracks time.
2. **Plate Number:** Type in the MicroPlate number.
3. **Sample Number:** Type in the sample ID number.

The screenshot shows the MicroLog 4.20.05 software interface. The 'Plate Info' section contains the following fields:

- Current Time: Apr 15 2007 21:42
- Plate Number: 1
- Sample Number: (empty)
- Plate Type: GN2 (dropdown menu)
- Strain Type: GN2 (dropdown menu)
- Strain Name: AN
- Strain Number: YT
- Incubation Time: FF
- Other: (empty)

The 'Pos/Neg' section shows a grid of 12 columns and 8 rows (A-H) with circles indicating results. The bottom section shows a table with columns ML4, NAME, PROB, SIM, DIST, and TYPE, with rows 1, 2, 3, and Other.

SAMPLE IDENTIFIER SECTION OF DATA WINDOW

4. **Plate Type:** Use the drop-down list to select the type of MicroPlate you used.
5. **Strain Type:** Use the drop-down list to select the strain type you are working with.
 - For an initial bacterial query use one of the following: GN-NENT, GN-ENT, GN-FAS, GP-COCCUS, GP-ROD, GP-ROD SB, AN-ALL (use the NON-APPLICABLE default for other MicroPlate types.)
 - Secondly, use GN-ALL or GP-ALL to screen organisms of questionable type and morphology as required.
 - If an ID is obtained using GN-ALL or GP-ALL be sure to check that the correct protocol for that ID was used.
 - For FF query Use the FF-AIR or FF-FOOD database for FF MicroPlates, depending on the source of the organism to be identified. These databases also contain the most commonly isolated species from the other databases (Yeast, Aspergillus, Colletotrichum, Fusarium, Penicillium, Trichoderma). The other

databases are intended for use by researchers working on specific genera, not for routine food and air isolates.

6. **Strain Name, Number:** If you have additional information, enter it in these optional fields.
7. **Incubation Time:** Use the drop-down list to select the hours that the MicroPlate has incubated.

Note: A 240-hour incubation time is listed for the FF database, however, there is no database for 240 hours. That listing is a research tool.

8. **Other:** If you have additional information, enter it in this optional field (colony morphology, culture conditions, etc.). An optional comment can be entered in this field.

Remember:
Light purple reactions are considered positive as long as the color is noticeable when compared to the A1 reference well.

Entering reactions manually


Entering reactions in manual mode requires visual assessment of MicroPlate reactions.

Table 5-1 will guide your visual assessment:

| Color Density | How To Assess It |
|---|------------------|
| Same as A1 well (usually colorless) | Negative |
| Noticeable purple color (more than A1 well) | Positive |
| <u>Extremely</u> faint color, Not sure | Borderline |

TABLE 5-1. VISUALLY ASSESSING WELL COLOR DENSITY

1. From the **Data** window, click **Read Next** (directly under sample information). The **Manual Entry of Data** window appears.



MANUAL ENTRY OF DATA WINDOW (TO ENTER MICROPLATE REACTIONS)

2. Enter reactions into the on-screen MicroPlate, using any of the methods shown below.
3. Click **OK** to read the Plate
4. The **Data** screen appears and displays the read result.

SELECT YOUR PREFERRED WAY TO MANUALLY ENTER REACTIONS

Mouse Keypad or Keyboard Method

Use your mouse to maneuver around the on-screen keypad or use the computer keyboard to enter reactions on the Manual Entry of Data window.

| | | |
|--------------------|---|----------------------------------|
| "1" key | ➔ | negative |
| "2" key | ➔ | borderline |
| "3" key | ➔ | positive |
| up and down arrows | ➔ | use to navigate around the wells |
| "home" key | ➔ | returns cursor to A1 well |
| "end" key | ➔ | moves cursor to H12 well |

Mouse Cursor Method

Select the radio button you want, then click your cursor in desired wells. With this method, you can also click and drag across areas to change large blocks of wells.

| | | |
|--------------|---|-----------------------------------|
| ● Positive | ➔ | makes selected well(s) positive |
| ● Borderline | ➔ | makes selected well(s) borderline |
| ● Negative | ➔ | makes selected well(s) negative |

All or Nothing Method


One click to make the whole MicroPlate positive or negative, then use radio buttons and your cursor to make alterations.

| | | |
|-------------|---|--------------------------|
| ■ - All Neg | ➔ | makes all wells negative |
| ■ + All Pos | ➔ | makes all wells positive |

Reading Plates Using the Plate Reader

Caution!
The Detailed Reader Setup window contains troubleshooting tools. Don't change anything without consulting Biolog Technical Service.

1. Use the **Data Input Mode** drop-down list to select **Reader** (to use the reader).
2. Select the correct MicroStation reader. **Molecular Devices** or **Bio-Tek**.
3. The **Initialize Reader** and **Detailed Reader Setup** selection bars appear. Use the **Com Port** drop-down list to assign communications (com) port 1 thru 5.



SET UP WINDOW (WITH READER MODE SELECTED)

4. Click the **Initialize Reader** selection bar. The **Reader Ready** field changes from **No** (red) to **Yes** (green).

Note: If the Reader Ready field does not change to Yes (green) after 1 minute, make sure the reader is connected properly and turned on. See Section 11 if you have trouble initializing the reader.

5. Select **Yes** or **No** in the **Use Worksheet** drop-down list. For information on using a Worksheet see Section 5, Page 11-14.
6. Select **Yes** or **No** in the **Write Data to Data File** drop-down list. For information on Saving Files see Section 5, page 15-16.

Entering Reactions Using the Reader

** For FF strain type query see section 4 page 5.*

1. Make sure the reader is initialized (on the Set Up window).
2. Enter sample identifiers.

Entering sample identifiers

-Current Time: Your computer's clock automatically tracks time.

-Plate Number: Type in the MicroPlate number.

-Sample Number: Type in the sample ID number.

MicroLog 4.20.05

Welcome Setup Data Data File Exit

Plate Info

Current Time: Apr 15 2007 21:42

Plate Number: 1

Sample Number:

Plate Type: GN2

Strain Type: GN2

Strain Name: AN

Strain Number: YT

Incubation Time: FF

Other: MT, ECO, OTHER

Pos/Neg

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| A | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| B | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| C | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| D | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| E | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| F | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| G | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |
| H | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> | <input type="radio"/> |

Re-Read Read Next Print Save View Database Cluster Menu

No ID Yet

| ML4 | NAME | PROB | SIM | DIST | TYPE |
|-------|------|------|-----|------|------|
| >1 | | | | | |
| 2 | | | | | |
| 3 | | | | | |
| Other | | | | | |

SAMPLE IDENTIFIER SECTION OF DATA WINDOW

-Plate Type: Use the drop-down list to select the type of MicroPlate you used.

-Strain Type: Use the drop-down list to select the strain type you are working with.

- For an initial bacterial query use one of the following: GN-NENT, GN-ENT, GN-FAS, GP-COCCUS, GP-ROD, GP-ROD SB, AN-ALL (use the NON-APPLICABLE default for other MicroPlate types.)
- Secondarily, use GN-ALL or GP-ALL to screen organisms of questionable type and morphology as required.
- If an ID is obtained using GN-ALL or GP-ALL be sure to check that the correct protocol for that ID was used.
- For FF query Use the FF-AIR or FF-FOOD database for FF MicroPlates, depending on the source of the organism to be identified. These databases also contain the most commonly isolated species from the other databases (Yeast, Aspergillus, Colletotrichum,

Fusarium, Penicillium, Trichoderma). The other databases are intended for use by researchers working on specific genera, not for routine food and air isolates.

-Strain Name, Number: If you have additional information, enter it in these optional fields.

-Incubation Time: Use the drop-down list to select the hours that the MicroPlate has incubated.

Note: A 240-hour incubation time is listed for the FF database, however, there is no database for 240 hours. That listing is a research tool.

-Other: If you have additional information, enter it in this optional field (colony morphology, culture conditions, etc.).

-Comment: An optional comment can be entered in this field.

*Check your
MicroPlate
before
loading.
Make sure
the reaction
pattern is
well
developed.*


3. Wipe the bottom of the MicroPlate to remove condensation and fingerprints.
4. Place the MicroPlate on the reader tray with the A1 well at the top left-hand corner.
5. Gently push the MicroPlate down until it snaps into a level position in the reader tray.
6. Remove the MicroPlate lid.

Note: For FF, leave the lid on. Replace it with a new lid if condensation is present.

7. Click **Read Next** on the **Data** window.
8. The reader calibrates itself, pulls the MicroPlate inside, takes a reading, displays the result on the **Data** window, and ejects the MicroPlate. Briefly do a visual check of the reaction with the +/-/- called by the auto-thresholding function (see pg 5.20 if adjustment is required).


Reading from a Saved File

1. Use the **Data Input Mode** drop-down list to select **File** (to view reactions from a saved file)
2. The **Input Data File Name** field appears. Click in the empty field and move through the listed folders to find the file you want to read and click on it. This will bring the file name up on the **Set up** screen.



SET UP WINDOW (WITH FILE MODE SELECTED)

3. In the **Write Data to Data File** drop-down list
 - a. Select **No**.
4. Click the **Data** tab.
5. Click **Select Read**
6. A list of saved data appears. Select the MicroPlate you want to read and click **OK**.



C:\Program Files\Biolog\ML3_42_05\datafiles\GNDemo.D4C

Data file = C:\Program Files\Biolog\ML3_42_05\datafiles\GNDemo.D4C

| # | OK | Not OK | Created | Date | Sample | by | Notes | Strain Type | Strain Name | Strain # |
|----|----|--------|-------------------|------|--------|-----|----------------|-----------------------|-------------|----------|
| 16 | OK | | Apr 16 1990 11:17 | | 24 HR | GN2 | ON-ENT | ESC. COL | 5574 | |
| 17 | OK | | Apr 15 1990 11:17 | | 4 HR | GN2 | ON-ENT | OT BRA | 13651 | |
| 18 | OK | | Apr 15 1990 11:17 | | 4 HR | GN2 | ON-ENT | ENT CLO | 11534 | |
| 19 | OK | | Apr 16 1990 11:14 | | 4 HR | GN2 | ON-ENT | OT MUR | 13656 | |
| 20 | OK | | Apr 17 1990 01:43 | | 24 HR | GN2 | ON-ENT | ENT ADO BO S (BPAV'S) | 11530 | |
| 21 | OK | | Apr 17 1990 01:43 | | 24 HR | GN2 | ON-ENT | ENT CLO | 11532 | |
| 22 | OK | | Apr 17 1990 01:43 | | 24 HR | GN2 | ON-ENT | KLUB PNE | 4145 | |
| 23 | OK | | May 14 1990 11:52 | | 4 HR | GN2 | ON-FAS OXI-HAE | ACT | 058 | |
| 24 | OK | | May 14 1990 11:52 | | 4 HR | GN2 | ON-FAS OXI-CAP | OCHSRU | 5618 | |
| 25 | OK | | May 15 1990 01:10 | | 24 HR | GN2 | ON-FAS OXI-HAE | INF | 5641 | |
| 26 | OK | | May 15 1990 01:10 | | 24 HR | GN2 | ON-FAS OXI-HAE | APH | 2704 | |
| 27 | OK | | May 15 1990 11:53 | | 4 HR | GN2 | ON-FAS OXI-HAE | ACT | 8847 | |
| 28 | OK | | May 15 1990 11:53 | | 4 HR | GN2 | ON-FAS OXI-HAE | HAE | 6755 | |
| 29 | OK | | May 15 1990 11:53 | | 4 HR | GN2 | ON-FAS OXI-HAE | INF | 813 | |
| 30 | OK | | May 17 1990 11:15 | | 24 HR | GN2 | ON-FAS OXI-CAP | OCHSRU | 7228 | |

Page 2/2

File View ON

OK

LIST OF SAVED DATA

7. The Data page will show the record.
 - a. Click **Print** to print the record.
 - b. Click **Read Next** to view the next record.
 - c. Click **Select Read** to select another record.
8. Repeat steps 6 thru 7 to view all records.

Setting Up a Worksheet

Worksheets offer an efficient way to log in and read multiple MicroPlates. By entering multiple sample identifiers and MicroPlate reactions, you simply read MicroPlates sequentially as the worksheet indicates, without stopping to enter more information. If you're working with over five MicroPlates, using a worksheet is faster. If you're reading ten MicroPlates or more, worksheets will save significant time.

Worksheets also allow you to:


- Have one shift set up a worksheet, and the second shift read the MicroPlates
- Avoid entering data fields more than once if you plan to read plates at multiple incubation time points

What's a Worksheet?
A worksheet is a spreadsheet within MicroStation/MicroLog software. Use it to log in and read multiple MicroPlates efficiently.


Printing worksheets gives a record of batch activity. You can use a worksheet whether in manual or reader mode. Worksheet use is optional.

Note: If the program is running in Restricted mode, you have access to this feature only if you have "Log-In" privileges.

1. On the **Set Up** window, select **Yes** on the **Use Worksheet** drop-down list.
2. Click on the **Worksheet Name** field. A **Worksheet File** dialog box appears. Type the desired file name (W4C extension) or choose a previously created worksheet to edit. (Note: Save worksheets that have been created in the "DataFile" folder, or create a separate Worksheet folder. This will assist you in recalling the files at a later date.)



3. Click the **Worksheet Entry** selection bar. The **Worksheet** window appears.



Entering sample identifiers

4. Click **Add Entry**. The **Plate Information** window appears.

Plate Information

C:\Program Files\Biolog\ML3_42_05\datafiles\test.W4C

Plate Number: 1

Sample Number:

Plate Type: GN2

Strain Type: GN-NENT OXI+

Strain Name:

Strain Number:

Other:

Buttons: Delete, Save, Clear, Save, Keep, Save, Done, Cancel

PLATE INFORMATION WINDOW (TO ENTER DATA INTO A WORKSHEET)


5. For each MicroPlate, enter the following information:

- **Plate Number:** Default begins with 1 (changeable)
- **Sample Number:** Type in the sample ID number (optional field).
- **Plate Type:** Use the drop-down list to select the type of MicroPlate you used.
- **Strain Type:**
 - For an initial bacterial query use one of the following: GN-NENT, GN-ENT, GN-FAS, GP-COCCUS, GP-ROD, GP-ROD SB, AN-ALL (use the NON-APPLICABLE default for other MicroPlate types.)
 - Secondly, use GN-ALL or GP-ALL to screen organisms of questionable type and morphology as required.
 - If an ID is obtained using GN-ALL or GP-ALL be sure to check that the correct protocol for that ID was used.
 - Use the FF-AIR or FF-FOOD database for FF MicroPlates, depending on the source of the organism to be identified. These databases also contain the most commonly isolated species from the other databases (Yeast, Aspergillus, Colletotrichum, Fusarium, Penicillium, Trichoderma). The other databases are intended for use by researchers working on specific genera, not for routine food and air isolates. The YT MicroPlates does not have a list of strain types from which to choose.

- **Strain Name, Number, Other:** Enter other information in these optional fields. (In the “Other” field, you can enter colony morphology, culture conditions, etc.).
6. After each entry, click the selection bars as desired:
- **Save, Clear** – to save data, to advance MicroPlate numbers incrementally, to clear all fields for entering information on the next MicroPlate
 - **Save, Keep** – to save data, to advance MicroPlate numbers incrementally, to keep current entries in all fields
 - **Save, Done** – to return to the worksheet
7. When you’re finished setting up the worksheet, click **OK**.

Editing a worksheet

1. To edit a worksheet, right click on the row you want to edit. The **Plate Information** window reappears. Edit fields as desired.



WORKSHEET (WITH SAMPLE ENTRIES)

2. Sort entries (Plate Number) by clicking **Sort Worksheet**.
3. To delete an entry click on the row and then click **Delete Record**.

Setting Up to Save Files

Note: If you are planning to compile your own database, see Section 7 before saving files.


If you plan to compile your own database, see Section 7 before saving files.

1. Click the **Write Data To Data File** drop-down list on the **Data** window to select **Yes** (if you want the software to save files) or **No** (if you don't want to save files).

*Note: This feature gives a default answer of **Yes** in the **Write Data To Data File** and the **Auto Write to Data File** fields when the program is running in Restricted mode. "Set Up" privileges are required to change these defaults.*

2. If you selected No, the two fields below will close. If you selected **Yes**, two new fields remain. Click on the **Output Data File Name** field. A **Save As** dialog box appears. Enter the desired file name (D4C extension) or choose a previously created data file.
3. Use the **Auto Write To Data File** drop-down list to select **Yes** (if you want to automatically save your results) or **No** (if you don't want to automatically save your results).

Note: For easier editing and compiling of saved data in the future, separate different plate types into individual data files.



CLOSE-UP OF OUTPUT PORTION OF SET UP WINDOW


*Note: If you are working in Restricted mode, new entries are saved as **Read Only** files. Original files cannot be edited. Make a backup for editing.*

Using Worksheets to Read Multiple MicroPlates

To efficiently read a batch of MicroPlates, use a worksheet. Make sure you're in the desired mode (reader or manual) and that you've selected **Yes** for worksheet use on the **Set Up** window. Select the desired file name.

Reading multiple MicroPlates in manual mode

1. Select **Manual** mode on the **Set Up** window.
2. Click the **Data** tab on the menu bar. The **Data** window appears.



DATA WINDOW

3. Use the **Incubation Time** drop-down list to select the desired incubation time.
4. To read from an open Worksheet:
 - Click **Select Read** to go to the desired worksheet.
 - Click on the sample number you want to read (an arrow will appear next to the sample number you select).
 - Click **OK** to read the MicroPlate. The **Manual Entry of Data** window appears.
 - If necessary, click **Cancel** to return to the worksheet.
5. To read MicroPlates sequentially:

- Click **Read Next** to sequentially read the next MicroPlate on the worksheet without returning to the **Worksheet** window.
- A **Confirmation** window appears. Click **OK** if it displays the MicroPlate you wish to read.
- The **Manual Entry of Data** window appears.
- If necessary, click **Cancel** to return to the **Data** window.

Assessing reactions

Entering reactions in manual mode requires visual assessment of MicroPlate reactions. SEE Section 5, page 5 for more information on reading reactions manually.

Table 5-1 will guide your visual assessment:


| Color Density | How To Assess It |
|---|------------------|
| Same as A1 well (usually colorless) | Negative |
| Noticeable purple color (more than A1 well) | Positive |
| <u>Extremely</u> faint color. Not sure | Borderline |

TABLE 5-1. VISUALLY ASSESSING WELL COLOR DENSITY

Remember:
Light purple reactions are considered positive as long as the color is noticeable when compared to the A1 reference well.

Entering reactions

1. Enter reactions into the on-screen MicroPlate using any of the methods shown on page 5.6.
2. Click **OK**.




MANUAL ENTRY OF DATA WINDOW (TO ENTER PLATE REACTIONS)

3. Click **Print** to print report (if you did not select Automatic Printout).

4. Click **Save** to save data (if you did not select Auto Write to Data file.)
5. Repeat Steps 1-4 for all MicroPlates.

Reading multiple MicroPlates in reader mode

1. Select **Reader** on the **Set Up** window. SEE Section 5, page 7 for further information on using the Reader.



SET UP WINDOW (TO SELECT READER MODE AND INITIALIZE READER)

2. Click **Initialize Reader** to set up the reader. Make sure that red **No** changes to a green **Yes**.
3. Click the **Data** tab on the menu bar. The **Data** window appears.
4. Use the **Incubation Time** drop-down list to select the number of hours the MicroPlate has incubated.
5. Place the correct MicroPlate in the reader.
6. To work from an open worksheet:
 - Click **Select Read** to go to the selected worksheet.
 - Click on the sample number you want to read. An arrow will appear next to the sample you select.
 - Click **OK**. The reader calibrates itself, pulls the MicroPlate inside, takes a reading, displays the result on the **Data** window, and ejects the MicroPlate.
 - If necessary, click **Cancel** to return to the worksheet.
7. To read plates sequentially:

- Click **Read Next** to sequentially read the next MicroPlate on the worksheet without returning to the Worksheet window.
 - A **Confirmation** window appears. Click **OK** if it displays the MicroPlate you want to read. The reader calibrates itself, pulls the MicroPlate inside, takes a reading, displays the result on the **Data** window, and ejects the MicroPlate.
 - If necessary, click **Cancel** to return to the **Data** window.
8. Click **Print** to print the report (if you did not select Automatic Printout).
 9. Click **Save** to save data (if you did not select Auto Write to Data file.) A **Confirmation** window appears.
 10. Repeat Steps 5-9 for all MicroPlates.


Adjusting Thresholds Manually

The software allows you to view histograms of MicroPlate patterns. In most cases, use of this feature should not be necessary to obtain accurate results. However, this advanced technique can be useful if you feel that the automatic reading of your MicroPlate does not match your visual interpretation. You can adjust either or both of the thresholds to override the automatic reading and obtain what you feel is the correct result. You can only adjust thresholds in Reader mode and for only GN, GP and AN plates

Caution!

Manually adjusting thresholds will change SIM and DIST values. These changes override automatic thresholds. This is an advanced technique. Call Biolog Technical Service if you need help.

1. Click on the small **DW Data** tab at the upper right edge of your on-screen Microplate (on the **Data** window). You'll see an on-screen MicroPlate of your sample, showing optical densities.
2. Click the **DW Histogram** tab. A histogram of the current ID selection appears.
3. The histogram is divided into three lateral sections, which represent color distribution within the MicroPlate. Wells with low and high color density are on the left and right sides of the histogram. The middle section (between the two vertical lines) represents wells with borderline reactions.




HISTOGRAM (ON LEFT) AND OPTICAL DENSITIES (ON RIGHT)

- Click the **Manual** radio button to make adjustments.
- Use the horizontal scroll bars under the histogram to make adjustments to the thresholds. The top scroll bar adjusts the negative threshold (white line); the bottom scroll bar adjusts the positive threshold (red line). Green lines represent the automatic threshold levels.
- Continue until you feel the adjusted threshold corresponds to your visual read of the MicroPlate.
- Click **Apply** button to update thresholds.

Note: FF MicroPlate types will have Pos/Neg, Data, and Histogram tabs for both color and turbidity.


Choosing a Database to Search

The software allows you to select which database you want to search. If you have created User (custom) databases you can search the MicroLog databases, the User database or both. For information on creating a User database, see Section 7. MicroLog databases may not be selected as a User database.



1. Click the **Database** tab on the **Set Up** window. The **Database** window appears.

You'll see a Data Management Menu bar on the Database window. See Section 7 for data management instructions.



DATABASE WINDOW (TO SELECT DATABASE TO SEARCH)

2. Use the drop-down list to select the database you want. If you select a User or MicroLog/User database, the **User Progressive Database** and **User End-Point Database** (YT & FF) fields appear. Click on the appropriate field and a **Database File** dialog box appears.
3. Select the desired database by choosing the file name. For End-Point databases, the file name will include the hour (e.g., for a 24 hour read, select "demo24.eid").
4. Click **OK**.

Note: When you're conducting a MicroLog/User database search, a "U" will appear on the Species ID ranked list on the Data window to denote species in the User database. Also, you will see a ML/user designation in the upper left corner of the organism choices.

Exiting

When you're finished, click the **Exit** tab on the menu bar. Click **Click Here to Exit**.

Note: If you change your mind after clicking the **Exit** menu tab, you can go back to your work by clicking any tab on the menu bar.

6. Interpreting Results

In this section:

- ➔ Reading Printouts
- ➔ Assessing the Accuracy of your ID
- ➔ Pinning Down an Uncertain ID
- ➔ Species Comparison
- ➔ Searching and Assessing Species Information
- ➔ Using Macroscopic and Microscopic Pictures
- ➔ Using the Other Row
- ➔ Accessing Raw OD values

Check the bottom section of the **Data** window. The software identifies your microbe in the form of a ranked list. Entry #1 is the best match selected from the database. The other identifications are ranked in descending order of pattern match. Scroll down to see the whole list.

| SPECIES ID : CITROBACTER FARMERI | | | | | |
|----------------------------------|--------------------------|------|-------|------|--------|
| | NAME | PROB | SIM | DIST | TYPE |
| => 1 | CITROBACTER FARMERI | 100% | 1.000 | 0.00 | GN-ENT |
| 2 | CITROBACTER AMALONATICUS | 0% | 0.000 | 5.26 | GN-ENT |
| 3 | CITROBACTER FREUNDII | 0% | 0.000 | 6.54 | GN-ENT |
| Other | - | - | - | - | - |

CLOSE-UP OF RESULTS SECTION OF DATA WINDOW

The ID box

Final identification results are displayed in the ID box at the top of the results area of the **Data** window. The software considers the possibility that even the #1 ranked species may not be a good enough match. In this case, you will see “No ID” in the ID box.

Reading Printouts

Every time you read a MicroPlate, the printer will issue a hard copy report. The information on these reports varies, depending on whether you are working in manual or reader mode. See Appendices for plate data statistics (A-1), and sample printouts and explanation of entries (A-5).

Assessing the Accuracy of Your ID

If the answer to the following three questions is yes, you can feel confident that ID #1 is accurate.

- Are the top-ranked ID choices on the list all the same (or closely related) genera?
- Check the SIM (similarity index value) rating of ID #1. Is it greater than .5?
- Check the DIST (distance) rating of ID #1. Is it greater than two distance points away from the 2nd choice.

Pinning Down an Uncertain ID

If the top-rated choices show any of the following results, you may want to do some additional probing:

*Understanding
the language of
microbe
identification:*

% PROB –
allows you to
compare our IDs
to other systems
that use this type
of calculation
SIM – tells you
how good each
match is
DIST – tells you
the number of
mismatches
between your
MicroPlate
results and the
database pattern
for that species

- The top-rated ID choices on the list are random and unrelated genera
- For GN, GP, & AN database (24 hour results):
 - SIM is very near 0.5
 - DIST is greater than 5.0
 - DIST of the first and second choices are nearly equal (less than two distance points apart)
 - There are many variable reactions (more than 15)
- For YT database (≥ 48 hour results):
 - SIM is near 0.5
 - DIST is greater than 5.0
 - DIST of the first and second choices are nearly equal
 - There are many variable reactions (more than 15)
- For FF database (all time points):
 - SIM is near cutoff values for each day (1-4)
 - DIST is greater than 5.0
 - DIST of the first, second and third choices are nearly equal
 - There are many variable reactions (more than 30)

Note: For FF MicroPlate identifications, verification of the culture morphology is necessary. Refer to the software database photo library or mycology textbooks as references.

For additional assistance in interpreting results, contact Biolog Technical Service.

ID Rules of Thumb:

SIM = 1.0 ➔ a perfect match

SIM = 0.0 ➔ no match

The closer to 1.0, the better the match.

1. The chart below will help you clarify your identification by using SIM and DIST values:

USING SIM AND DIST TO ASSESS IDENTIFICATIONS

Assess SIM for ID #1

| | |
|-------------------------------------|--|
| Gram-negative aerobes | ➔ must be ≥ 0.5 at 16-24 hours incubation |
| Gram-positive aerobes | ➔ must be ≥ 0.75 at 4-6 hours incubation |
| Anaerobic microbes | ➔ must be ≥ 0.5 at 20-24 hours incubation |
| Yeasts | ➔ must be ≥ 0.75 at 24 hours incubation ➔ must be ≥ 0.5 at 48 or 72 hours incubation |
| Filamentous Fungi (& select yeasts) | ➔ must be ≥ 0.90 at 24 hours incubation ➔ must be ≥ 0.70 at 48 hours incubation ➔ must be ≥ 0.65 at 72 hours incubation ➔ must be ≥ 0.60 at 96 hours incubation |

Assess SIM for top IDs

If SIMs for top IDs are $<$ the SIM required **AND** all belong to the same genus **AND** if their total is \geq the SIM required, you can be confident that the microbe is in that genus. MicroLog will give you a genus-level ID rather than a species-level ID.

Assess DIST for top IDs

| | |
|--|--|
| If SIM is OK, but DIST between ID#1 and ID#2 are very close | ➔ Good match to both species (some are quite close) |
| If you remain unconfident of the ID, try to differentiate organisms by | ➔ Gram stain, colony morphology, oxidase and/or catalase tests, additional tests (see <i>Bergey's Manual</i>) |

What Should I Do About "No Identification"?

Check the on-screen + and – signs for that MicroPlate. You can use these mismatch indicators to assess what happened.

2. If you get a "No Identification" result, the % Prob value will not be displayed. When this happens, check for black + and – signs in the on-screen MicroPlate.

- If there are both + and – random mismatches, this is probably a true mismatch. If so, this is probably a species not in Biolog's database. You can add it to your own user-created database. See Section 8.
- If the mismatches are all + or all –, you may have made a testing error. Refer to Table 6-1, and see Section 10.

Table 6-1 will assist you in figuring out the cause of these mismatches:

| Mismatch Type | What It Might Mean |
|---|--|
| All + (your pattern is giving fewer positive reactions than the species you're comparing it to) | Under-inoculation |
| | A1 well is overfilled, contains clumps, or is cloudy |
| | Cells were mishandled, too old, cultured on wrong medium, suspended in wrong inoculating fluid, incubated at wrong temperature, etc. |
| All – (your pattern is giving more positive reactions than the species you're comparing it to) | Over-inoculation (especially with enterics) |
| | A1 well is underfilled |
| | Contamination (mixed culture) |

TABLE 6-1. ASSESSING MISMATCHES

FF database

Perform a visual macroscopic and microscopic assessment of your isolate compared to the first three ID ranked species. Refer to MicroLog software photo library or mycology textbooks as references.

Note: The A1 well is not used as a negative control in the FF database, therefore this well may not appear colorless.

FF ID Workflow

1. Read daily (at 24, 48, 72 hours) until ID is called.
2. After an ID is called, verify the first three choice macro/microscopically
3. If no ID, re-incubate the MicroPlate
4. Read at 96 hours. If an ID is called, verify the first three choices macro/microscopically
5. If no ID, launch additional investigation. Check first three choices macro/microscopically.

*When you're
interpreting
results:*

- 1. Trust your
eye over the
reader.*
- 2. When in
doubt,
repeat the
test, using
stringent
technique.*

GN, GP, and AN databases

1. Perform a visual assessment by checking the A1 well.
 - If it contains a significant purple color, this color may be obscuring the pattern. Re-examine your general setup procedures.
 - The strain may need to be retested with thioglycolate added to the inoculating fluid.
 - With some rare gram positive species, salicylate must be added (in addition to thioglycolate) to eliminate the false positive purple color.
2. Check for borderline reactions. The system essentially disregards these when making identifications.
 - Review and verify your testing procedures. Did you use all optimal conditions?
 - If there are more than 15 borderline reactions, visually check the MicroPlate. If you used the reader, compare how it read the pattern in relation to how you read it by eye. If you disagree with the reader, click the **Set Up** tab on the menu bar. Select **Manual** mode and correct reactions manually. Or you can use the histogram function to adjust the thresholds (see Section 5). Re-identify based on the modified pattern or incubate the MicroPlate for a longer time (do not exceed time limits).
3. Check for a message in the ID box. When appropriate, the ID box will offer additional information to guide your identification. For example:
 - ID of *Salmonella*, *Shigella*, or *E. coli* O157, ID box message says "Confirm by Serology"
 - ID of *Listeria* species or *Neisseria gonorrhoeae*, ID box message says "Confirmation of ID is recommended."

Species Comparison


MicroLog software has several features that can help with species identification. You can compare any pattern with the expected pattern of any species or taxa in the database. You can also display expected species patterns visually.

To compare your #1 ID with others on the list of runners-up:

*What are those +
and – signs?*

*These are signs
showing
mismatches
between your
sample and the
database record.
+ ≥ 80% pos
– ≤ 20% pos*

1. After a pattern is entered, the on-screen MicroPlate automatically displays the major mismatches between the #1 ranked species and your MicroPlate.
2. To compare your ID with other species on the list, click on the species you want to view. A small arrow appears next to the selected organism and the major mismatches are displayed.



PLUS (+) AND MINUS (–) MISMATCHES

3. To compare against a species that is not in the 10 top-ranking, click on the **Other** row. See page 13

Searching and Assessing Species Information


MicroLog software gives you an opportunity to find out more about the range of reaction patterns of species in the MicroLog database (or any record in a user database). There are three possible ways to access the species information feature, all from the **Data** window:

- Using the View Database button on the selection bar.
- Using the ID box
- Using the “Other” row


Using the View Database Button

You can get advanced species information without running identification, for single or multiple records in any database.

1. Click the **View Database** selection bar on the **Data** window. The **Database List** appears. This is a species list from the database, based on the strain type you last selected on the Data window.
2. Click the **Database Options** tab to select the database you want to view, based on strain type.
3. The **Database List** contains 15 records per page. Click **Page Up** or **Page Down** or use the scroll bar to move through the list.
4. Click on the row containing the species you want.
5. This will take you to the **End-Point** and **Progressive** data pages.



DATABASE OPTIONS WINDOW



DATABASE LIST WINDOW

Endpoint data analysis

For each well, the number displayed represents the percentage of strains of the displayed species that were positive at the selected incubation time.

1. In the ID box in the **Data** window, right-click on the species you want to investigate. The **End-Point** window in the **Database Data** window appears.

NOTE: A1 value may be greater than 0 if the organism tends to give false positive reactions.

Database Data

Achromobacter xylosoxidans ss xylosoxidans
GN-NENT OXI+

End-Point Progressive

Percent of Times Wells Are Positive At End-Point

| Color | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|-------|-----|-----|----|----|----|-----|----|----|----|-----|----|-----|
| A | 2 | 0 | 6 | 18 | 26 | 40 | 0 | 0 | 0 | 25 | 0 | 3 |
| B | 0 | 1 | 0 | 0 | 0 | 73 | 0 | 0 | 0 | 0 | 0 | 7 |
| C | 0 | 1 | 2 | 0 | 0 | 0 | 0 | 0 | 3 | 0 | 99 | 100 |
| D | 67 | 91 | 85 | 95 | 59 | 0 | 55 | 0 | 0 | 98 | 98 | 9 |
| E | 31 | 80 | 76 | 74 | 69 | 100 | 0 | 86 | 0 | 39 | 24 | 100 |
| F | 100 | 100 | 0 | 53 | 59 | 88 | 83 | 80 | 94 | 100 | 66 | 97 |
| G | 25 | 0 | 92 | 6 | 92 | 95 | 42 | 9 | 86 | 50 | 0 | 3 |
| H | 81 | 0 | 0 | 0 | 44 | 0 | 0 | 0 | 0 | 5 | 0 | 0 |

Print Data Close

DATABASE DATA WINDOWS (END-POINT WINDOW)

YT & FF
MicroPlates are limited to endpoint data analysis only

Database Data

Aspergillus fumigatus Fresen.
ASP

End-Point Macro Pictures Micro Pictures

Percent Of Times Wells Are Positive After 3 Days

| Color | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 100 | 0 | 100 | 0 | 19 | 84 | 6 | 94 | 34 | 100 | 91 |
| B | 0 | 56 | 100 | 91 | 81 | 0 | 25 | 100 | 88 | 100 | 0 | 91 |
| C | 100 | 0 | 100 | 100 | 100 | 88 | 100 | 12 | 6 | 69 | 88 | 94 |
| D | 100 | 100 | 72 | 47 | 6 | 19 | 0 | 81 | 91 | 0 | 91 | 75 |
| E | 84 | 100 | 9 | 100 | 100 | 100 | 100 | 9 | 91 | 81 | 75 | 94 |
| F | 94 | 6 | 56 | 31 | 28 | 100 | 88 | 91 | 59 | 69 | 72 | 100 |
| G | 91 | 0 | 63 | 91 | 88 | 19 | 91 | 100 | 100 | 100 | 100 | 100 |
| H | 59 | 100 | 91 | 100 | 100 | 100 | 100 | 100 | 91 | 50 | 88 | 88 |

| Turbidity | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|-----------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 100 | 0 | 100 | 0 | 38 | 100 | 28 | 100 | 66 | 100 | 100 |
| B | 56 | 100 | 100 | 88 | 100 | 9 | 75 | 100 | 100 | 100 | 16 | 100 |
| C | 72 | 0 | 100 | 100 | 100 | 94 | 100 | 34 | 28 | 100 | 100 | 100 |
| D | 100 | 100 | 94 | 94 | 0 | 50 | 0 | 100 | 100 | 19 | 100 | 94 |
| E | 100 | 100 | 6 | 100 | 100 | 100 | 100 | 12 | 100 | 100 | 81 | 100 |
| F | 100 | 12 | 66 | 59 | 16 | 91 | 72 | 56 | 66 | 66 | 69 | 100 |
| G | 81 | 81 | 69 | 75 | 100 | 25 | 84 | 100 | 100 | 100 | 100 | 100 |
| H | 53 | 100 | 88 | 100 | 100 | 100 | 88 | 84 | 81 | 56 | 22 | 22 |

Print Data

- The numbers appearing in the on-screen MicroPlate show (for each well) the average “percent positive” value, considering all strains of the species that were tested. This is a simple way to view species data.
- Click the **Progressive** tab to view more sophisticated bacterial species data, namely, the range of reactions you can get for this species as the pattern develops in a progressive fashion.

Progressive data
analysis applies to bacterial identifications only

On all progressive data windows:

Click on a specific well to select it.

Find well information in "Key" box.

Click "Set to Current Data Value" to return all settings to number of positive reactions your current MicroPlate has.


Progressive data analysis

MicroLog software uses a new method of pattern matching called Progressive ID (PID). The PID method takes into account the progressive nature by which purple color forms in the wells of the MicroPlate. What typically happens after cells are inoculated is that some wells turn purple quickly and strongly, resulting in dark purple wells. These are usually the preferred carbon source for the microbe. After additional incubation, other wells also turn purple, sometimes forming weaker purple wells.

For example, if a MicroPlate has 23 purple wells, PID will compare the pattern to the other species in the database *at their stage of development when they also had 23 purple wells*. Each species in the database responds to the search, as if to say "This is my pattern when I had 23 positive reactions." PID finds the best match among the species in the database. The mathematics of matching the species to the pattern is otherwise identical to that used in the previous method of End-point ID (EID).

Reaction range

This display shows the entire metabolic progression of the species in a single image. The boxes are entirely purple when the well is positive 100% of the time. A teal color indicates that the well is positive from 20% to 80% of the time. Empty boxes indicate that the substrate is used < 20% of the time. Each box contains a black line that shows how many positive reactions are on the given MicroPlate. This black line is a gauge of progressive identification. The top right corner of the window contains a key showing the current reaction of a specific well and how often it is positive.



REACTION RANGE WINDOW


Slide the scroll bar left and right to represent various levels of positive reactions. As the vertical line moves across each well of the on-screen MicroPlate, it indicates where in the histogram a chosen positive well value falls.

- The purple area under the graph (histograms) represents when a well is positive for that species.
- A box that is completely purple is always positive for that species. This is a preferred carbon source for that species.
- A box that has limited levels of purple is only positive when there are relatively high numbers of positive wells on the MicroPlate.
- Teal represents borderline well reactions.
- You can click on any well and it will appear as the example shown on the key at the top.

Specific pattern

This display shows the numerical pattern for the species at a specific level of positive reactions. This is a more explicit display of the **Reaction Range** display. For the number of positive reactions (indicated below the display and adjusted by using the scroll bar), this is the percentage of times that a specific well is positive.

- A well that is 100 when the scroll bar is far to the left is a well that is always positive for that species, even when it has few positive reactions.
- A well that is only 100 when the scroll bar is to the far right is a well that is only positive when the species has many positive reactions.
- A well that contains a low number or zero is a well that is rarely or never positive, even when the species has expressed its full metabolic potential.




SPECIFIC PATTERN WINDOW

Vs. current MicroPlate

This display shows a graphical representation of the pattern for the species at a specific level of positive reactions versus the MicroPlate currently on the **Data** window. The circle at the left of each box represents the reactions in the current MicroPlate. The bar at the right of each box represents the percent of times the species is positive at the given level of positive reactions.

- Slide the scroll bar left and right to see the animation of the metabolic progression for that species.
- Purple indicates positive wells.
- Teal represents borderline well reactions.
- Empty boxes are always negative.
- The display shows major mismatches between the current MicroPlate and the species, at the current level of positive reactions.
- If you set the number of positive reactions to be equal to the number of positive reactions in the current MicroPlate, you can see how the ID search treats the match between your MicroPlate and the displayed species.



VS. CURRENT MICROPLATE WINDOW


Using Macroscopic and Microscopic Pictures

You can access macroscopic and microscopic species of selected fungi in the MicroLog photo library of the FF database.


1. From the **Data** window, click the **View Database** selection bar. The **Database List** appears. This is a species list from the database based on strain type last selected on the Data screen.
2. Click the **Database Options** tab and select the database you want to view, based on strain type.
3. Click back to the **Database List**. Click **Page Up** or **Page Down** or use the scroll bar to move through the list.
4. Click on the row containing the species you want.
5. The **End-Point**, **Micro Pictures**, and **Macro Pictures** tabs appear.
6. Click on the **Macro** or **Micro** tab to access the photo library.



MACRO AND MICRO TABS ON DATABASE DATA WINDOW



MACROSCOPIC PICTURE



MICROSCOPIC PICTURE

Using the Other Row


You can use the **Other** row at the bottom of the **Data** window to get additional species information.

1. Double-click anywhere on the **Other** row at the bottom of the ID box (watch for the small arrow).

2. The **Database List** appears. Click on the row containing the species you want to highlight it.
3. Click OK to return to the **Data** window. Note that your selection now appears in the **Other** row. In addition, positive (+) and negative (–) mismatches appear in the on-screen MicroPlate, as does the DIST value in the ID box. This performs a comparison between the database record and your MicroPlate ID.

Note: If you are searching both the MicroLog and User databases, the program will request that you designate which database you wish to compare to before the Database List appears.

4. Click on the **Other** row and press Enter or Backspace on your keyboard to clear the entry.



DATA WINDOW WITH “OTHER” FIELD FILLED (NOTE DIST AND +/- ENTRIES)

Accessing Raw OD Values

This function allows you to access the raw OD values for the wavelengths read.

1. Click the **Raw Data** tab on the **Data** window. This tab is located below the reaction section of the **Data** window. The raw OD values for the wavelengths read will be displayed.
2. To return to the previous data view, click the **Data** tab.

Accessing Data Values

This function allows you to access the Data values for the wavelengths read. This is useful to view both sets of data values for the FF results on one screen.

1. Click the **Data Stats** tab on the **Data** window. This tab is located below the reaction section of the **Data** window. The Data values for the wavelengths read will be displayed.
2. To return to the previous data view, click the **Data** tab.


7. Data Management and Compiling Databases

In this section:

- ➔ BackUp a Data File
- ➔ Combining Data Files
- ➔ Viewing and Editing Files
- ➔ Compiling Data Files
- ➔ Exporting ASCII Data
- ➔ Cluster Analysis

The heart of the software is a library of over 1,900 microbes, distributed over a number of organism databases. You can replace or augment this library with other species or strains. The user can compile an ID File (customized database) from a Data File and perform identifications using either our database or your own customized database. The software also provides functionality for managing Data Files, which is also useful in preparing a Data File for compilation. These functions include backup (copying), combining, editing, printing, exporting to other applications (such as Microsoft Excel® or Microsoft Access™).

The software also provides dendrograms and Multi-Dimensional Scaling diagrams, useful for cluster analysis.



DATA FILE DATA MANAGEMENT MENU PAGE

Data File Management under Restricted Access mode

If the software is running under Unrestricted Access mode, then any user may access the Data File Management page. Under Restricted Access mode, only users with at least the View/Print privileges may have access. If a user with the Log-In and Set Up privileges only, then a warning message box pops up stating that the user: “Must have the View/Print privilege to navigate to the Data File Management page”.

Users with the View/Print privilege may access the following features:

- Backup Data File

- View and Print functions (from the **View/Edit Data File** button)
- Export to ASCII File

Users with Edit Data File privileges have access to advanced features, including:


- Add File to Data
- Editing functions (from the **View/Edit Data File** button)
- Compile/View Data File

Backup Data File (Copying a Data File)

Before altering a Data File in any way, it is good practice to make a back-up copy of the Data File. Frozen files cannot be edited; therefore, a backup copy must be created on which to perform any editing. The backup copy of the Data file keeps track of the file from which a record was copied (parent file). This creates an electronic trail of the record.

*If desired,
view the
destination
/edited file
using the
View/Edit
feature to
confirm
data file
modifications.*

1. Click the **Data File** tab on the menu bar to view the **Data Management** page. Click on the **Data Management File** field and a **Data File** dialog box appears which lists available Data files. To select the file, click on the name of the file.
2. Click **Open**.
3. Click the **Backup Data File** button to view the BackUp Data File save dialog.
4. Enter the new desired file name (backup copy) into the **File name:** field. You can stay in the same directory as the original or navigate to another directory to save the file.




BACKUP DATA FILE SAVE DIALOG

5. Click **Save**.

Combining Data Files

1. Click on the **Data Management File** field and a **Data File** dialog box appears which lists available Data files. Click on the name of a non-Frozen Data File (non-original data file or backup data file) to select that file as the destination file.
2. Click **Open**.
3. Click the **Add File to Data File** button to view the **Add Single Files** dialog. *Note: Destination file is listed in the Old File section.*




4. Click in the **Source File** section to view the **Find New Data File** open dialog.
5. Click on the desired data file name (file to add).
6. Click **Open**. *Note: Last selected file is listed in the New File section.*
7. Click the **Add Source to Target** button. The file you've added will be listed in the 'Files Added This Session' section. The file listed in the **Target File** box now also contains a copy of the file listed in the **Source File** box. *(The Source File box still contains the name of the added file until another file is selected.)*
8. Follow steps 4 thru 7 to add additional files.

***Note:** As a safety feature(to avoid spurious duplication of records), the software will not allow you to add a file more than once in one session. However, if the user closes and re-opens this dialog, then the software does not prevent the user from adding a file that may have already been added.*

9. To freeze this file after you have added to it, click **Freeze Data Records** at the bottom of the **Data Management** window.


Viewing and Editing Data Files

1. Click on the **Data Management File** field and a **Data File** dialog box appears which lists available Data files. Click on the name of a Data File to select that file (non-original data file or backup data file).
2. Click **Open**.
3. Click the **View/Edit Data File** button to view the **Data File Edit** dialog. *Note: Dialog contains a spreadsheet, in which each record corresponds to a read of a plate. The dialog is initially in the View Only mode.*



VIEW/EDIT DATA FILE DIALOG

4. You will be in **View Only** mode. You can examine the records, but not edit them even if the file is not frozen. *If desired, navigate to other columns or rows by using the horizontal or vertical scroll bar or adjust the column widths by carefully positioning the cursor between column headings until the cursor icon changes and, subsequently, dragging the cursor.*
5. Select the **Print** tab to access the print options.




CLOSE-UP OF PRINT OPTIONS IN DATA FILE EDIT DIALOG

OK, Out, Hold, Atypical
Sorting labels may be applied to individual identification records. This sorting facilitates the task of organizing records within in a datafile for the compiling user databases.

- Click desired radio button to select ‘**Text**’ or ‘**Graphics**’ as the **Print Mode**.
- Click desired radio button to select a **Print Selection** option: ‘**Print Record**’ allows printing of the selected record, which is the record with an arrow icon in the row heading. ‘**Print Page**’ allows printing of all the records that appear on the spreadsheet. ‘**Print All**’ allows printing of all the records in the Data File.
- Click the **Print Fields** button to print the ‘information’ fields that are in the spreadsheet.
- Click the **Print Results** button to print each identification report printout.

6. Select the **Edit In/Out** tab to label a record (used in the Data File compilation feature, Sec. 7 page 7). *Users without the Data File Edit privilege are prevented from accessing this page.*




CLOSE-UP OF DATA FILE COMPILATION OPTIONS IN DATA FILE EDIT DIALOG

- Click desired radio button to select a ‘**Prompt Before Edit**’ option. Note that the ‘**Prompt**’ option will prompt the user for an OK button click before applying any edit.
- Click desired radio button to select an **Edit Mode** option:


| | |
|-----------------|---|
| OK | (to include a record in compilation) |
| Out | (to exclude a record from compilation, while indicating that it will never become eligible in the future) |
| Hold | (to exclude a record, while indicating that it may become eligible in the future) |
| Atypical | (to exclude a record, while indicating that it may become eligible on a contingency basis) |
- Once you have selected to **Edit Mode**, click on the row to label the record in the spreadsheet. Note that the same edit may rapidly be applied to many records by clicking many records in succession.

7. Select the **Edit Fields** tab to access options that allow rapid editing of many records. *Note that users without the Data File Edit privilege are prevented from accessing this page.*



CLOSE-UP OF EDIT FIELDS OPTIONS IN DATA FILE EDIT DIALOG

- Click radio buttons to select a field for editing from among the **Data Field** options.
 - Enter the value of the selected field into the Edit box.
 - Once you have selected the **Data Field** and entered the value (new), click on the row on the spreadsheet to edit (replace value). *Note that the same edit may rapidly be applied to many records by clicking many records in succession.*
8. Select the **Edit Values** tab to access options that allow extensive editing of a single record. *Note that users without the Data File Edit privilege are prevented from accessing this page.*
 - Click on a row in the spreadsheet to view a **Plate Information** dialog for the record corresponding to the selected row.
 - Enter a value into or change a value in each Edit box, as desired. Select a different value in each Pull Down box, as desired.
 - Click Save, Done or Cancel.
 9. Select the **Copy** tab to access options that allow copying of selected records into another Data File. *Note that users without the Data File Edit privilege are prevented from accessing this page.*



CLOSE-UP OF COPY OPTIONS IN DATA FILE EDIT DIALOG

- Click desired radio button to select a **Mark Record** option:

Remember!

When building your own databases, strain names must be consistent when you group patterns. Using inconsistent strain names will result in incorrect database groupings.

Mark (to select a record for copying)
Clear (to de-select a record for copying)

- Once you have selected the **Mark Record** option, click on the row to label the record in the spreadsheet. Note that the same edit may rapidly be applied to many records by clicking many records in succession.
- To mark or clear a large numbers of records simultaneously, click one of the following-buttons:

Mark All (to select all the records in the Data File)

Mark OK's (to select for copying only those records that are also 'OK' for compilation)

Clear All (to de-select all records)

- Click the **Write Marked Recs to New File** button to view the **Save Marked Data to this File** save dialog.
- Enter the desired file name of the copy into the **File name:** field. Click the **Save** button.

10. Click **Freeze Data Records** if you've finished editing the file and want to save it as a "**Read Only**" file. This will create a file that is non-editable until another backup copy is made (Backup Section 7, page 2).

Compiling Data Files

The software also allows you to compile data files and produce databases.

Producing databases

Be sure to save GN, GP, AN, YT, and FF MicroPlate-type readings in separate files for database compilation. Certain features of the compilation menu are specifically linked to MicroPlate type.

Choose the existing file to compile by clicking on the **Data Management File** field on the **Data Management** window. A **Data File** dialog box appears.


1. Click on the desired file name.
2. Click **Open**.
3. Click the **Compile/View Data File** selection bar on the **Data Management** window. The **Compile/Index/Data** tabs appear, defaulting to the **Compile** window.

4. In the **Data To Include** section, enable check marks in the desired box(es). (See pages 7.4-7.6 for the Edit In/Out process).

| | |
|-----------------|--|
| OK | (to include records designated OK) |
| Hold | (to include records designated Hold) |
| Atypical | (to include records designated Atypical) |
| Out | (to include records designated Out) |


6. In the **What Processing To Do** section:

- Enable the check mark next to **Make Index Files** (this must be enabled to use the rest of this feature below).
- Enable **Make ID Files** to make your own database.
- If you are compiling files with GN, GP, or AN plate types, the window below opens.




COMPILE/INDEX/DATA WINDOWS (FOR GN/GP/AN PLATE TYPES)

- If the file you are compiling contains YT or FF MicroPlate data, you must specify the incubation time for the data you are compiling. If the file contains several incubation times, select **All**. The software will compile the entries by strain name and automatically make an endpoint database for each incubation time. For further information, follow the on-screen instructions.




COMPILE/INDEX/DATA WINDOWS (FOR YT/FF PLATE TYPES)

- Click **Go** to compile database.
 - Indication bars at the bottom of the window show compilation progress.
 - The Progressive database file will have the same name as your data file, with a ".KID" extension (GN,GP, and AN files) and the Endpoint database file, with a "EID" extension (YT and FF files). The user databases are saved as read only files.
7. Click the **Index** tab on the **Compile/View Data File** window to show the **Index** window.
- You will be able to view a data file list by **Sample Number**, **Strain Name**, or **Strain Number**. Click the desired radio button and the list will change accordingly.
 - Click **Page Up** or **Page Down** or use the scroll bar to view the entire list.
 - To view an entry more closely, click on a row, then click the **Data** tab at the top of the page.




INDEX WINDOW OF COMPILE/VIEW DATABASE (WITH SAMPLE ENTRIES)

- The data spreadsheet will show you all information for that sample. You can perform all functions (such as editing and printing) from the **View/Edit Data File** menu.
- If you edit your data file after compilation, you must recompile the data before using the database to identify organisms.
- Click **OK** to close.




DATA WINDOW OF COMPILE/VIEW DATABASE

- Click on the **Setup** tab to view the User Database you created.
 - Click on the **Database** tab.
 - Click on the **Database to Search** pull-down list.



DATABASE WINDOW

- Select **User**.
- If the database contains GN/GP/AN plate types, click on the **User Progressive Database** field.
- If the database contains YT/FF plate types, click on the **User End-point Database** field.
- A **Progressive** or **Endpoint Database File** dialog box appears.



DATABASE FILE DIALOG BOX

- Select the user database you want.
- Click **Save**.
- Freeze and/or comment on the database file when you finish editing.
- Click on the **Data** tab.


- Select the appropriate Plate Type, Strain Type, and Incubation time for the User Database.
- Click on the **View Database** bar under the picture of the MicroPlate.
- A list of organisms in your User Database appears similar to a MicroStation/MicroLog database list.

User database development
To evaluate various database versions during development, the Re-ID function allows you to perform ID test evaluations on test backup data files.

!Note: Use of the compile tab re-ID feature is not 21CFR part 11 compliant.


Compile tab Re-ID feature

1. Choose the existing backup file for User database development and analysis to Re-ID by clicking on the **Data Management File** field on the **Data Management** window. A **Data File** dialog box appears.
2. Click on the desired file name.
3. Click **Open**.
4. Click **Compile/View Data File**.
5. In the **What Processing To Do** section:
 - Enable the check mark next to **Make Index File** (this must be enabled to use the rest of this feature).
 - Enable **Re-ID Data** identify your data against a database you have created or re-threshold data collected in an earlier software release to the thresholds in the current release.
 - Click on the **Re-ID Database** field to select the appropriate database. The software will re-identify the data file based on this entry.
 - If the file to **Re-ID** contains GN, GP or AN MicroPlate types, the appropriate progressive database is selected.



SELECT RE-ID DATABASE WINDOW

- If the file to Re-ID contains YT or FF plate data the with only one incubation time, choose the corresponding endpoint database file for that time in the **Re-ID Database** field. For example, if the FF file is from 24-hour incubation time in Step 1, choose 1 Day. For the **Re-ID Database**, select FF60124.eid (24 Hour).



SELECT RE-ID DATABASE WINDOW

- If the YT/FF file contains several incubation times, select **All** at Step 1. The software will compile and Re-ID the entries by strain name and incubation time. You do not need to select a database for the **Re-ID Database** field.
- Click **Go** to Re-ID the file.
- Indication bars at the bottom of the window shows the Re-ID progress.
- Click **OK** when the process is finished.

Printing data headers, species counts and dendrograms

These features are automatically printed once they are compiled. The software does not save a permanent record.


1. Choose the existing file to compile by clicking on the **Data Management File** field on the **Data Management** window. A **Data File** dialog box appears.
2. Click on the desired file name.
3. Click **Open**.
4. In the **What Processing To Do** section:
 - Enable the check mark next to **Make Index File** (this must be enabled to use the rest of this feature).

- Enable **Print Ranked Data Headers** if you want to print a list of all the information in the data file as entered in a spreadsheet row, ranked by MicroPlate type and species (see Appendices for sample printout).
- Enable **Print Species Counts** if you want to print by strain type, one row per species, with species name and number of MicroPlates for that species (this will give you a species count) (see Appendices for sample printout).
- Enable **Print Species/Plate Dendrograms** if you want to generate a dendrogram printout, with each dendrogram entry equaling one MicroPlate. The same species must have at least three entries in the file to make a dendrogram. Each dendrogram is determined by the strain name (see Appendix 5 for sample printouts). The dendrograms display linear relationships.

Entering Reactions from a Saved File for Re-ID


*Re-ID using
input mode
features
You must use
this method to
re-ID previously
run data against
your user
database and
save as a new
file to be 21
CFR part 11
compliant.*

1. Use the **Data Input Mode** drop-down list to select **File** (to re-ID saved data records from a saved file)
2. The **Input Data File Name** field appears. Click in the empty field and move through the listed folders to find the file you want to read and click on it. This will bring the file name up on the **Set up** screen.



SET UP WINDOW (WITH FILE MODE SELECTED)

3. In the **Write Data to Data File** drop-down list
 - a. Select **Yes** to save the generated data to another file.
 - b. Click on the **Output Data File Name** field. A **Save As** dialog box appears. Enter the desired file name (D4C extension).
4. Click the **Data** tab.
5. Click **Select Read**
6. A list of saved data appears. Select the MicroPlate you want to read and click **OK**.



| # | OK | Created | Date | Sample | V | Plate | Strain Type | Strain Name | Strain # |
|----|----|-------------------|-------|--------|--------|---------------------|-------------|-------------|----------|
| 16 | OK | Apr 16 1998 11#17 | 24 HR | ON2 | ON-ENT | ESC. COL | | 5574 | |
| 17 | OK | Apr 15 1998 11#7 | 4 HR | ON2 | ON-ENT | OIT BRA | | 13651 | |
| 18 | OK | Apr 15 1998 11#8 | 4 HR | ON2 | ON-ENT | ENT CLO | | 11534 | |
| 19 | OK | Apr 15 1998 11#4 | 4 HR | ON2 | ON-ENT | OIT MUR | | 13656 | |
| 20 | OK | Apr 17 1998 01#3 | 24 HR | ON2 | ON-ENT | ENT ADO RO 6 (ERRV) | | 11520 | |
| 21 | OK | Apr 17 1998 01#2 | 24 HR | ON2 | ON-ENT | ENT CLO | | 11532 | |
| 22 | OK | Apr 17 1998 01#22 | 24 HR | ON2 | ON-ENT | KLB.FNE | | 4145 | |
| 23 | OK | May 14 1998 11#52 | 4 HR | ON2 | ON-FAS | OXI.HAE.ACT | | 958 | |
| 24 | OK | May 14 1998 11#52 | 4 HR | ON2 | ON-FAS | OXI.CAP.OCHSRU | | 5618 | |
| 25 | OK | May 15 1998 01#70 | 24 HR | ON2 | ON-FAS | OXI.HAE.INF | | 5641 | |
| 26 | OK | May 15 1998 01#60 | 24 HR | ON2 | ON-FAS | OXI.HAE.APH | | 2704 | |
| 27 | OK | May 15 1998 11#53 | 4 HR | ON2 | ON-FAS | OXI.HAE.ACT | | 8847 | |
| 28 | OK | May 15 1998 11#15 | 4 HR | ON2 | ON-FAS | OXI.HAE.HAE | | 8755 | |
| 29 | OK | May 15 1998 11#35 | 4 HR | ON2 | ON-FAS | OXI.HAE.INF | | 813 | |
| 30 | OK | May 17 1998 11#15 | 24 HR | ON2 | ON-FAS | OXI.CAP.OCHSRU | | 7228 | |

LIST OF SAVED DATA

7. The Data page will show the MicroLog ID record.
 - a. Click the **Setup** tab.
 - b. Click the **Database** tab.
 - i. Select **User**.
 - ii. Enter database file name, as required for User database search options.
8. Click the **Data** tab.
 - a. The re-ID based on the MicroLog database is in view.
 - b. Click the Setup tab.
 - c. Select Manual input mode.
 - d. Click the Data tab.
 - e. Click on Re-read button
 - f. Manual entry display appears with the reactions populated.
 - g. Click OK to re-ID.
 - h. The Data screen is displayed with the re-ID result.
 - i. Click **Save** to save the ID generated for the file record entered to your Output file.


9. Click the **Setup** tab.
10. Use the **Data Input Mode** drop-down list to select **File**
11. Click the **Data** tab.
12. Click **Read Next** (to select the next record) or click **Select Read** (to choose another record) and click **OK** after the selection. *The program defaults to the MicroLog database in file input mode.*
13. Repeat steps 7 thru 12 until you have re-ID'd all of the records that you wish. Click **Save** to save the ID generated for the file record entered to your Output file.

Exporting ASCII Data

The software enables you to export data into an ASCII file, and export data to Excel and/or Microsoft Word.


Exporting data into an ASCII file

1. Choose an existing file to export to by clicking on the **Data Management File** field on the **Data Management** window. A **Data File** dialog box appears.
2. Click on the desired file name.
3. Click **Open**.
4. Click the **Export to ASCII File** selection bar on the **Data Management** window. The **Export to ASCII File** window appears. The file you chose appears in the **Data File To Give Data** field.




ASCII EXPORT WINDOW

5. Click on the **ASCII File to Get Data** field. An **ASCII Export File** dialog box appears.



ASCII EXPORT FILE DIALOG BOX

6. Use the **Save as type** pull-down list to select the format of the ASCII file you want to import.
 - **Text File** shows a list with a .prn or .txt extension.
 - **Comma Delimited** shows a list with a .csv extension.
 - **All Files** shows a list of all files
7. Type in the desired name for the ASCII file where you want to export the data or click on a previously created ASCII file to append to it. Be sure to use the correct file extension.
8. Click **Save**. The **ASCII Files** window reappears.
9. **Choose Type of Format**, click the radio button designating the format you desire for your ASCII file.
 - One-Line, Comma Delimited (see Section 7, page 17)
 - Line Delimited, Data 8x12 Space Delimited (see Section 7, pages 18-19)
 - Fixed 8x12 O.D. Format
 - Comma Delimited Files for Importing into MicroSoft's ACCESS database software
10. Click the **Format** tab or **Go To Format Page**. The **Export to ASCII File** window appears. The choices on the **Data Fields** list depend on the MicroPlate in the file you've chosen. For a definition of each of the **Data Fields** (refer to Section 7, page 20 lists).



11. If necessary, click **Clear All Fields** to clear fields from the right list.
12. Click **Add All Fields** to select all data fields for ASCII File Fields.
13. Click on the left (data) field you want, then click on the next empty row on the right (ASCII) field. This information will be exported from the ASCII file. Add all the desired information.
Note: Clicking on a previously used row inserts a field.
14. Click **Remove Field** and then the **ASCII File Fields** field to remove a selected Data Field.
15. Click **OK**. The **Files** window reappears.
16. Click **Convert Data File → ASCII File**.
17. Click **OK** to exit.
18. Minimize the MicroStation/MicroLog software.
19. Open the software where you wish to view the exported data.

Exporting to a .csv file for EXCEL

1. Choose **Comma Delimited** as the **Save as Type** at Step 6 in the procedure above.
2. Choose **One Line, Comma Delimited** at Step 9 above. All fields selected for export will be in one line across the columns of the EXCEL spreadsheet. Data Field titles will also be exported.

3. Open EXCEL after you have exported the desired data and minimized the MicroStation/MicroLog software.
4. Click on **File → Open**. Under **Files of type**, choose **Text Files**.
5. Move through the folders to the C: drive. Open the Biolog folder and look for the file where the data was exported

Note: The file extension may not appear, but you will see the EXCEL symbol on the icon next to the file.

6. Click on your file name. The file will open and data will appear with Data Field Headers.
7. Select **Format** from the toolbar.
8. Select **Column**.
9. Select **Auto Fit Selection**.

Note: The instructions for opening the file in EXCEL may vary slightly, depending on the version you're using.

Exporting as a .txt or .prn file for EXCEL

Note: Use this to export reaction results in an 8x12 format similar to software printouts.

1. Choose **Text File** as the **Save as type** at Step 6 in Section 7, page 15.
 2. Choose **Line delimited, Data 8x12 Space Delimited** at Step 9 in Section 7, page 16. All fields selected for export will be on separate lines on the EXCEL spreadsheet. Data Field titles will not be exported.
 3. Open EXCEL after you have exported the desired data and minimized the MicroStation software.
 4. Click on **File → Open**. Under **Files of type**: choose **Text Files**.
 5. Move up through the folders to the C: drive. Open the Biolog folder and look for the file where the data was exported
- Note: The file extension may not appear, but you will see a text icon next to the file.*
6. Click on your file. A **Text Import Wizard** box appears.

7. Choose file type as **Delimited**.
8. Click **Next**.
9. Under **Delimiters**, deselect **Tab** and select **Space**.
10. Click **Next**.
11. Click **Finish**.
12. Your file will open with all the Data Fields on separate lines. If you chose to export any +/- values or OD values they will be in an 8x12 format.

*Note: The instructions for opening the file in EXCEL may vary slightly, depending on the version you're using. There are **no** Data Field Headers for each field exported.*

Exporting as a .txt or .prn file for Microsoft Word

Note: Use this option to export reaction results in an 8x12 format similar to software printouts.

1. Choose **Text File** as the **Save as type** at Step 6 in Section 7, page 15.
2. Choose **Line Delimited, Data 8x12 Space Delimited** at Step 9 in Section 7, page 16. All fields selected for export will be on separate lines on the EXCEL spreadsheet. Data Field titles will not be exported.
3. Open Microsoft Word after you have exported the desired data and minimized the MicroStation/MicroLog software.
4. Click on **File → Open**. Under **Files of type:** choose **Text Files**.
5. Move up through the folders to the C: drive. Open the Biolog folder and look for the file where the data was exported.


Note: The file extension may not appear, but you will see a text icon next to the file.

6. Click on your file name. A **File Conversion** box appears.
7. Choose **Plain text** as the encoding loading this file.
8. Click **OK**.

9. Your file will open with all the Data Fields on separate lines. If you chose to export any +/- values or OD values they will be in an 8x12 format.

Exporting to Microsoft Access

You may also wish to export data to Microsoft Access. This can be done by clicking the **Export to ASCII File** button on the **Data File** window. In the Step 2 box mark the comma delimited for Access format. For further information on exporting into Access, please see the Help tab on the **Export to ASCII File** window (see below).



ASCII Export Window Entries

- Date: Date data was saved
- Plate Type: GN, GP, YT, FF, etc. required and selected from saved file
- Plate Number: number you designated for plate information
- Incubation Time: time you designated depending on MicroPlate type
- Sample Number: number you designated for plate information
- Strain Type: GN-ENT, or GN-NENT, etc., determined by the MicroPlate type
- Strain Name: name you designated for plate information
- Strain Number: number you designated for plate information
- Other: information you add to further describe sample
- “+/-“ Reactions: results from reading MicroPlate displayed as:
 - “+” = positive
 - “-“ = negative
 - “b” = borderline
- 1/0.5/0 Reactions: OD display of results if reactions read manually
 - 1 = positive
 - 0 = negative
 - 0.5 = borderline
- Dual Wavelength Data (DWD): O.D. data on printouts if GN, GP, AN, MT or Eco plate is read from a MicroStation Reader = $[(590\text{nm} - 750\text{nm})_x - (590\text{nm} - 750\text{nm})_{A1}] \times 1000$, where x=any well
- Dual Wavelength O.D. (DW): OD difference data if GN, GP, AN, MT, or Eco plate read from a MicroStation Reader = $(590\text{nm} - 750\text{nm})_x$, where x=any well
- 750nm O.D.: Raw O.D. read at 750 nm from MicroStation Reader
- 590nm O.D.: Raw O.D. read at 590 nm from MicroStation Reader
- 590nm %: O.D. data on printouts if YT, SFP, or SFP plate read from a MicroStation Reader=
- 590nm % = $(590\text{nm}_x - A1) \times 100$, where x=any well
- 590nm Filter Position: position 5
- 750nm Filter Position: position 6
- Data Generation Mode: Manual or Reader
- ID Result: Message if there is NO ID or Genus species name if ID obtained from reading MicroPlate
- ID #1 Name: Genus species name of organism listed as first choice on printout of results
- ID #1 SIM: SIM of organism listed as first choice on printout of results
- ID # 1 DIS: DIS of organism listed as first choice on printout of results
- ID #2 Name: Genus species name of organism listed as second choice on printout of results
- ID #2 SIM: SIM of organism listed as second choice on printout of results
- ID # 2 DIS: DIS of organism listed as second choice on printout of results
- ID #3 Name: Genus species name of organism listed as third choice on printout of results
- ID #3 SIM: SIM of organism listed as third choice on printout of results
- ID # 3 DIS: DIS of organism listed as third choice on printout of results

ASCII Export Window Entries (Unique Entries for FF)

- Color Data: O.D. Data on FF printouts for color data =
- $[(490\text{nm}_x - 490\text{nm}_{A1}) - (750\text{nm}_x - 750\text{nm}_{A1})] \times 1000$, where x=any well.
- Turbidity Data: O.D. Data on FF printouts for turbidity data =
- $(750\text{nm}_x - 750\text{nm}_{A1}) \times 1600$, where x=any well.
- Dual Wavelength O.D. (DW): OD readings at 490 nm and 750 nm minus the A-1 well:
- $(490\text{nm}_x - 450\text{nm}_{A1})$, where x=any well.
- $(750\text{nm}_x - 750\text{nm}_{A1})$, where x=any well.
- 750nm O.D.: Raw O.D. read at 750 nm from MicroStation Reader.
- 490nm O D: Raw O.D. read at 490 nm from MicroStation Reader.
- Color Filter Position: 490nm, position 3
- Turbidity Filter Position: 750nm, position 6

Cluster Analysis

Cluster analysis uses a mathematical approach to display the relative similarities and differences between strains or species based on pattern matching. It offers a visual method for examining groups using linear dendrograms and two- and three-dimensional cluster diagrams. You can generate cluster analyses on groups of entries in Biolog's database or your own database.

The math behind cluster analysis

Dendrograms are generated using a modified UPGMA analysis developed by Biolog (algorithm for constructing rooted phylogenetic trees). The algorithm uses DIST values to generate the branching structure of the dendrogram. A DIST value of 14 is usually used as a cutting level to separate genera. This value has been empirically determined to be an appropriate one for most genera of gram-negative bacteria.

Once the software generates a preliminary dendrogram, it makes a second pass, rotating all branch points to their most preferred position. This best maintains proximity of related strains or species. If the distance between strains or species adjacent to cut points is less than 14, they are joined with a dashed line that indicates their actual distance.

By mathematically removing the linear constraints of the dendrogram, the software can generate two- and three-dimensional diagrams of DIST relationships.

All of these cluster diagrams are only approximations of true relationships. With 95 carbon source tests, it would take 95-dimensional space to provide fully detailed clustering. As a practical matter, the linear dendrograms and two- and three-dimensional diagrams give very satisfactory representations.


Using cluster analysis with IDs

1. Once you get a species ID, click the **Cluster Menu** selection bar on the **Data** window.

2. The screen shows a dendrogram of your ID.

Dendrogram numbering:


→ The number in () is one less than the corresponding number found in the data file record.
→ The numbering here begins with 0
→ Data file numbering begins with 1




CLUSTER MENU (TEXT)

3. Click on the black line below or above the dendrogram to scroll.


4. Click the **Cluster Menu (Graphics)** tab to view one-dimensional, two-dimensional, and three-dimensional cluster analyses.



ONE-DIMENSIONAL (1-D) CLUSTER ANALYSIS



TWO-DIMENSIONAL (2-D) CLUSTER ANALYSIS



THREE-DIMENSIONAL (3-D) CLUSTER ANALYSIS

5. Design your cluster according to your preferences.

- Select **B/W** for a black and white representation.
- Select **Copy** to paste your cluster into a text document (print screen).

Note: Copy to Paint and use the cut and paste feature to place a “cluster only” picture into a Word document.

- Select **Color** to choose the color scheme you want.
- Select **Big** to enlarge your species.
- Select **Small** to shrink your cluster.

- Select **In** to zoom in on cluster.
 - Select **Out** to zoom out on cluster.
 - Select **Distance** to designate the distances between organisms with lines and to note the distances with a key in the lower right corner of the cluster.
 - Select **Groups** to encircle the entries that are in the same groups in the dendrogram.
 - Select **Lines** in the three-dimensional cluster to draw a line from each organism to the base.
 - Select to the left of **Live** to enable the spin of the three-dimensional cluster left.
 - Select to the right of **Live** to enable the spin of the three-dimensional cluster right.
 - Click on **Spin, Roll, and /or Rotate, Fast or Slow**, to capture the best view of your cluster.
6. Click **Print** to print any diagram. (Note: “Live” feature must be disabled to print cluster diagram)

Note: When searching your own database, you can view and print dendrograms, 2-D, and 3-D cluster diagrams by following the instructions above in “Using cluster analysis with IDs.”

Using cluster analysis with collected data

To obtain a dendrogram displaying the spatial relationships between your own organisms in a saved data file, refer to “Printing Data Headers, Species Counts, and Dendrograms” beginning on page 7-13.

1. From the **Compile** window, enable the desired box:
 - Make Index Files
 - **Print Species/Plate Dendrograms** (to generate one dendrogram for each MicroPlate)
2. When the data file compiles, printouts of your dendrograms will display linear relationships.

8. Administration and Security

In this section:

- ➔ **What Is Restricted Access Mode?**
- ➔ **First Log-In and Setting Up an Administrator**
- ➔ **Administrator Options**
- ➔ **Viewing the Log-In Log**
- ➔ **Changing or Lost Passwords**
- ➔ **Unauthorized Log-In Attempts**

The software runs in Restricted or Unrestricted Access Mode, a term that may be familiar from experience with the OmniLog® System. Let's review what Restricted Access Mode does and why it is important.

What is Restricted Access Mode?

Restricted Access Mode requires that only users with Administrator privileges have the ability to oversee and control access to the software.

The Administrator will:

- Assign Usernames and passwords for those who will use the system.
- Assign access privileges to each user.
- Oversee the security of the system

Restricted Access Mode requires all users to Log-In with a Username and password when entering the software or changing users. The software maintains files of the **User List**, **Log-In Log** and **Log-In Log Archive** to keep track of registered users, access privileges, and Log-In/out activity in the system.

All of these files are encrypted. The Log-In Log Archive files can be placed by the Administrator in the computer/network location of their choosing. This allows the Administrator to place the files in a location that has secured access.

Why Restricted Access Mode?

Restricted Access Mode ensures that data integrity and security controls are implemented in accordance with the guidelines of 21 CFR Part 11. It assists in compliance with federal Current Good Manufacturing Practices (cGMP) by ensuring the integrity of software use and electronic files.

Security Features:

- User List and Log-in Log
- Notification of failed Log-In attempts

**21 CFR Part 11:
Only restricted
data should be
imported.
(global setting)**

- User password change
- User privileges
- Session time-out
- Password expiration
- Log-in archive

Electronic record integrity:

- Limiting access
- Original record integrity
- Documentation of changes (by whom and when)
- Audit trails

What if I don't work in a 21 CFR compliant environment?

No problem!

- It is easy for the software Administrator to assign full access privileges to anyone in your organization who will be using the software.
 - This will give everyone the ability to use all of the system's features, with the exception of the ability to access the Administration Options button.
 - All records will contain the Username of the individual users

Other options are:

- Simply to use the Administrator ID and password for all users, thereby giving everyone equal access.
- Select to run the software in Unrestricted mode.

First Log-In and Setting up an Administrator

Please refer to *Section 2, page 2* for full instructions on how to set up an Administrator and Log-In to the software for the first time. Ideally, this will be done soon after the software is installed.

Administrator Options

The Administrator (or a user with Administrator access rights) uses the Administration Options button, located on the Welcome window, to register new users and assign access privileges, as well as manage security controls. This button is only available when the software Administrator or a user with Administration access is Logged In; in all other cases it will not appear on the Welcome window.

Important!


- If you test Password Expiration Period or Time Out Period features using system clock changes, the software must be closed when system clock changes are made.
- If not the Administration screen may appear with each Administrator Log-in until 100 Log-In/Log-Outs occur.

Remember:

- A user name must be at least 1 character in length.
- The password must be at least 6 characters in length, contain at least one number, and is case sensitive.

Options Tab

Once you have selected the Administration Options buttons, the **Administration** window will open. By default, it will appear with the **Options** tab selected. Make your selections from the available options and then click the **Save and Close** button.



- **Destination Directory for Log-In-Log files**
Shows the default directory where Log-In Log Archive files are saved. Select the desired directory to place in the computer/network location of your choosing. You should select a secure location (for example, in a secure server or a password protected file folder).
- **Password Expiration Period**
Use the pull down menu to select either **Three Months** or **One Minute**. The default is three months. This requires all users to select new passwords after 3 months. Select one minute to expedite validation testing only.
- **Time-Out Period**
Use the pull down menu to select either **Fifteen Minutes** or **10 seconds**. The default is fifteen minutes. Select **10 seconds** to expedite validation testing only. Please refer to *Timed Log Out* in *Section 3, Pg. 5* for more information.
- **Restricted Access Mode**
Click on the checkbox to change the mode (make sure there's a check in the box if you want Restricted Access mode).

Creating a User List

From the **Administration** window, select the **User List** tab. Go to this list to add new users and assign levels of software access to those users, or to make changes to their level of access. User names can never be changed or deleted.

Once the administrator adds a new user name to the user list, that name can never be deleted or changed.

Adding new users and Assigning Privileges

1. Click the **User List** tab. A numbered list will appear, showing all users registered to that point, starting with the original Software Administrator (in row 1).
2. To add a new user, click in the next blank field in the **Username** column. Enter the new user name.

| Options | User List | Log-In Log | | | | | | | |
|---------|---------------------|-------------------|----------|------------------------|--------|--------|-------|--------|-------|
| | Last Log-In Attempt | Last Log-Out | Username | Assigned Password | Log-In | Set-Up | View/ | Edit D | Admin |
| 1 | Aug 23 2004 20:06 | Aug 23 2004 20:07 | biolog | Original Administrator | Yes | Yes | Yes | Yes | Yes |
| 2 | Aug 23 2004 20:05 | Aug 23 2004 20:06 | biolog2 | microbe2 | Yes | Yes | Yes | Yes | No |
| 3 | Aug 23 2004 20:08 | Aug 23 2004 20:09 | biolog3 | microbe4 | Yes | Yes | No | No | No |
| 4 | Aug 23 2004 20:09 | Aug 23 2004 20:09 | biolog4 | microbe6 | Yes | Yes | No | No | No |
| 5 | Aug 23 2004 20:09 | Aug 23 2004 20:10 | biolog5 | microbe8 | Yes | No | No | No | No |
| 6 | never | never | | | No | No | No | No | No |

3. Click in the blank **Assigned Password** field next to that new user name. Enter a temporary password for that new user.
4. Click in each Privilege box to the right, toggling between **Yes** and **No** to assign or deny specific access levels to that user.
5. Click the **Save and Close** button when you are finished.
6. Give the Username and Password to the person you have registered, and refer them to *Logging In and Out* (starting on Pg. 4, Section 3), if they need help logging into and out of the software. Remember that the password you have chosen is only temporary; the user will be prompted to enter a new password the first time they Log-In to the system.

Each new user must be assigned a level of software access by the Administrator. Consult Table 9-1 for a listing of these privileges.

Log-In Privileges

| Privilege | What It Allows | What It Does Not Allow |
|------------|---|---|
| Log-In | User will be able to Log-In and out of MicroLog software using a password, and will be able to use the Set Up Input and Data features (to enter MicroPlate information, read MicroPlates, and print results. All Data will be saved automatically). | User cannot use Data Management functions |
| Set-Up | User will be able change screen colors and fonts, will be able to access the Detailed Reader Setup page (for troubleshooting), and can disable the Auto Save feature. | User cannot use Data Management functions |
| View/Print | User will be able to view or print data files | User cannot edit data files |
| Edit | User will be able to use all Data Management features, including editing data files and compiling databases. | User cannot perform any Administrator functions |
| Admin | User will have complete access to all aspects of the software (just like the software administrator). | User cannot delete or change user names |

TABLE 8-1: USER ACCESS PRIVILEGES, IN DESCENDING ORDER OF LOWEST TO HIGHEST LEVELS OF ACCESS.

Viewing the Log-In Log

View the Log-In Log by selecting the **Log-In Log** tab on the **Administration** window. This feature of the software keeps meticulous track of the software's use in descending order of date (with the most recent date and log-in first). This log is non-editable.

The Log-In Log records the past 100 log-in attempts. For each log-in attempt, the Log-In Log records the username, date and time logged-in, whether the username is registered, and the date and time the person logged-out. The log also records which privileges registered users have, and shows the software is in Restricted Access Mode. This creates an audit trail. As log-in records in excess of 100 drop from the list, they are saved to a "Read-Only" file (Logins/LoginLogXXXXX.csv¹). When this saved Log-In Log archive file contains 100 records, subsequent records will be saved to a new Log-In Log archive file with a different date stamp¹.

¹ XXXXX is a five-digit date stamp

| Options User List Log-In Log | | | | | | | | | | |
|----------------------------------|----------|------------------|----------|-------------------|-------------------|--------|---------|--------|----------|-------|
| | Restrict | Entered Username | Register | Attempt Date | Log-Out Time | Log-In | Options | View/F | Edit Pri | Admin |
| 1 | Yes | biolog | Yes | Aug 23 2004 20:10 | never | Yes | Yes | Yes | Yes | Yes |
| 2 | Yes | biolog5 | Yes | Aug 23 2004 20:09 | Aug 23 2004 20:10 | Yes | No | No | No | No |
| 3 | Yes | biolog4 | Yes | Aug 23 2004 20:09 | Aug 23 2004 20:09 | Yes | Yes | No | No | No |
| 4 | Yes | biolog3 | Yes | Aug 23 2004 20:08 | Aug 23 2004 20:09 | Yes | Yes | Yes | No | No |
| 5 | Yes | biolog | Yes | Aug 23 2004 20:06 | Aug 23 2004 20:07 | Yes | Yes | Yes | Yes | Yes |
| 6 | Yes | biolog2 | Yes | Aug 23 2004 20:05 | Aug 23 2004 20:06 | Yes | Yes | Yes | Yes | No |
| 7 | Yes | biolog | Yes | Aug 23 2004 20:04 | Aug 23 2004 20:05 | Yes | Yes | Yes | Yes | Yes |
| 8 | No | | No | never | never | No | No | No | No | No |
| 9 | No | | No | never | never | No | No | No | No | No |
| 10 | No | | No | never | never | No | No | No | No | No |
| 11 | No | | No | never | never | No | No | No | No | No |
| 12 | No | | No | never | never | No | No | No | No | No |
| 13 | No | | No | never | never | No | No | No | No | No |
| 14 | No | | No | never | never | No | No | No | No | No |
| 15 | No | | No | never | never | No | No | No | No | No |

THE LOG-IN LOG

Table 8-2 defines each of the columns in the Log-In Log and describes what they mean.

| Interpreting the Log-In Log | |
|-----------------------------|---|
| Column Name | Information Given |
| Restricted | If software was in Restricted Access mode when user logged-in, this entry will say Yes ; if software was operating in Unrestricted Access mode, this entry will say No . |
| Entered Username | The user name of each person who logged in is listed here. |
| Registered | If this person was an approved user, this entry will say Yes ; if not, this entry will say No . |
| Attempt Date | This is the exact date and time the user logged in (or attempted to do so). |
| Log-Out Time | This is the exact date and time the user logged out. The entry here will read Never in the event of a failed log-in and when the software administrator is currently using the software. |
| User Privileges | |
| Log-In | Yes if access was given during that log-in period; No if not |
| Option | Yes if access was given during that log-in period; No if not |
| View/Print | Yes if access was given during that log-in period; No if not |
| Edit | Yes if access was given during that log-in period; No if not |
| Admin | Yes if access was given during that log-in period; No if not |

TABLE 8-2: INFORMATION IN THE LOG-IN LOG.


Changing or Lost Passwords

The software requires that all users **Log-In** with a password to enter the software. This ensures a secure environment by controlling access to the software. A user can change his or her password at any time.

Remember that the password must be at least 6 characters in length, contain at least one number, and is case sensitive.

Changing your Password

1. **Log-In** to the software.
2. At the **Welcome** window, select the **Change Password** button.
3. The Change Password dialog box appears.
4. Enter the password that was assigned to you, then enter a password chosen by you and reenter it.
5. Click **OK**.



CHANGE PASSWORD DIALOG WINDOW

NOTE: The software will not allow the use of the last four passwords established for a given user name.

Lost or Revoked Password

Occasionally, a user might simply forget or misplace their password. The user may have attempted entry 5 times with the incorrect password. If the user is NOT the only Administrator, this situation is easily remedied. Any Software Administrator can easily go to the **User List** under **Administration** and assign the person a new password. Make sure to re-assign the user's privileges as well, or the new password will not work.

Logging In

When the system is operating in Restricted Access Mode, all users must Log-In when starting up the software or changing users.

1. When starting up the software, the software will automatically show a **Password Dialog** window. Enter your user name and password, and then click **OK**. After a few seconds, the **Welcome** window will appear.



PASSWORD DIALOG WINDOW

2. If the software is already open and it is necessary to change users, go to the **Welcome** window and click **Change Users**. A **Password Dialog** window will appear. Enter the new user's name and password.

Unauthorized log-in attempts

If someone enters an incorrect user name and/or password, the software will allow five attempts to enter the information correctly. After five attempts, the software will automatically close.

Subsequent attempts to open the software and Log-In using an incorrect user name or password will result in a warning tone and screen message noting that "Unauthorized Access Has Been Attempted."

The message will remain until the Administrator clicks on the message while in the **Administration** screens.

9. Technical Notes

In this section:

- ➔Materials List
- ➔Media Preparation
- ➔Specialized MicroPlates

This section gives specific information regarding consumable materials you'll need to perform microbial identifications, as well as instructions for how to prepare isolation media and inoculating fluid.

Materials List

| | Item | Description |
|------------------------|------------------------|--|
| Media | BUG + B | Growth medium for gram-negative and gram-positive bacteria |
| | BUG + M | Growth medium for gram-positive spore-forming rods |
| | BUG | Growth medium for agricultural bacteria |
| | BUA | Growth medium for anaerobes |
| | BUY | Growth medium for yeasts |
| | CHOC | Growth medium for fastidious gram-negative bacteria |
| | 2% ME | Growth medium for filamentous fungi (and select yeast species) |
| Inoculating Fluid | GN/GP-IF | Inoculating fluid for gram-negative and gram-positive bacteria |
| | AN-IF | Inoculating fluid for anaerobic bacteria |
| | FF-IF | Inoculating fluid for filamentous fungi (and select yeast species) |
| | Sterile water | Inoculating fluid for yeasts |
| | Salicylate | Inoculating fluid additive for some microbes |
| | Thioglycolate | Inoculating fluid additive for some microbes (available in droppers) |
| Miscellaneous Supplies | Streakerz | Sterile 6" pointed streaking sticks |
| | LongSwabs™ | Sterile 7" cotton-tipped swabs |
| | Transfer pipettes | Sterile 6" disposable pipettes, graduated tip, 5 ml |
| | Reagent reservoirs | Sterile reservoirs |
| | Test tube rack | Autoclavable, holds 40 tubes, 20 mm diameter |
| | Pipette tips | Sterile, racked, for repeating multichannel pipettor |
| | Pipette tips, filtered | Sterile, filtered, racked, for repeating multichannel pipettor |
| | MicroPlate lids | Gamma irradiated MicroPlate lids |
| MicroPlates | GN2 MicroPlates | MicroPlates for gram-negative bacteria |
| | GP2 MicroPlates | MicroPlates for gram-positive bacteria |
| | AN MicroPlates | MicroPlates for anaerobic bacteria |
| | YT MicroPlates | MicroPlate for yeasts |
| | FF MicroPlates | MicroPlates for filamentous fungi (and select yeast species) |

TABLE 9-1. CONSUMABLE MATERIALS AND EQUIPMENT

Media Preparation

Caution!

If you're making your own media, follow instructions exactly. Take care to avoid contamination.

Biolog can provide you with all the pre-made media you will need for microbe identification. However, you can also prepare your own media, using Biolog's formulas.

Making Biolog Universal Growth Agar + blood (BUG + B)

1. Mix the following in a 2-3 liter container:
 - 57 grams BUG Agar
 - 950 ml purified, distilled, or deionized water
2. Gently boil mixture while stirring to dissolve the agar and other components.
3. Cool an aliquot and measure the pH. Adjust pH with NaOH or HCl to achieve a final pH of 7.3 ± 0.1 at 25°C .
4. Sterilize by autoclaving at 15 pounds of pressure and 121°C for 15 minutes.
5. Cool to $45\text{-}50^{\circ}\text{C}$.
6. Add 50 ml sterile fresh defibrinated sheep's blood just prior to dispensing and mix gently. (Use good quality blood with a hematocrit of at least 40%).
7. Dispense into sterile petri dishes.

Making Biolog Universal Growth Agar + maltose (BUG + M)

1. Follow steps 1-5 above for BUG + B, except use 990 ml purified water.
2. Add 10 ml sterile solution of 25% maltose and mix thoroughly.
3. Dispense into sterile petri dishes.

Making Biolog Universal Growth Agar (BUG)

1. Mix the following in a 2-3 liter flask:
 - 57 grams BUG Agar
 - 1,000 ml purified, distilled, or deionized water
2. Gently boil while stirring to dissolve the agar and other components.
3. Cool an aliquot and measure the pH. Adjust pH with NaOH or HCl to achieve a final pH of 7.3 ± 0.1 .
4. Sterilize by autoclaving at 15 pounds of pressure, at 121°C , for 15 minutes.
5. Cool to $45\text{-}50^{\circ}\text{C}$.
6. Dispense into sterile petri dishes.

Making Biolog Universal Anaerobe Agar + blood (BUA + B)

1. Mix the following in a 2-3 liter container:
 - 51.7 grams BUA Agar
 - 950 ml purified, distilled, or deionized water
2. While flushing with oxygen-free nitrogen gas, gently boil mixture while stirring to dissolve the agar and other components.
3. Cool an aliquot and measure the pH. Adjust pH with NaOH or HCl to achieve a final pH of 7.2 ± 0.1 at 25°C .
4. Sterilize by autoclaving at 15 pounds of pressure and 121°C for 15 minutes. Make sure the bottle is tightly capped to prevent the entry of oxygen.
5. Cool to $45\text{-}50^{\circ}\text{C}$ under a stream of oxygen-free nitrogen gas.
6. Add 50 ml sterile fresh defibrinated sheep's blood just prior to dispensing and mix gently. (Use good quality blood with a hematocrit of at least 40%).
7. Dispense into sterile petri dishes in an anaerobic chamber.

Making Biolog Universal Yeast Agar (BUY)

1. Mix the following in a 2-3 liter container:
 - 60 grams BUY Agar
 - 1,000 ml purified, distilled, or deionized water
2. Gently boil while stirring to dissolve the agar and other components.
3. Cool an aliquot and measure the pH. Adjust pH with NaOH or HCl to achieve a final pH of 5.6 ± 0.4 at 25° C.
4. Sterilize by autoclaving at 15 pounds of pressure and 121° C for 15 minutes.
5. Cool to 45-50° C.
6. Dispense into sterile petri dishes.

Making 2% Malt Extract Agar

1. Mix the following in a 2-3 liter container:
 - 20 grams Oxoid Malt Extract (LP0039B)
 - 18 grams Bacteriological Grade Agar
 - 1,000 ml purified, distilled, or deionized water
2. Gently boil while stirring to dissolve the agar and other components.
3. Cool an aliquot and measure the pH. Adjust pH with NaOH or HCl to achieve a final pH of 5.5 ± 0.2 at 25° C.
4. Sterilize by autoclaving at 15 pounds of pressure and 121° C for 15 minutes.
5. Cool to 45-50° C.
6. Dispense into sterile petri dishes.

General hints for culture media preparation



Keep dehydrated media powder in original bottles with lids tightly closed to avoid water absorption and deterioration.



Use clean dry glassware that has been rinsed free of all soap residue.



Add water to the vessel first, then weigh agar powder and add it to the vessel. Mix to obtain an even suspension. Do NOT fill the vessel more than two-thirds full (to avoid boiling over during heating and autoclaving).



Heat agar gently, with constant stirring.

Making inoculating fluid (GN/GP-IF)

Biolog supplies pre-made inoculating fluid, with the following composition:

- 0.40% sodium chloride (NaCl)
- 0.03% Pluronic F-68 (e.g., Sigma P7061)
- 0.02% Gellan Gum (e.g., Phytigel™, Sigma P8169)

However, if you wish to prepare your own, using Biolog's formula.

1. Add 0.2 grams Gellan Gum to 1 liter of H₂O.
2. Heat to boiling, stirring constantly, until the Gellan Gum is completely dissolved.
3. Turn off the heat but continue stirring.
4. Add 4 grams NaCl and continue stirring until the NaCl is completely dissolved.
5. Add 0.3 grams of Pluronic F-68. Continue stirring until the Pluronic F-68 is completely dissolved.
6. Dispense the volume appropriate to obtain 19 ml post-autoclaving into 20 x 150 mm tubes.
7. Sterilize by autoclaving at 15 pounds of pressure and 121° C for 30 minutes.

Note: Gellan Gum concentration of 0.01-0.02% is acceptable.

Making inoculating fluid (FF-IF)

Biolog supplies pre-made inoculating fluid, with the following composition:

- 0.03% Tween 40 (e.g., Sigma P1504)
- 0.25% Gellan Gum (e.g., Phytigel, Sigma P8169)

However, if you wish to prepare your own, using Biolog's formula.

1. Heat to boiling 1 liter of H₂O.
2. Add 2.5 grams Gellan Gum and 0.3 grams Tween 40. Stir.
3. Turn off the heat but continue stirring until all components are completely dissolved and the solution becomes clear.
4. Dispense the volume appropriate to obtain 16 ml post-autoclaving into 20 x 150 mm tubes.
5. Sterilize by autoclaving at 15 pounds of pressure and 121° C for 30 minutes.

Specialized MicroPlates

Biolog makes a series of additional MicroPlates for special applications. While there is currently no database for these, you can use them to create your own database with the software.

Table 9-2 lists these specialized MicroPlates and their uses. Contact Biolog or your Biolog Distribution Partner for more information.

| Specialized MicroPlate | Use |
|-----------------------------|--|
| SF-N2 and SF-P2 MicroPlates | For testing of Sporulating and Filamentous microorganisms (such as actinomycetes and fungi) |
| MT2 MicroPlate | Contains same nutrient base and color chemistry as GN MicroPlate, but without carbon sources. Add your own carbon sources to study the metabolic capabilities of any microbe of your choice. |
| EcoPlate | Contains the same color chemistry as the GN2 plate, but has three sets of 31 the identical carbon sources. <i>Note: the A1 well is used as the threshold algorithm control. Running different samples on one plate may affect the positive/negative/borderline calls.</i> |

TABLE 9-2. BIOLOG SPECIALIZED MICROPLATES

10. Frequently Asked Questions

In this section:

➔ Answers to
Common
Questions

Q *How I should I store my MicroPlates?*

A Unpack MicroPlates as soon as you receive them and keep them refrigerated until use. Biolog MicroPlates should be kept refrigerated (not frozen). MicroPlates are fairly stable at room temperature; we suggest refrigeration because they must be kept cold for maximum shelf life. At moderately warm temperatures, they slowly begin to deteriorate. Visually examine MicroPlates before using them. You may see faint yellow-brown or pink shades in wells when they are dry (this is OK). However, if wells show significant color immediately after inoculating the MicroPlate with the cell suspension, they have most likely been heat damaged.

Q *Do I have to handle frozen or lyophilized cultures in a special way?*

A Yes. Subculture frozen or lyophilized cultures two to three times before testing.

Q *Our lab would really like to use our own media. We know that our strains will grow on it. Is that all right?*

A No. However tempting it may be to use nonrecommended media, using anything other than Biolog-recommended media is a mistake. The carbon source tests in our MicroPlates are configured for precise metabolic reactions. The metabolic state of an organism has a profound influence on the resulting MicroPlate pattern that is the basis for identification. Many species will give fewer positive reactions when grown on non-recommended media. Use Biolog-recommended media only!

Q *Why is there a defined range for the turbidity of inocula?*

A It is essential to prepare inocula in a consistent and precise manner. The Biolog redox test chemistry is sensitive to oxygen concentration, which is determined by cell density. Use Biolog turbidity standards and calibrate your turbidimeter.

Q *Are Biolog's Turbidity standards comparable to McFarland standards?*

A Biolog's Turbidity standards have different set points and are not comparable to McFarland standards. Biolog's standards are set for cell densities at which the Biolog test performs optimally. Use of Biolog standards is necessary to achieve consistent and accurate results.

Q *Can I use one incubator for all my strains and set it at 32 °-33 °c instead of having two incubators at 30 °c and 35 °-37 °c respectively?*

A Organisms grow better depending on the temperature of their environment. Biolog has chosen the temperatures at which the organisms in our database grow optimally. Using different temperatures than those recommended will decrease the performance of your Biolog ID system.

Q *We read our MicroPlates visually. Occasionally we are unsure whether a reaction is positive or negative. How should we enter those wells?*

A Light purple reactions are considered positive as long as the color is noticeable when compared to the A1 reference well. However, if you're still unsure, enter these reactions as "borderline." If, for example, you enter 4 out of 95 tests as borderline, the software will ignore the 4 borderline tests and base its identification on the remaining 91 tests. This will give you far more accurate results than guessing wrong about reactions you're not sure of.

Q *What should I do if I get different ID results at 4 hours and at 24 hours?*

A This is most often the result of having a mixed culture. See Section 4 for recommended procedures and examine your culture carefully to see if it is mixed. Also, consider the similarity value at one incubation time vs. the other. Environmental strains will often grow a little differently than the lab strains used to make the databases. Some strains may have a high similarity in our 4 hour database, while others match our 24 hour database more closely. Compare your similarity values to see which is a better match.

Q *Are there any genera where the taxonomy is still undergoing change?*

A Taxonomists have not yet agreed upon how all species should be delineated and microbiology is an ever-evolving science. We update and expand our software library on a continual basis, but at this point there are certain genera still undergoing revision. These include *Bacillus*, *Corynebacterium*, *Enterobacter*, *Pseudomonas*, *Aeromonas*, *Vibrio*, *Acinetobacter*, *Moraxella*, *Clostridium*, and *Candida*. www.dsmz.de/bactnom is a good site for current bacterial nomenclature. You can also link to this site through Biolog's web page. www.cbs.knaw.nl is a good site for current yeast and filamentous fungi nomenclature and information.

Q *Our lab is interested in your quality control and validation procedures. Where do we get this information?*

A Material Safety Data Sheets, Quality Control organism information, and Certificates of Performance are available on the website. Validation procedures, reagents kits and bacterial strains are available from Biolog.

11. Troubleshooting

In this section:

- Gram Stain Identification
- Culturing Microorganisms
- Preparing Inocula
- Inoculating MicroPlates
- Incubating MicroPlates
- MicroStation Reader

Carefully following the instructions in this guide will greatly minimize problems. Occasionally, however, you may get stuck or encounter difficulties. This section addresses the symptoms, causes, and solutions to those occasional problems. If you still can't figure out the cause of the problem, call Biolog Technical Service. We're always glad to help.

Additional help can be found on Biolog's website at www.biolog.com.

Gram Stain Identification

Since Gram stain readings start the chain of steps leading to a MicroStation/ MicroLog identification, it's essential to perform them according to standard lab protocol and to interpret results correctly.



View the smear with a light microscope, using the oil-immersion objective. Gram-positive bacteria appear blue or violet; gram-negative bacteria appear pink or red.



Gram-negative bacteria may appear gram-positive if the smear is too thick and decolorization is incomplete.



Gram-positive bacteria may appear gram-negative if the smear is over-decolorized. This can also occur if the culture is not fresh and has reached stationary phase (some *Bacillus* species are gram positive for only a few divisions after spore germination). In addition, gram-positive bacteria will appear gram-negative if the integrity of their cell walls is damaged.



To prevent misidentification, prepare light smears of young, actively growing cultures. Use known gram-positive and gram-negative controls. When you're determining organism morphology, perform Gram stains from liquid medium. Solid agar can affect the appearance of the organism.

| Symptom | Cause | Solution |
|---------------------------------|---|---|
| Gram variability | Smear is too thick | Make a single-cell-layer smear. Run positive and negative controls. |
| All smears appear gram negative | Overwarming of slides | Air dry slides completely before warming. Don't use a flame to fix slides. |
| | Using old cultures | Prepare slides from fresh log-phase growth cultures. |
| | Over-decolorizing | Repeat Gram stain, using specified timing. Run positive and negative controls. |
| Morphology unclear | Growth on agar media affects morphological appearance | Do a Gram stain from broth culture. |
| Too-intense color | Under-decolorizing | Use decolorizer made by same company that makes your stains. Prepare decolorizer or acetone/alcohol at several concentrations. Test solutions at 5, 10, and 15 seconds with known positive, negative, and weakly positive cultures. Run positive and negative controls. |

Additional differentiation techniques

Make sure you're using proper Gram stain procedure. If you still get ambiguous results, use the following table as a guide:

| Test | Reactions |
|---|--|
| Vancomycin sensitivity (30 µg disks) | Gram-negative bacteria → resistant Gram-positive bacteria → most are sensitive |
| Growth on MacConkey and CNA plates | Gram-negative bacteria → MacConkey +, CNA – Gram-positive bacteria → MacConkey –, CNA + |
| KOH test (3%)* | Gram-negative bacteria → suspension becomes thick and stringy |
| L-alanine aminopeptidase activity (use commercially available test strips impregnated with a colorimetric substrate for this enzyme)* | Gram-negative bacteria → + Gram-positive bacteria → – |

*For test procedures, see *Bailey & Scott's Diagnostic Microbiology*, Baron E.J. and Finegold, S.M., C.V. Mosby Co., 1990, pages 102-103.

Culturing Microorganisms

| Symptom | Cause | Solution |
|---|---|--|
| Poor overall growth | Using nonrecommended media | Use Biolog-recommended media. |
| Isolate takes several days to form a colony | Some environmental organisms take several days to become visible on growth plate, at which point culture may be too old for successful ID | Subculture one or two passages in broth medium. Set up two or three agar plates to obtain sufficient growth. |
| Bacterium will not grow on BUG + B or forms pinpoint colonies | Fastidious gram-negative bacteria need special culture conditions | Use Chocolate agar and incubate with elevated CO ₂ at 35-37° C. |
| | Some environmental bacteria may be oligotrophic or temporarily in an oligotrophic state (i.e., they will only grow on low nutrient media such as R2A) | Try to subculture the bacterium from R2A and see if they will gradually adapt to growing on BUG + B. If not, it is a species that is not included in our MicroStation/MicroLog database. Agricultural bacteria may be grown on BUG without blood. |
| Slow growth and/or weak pattern formation | Sub-optimal growth temperature and/or humidity and/or atmosphere | Use specified incubating temperatures (see Section 4 and Appendices). If your unknown organism came from a cold or warm environment, grow it first at its environmental temperature, then try to grow it at 35-37° C, 30° C, or 26° C. Add a pan of water to provide humidity in your incubator. Incubate the MicroPlate at the same temperature as the growth plate. |
| Mixed growth on agar plates | Sticky bacterial surfaces | Transfer a colony into a few ml of GN/GP-If or sterile water, vortex for several seconds, and streak out from the cell suspension onto a medium that aids the detection of subtle differences in colony morphology. Take a single colony and streak for isolation on a separate plate. Repeat as necessary until you have a pure culture. |
| | Filamentous fungi mixed with bacteria | Cultures can be streaked on media containing antibacterial antibiotics. If this fails, use 1 drop of saline suspension of the mixed culture to inoculate each of four tubes of Sabouraud broth (1.0 ml). Add in series either 1,2,3 or 4 drops of concentrated HCl. Streak out from each tube. |

Preparing Inocula

| Symptom | Cause | Solution |
|--|--|--|
| Incorrect turbidity | Turbidimeter out of calibration | Calibrate accurately. Use correct standards. Make sure standards have not expired. |
| | Tubes not properly blanked | Once you blank a tube, do not rotate it while making suspension. Blank each suspension tube (tubes are not optically precise and can vary tube-to-tube and with rotation). |
| | Significant scratches or smudges on tube | Use suspension tubes only once, then discard. Visually inspect tubes before using. Wipe outside of tubes with a tissue before placing in turbidimeter. |
| Inaccurate %T measurements | Nonuniform suspension of cells | Use 7" swab to mix cells all the way to the bottom of the tube. The light path that looks through the tube is only viewing the bottom of the tube. |
| Trouble making a homogenous suspension | Mucoid or clumpy cells | For Spore Forming Gram positive Rods, Follow the recommendations on pages 4-10 to 4-12. |
| | | For Non-Spore formers: 1. Roll swab over young colonies in 3 rd and 4 th quadrants. 2. Mash colonies against side of a sterile tube (above meniscus). 3. Add several ml of inoculating fluid and wash the sides of the tube with a cotton swab. 4. Mix well, then examine for clumps. If clumps are present, allow them to settle. 5. Transfer this concentrated inoculum into a fresh Inoculating Fluid tube until the recommended turbidity is reached. * Be careful not to transfer any colony clumps into the Inoculating fluid. |
| | | 1. For difficult organisms such as <i>Gornia</i> or <i>Tsukamurella</i> , suspend colonies as a dense suspension in a small amount of inoculating fluid as above. 2. Incubate as needed at 35° C (do not exceed 10 minutes). Vortex (if clumps still present, pull off supernatant and repeat). 3. Bring supernatant volume up to 18 ml and mix until homogenous. |
| Too few positive reactions in MicroPlates | Detergent traces in tube | Use suspension tubes only once, then discard. |
| | Use of "bad" water | Use high quality purified water without preservatives. |

Inoculating MicroPlates

| Symptom | Cause | Solution |
|---------------------------|---------------------------|---|
| Pipetting problems | Wells inoculated unevenly | <p>Dispense first aliquot from pipettor back into reservoir.</p> <p>Depress lever smoothly and completely.</p> <p>Make sure tips are lined up evenly and tightly seated.</p> <p>Avoid pipetting or trapping air bubbles.</p> <p>Visually check wells after inoculating.</p> |
| | Incorrect volumes used | 100 µl for anaerobic bacteria, yeasts and fungal. |
| | | 150 µl for aerobic bacteria. |

Incubating MicroPlates

| Symptom | Cause | Solution |
|---|---|--|
| Poor growth or poor pattern formation in MicroPlates | Wrong incubation conditions | Incubate according to specified temperature and conditions. See Section 4 and Appendices. |
| | Wrong MicroPlate | Incubate the MicroPlate at the same temperature as the growth plate. |
| | Organism not in Database | Make sure incubator humidity is sufficient. |
| Wells at the four corners or edges of the MicroPlate have lower liquid levels after incubation | Incubator is too dry | Add a pan of water to humidify incubator, or incubate MicroPlates in a plastic container humidified with moist paper towels. |
| | Filamentous fungi not incubated in a closed container | Make sure you have the FF MicroPlates in a closed container or in plastic bags on trays. Do not add additional moisture when incubating filamentous fungi. |

MicroStation Reader

| Symptom | Cause | Solution |
|---|---|--|
| Erratic or inaccurate reading | Moisture, scratches, or smudges on MicroPlate | Remove MicroPlate lid before reading. Wipe bottom of MicroPlate before putting into reader. |
| | Many borderline reactions on MicroPlate | Review all sample preparation procedures for correctness and accuracy. |
| | Reading artifacts | Always double-check results by eye. Verify that your visual reading coincides with reader results. Use the Histogram software function to manually adjust threshold if necessary. |
| Software won't communicate with or initialize reader | Wrong com port selected | Choose correct com port. |
| | Loose cable connection | Turn reader off. Unplug cable, then plug it back in. Turn reader on and try again. |
| | Unknown cause | Turn reader off and exit program. Restart program, turn reader on, and try again. The reader has its own error messages, which should be self-explanatory. Call Biolog Technical Service if you need further assistance. |
| | Reader Bulb is Out/Low | Use the System Check test in the Bio-Tek reader manual. (Use Air reference test in the Molecular Devices reader manual.) |
| MicroPlate won't go into reader | MicroPlate mispositioned | Make sure MicroPlate snaps into place and is seated levelly. Make sure MicroPlate lid is off before reading. Remove MicroPlate and reposition in reader tray. Make sure A1 well is at the top left. Verify that software is communicating with reader. |

12. Glossary

In this section:
➔ **Definitions of
Terms**

Anamorph

Asexual state for yeasts and filamentous fungi.

Ascospore

Sexual spore formed within an ascus in fungi of the division Ascomycota. There are usually 4 or 8 ascospores per ascus (sometimes 2 or multiples of 4).

Ascus

Cell in the shape of a sac (round, club-shaped or cylindrical) within which ascospores are produced.

Aseptic technique

Standard lab procedures used to prevent contamination.

Basidiospore

A sexual spore formed on the upper surface of a basidium in fungi of the division Basidiomycota. There are usually 4 basidiospores, sometimes 2, and rarely multiples of 4.

Basidium

A club shaped cell upon which basidiospores form.

Carbon source utilization

Basic process used to identify microbes based on the chemicals they can utilize.

Catalase test

Additional test used to characterize gram-positive bacteria.

Cluster analysis

Mathematical depiction of the relationship between closely related microbe identification patterns. Shows dendrograms and two- and three-dimensional diagrams.

Conidium

A unicellular or multicellular fungal element specialized to detach from the mycelium and disseminate, thus serving as an asexual reproductive structure.

Dendrogram

Cluster diagram in the form of a branching tree.

Endpoint ID (EID)

Biolog, Inc. developed pattern matching method which considers the daily endpoint determination.

Enteric

Gram-negative bacteria belonging to the group *Enterobacteriaceae*.

Freeze

The act of converting a database into a Read Only and non-editable format.

Gram negative

Bacteria showing typical pink or red reaction on Gram stains.

Gram positive

Bacteria show typical blue or violet reaction on Gram stains.

Histogram

Visual representation of MicroPlate color distribution and thresholds.

Hypha

Septate or aseptate filaments of a fungus.

Inocula

Cell suspension used to inoculate MicroPlates.

Inoculating fluid

Fluid used to prepare inocula.

Maltose

A sugar added to BUG Agar. Required for culturing spore-forming gram-positive rods (e.g., *Bacillus* species).

Manual mode

Software mode used to record MicroPlate reactions visually.

MicroStation/MicroLog software

Windows-based program for microbe identification.

MicroPlate

Plate with 95 carbon source utilization tests (one in each well), with A1 well as control. The MicroPlate for yeasts has 94 tests and two control wells (A1 and D1).

MicroPlate reactions

Positive, negative, and borderline color responses used to identify microbes.

MicroStation System

Comprehensive system for microbe identification, including computer, MicroStation/MicroLog software, MicroStation reader, pipettor, turbidimeter, and ancillary materials.

Non-enteric

Gram-negative bacteria not belonging to the group *Enterobacteriaceae*.

Oxidase test

Additional test used to characterize gram-negative bacteria.

Pattern

Color responses in MicroPlates.

Pleomorphic

Having various distinct forms or shapes exhibited by a single strain or species

Progressive ID (PID)

Biolog, Inc. developed pattern matching method which considers the progressive sequence in which purple wells are formed.

Pure culture

Culture containing only one microbe species.

Reader mode

Software mode used to register MicroPlate reactions from MicroStation Reader.

Read Only File

The MicroStation/MicroLog Program will prevent editing of designated Read Only files

Restricted Access mode

Mode set by system administrator to prevent unregistered users from using system; assign user names, passwords and privileges to each user; create an audit trail, and freeze data files to maintain data integrity.

Salicylate

Anticapsulate agent required as an addition to inoculating fluid for a few gram-positive species.

Sporangiospore

An asexual spore produced by cytoplasmic cleavage within a sporangium.

Spore, fungal

A fungal propagative element, produced either as part of a sexual process (ascospore, basidiospore, zygosporium), or by an asexual reproductive process involving a process of cytoplasmic cleavage (sporangiosporium).

Sterile, fungi

Used to describe a fungal culture which produces no spores or conidia.

Teleomorph

Sexual state for yeasts and filamentous fungi.

Thioglycolate

Anticapsulate agent required as an addition to inoculating fluid for gram-negative enteric bacteria, gram-negative fastidious bacteria, and gram-positive cocci and rods.

Thresholds

The optical boundaries between negative, borderline, and positive reactions.

TSI slant

Additional test used to characterize gram-negative bacteria.

Turbidity

Measurement of cloudiness, which is indicative of inocula cell densities.

Unrestricted Access mode

Mode set by system administrator to allow any user to use all software functions.

Zygosporium

A sexual spore formed from the fusion of two similar cells called gametangia; this spore is characteristic of fungi in the division Zygomycota.

13. Appendices


In this section:

→ Plate Data
Statistics
→ Setup
Flowcharts
→ Database
Species Lists and
Their
Characteristics
→ Program
Printouts

Appendix 1: Data Statistics For Various Plate Types

| Plate Type | General Read Statistics | ID specific limitations: | |
|--------------------|---|--|--|
| | | Minimum Number of positives to call an ID | Maximum number of borderline Reactions |
| GN, GP, MT or Eco: | -Dual wavelength color statistic (590, 750) -A1-zeroed OD 750 turbidity statistic -Progressive identification algorithm -0.75, 0.5 similarity cut-offs for 4-6 hour, and 24 hours | 4-6 hour read 3 24 hour read 3 | 25 Borderline Reactions |
| AN: | -Dual wavelength color statistic (590, 750) -A1-zeroed OD 750 turbidity statistic -Progressive identification algorithm -0.5 similarity cut-offs for 20-24 hours | 20-24 hour read 7 | 25 Borderline Reactions |
| YT or Other: | -Single wavelength color statistic (590) -End Point identification algorithm -0.75, 0.5, 0.5 similarity cut-offs for 24 hours, 48 hours, and 72 hours | 24 hour read 3 48 hour read 3 72 hour read 3 | 25 Borderline Reactions |
| FF: | -Dual wavelength color statistic (490, 750) -A1-zeroed OD 750 turbidity statistic -End-point identification algorithm -0.9, 0.7, 0.65, 0.6, 0.6 similarity cut-offs for 1,2,3,4, and 7 days (7 day reading is present, however not included in the identification procedure) | 24 hour 30 48 hour 25 72 hour 20 96 hour 15 168 hour 5 | 50 Borderline Reactions |

Appendix 2: Release 4.2 Setup Flowchart (non-FF)



^aWith the exception of thermophiles, all species are either 26°C, 30° or 35°-37°C.

^bRequiring CHOC or CO₂ for growth or forming <1 mm colonies on BUG+B.

^cAgricultural species may be grown on BUG without blood.

^dAir, unless organism requires CO₂ for growth.

^eAN MicroPlate must be incubated in a hydrogen-free anaerobic atmosphere

^fIf control well (A1) is positive, add thioglycolate to inoculating fluid.

^gIf control well (A1) is positive, add 1 ml salicylate along with thioglycolate to inoculating fluid.

ABBREVIATIONS

T= Thioglycolate

Oxi = Oxidase; TSI = Triple Sugar Iron Slant; A = Acid, A^w = Weak Acid, K = Alkaline

Appendix 3: FF Setup Flowchart

Sample Preparation Process for FF MicroPlate

| Initial culture medium | 2% ME | |
|---|-------------------|----------------|
| | | |
| Wet prep results | Hyphal elements | Yeast cells |
| | | |
| Microbe type | Filamentous Fungi | Yeast |
| Culture medium | 2% ME | 2% ME |
| Temperature | 26° C | 26° C |
| Atmosphere | Air | Air |
| Culture Incubation time | 5 - 10 days | 2 days |
| Inoculating fluid | FF-IF | FF-IF |
| Inoculum turbidity/ Turbidity standard | 75% T FF | 75% T FF |
| MicroPlate type/ μl per well | FF 100 | FF 100 |
| Incubation time (hours) | 24, 48, 72, 96 | 24, 48, 72, 96 |

**Note: We recommend that filamentous fungi be manipulated in a biological safety cabinet.*

ABBREVIATIONS

| | | |
|-------|---|-----------------------|
| 2% ME | = | 2% Malt Extract Agar |
| FF | = | Filamentous Fungi |
| %T | = | percent transmittance |
| IF | = | inoculating fluid |

Notes:

Appendix 4: Database Species Lists and Their Characteristics

Gram-Negative Aerobic Bacteria

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|----|---|---------|------|------|--------|----------------------|-------|
| 1 | "Achromobacter cholinophagum" | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 2 | Achromobacter xylosoxidans ss denitrificans | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 3 | Achromobacter xylosoxidans ss xylosoxidans | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 4 | Acidovorax avenae ss avenae | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 5 | Acidovorax avenae ss cattleyae | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 6 | Acidovorax avenae ss citruli | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 7 | Acidovorax delafieldii | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 8 | Acidovorax facilis | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 9 | Acidovorax konjaci | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 10 | Acidovorax temperans | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 11 | Acinetobacter baumannii/genospecies 2 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 12 | Acinetobacter calcoaceticus bv alc | GN-NENT | O- | N | BUG+B | Air | 30 |
| 13 | Acinetobacter calcoaceticus/genospecies 1 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 14 | Acinetobacter calcoaceticus/genospecies 3 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 15 | Acinetobacter calcoaceticus/genospecies 13 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 16 | Acinetobacter genospecies 6 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 17 | Acinetobacter genospecies 10 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 18 | Acinetobacter genospecies 11 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 19 | Acinetobacter genospecies 14 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 20 | Acinetobacter genospecies 15 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 21 | Acinetobacter haemolyticus/genospecies 4 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 22 | Acinetobacter johnsonii/genospecies 7 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 23 | Acinetobacter junii/genospecies 5 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 24 | Acinetobacter lwoffii/genospecies 8/9 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 25 | Acinetobacter radioresistens/genospecies 12 | GN-NENT | O- | N | BUG+B | Air | 30 |
| 26 | Actinobacillus capsulatus | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 27 | Actinobacillus equuli | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 28 | Actinobacillus hominis | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 29 | Actinobacillus indolicus | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 30 | Actinobacillus lignieresii | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 31 | Actinobacillus minor | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 32 | Actinobacillus muris | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 33 | Actinobacillus pleuropneumoniae | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 34 | Actinobacillus porcinus | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 35 | Actinobacillus rossii | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 36 | Actinobacillus salpingitidis | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 37 | Actinobacillus seminis | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 38 | Actinobacillus suis | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 39 | Actinobacillus ureae | GN-FAS | O+/- | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 40 | Aeromonas allosaccharophila | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 41 | Aeromonas caviae DNA group 4 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 42 | Aeromonas encheleia | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 43 | Aeromonas enteropelogenes | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 44 | Aeromonas eucrenophila DNA group 6 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 45 | Aeromonas hydrophila DNA group 1 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 46 | Aeromonas hydrophila-like DNA group 2 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 47 | Aeromonas hydrophila-like DNA group 3 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 48 | Aeromonas ichthiosmia | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 49 | Aeromonas jandaei DNA group 9 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 50 | Aeromonas media DNA group 5b | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 51 | Aeromonas media-like DNA group 5a | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 52 | Aeromonas salmonicida ss achromogenes | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 53 | Aeromonas salmonicida ss masoucida | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 54 | Aeromonas salmonicida ss salmonicida | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 55 | Aeromonas schubertii DNA group 12 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 56 | Aeromonas sobria DNA group 7 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 57 | Aeromonas trota DNA group 13 | GN-NENT | O+ | N | BUG+B | Air | 30 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|---------|------|------|-----------|---------|-------|
| 58 | <i>Aeromonas veronii</i> DNA group 10 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 59 | <i>Aeromonas veronii/sobria</i> DNA group 8 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 60 | <i>Alcaligenes faecalis</i> ss faecalis | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 61 | <i>Alysiella filiformis</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 62 | <i>Aminobacter aminovorans</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 63 | <i>Ancylobacter aquaticus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 64 | <i>Aquaspirillum autotrophicum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 65 | <i>Aquaspirillum dispar</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 66 | <i>Aquaspirillum metamorphum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 67 | <i>Aquaspirillum peregrinum</i> ss integrum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 68 | <i>Aquaspirillum peregrinum</i> ss peregrinum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 69 | <i>Aquaspirillum putridiconchylum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 70 | <i>Bergeyella zoohelcum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 71 | <i>Bordetella avium</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 72 | <i>Bordetella bronchiseptica</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 73 | <i>Bordetella hinzii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 74 | <i>Bordetella holmesii</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 75 | <i>Bordetella parapertussis</i> | GN-NENT | O- | N | BUG+B | Air | 35-37 |
| 76 | <i>Bordetella pertussis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 77 | <i>Bordetella trematum</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 78 | <i>Bordetella-like species</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 79 | <i>Brenneria rubrifaciens</i> | GN-ENT | O- | N | BUG+B† | Air | 30 |
| 80 | <i>Brenneria salicis</i> | GN-ENT | O- | N | BUG+B† | Air | 30 |
| 81 | <i>Brevundimonas diminuta</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 82 | <i>Brevundimonas vesicularis</i> | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 83 | <i>Brucella melitensis</i> (abortus) | GN-FAS | O+ | N | CHOC** | 6.5%CO2 | 35-37 |
| 84 | <i>Brucella melitensis</i> (abortus, melitensis) | GN-FAS | O+ | N | CHOC** | 6.5%CO2 | 35-37 |
| 85 | <i>Brucella melitensis</i> (melitensis) | GN-FAS | O+ | N | CHOC** | 6.5%CO2 | 35-37 |
| 86 | <i>Brucella melitensis</i> (suis, canis) | GN-FAS | O+ | N | CHOC** | 6.5%CO2 | 35-37 |
| 87 | <i>Brucella melitensis</i> (suis) | GN-FAS | O+ | N | CHOC** | 6.5%CO2 | 35-37 |
| 88 | <i>Budvicia aquatica</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 89 | <i>Burkholderia andropogonis</i> | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 90 | <i>Burkholderia caryophylli</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 91 | <i>Burkholderia cepacia</i> | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 92 | <i>Burkholderia gladioli</i> | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 93 | <i>Burkholderia glathei</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 94 | <i>Burkholderia glumae</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 95 | <i>Burkholderia mallei</i> | GN-NENT | O+/- | N | BUG+B†*** | Air | 30 |
| 96 | <i>Burkholderia multivorans</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 97 | <i>Burkholderia phenazinium</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 98 | <i>Burkholderia plantarii</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 99 | <i>Burkholderia pseudomallei</i> | GN-NENT | O+ | N | BUG+B†*** | Air | 30 |
| 100 | <i>Burkholderia pyrrocinia</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 101 | <i>Burkholderia vietnamiensis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 102 | <i>Buttiauxella agrestis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 103 | <i>Buttiauxella brennerae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 104 | <i>Buttiauxella ferragutiae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 105 | <i>Buttiauxella gaviniae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 106 | <i>Buttiauxella izardii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 107 | <i>Buttiauxella noackiae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 108 | <i>Buttiauxella warmboldiae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 109 | CDC group DF-3 (Capnocytophaga) | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 110 | CDC group EF-4 (Neisseria) | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 111 | CDC group EO-2 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 112 | CDC group II-E subgroup A | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 113 | CDC group II-E subgroup B | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 114 | CDC group II-H | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 115 | <i>Capnocytophaga canimorsus</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 116 | <i>Capnocytophaga cynodegmi</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 117 | <i>Capnocytophaga gingivalis</i> | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 118 | <i>Capnocytophaga granulosa</i> | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 119 | <i>Capnocytophaga haemolytica</i> | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 120 | <i>Capnocytophaga ochracea/sputigena</i> | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 121 | <i>Cardiobacterium hominis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 122 | <i>Cedecea davisae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 123 | <i>Cedecea lapagei</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 124 | <i>Cedecea neteri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 125 | <i>Chromobacterium violaceum</i> | GN-NENT | O+/- | N | BUG+B | Air | 30 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|---------|------|------|--------|---------|-------|
| 126 | <i>Chryseobacterium balustinum</i> | GN-NENT | O+ | Y | BUG+B | Air | 30 |
| 127 | <i>Chryseobacterium gleum</i> /indologenes | GN-NENT | O+ | Y | BUG+B | Air | 30 |
| 128 | <i>Chryseobacterium indoltheticum</i> | GN-NENT | O+ | Y | BUG+B | Air | 30 |
| 129 | <i>Chryseobacterium meningosepticum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 130 | <i>Chryseobacterium scophthalmum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 131 | <i>Chryseomonas luteola</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 132 | <i>Citrobacter amalonaticus</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 133 | <i>Citrobacter braakii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 134 | <i>Citrobacter farmeri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 135 | <i>Citrobacter freundii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 136 | <i>Citrobacter gillenii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 137 | <i>Citrobacter koseri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 138 | <i>Citrobacter murlinae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 139 | <i>Citrobacter rodentium</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 140 | <i>Citrobacter sedlakii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 141 | <i>Citrobacter werkmanii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 142 | <i>Citrobacter youngae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 143 | <i>Comamonas terrigena</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 144 | <i>Comamonas testosteroni</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 145 | <i>Cytophaga fermentans</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 146 | <i>Delftia acidovorans</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 147 | <i>Edwardsiella hoshinae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 148 | <i>Edwardsiella ictaluri</i> | GN-ENT | O- | Y | BUG+B | Air | 30 |
| 149 | <i>Edwardsiella tarda</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 150 | <i>Eikenella corrodens</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 151 | <i>Empedobacter brevis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 152 | <i>Enterobacter aerogenes</i> (<i>Klebsiella mobilis</i>) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 153 | " <i>Enterobacter agglomerans</i> " bgp 2 | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 154 | " <i>Enterobacter agglomerans</i> " bgp 3 | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 155 | " <i>Enterobacter agglomerans</i> " bgp 4 (<i>Pantoea</i>) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 156 | " <i>Enterobacter agglomerans</i> " bgp 5 | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 157 | " <i>Enterobacter agglomerans</i> " bgp 6 (<i>Pectobacterium</i>) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 158 | " <i>Enterobacter agglomerans</i> " bgp 7 (<i>Pantoea</i>) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 159 | <i>Enterobacter amnigenus</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 160 | <i>Enterobacter asburiae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 161 | <i>Enterobacter cancerogenus</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 162 | <i>Enterobacter cloacae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 163 | <i>Enterobacter gergoviae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 164 | <i>Enterobacter hormaechei</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 165 | <i>Enterobacter intermedius</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 166 | <i>Enterobacter nimipressuralis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 167 | <i>Enterobacter sakazakii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 168 | <i>Erwinia amylovora</i> | GN-ENT | O- | Y | BUG+B† | Air | 30 |
| 169 | <i>Escherichia blattae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 170 | <i>Escherichia coli</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 171 | <i>Escherichia coli</i> (USP5-7085) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 172 | <i>Escherichia coli</i> inactive | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 173 | <i>Escherichia coli</i> O157:H7 | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 174 | <i>Escherichia fergusonii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 175 | <i>Escherichia hermannii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 176 | <i>Escherichia vulneris</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 177 | <i>Ewingella americana</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 178 | <i>Flavimonas oryzae</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 179 | <i>Flavobacterium ferrugineum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 180 | <i>Flavobacterium flevense</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 181 | <i>Flavobacterium hydati</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 182 | <i>Flavobacterium johnsoniae</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 183 | <i>Flavobacterium mizutaii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 184 | <i>Flavobacterium mizutaii</i> -like (CDC group II-I) | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 185 | <i>Flavobacterium resinovorum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 186 | " <i>Flavobacterium tirrenicum</i> " (<i>Chryseobacterium</i>) | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 187 | <i>Francisella philomiragia</i> | GN-FAS | O-/+ | N | CHOC | 6.5%CO2 | 35-37 |
| 188 | <i>Francisella tularensis</i> | GN-FAS | O- | N | CHOC** | 6.5%CO2 | 35-37 |
| 189 | <i>Gilardi unnamed rod group 1</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 190 | <i>Haemophilus actinomycetemcomitans</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 191 | <i>Haemophilus aegyptius</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 192 | <i>Haemophilus aphrophilus</i> | GN-FAS | O- | Y | CHOC | 6.5%CO2 | 35-37 |
| 193 | <i>Haemophilus avium</i> (<i>Pasteurella</i>) | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|---------|------|------|--------|----------|-------|
| 194 | <i>Haemophilus ducreyi</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 195 | <i>Haemophilus haemoglobinophilus</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 196 | <i>Haemophilus haemolyticus</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 197 | <i>Haemophilus influenzae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 198 | <i>Haemophilus paracuniculus</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 199 | <i>Haemophilus paragallinarum</i> | GN-FAS | O- | Y | CHOC | 6.5% CO2 | 35-37 |
| 200 | <i>Haemophilus parahaemolyticus</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 201 | <i>Haemophilus parainfluenzae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 202 | <i>Haemophilus paraphrohaemolyticus</i> | GN-FAS | O+/- | Y | CHOC | 6.5% CO2 | 35-37 |
| 203 | <i>Haemophilus paraphrophilus</i> | GN-FAS | O+/- | Y | CHOC | 6.5% CO2 | 35-37 |
| 204 | <i>Haemophilus parasuis</i> | GN-FAS | O- | Y | CHOC | 6.5% CO2 | 35-37 |
| 205 | <i>Haemophilus segnis</i> | GN-FAS | O- | Y | CHOC | 6.5% CO2 | 35-37 |
| 206 | <i>Haemophilus somnus</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 207 | <i>Hafnia alvei</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 208 | <i>Herbaspirillum rubrisubalbicans</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 209 | <i>Herbaspirillum seropedicae</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 210 | <i>Hydrogenophaga flava</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 211 | <i>Hydrogenophaga palleronii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 212 | <i>Hydrogenophaga pseudoflava</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 213 | <i>Hydrogenophaga taeniospiralis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 214 | <i>Iodobacter fluviatilis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 215 | <i>Janthinobacterium lividum</i> | GN-NENT | O+/- | N | BUG+B | Air | 26 |
| 216 | <i>Kingella denitrificans</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 217 | <i>Kingella kingae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 218 | <i>Kingella oralis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 219 | <i>Klebsiella oxytoca</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 220 | <i>Klebsiella pneumoniae</i> ss <i>ozaenae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 221 | <i>Klebsiella pneumoniae</i> ss <i>pneumoniae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 222 | <i>Klebsiella pneumoniae</i> ss <i>rhinoscleromatis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 223 | <i>Kluyvera ascorbata</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 224 | <i>Kluyvera cochleae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 225 | <i>Kluyvera cryocrescens</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 226 | <i>Kluyvera georgiana</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 227 | <i>Lampropedia hyalina</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 228 | <i>Leclercia adecarboxylata</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 229 | <i>Leminorella grimonitii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 230 | <i>Leminorella richardii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 231 | <i>Listonella anguillarum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 232 | <i>Listonella pelagia</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 233 | <i>Mannheimia granulomatis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 234 | <i>Mannheimia haemolytica</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 235 | <i>Moellerella wisconsensis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 236 | <i>Moraxella bovis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 237 | <i>Moraxella canis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 238 | <i>Moraxella caprae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 239 | <i>Moraxella catarrhalis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 240 | <i>Moraxella caviae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 241 | <i>Moraxella cuniculi</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 242 | <i>Moraxella equi</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 243 | <i>Moraxella lacunata</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 244 | <i>Moraxella nonliquefaciens</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 245 | <i>Moraxella osloensis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 246 | <i>Moraxella ovis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 247 | <i>Morganella morganii</i> ss <i>morganii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 248 | <i>Myroides odoratimimus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 249 | <i>Myroides odoratus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 250 | <i>Neisseria animalis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 251 | <i>Neisseria canis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 252 | <i>Neisseria cinerea</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 253 | <i>Neisseria denitrificans</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 254 | <i>Neisseria elongata</i> ss <i>elongata</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 255 | <i>Neisseria flava</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 256 | <i>Neisseria flavescens</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 257 | <i>Neisseria gonorrhoeae</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 258 | <i>Neisseria lactamica</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 259 | <i>Neisseria meningitidis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 260 | <i>Neisseria mucosa</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 261 | <i>Neisseria perflava</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|---------|------|------|--------|---------|-------|
| 262 | <i>Neisseria sicca</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 263 | <i>Neisseria subflava</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 264 | <i>Neisseria weaveri</i> | GN-FAS | O+ | Y | CHOC | 6.5%CO2 | 35-37 |
| 265 | <i>Obesumbacterium proteus</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 266 | <i>Obesumbacterium proteus</i> biogroup 2 | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 267 | <i>Oceanomonas doudoroffii</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 268 | <i>Ochrobactrum anthropi</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 269 | <i>Oligella ureolytica</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 270 | <i>Oligella urethralis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 271 | <i>Ornithobacterium rhinotracheale</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 272 | <i>Pandoraea norimbergensis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 273 | <i>Pantoea agglomerans</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 274 | <i>Pantoea citrea</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 275 | <i>Pantoea dispersa</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 276 | <i>Pantoea punctata</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 277 | <i>Pantoea stewartii</i> ss <i>stewartii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 278 | <i>Pantoea terreia</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 279 | " <i>Pasteurella</i> " <i>aerogenes</i> | GN-NENT | O+/- | N | BUG+B | Air | 35-37 |
| 280 | <i>Pasteurella anatis</i> | GN-NENT | O+/- | N | BUG+B | Air | 35-37 |
| 281 | <i>Pasteurella bettyae</i> | GN-NENT | O- | N | BUG+B | Air | 35-37 |
| 282 | <i>Pasteurella caballi</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 283 | <i>Pasteurella canis/stomatis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 284 | <i>Pasteurella dagmatis</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 285 | <i>Pasteurella gallinarum</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 286 | <i>Pasteurella langaaensis</i> | GN-NENT | O+/- | N | BUG+B | Air | 35-37 |
| 287 | <i>Pasteurella lymphangitidis</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 288 | <i>Pasteurella mairii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 289 | <i>Pasteurella multocida</i> ss <i>multocida</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 290 | <i>Pasteurella pneumotropica</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 291 | " <i>Pasteurella</i> " <i>testudinis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 292 | <i>Pasteurella trehalosi</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 293 | <i>Pasteurella volantium</i> | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 294 | <i>Paucimonas lemoignei</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 295 | <i>Pectobacterium carotovorum</i> ss <i>atrosepticum</i> | GN-ENT | O- | Y | BUG+B† | Air | 35-37 |
| 296 | <i>Pectobacterium carotovorum</i> ss <i>betavascolorum</i> | GN-ENT | O- | Y | BUG+B† | Air | 35-37 |
| 297 | <i>Pectobacterium carotovorum</i> ss <i>carotovorum</i> | GN-ENT | O- | Y | BUG+B† | Air | 35-37 |
| 298 | <i>Pectobacterium chrysanthemi</i> | GN-ENT | O- | Y | BUG+B† | Air | 35-37 |
| 299 | <i>Pectobacterium cypripedii</i> | GN-ENT | O- | Y | BUG+B† | Air | 35-37 |
| 300 | <i>Pedobacter heparinus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 301 | <i>Photobacterium angustum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 302 | <i>Photobacterium damsela</i> ss <i>damselae</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 303 | <i>Photobacterium fischeri</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 304 | <i>Photobacterium leiognathi</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 305 | <i>Photobacterium luminescens</i> ss <i>luminescens</i> | GN-ENT | O- | N | BUG+B | Air | 30 |
| 306 | <i>Phyllobacterium myrsinacearum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 307 | <i>Phyllobacterium rubiacearum</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 308 | <i>Plesiomonas shigelloides</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 309 | <i>Pragia fontium</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 310 | <i>Proteus mirabilis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 311 | <i>Proteus myxofaciens</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 312 | <i>Proteus penneri/vulgaris</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 313 | <i>Providencia alcalifaciens</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 314 | <i>Providencia heimbachae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 315 | <i>Providencia rettgeri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 316 | <i>Providencia rustigianii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 317 | <i>Providencia stuartii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 318 | <i>Pseudomonas aeruginosa</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 319 | <i>Pseudomonas agarici</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 320 | <i>Pseudomonas alcaligenes</i> | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 321 | <i>Pseudomonas asplenii</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 322 | <i>Pseudomonas aurantiaca</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 323 | " <i>Pseudomonas</i> " <i>bathycetes</i> " | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 324 | " <i>Pseudomonas</i> " <i>boreopolis</i> (<i>Stenotrophomonas</i>) | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 325 | <i>Pseudomonas caricapapayae</i> (<i>syringae</i>) | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 326 | <i>Pseudomonas chlororaphis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 327 | <i>Pseudomonas cichorii</i> (<i>syringae</i>) | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 328 | " <i>Pseudomonas</i> " <i>cissicola</i> (<i>Xanthomonas</i> -like) | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 329 | <i>Pseudomonas citronellolis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|---------|------|------|--------|-----|-------|
| 330 | <i>Pseudomonas corrugata</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 331 | " <i>Pseudomonas floridana</i> " (Burkholderia-like) | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 332 | <i>Pseudomonas fluorescens</i> | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 333 | <i>Pseudomonas fluorescens</i> biotype A | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 334 | <i>Pseudomonas fluorescens</i> biotype C | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 335 | <i>Pseudomonas fluorescens</i> biotype F | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 336 | <i>Pseudomonas fluorescens</i> biotype G | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 337 | <i>Pseudomonas fragi</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 338 | <i>Pseudomonas fulva</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 339 | <i>Pseudomonas fuscovaginae</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 340 | " <i>Pseudomonas</i> group 2" (Burkholderia-like) | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 341 | " <i>Pseudomonas</i> " <i>huttiensis</i> (Burkholderia-like) | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 342 | <i>Pseudomonas lundensis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 343 | <i>Pseudomonas</i> "maculicola" | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 344 | <i>Pseudomonas marginalis</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 345 | <i>Pseudomonas mendocina</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 346 | " <i>Pseudomonas mephitica</i> " | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 347 | <i>Pseudomonas mucidolens</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 348 | <i>Pseudomonas nitroreducens/azelaica</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 349 | <i>Pseudomonas oleovorans</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 350 | <i>Pseudomonas pertucinogena</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 351 | <i>Pseudomonas pseudoalcaligenes</i> ss <i>pseudoalcaligenes</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 352 | <i>Pseudomonas putida</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 353 | <i>Pseudomonas putida</i> biotype A | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 354 | <i>Pseudomonas putida</i> biotype B | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 355 | <i>Pseudomonas resinovorans</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 356 | <i>Pseudomonas savastanoi</i> pv <i>fraxini</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 357 | <i>Pseudomonas savastanoi</i> pv <i>glycinea</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 358 | <i>Pseudomonas savastanoi</i> pv <i>nerii</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 359 | <i>Pseudomonas spinosa</i> (Burkholderia) | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 360 | <i>Pseudomonas straminea</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 361 | <i>Pseudomonas stutzeri</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 362 | <i>Pseudomonas synxantha</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 363 | <i>Pseudomonas syringae</i> pv <i>aceris</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 364 | <i>Pseudomonas syringae</i> pv <i>antirrhini</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 365 | <i>Pseudomonas syringae</i> pv <i>apii</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 366 | <i>Pseudomonas syringae</i> pv <i>aptata</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 367 | <i>Pseudomonas syringae</i> pv <i>atrofaciens</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 368 | <i>Pseudomonas syringae</i> pv <i>coronafaciens</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 369 | <i>Pseudomonas syringae</i> pv <i>cunninghamiae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 370 | <i>Pseudomonas syringae</i> pv <i>delphinii</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 371 | <i>Pseudomonas syringae</i> pv <i>erobotryae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 372 | <i>Pseudomonas syringae</i> pv <i>glycinea</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 373 | <i>Pseudomonas syringae</i> pv <i>helianthi</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 374 | <i>Pseudomonas syringae</i> pv <i>lachrymans</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 375 | <i>Pseudomonas syringae</i> pv <i>mori</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 376 | <i>Pseudomonas syringae</i> pv <i>myricae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 377 | <i>Pseudomonas syringae</i> pv <i>oryzae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 378 | <i>Pseudomonas syringae</i> pv <i>papulans</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 379 | <i>Pseudomonas syringae</i> pv <i>persicae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 380 | <i>Pseudomonas syringae</i> pv <i>phaseolicola</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 381 | <i>Pseudomonas syringae</i> pv <i>pisi</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 382 | <i>Pseudomonas syringae</i> pv <i>porri</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 383 | <i>Pseudomonas syringae</i> pv <i>primulae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 384 | <i>Pseudomonas syringae</i> pv <i>sesami</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 385 | <i>Pseudomonas syringae</i> pv <i>syringae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 386 | <i>Pseudomonas syringae</i> pv <i>tabaci</i> A | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 387 | <i>Pseudomonas syringae</i> pv <i>tabaci</i> B | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 388 | <i>Pseudomonas syringae</i> pv <i>tagetis</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 389 | <i>Pseudomonas syringae</i> pv <i>tomato</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 390 | <i>Pseudomonas syringae</i> pv <i>zinniae</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 391 | <i>Pseudomonas taetrolens</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 392 | <i>Pseudomonas tolaasii</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 393 | <i>Pseudomonas viridiflava</i> (syringae) | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 394 | <i>Pseudomonas viridilivida</i> | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 395 | <i>Psychrobacter immobilis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 396 | <i>Psychrobacter phenylpyruvicus</i> | GN-NENT | O+ | Y | BUG+B | Air | 35-37 |
| 397 | <i>Rahnella aquatilis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|---------|------|------|--------|----------|-------|
| 398 | Ralstonia eutropha | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 399 | Ralstonia pauca | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 400 | Ralstonia pickettii | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 401 | Ralstonia solanacearum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 402 | Raoultella planticola | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 403 | Raoultella planticola/ornithinolytica | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 404 | Raoultella terrigena | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 405 | Rhizobium like-Cystic Fibrosis | GN-NENT | O+ | N | BUG+B† | Air | 35-37 |
| 406 | Rhizobium radiobacter | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 407 | Rhizobium rhizogenes | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 408 | Rhizobium vitis | GN-NENT | O+/- | N | BUG+B† | Air | 30 |
| 409 | Riemerella anatipestifer | GN-NENT | O+ | N | BUG+B | Air | 35-37 |
| 410 | Roseomonas cervicalis | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 411 | Roseomonas fauriae | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 412 | Roseomonas genomospecies 4 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 413 | Roseomonas genomospecies 5 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 414 | Roseomonas genomospecies 6 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 415 | Roseomonas gilardii | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 416 | Salmonella gp 1 (choleraesuis) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 417 | Salmonella gp 1 (choleraesuis) st choleraesuis | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 418 | Salmonella gp 1 (choleraesuis) st gallinarum | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 419 | Salmonella gp 1 (choleraesuis) st paratyphi A | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 420 | Salmonella gp 1 (choleraesuis) st pullorum | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 421 | Salmonella gp 1 (choleraesuis) st typhi | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 422 | Salmonella gp 1 (choleraesuis) st typhimurium | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 423 | Salmonella gp 3a (arizonae) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 424 | Salmonella gp 3b (diarizonae) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 425 | Salmonella gp 4 (houtenae) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 426 | Salmonella gp 5 (bongori) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 427 | Salmonella gp 6 (indica) | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 428 | Serpens flexibilis | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 429 | Serratia entomophila | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 430 | Serratia ficaria | GN-ENT | O-/v | Y | BUG+B | Air | 35-37 |
| 431 | Serratia fonticola | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 432 | Serratia liquefaciens/grimesii | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 433 | Serratia marcescens | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 434 | Serratia odorifera | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 435 | Serratia plymuthica | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 436 | Serratia proteamaculans ss proteamaculans | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 437 | Serratia rubidaea | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 438 | Shewanella algae | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 439 | Shewanella putrefaciens A | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 440 | Shewanella putrefaciens B | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 441 | Shigella boydii | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 442 | Shigella dysenteriae | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 443 | Shigella flexneri | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 444 | Shigella sonnei | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 445 | Simonsiella crassa | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 446 | Simonsiella muelleri | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 447 | Simonsiella steedae | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 448 | Sinorhizobium meliloti | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 449 | Sphingobacterium multivorum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 450 | Sphingobacterium multivorum-like | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 451 | Sphingobacterium spiritivorum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 452 | Sphingobacterium thalpophilum | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 453 | Sphingomonas adhaesiva | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 454 | Sphingomonas capsulata | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 455 | Sphingomonas echinoides | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 456 | Sphingomonas macrogoltabidus | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 457 | Sphingomonas parapaucimobilis | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 458 | Sphingomonas paucimobilis A | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 459 | Sphingomonas paucimobilis B | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 460 | Sphingomonas sanguinis | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 461 | Sphingomonas terrae | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 462 | Sphingomonas yanoikuyae | GN-NENT | O+/- | N | BUG+B | Air | 30 |
| 463 | Stenotrophomonas maltophilia | GN-NENT | O- | N | BUG+B† | Air | 30 |
| 464 | Suttonella indologenes | GN-FAS | O+ | Y | CHOC | 6.5% CO2 | 35-37 |
| 465 | Tatumella ptyseos | GN-ENT | O- | Y | BUG+B | Air | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|---------|------|------|---------|----------------------|-------|
| 466 | <i>Taylorella equigenitalis</i> | GN-FAS | O+ | Y | CHOC | 6.5% CO ₂ | 35-37 |
| 467 | <i>Trabulsiella guamensis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 468 | <i>Variovorax paradoxus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 469 | <i>Vibrio aestuarianus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 470 | <i>Vibrio alginolyticus</i> | GN-NENT | O+ | N | BUG+B† | Air | 30 |
| 471 | <i>Vibrio campbelli</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 472 | <i>Vibrio carchariae</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 473 | <i>Vibrio cholerae</i> O1 (ATCC 25870) | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 474 | <i>Vibrio cholerae</i> O1/0139 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 475 | <i>Vibrio cholerae</i> non O1 | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 476 | <i>Vibrio cincinnatiensis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 477 | <i>Vibrio diazotrophicus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 478 | <i>Vibrio fluvialis</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 479 | <i>Vibrio furnissii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 480 | <i>Vibrio harveyi</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 481 | <i>Vibrio mediterranei</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 482 | <i>Vibrio metschnikovii</i> | GN-NENT | O- | N | BUG+B | Air | 30 |
| 483 | <i>Vibrio mimicus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 484 | <i>Vibrio natriegens</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 485 | <i>Vibrio ordalii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 486 | <i>Vibrio parahaemolyticus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 487 | <i>Vibrio proteolyticus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 488 | <i>Vibrio splendidus</i> | GN-NENT | O+ | N | BUG+B | Air | 26 |
| 489 | <i>Vibrio tubiashii</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 490 | <i>Vibrio vulnificus</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 491 | <i>Vogesella indigofera</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 492 | <i>Weeksella virosa</i> | GN-NENT | O+ | N | BUG+B | Air | 30 |
| 493 | <i>Xanthomonas albilineans</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 494 | <i>Xanthomonas campestris</i> pv <i>begoniae</i> A | GN-NENT | O- | N | BUG | Air | 30 |
| 495 | <i>Xanthomonas campestris</i> pv <i>begoniae</i> B | GN-NENT | O- | N | BUG | Air | 30 |
| 496 | <i>Xanthomonas campestris</i> pv <i>campestris</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 497 | <i>Xanthomonas campestris</i> pv <i>carotae</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 498 | <i>Xanthomonas campestris</i> pv <i>dieffenbachiae</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 499 | <i>Xanthomonas campestris</i> pv <i>hyacinthi</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 500 | <i>Xanthomonas campestris</i> pv <i>juglandis</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 501 | <i>Xanthomonas campestris</i> pv <i>malvacearum</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 502 | <i>Xanthomonas campestris</i> pv <i>nigromaculans</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 503 | <i>Xanthomonas campestris</i> pv <i>pelargonii</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 504 | <i>Xanthomonas campestris</i> pv <i>phaseoli</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 505 | <i>Xanthomonas campestris</i> pv <i>poinsetticola</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 506 | <i>Xanthomonas campestris</i> pv <i>raphani</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 507 | <i>Xanthomonas campestris</i> pv <i>tardicrescens</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 508 | <i>Xanthomonas campestris</i> pv <i>translucens</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 509 | <i>Xanthomonas campestris</i> pv <i>vesicatoria</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 510 | <i>Xanthomonas oryzae</i> pv <i>oryzicola</i> | GN-NENT | O- | N | BUG | Air | 30 |
| 511 | <i>Xenorhabdus bovienii</i> | GN-ENT | O- | N | BUG+B | Air | 30 |
| 512 | <i>Xenorhabdus nematophila</i> | GN-ENT | O- | N | BUG+B | Air | 30 |
| 513 | <i>Yersinia aldovae</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 514 | <i>Yersinia bercovieri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 515 | <i>Yersinia enterocolitica</i> ss <i>enterocolitica</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 516 | <i>Yersinia fredericksonii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 517 | <i>Yersinia intermedia</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 518 | <i>Yersinia kristensenii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 519 | <i>Yersinia mollaretii</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 520 | <i>Yersinia pestis</i> | GN-ENT | O- | N | BUG+B** | Air | 35-37 |
| 521 | <i>Yersinia pseudotuberculosis</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 522 | <i>Yersinia rohdei</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 523 | <i>Yersinia ruckeri</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |
| 524 | <i>Yokenella regensburgei</i> | GN-ENT | O- | Y | BUG+B | Air | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Gram-Positive Aerobic Bacteria

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|----|--|-----------|-------|------|---------|----------------------|-------|
| 1 | <i>Actinomyces bovis</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 2 | <i>Actinomyces canis</i> | GP-ROD | C+/- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 3 | <i>Actinomyces hordeovulneris</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 4 | <i>Actinomyces hyovaginalis</i> | GP-ROD | C-/W+ | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 5 | <i>Actinomyces naeslundii</i> | GP-ROD | C-/W+ | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 6 | <i>Actinomyces neuii ss anitratus</i> | GP-ROD | C+/- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 7 | <i>Actinomyces neuii ss neuii</i> | GP-ROD | C+/- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 8 | <i>Actinomyces odontolyticus</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 9 | <i>Actinomyces radingae/turicensis (CDC.E)</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 10 | <i>Actinomyces viscosus</i> | GP-ROD | C+ | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 11 | <i>Aerococcus christensenii</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 12 | <i>Aerococcus urinae</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 13 | <i>Aerococcus viridans</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 14 | <i>Alloioicoccus otitis</i> | GP-COCCUS | C-/W+ | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 15 | <i>Arcanobacterium bernardiae (CDC.2)</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 16 | <i>Arcanobacterium haemolyticum</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 17 | <i>Arcanobacterium pyogenes</i> | GP-ROD | C- | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 18 | <i>Arthrobacter cummingsii</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 19 | <i>Arthrobacter histidinolovorans</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 20 | <i>Arthrobacter ilicis</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 21 | <i>Arthrobacter woluwensis</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 22 | <i>Aureobacterium resistens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 23 | <i>Bacillus alcalophilus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 24 | <i>Bacillus amyloliquefaciens A</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 25 | <i>Bacillus amyloliquefaciens B</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 26 | <i>Bacillus anthracis subgroup A</i> | GP-ROD SB | C-/W+ | N | BUG+M** | Air | 30 |
| 27 | <i>Bacillus anthracis subgroup B</i> | GP-ROD SB | C-/W+ | N | BUG+M** | Air | 30 |
| 28 | <i>Bacillus anthracis subgroup C</i> | GP-ROD SB | C-/W+ | N | BUG+M** | Air | 30 |
| 29 | <i>Bacillus anthracis subgroup D</i> | GP-ROD SB | C-/W+ | N | BUG+M** | Air | 30 |
| 30 | <i>Bacillus badius</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 31 | <i>Bacillus cereus/thuringiensis A</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 32 | <i>Bacillus cereus/thuringiensis B</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 33 | <i>Bacillus cereus/thuringiensis C</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 34 | <i>Bacillus circulans</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 35 | <i>Bacillus coagulans</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 36 | <i>Bacillus fastidiosus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 37 | <i>Bacillus firmus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 38 | <i>Bacillus halodurans</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 39 | <i>Bacillus laevolacticus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 40 | <i>Bacillus licheniformis</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 41 | <i>Bacillus maroccanus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 42 | <i>Bacillus megaterium A</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 43 | <i>Bacillus megaterium B</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 44 | <i>Bacillus mycoides</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 45 | <i>Bacillus psychrosaccharolyticus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 46 | <i>Bacillus pumilus A</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 47 | <i>Bacillus pumilus B</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 48 | <i>Bacillus pumilus C</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 49 | <i>Bacillus racemilacticus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 50 | <i>Bacillus sphaericus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 51 | <i>Bacillus subtilis A</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 52 | <i>Bacillus subtilis B</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 53 | <i>Bacillus subtilis C</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 54 | <i>Brevibacillus brevis</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 55 | <i>Brevibacterium casei</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 56 | <i>Brevibacterium epidermidis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 57 | <i>Brevibacterium linens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 58 | <i>Brevibacterium liquifaciens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 59 | <i>Brevibacterium mcbrellneri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 60 | <i>Brevibacterium otitidis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 61 | <i>Brochothrix campestris</i> | GP-ROD | C+ | Y | BUG+B | 6.5% CO ₂ | 35-37 |
| 62 | <i>Brochothrix thermosphacta</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 63 | <i>Carnobacterium alterfunditum</i> | GP-ROD | C- | Y | BUG+B | Air | 26 |
| 64 | <i>Carnobacterium divergens</i> | GP-ROD | C- | Y | BUG+B | Air | 26 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|-----------|------|------|--------|---------------------|-------|
| 65 | <i>Carnobacterium gallinarum</i> | GP-ROD | C- | Y | BUG+B | Air | 26 |
| 66 | <i>Carnobacterium mobile</i> | GP-ROD | C- | Y | BUG+B | Air | 26 |
| 67 | <i>Carnobacterium piscicola</i> | GP-ROD | C- | Y | BUG+B | Air | 26 |
| 68 | <i>Cellulomonas biazotea</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 69 | <i>Cellulomonas cellasea</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 70 | <i>Cellulomonas fimi</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 71 | <i>Cellulomonas flavigena</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 72 | <i>Cellulomonas gelida</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 73 | <i>Cellulomonas hominis</i> (CDC.A-3) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 74 | <i>Cellulomonas turbata</i> (Oerskovia turbata) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 75 | <i>Cellulomonas uda</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 76 | <i>Cellulosimicrobium cellulans</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 77 | <i>Clavibacter agropyri</i> (Corynebacterium) | GP-ROD | C+ | N | BUG | Air | 30 |
| 78 | <i>Clavibacter michiganensis</i> ss insidiosus | GP-ROD | C+ | N | BUG | Air | 30 |
| 79 | <i>Clavibacter michiganensis</i> ss michiganensis | GP-ROD | C+ | N | BUG | Air | 30 |
| 80 | <i>Clavibacter michiganensis</i> ss nebraskensis | GP-ROD | C+ | N | BUG | Air | 30 |
| 81 | <i>Clavibacter michiganensis</i> ss sepedonicus | GP-ROD | C+ | N | BUG | Air | 30 |
| 82 | <i>Clavibacter michiganensis</i> ss tessellarius | GP-ROD | C+ | N | BUG | Air | 30 |
| 83 | <i>Corynebacterium accolens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 84 | <i>Corynebacterium afermentans</i> ss afermentans(CDC.ANF-1) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 85 | <i>Corynebacterium afermentans</i> ss lipophilum | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 86 | " <i>Corynebacterium</i> " ammoniagenes (Brevibacterium-like) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 87 | <i>Corynebacterium amycolatum</i> (CDC.F-2) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 88 | <i>Corynebacterium argenteorotense</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 89 | <i>Corynebacterium auris</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 90 | <i>Corynebacterium bovis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 91 | <i>Corynebacterium callunae</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 92 | <i>Corynebacterium camporealensis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 93 | <i>Corynebacterium coyleae</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 94 | <i>Corynebacterium cystitidis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 95 | <i>Corynebacterium diphtheriae</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 96 | <i>Corynebacterium durum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 97 | <i>Corynebacterium falsenii</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 98 | <i>Corynebacterium flavescens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 99 | <i>Corynebacterium glucuronolyticum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 100 | <i>Corynebacterium glutamicum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 101 | <i>Corynebacterium imitans</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 102 | <i>Corynebacterium jeikeium</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 103 | <i>Corynebacterium kutscheri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 104 | <i>Corynebacterium lipophiloflavum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 105 | <i>Corynebacterium macginleyi</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 106 | <i>Corynebacterium mastitidis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 107 | <i>Corynebacterium matruchotii</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 108 | <i>Corynebacterium minutissimum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 109 | <i>Corynebacterium mucifaciens</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 110 | <i>Corynebacterium mycetoides</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 111 | " <i>Corynebacterium nitrilophilus</i> " | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 112 | <i>Corynebacterium pilosum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 113 | <i>Corynebacterium pseudodiphtheriticum/propinquum</i> (CDC.ANF-3) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 114 | <i>Corynebacterium pseudotuberculosis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 115 | <i>Corynebacterium renale</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 116 | <i>Corynebacterium riegeli</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 117 | <i>Corynebacterium seminale</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 118 | <i>Corynebacterium singulare</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 119 | <i>Corynebacterium</i> spp. (CDC.G) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 120 | <i>Corynebacterium striatum</i> (CDC.I-1) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 121 | <i>Corynebacterium thomassenii</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 122 | <i>Corynebacterium ulcerans</i> | GP-ROD | C+ | Y | BUG+B | 6.5%CO ₂ | 35-37 |
| 123 | <i>Corynebacterium urealyticum</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 124 | <i>Corynebacterium variabile</i> (Caseobacter polymorphus) | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 125 | <i>Corynebacterium vitaeruminis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 126 | " <i>Corynebacterium</i> " xerosis (GPC) | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 127 | <i>Curtobacterium albidum</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 128 | <i>Curtobacterium citreum</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 129 | <i>Curtobacterium flaccumfaciens</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 130 | <i>Curtobacterium luteum</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 131 | <i>Curtobacterium pusillum</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 132 | <i>Deinococcus grandis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|-----------|-------|------|--------|----------|-------|
| 133 | Deinococcus proteolyticus | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 134 | Deinococcus radiodurans | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 135 | Deinococcus radiophilus | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 136 | Deinococcus radiopugnans | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 137 | Dermabacter hominis | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 138 | Dermaococcus nishinomiyaensis | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 139 | Dolosicoccus paucivorans | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 140 | Dolosigranulum pigrum | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 141 | Enterococcus avium | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 142 | Enterococcus casseliflavus | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 143 | Enterococcus cecorum | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 144 | Enterococcus columbae | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 145 | Enterococcus dispar | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 146 | Enterococcus durans | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 147 | Enterococcus faecalis | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 148 | Enterococcus faecium | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 149 | Enterococcus flavescens | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 150 | Enterococcus gallinarum | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 151 | Enterococcus hirae | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 152 | Enterococcus malodoratus | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 153 | Enterococcus mundtii | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 154 | Enterococcus pseudoavium | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 155 | Enterococcus raffinosus | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 156 | Enterococcus saccharolyticus | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 157 | Enterococcus solitarius | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 158 | Enterococcus sulfureus | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 159 | Eremococcus coleocola | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 160 | Erysipelothrix rhusiopathiae/tonsillarum | GP-ROD | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 161 | Exiguobacterium acetylicum | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 162 | Gardnerella vaginalis | GP-ROD | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 163 | Gemella bergeri | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 164 | Gemella haemolysans/morbilorum | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 165 | Gemella palaticanis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 166 | Gemella sanguinis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 167 | Geobacillus stearothermophilus | GP-ROD SB | C-/W+ | N | BUG+M | Air | 55 |
| 168 | Geobacillus thermoglucosidasius | GP-ROD SB | C-/W+ | N | BUG+M | Air | 55 |
| 169 | Globicatella sanguinis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 170 | Gordonia aichiensis | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 171 | Gordonia bronchialis | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 172 | Gordonia rubropincta | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 173 | Gordonia sputi | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 174 | Gordonia terrae | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 175 | Helcococcus kunzii | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 176 | Ignavigranum rouffiae | GP-COCCUS | C+/- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 177 | Jonesia denitrificans | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 178 | Kocuria kristinae | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 179 | Kocuria rosea | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 180 | Kocuria varians | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 181 | Kurthia gibsonii | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 182 | Kurthia sibirica | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 183 | Kurthia zopfii | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 184 | Kytococcus sedentarius | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 185 | Lactococcus garvieae | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 186 | Lactococcus lactis ss cremoris | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 187 | Lactococcus lactis ss diacetylactis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 188 | Lactococcus lactis ss hordniae | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 189 | Lactococcus lactis ss lactis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 190 | Lactococcus plantarum | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 191 | Lactococcus raffinolactis | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 192 | Leifsonia aquatica | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 193 | Leuconostoc carnosum | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 194 | Leuconostoc citreum | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 195 | Leuconostoc fallax | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 196 | Leuconostoc gelidum | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 197 | Leuconostoc lactis | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 198 | Leuconostoc mesenteroides | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 199 | Leuconostoc mesenteroides ss dextranicum | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 200 | Leuconostoc mesenteroides ss mesenteroides | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|-----------|-------|------|--------|----------|-------|
| 201 | <i>Listeria grayi</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 202 | <i>Listeria innocua/monocytogenes/seeligeri/welshimeri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 203 | <i>Listeria ivanovii</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 204 | <i>Listeria ivanovii ss londoniensis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 205 | <i>Listeria monocytogenes/innocua/seeligeri/welshimeri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 206 | <i>Listeria seeligeri/innocua/monocytogenes/welshimeri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 207 | <i>Listeria welshimeri/innocua/monocytogenes/seeligeri</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 208 | <i>Macrococcus bovis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 209 | <i>Macrococcus caroulesicus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 210 | <i>Macrococcus caseolyticus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 211 | <i>Macrococcus equiperdus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 212 | <i>Microbacterium arborescens</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 213 | <i>Microbacterium esteraromaticum</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 214 | <i>Microbacterium flavescens</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 215 | <i>Microbacterium imperiale</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 216 | <i>Microbacterium lacticum</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 217 | <i>Microbacterium laevaniformans</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 218 | <i>Microbacterium liquefaciens</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 219 | <i>Microbacterium maritimum</i> | GP-ROD | C+/- | Y | BUG+B | Air | 30 |
| 220 | <i>Microbacterium saperdae</i> | GP-ROD | C+/- | Y | BUG+B | Air | 30 |
| 221 | <i>Microbacterium spp. (CDC.A-4)</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 222 | <i>Microbacterium spp. (CDC.A-5)</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 223 | <i>Microbacterium terregens</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 224 | <i>Microbacterium testaceum</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 225 | " <i>Micrococcus diversus</i> " | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 226 | <i>Micrococcus luteus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 227 | <i>Micrococcus luteus (ATCC 9341)</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 228 | <i>Micrococcus lylae</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 30 |
| 229 | <i>Paenibacillus azotofixans</i> | GP-ROD SB | C- | N | BUG+M | Air | 30 |
| 230 | <i>Paenibacillus larvae ss larvae</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 231 | <i>Paenibacillus macerans</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 232 | <i>Paenibacillus pabuli</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 233 | <i>Paenibacillus polymyxa</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 234 | <i>Paenibacillus popilliae</i> | GP-ROD SB | C- | N | BUG+M | Air | 30 |
| 235 | <i>Paenibacillus validus</i> | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |
| 236 | <i>Pediococcus acidilactici/parvulus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 237 | <i>Pediococcus dextrinicus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 238 | <i>Pediococcus pentosaceus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 239 | <i>Pediococcus urinaequi</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 240 | <i>Pediococcus urinaequi-like</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 241 | <i>Rathayibacter rathayi</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 242 | <i>Rathayibacter tritici</i> | GP-ROD | C+ | N | BUG | Air | 30 |
| 243 | <i>Rhodococcus australis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 244 | <i>Rhodococcus coprophilus</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 245 | <i>Rhodococcus equi</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 246 | <i>Rhodococcus erythropolis</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 247 | <i>Rhodococcus fascians</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 248 | <i>Rhodococcus globulus</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 249 | <i>Rhodococcus rhodnii</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 250 | <i>Rhodococcus rhodochrous</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 251 | <i>Rhodococcus ruber</i> | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 252 | <i>Rothia dentocariosa</i> | GP-ROD | C+/- | Y | BUG+B | Air | 35-37 |
| 253 | <i>Rothia mucilaginosa</i> | GP-COCCUS | C-/W+ | Y | BUG+B | Air | 35-37 |
| 254 | <i>Sanguibacter inulinus</i> | GP-ROD | C+ | Y | BUG+B | 6.5% CO2 | 35-37 |
| 255 | <i>Sanguibacter keddiei</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 256 | <i>Sanguibacter suarezi</i> | GP-ROD | C+ | Y | BUG+B | Air | 30 |
| 257 | <i>Staphylococcus arlettae</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 258 | <i>Staphylococcus aureus ss anaerobius</i> | GP-COCCUS | C- | Y | BUG+B | Air | 35-37 |
| 259 | <i>Staphylococcus aureus ss aureus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 260 | <i>Staphylococcus auricularis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 261 | <i>Staphylococcus capitis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 262 | <i>Staphylococcus caprae</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 263 | <i>Staphylococcus carnosus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 264 | <i>Staphylococcus chromogenes</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 265 | <i>Staphylococcus cohnii</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 266 | <i>Staphylococcus delphini</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 267 | <i>Staphylococcus epidermidis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 268 | <i>Staphylococcus equorum</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|-----------|------|------|--------|----------|-------|
| 269 | <i>Staphylococcus felis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 270 | <i>Staphylococcus gallinarum</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 271 | <i>Staphylococcus haemolyticus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 272 | <i>Staphylococcus hominis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 273 | <i>Staphylococcus hominis</i> ss novobiosepticus | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 274 | <i>Staphylococcus hyicus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 275 | <i>Staphylococcus intermedius</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 276 | <i>Staphylococcus kloosii</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 277 | <i>Staphylococcus lentus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 278 | <i>Staphylococcus lugdunensis</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 279 | <i>Staphylococcus lutrae</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 280 | <i>Staphylococcus muscae</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 281 | <i>Staphylococcus pasteurii</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 282 | <i>Staphylococcus piscifermentans</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 283 | <i>Staphylococcus pulvereri/vitulinus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 284 | <i>Staphylococcus saprophyticus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 285 | <i>Staphylococcus schleiferi</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 286 | <i>Staphylococcus sciuri</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 287 | <i>Staphylococcus sciuri</i> ss rodentium | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 288 | <i>Staphylococcus simulans</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 289 | <i>Staphylococcus warneri</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 290 | <i>Staphylococcus xylosus</i> | GP-COCCUS | C+ | Y | BUG+B | Air | 35-37 |
| 291 | <i>Streptococcus acidominimus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 292 | <i>Streptococcus agalactiae</i> A (GP B) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 293 | <i>Streptococcus agalactiae</i> B (GP B) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 294 | <i>Streptococcus alactolyticus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 295 | <i>Streptococcus anginosus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 296 | <i>Streptococcus bovis</i> (GP D) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 297 | <i>Streptococcus canis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 298 | <i>Streptococcus constellatus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 299 | <i>Streptococcus criceti</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 300 | <i>Streptococcus cristatus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 301 | <i>Streptococcus downei</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 302 | <i>Streptococcus dysgalactiae</i> ss dysgalactiae | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 303 | <i>Streptococcus dysgalactiae</i> ss equisimilis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 304 | <i>Streptococcus equi</i> ss equi | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 305 | <i>Streptococcus equi</i> ss zooepidemicus (GP C) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 306 | <i>Streptococcus equinus</i> (GP D) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 307 | <i>Streptococcus ferus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 308 | <i>Streptococcus gallolyticus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 309 | <i>Streptococcus gordonii</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 310 | <i>Streptococcus hyointestinalis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 311 | <i>Streptococcus hyovaginalis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 312 | <i>Streptococcus infantarius</i> ss coli | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 313 | <i>Streptococcus infantarius</i> ss infantarius | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 314 | <i>Streptococcus infantis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 315 | <i>Streptococcus iniae</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 316 | <i>Streptococcus intermedius</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 317 | <i>Streptococcus intestinalis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 318 | <i>Streptococcus macacae</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 319 | <i>Streptococcus macedonicus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 320 | <i>Streptococcus mitis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 321 | <i>Streptococcus mutans/ratti</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 322 | <i>Streptococcus oralis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 323 | <i>Streptococcus parasanguinis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 324 | <i>Streptococcus peroris</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 325 | <i>Streptococcus phocae</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 326 | <i>Streptococcus plurianimalium</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 327 | <i>Streptococcus pneumoniae</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 328 | <i>Streptococcus porcinus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 329 | <i>Streptococcus pyogenes</i> A (GP A) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 330 | <i>Streptococcus pyogenes</i> B (GP A) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 331 | <i>Streptococcus pyogenes</i> C (GP A) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 332 | <i>Streptococcus salivarius</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 333 | <i>Streptococcus sanguinis</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 334 | <i>Streptococcus sobrinus</i> | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 335 | <i>Streptococcus suis</i> (GP RST) | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 336 | <i>Streptococcus suis</i> serogroup 1/2 | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|-----------|-------|------|--------|----------|-------|
| 337 | Streptococcus suis serogroup 4/6 | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 338 | Streptococcus suis serogroup 5 | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 339 | Streptococcus suis serogroup 7 | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 340 | Streptococcus thoraltensis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 341 | Streptococcus uberis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 342 | Streptococcus vestibularis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 343 | Streptococcus waius | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 344 | Tetragenococcus halophilus | GP-COCCUS | C- | Y | BUG+B | Air | 30 |
| 345 | Tsukamurella inchoensis | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 346 | Tsukamurella paurometabola | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 347 | Turicella otitidis | GP-ROD | C+ | Y | BUG+B | Air | 35-37 |
| 348 | Vagococcus fluvialis | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 349 | Vagococcus lutrae | GP-COCCUS | C- | Y | BUG+B | 6.5% CO2 | 35-37 |
| 350 | Vagococcus salmoninarum | GP-COCCUS | C- | Y | BUG+B | Air | 26 |
| 351 | Virgibacillus pantothenicus | GP-ROD SB | C-/W+ | N | BUG+M | Air | 30 |

†Agricultural bacteria that may be grown on BUG without blood.

** Found in the Dangerous Pathogen database

Anaerobic Bacteria

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|----|--|----------------|------|------|--------|-----|-------|
| 1 | Abiotrophia defectiva | AN GP-COCOCCUS | N/A | N | BUA | ANA | 35-37 |
| 2 | Acetivibrio ethanolognens | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 3 | Acetoanaerobium noterae | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 4 | Acidaminococcus fermentans | AN GN-COCOCCUS | N/A | N | BUA | ANA | 35-37 |
| 5 | Actinobaculum schaalii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 6 | Actinobaculum suis* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 7 | Actinomyces bovis BGA* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 8 | Actinomyces bovis BGB* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 9 | Actinomyces europaeus* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 10 | Actinomyces graevenitzi* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 11 | Actinomyces howellii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 12 | Actinomyces hyovaginalis* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 13 | Actinomyces israelii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 14 | Actinomyces meyeri* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 15 | Actinomyces naeslundii/viscosus* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 16 | Actinomyces neuii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 17 | Actinomyces odontolyticus* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 18 | Actinomyces radingae* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 19 | Actinomyces slackii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 20 | Actinomyces turicensis* | AN GP-ROD | N/A | N | BUA | ANA | 35-3 |
| 21 | Anaerobiospirillum succiniproducens | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 22 | Anaeromusa acidaminophila | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 23 | Anaerorhabdus furcosa | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 24 | Anaerosinus glycerini | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 25 | Arcanobacterium pyogenes* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 26 | Atopobium fossor | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 27 | Atopobium minutum | AN GP-COCOCCUS | N/A | N | BUA | ANA | 35-37 |
| 28 | Atopobium parvulum | AN GP-COCOCCUS | N/A | N | BUA | ANA | 35-37 |
| 29 | Bacteroides caccae* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 30 | Bacteroides coagulans* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 31 | Bacteroides distasonis* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 32 | Bacteroides eggerthii* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 33 | Bacteroides forsythus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 34 | Bacteroides fragilis BGA* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 35 | Bacteroides fragilis BGB* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 36 | Bacteroides helcogenes* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 37 | Bacteroides merdae* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 38 | Bacteroides ovatus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 39 | Bacteroides pectinophilus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 40 | Bacteroides putredinis* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 41 | Bacteroides pyogenes* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 42 | Bacteroides splanchnicus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 43 | Bacteroides stercoris* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 44 | Bacteroides suis* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 45 | Bacteroides tectus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 46 | Bacteroides thetaiotaomicron* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 47 | Bacteroides uniformis* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 48 | Bacteroides ureolyticus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 49 | Bacteroides vulgatus* | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 50 | Bifidobacterium adolescentis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 51 | Bifidobacterium angulatum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 52 | Bifidobacterium animalis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 53 | Bifidobacterium asteroides | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 54 | Bifidobacterium bifidum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 55 | Bifidobacterium boum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 56 | Bifidobacterium breve | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 57 | Bifidobacterium catenulatum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 58 | Bifidobacterium choerinum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 59 | Bifidobacterium coryneforme | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 60 | Bifidobacterium cuniculi | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 61 | Bifidobacterium dentium | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 62 | Bifidobacterium gallicum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 63 | Bifidobacterium gallinarum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 64 | Bifidobacterium indicum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 65 | Bifidobacterium infantis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 66 | Bifidobacterium longum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |

*Indicates species that require a MicroStation Reader for identification.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|-----------|------|------|--------|-----|-------|
| 67 | <i>Bifidobacterium magnum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 68 | <i>Bifidobacterium merycicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 69 | <i>Bifidobacterium minimum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 70 | <i>Bifidobacterium pseudocatenulatum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 71 | <i>Bifidobacterium pseudolongum</i> ss <i>globosum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 72 | <i>Bifidobacterium pseudolongum</i> ss <i>pseudolongum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 73 | <i>Bifidobacterium pullorum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 74 | <i>Bifidobacterium ruminantium</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 75 | <i>Bifidobacterium saeculare</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 76 | <i>Bifidobacterium subtile</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 77 | <i>Bifidobacterium suis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 78 | <i>Bifidobacterium thermophilum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 79 | <i>Bilophila wadsworthia</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 80 | <i>Butyrivibrio fibrisolvens</i> * | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 81 | <i>Campylobacter curvus/gracilis</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 82 | <i>Campylobacter rectus</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 83 | <i>Catonella morbi</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 84 | <i>Centipeda periodontii</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 85 | <i>Clostridium absonum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 86 | <i>Clostridium acetobutylicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 87 | <i>Clostridium aerotolerans</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 88 | <i>Clostridium "aminobutyricum"</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 89 | <i>Clostridium aminophilum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 90 | <i>Clostridium aminovalericum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 91 | <i>Clostridium barati</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 92 | <i>Clostridium botulinum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 93 | <i>Clostridium cadaveris</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 94 | <i>Clostridium carnis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 95 | <i>Clostridium cellobioparum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 96 | <i>Clostridium cellulolyticum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 97 | <i>Clostridium clostridioforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 98 | <i>Clostridium coccooides</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 99 | <i>Clostridium cochlearium</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 100 | <i>Clostridium cocleatum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 101 | <i>Clostridium colinum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 102 | <i>Clostridium cylindrosporum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 103 | <i>Clostridium difficile</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 104 | <i>Clostridium fallax</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 105 | <i>Clostridium glycolicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 106 | <i>Clostridium hastiforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 107 | <i>Clostridium hydroxybenzoicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 108 | <i>Clostridium innocuum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 109 | <i>Clostridium intestinale</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 110 | <i>Clostridium irregularis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 111 | <i>Clostridium magnum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 112 | <i>Clostridium malenominatum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 113 | <i>Clostridium nexile</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 114 | <i>Clostridium novyi</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 115 | <i>Clostridium orbiscindens</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 116 | <i>Clostridium oroticum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 117 | <i>Clostridium papyrosolvens</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 118 | <i>Clostridium paraputrificum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 119 | <i>Clostridium pasteurianum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 120 | <i>Clostridium perfringens</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 121 | <i>Clostridium propionicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 122 | <i>Clostridium purinolyticum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 123 | <i>Clostridium putrificum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 124 | <i>Clostridium ramosum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 125 | <i>Clostridium rectum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 126 | <i>Clostridium saccharolyticum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 127 | <i>Clostridium sardiniensis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 128 | <i>Clostridium sartagoformum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 129 | <i>Clostridium scatologenes</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 130 | <i>Clostridium scindens</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 131 | <i>Clostridium septicum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 132 | <i>Clostridium sordelli</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 133 | <i>Clostridium sphenoides</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 134 | <i>Clostridium spiroforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 135 | <i>Clostridium sporogenes</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |

*Indicates species that require a MicroStation Reader for identification.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|--------------|------|------|--------|-----|-------|
| 136 | <i>Clostridium sporosphaeroides</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 137 | <i>Clostridium sticklandii</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 138 | <i>Clostridium subterminale</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 139 | <i>Clostridium symbiosum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 140 | <i>Clostridium tertium</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 141 | <i>Clostridium tetani</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 142 | <i>Clostridium tetanomorphum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 143 | <i>Clostridium tyrobutyricum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 144 | <i>Clostridium/Eubacterium species</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 145 | <i>Collinsella aerofaciens</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 146 | <i>Coprococcus catus</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 147 | <i>Dendrosporobacter quercicolus</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 148 | <i>Desulfobulbus elongatus</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 149 | <i>Desulfomonas pigra</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 150 | <i>Desulfomonile tiedjei</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 151 | <i>Desulfovibrio desulfuricans ss desulfuricans</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 152 | <i>Desulfovibrio fructosivorans</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 153 | <i>Desulfovibrio vulgaris ss vulgaris</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 154 | <i>Eggerthella lenta*</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 155 | <i>Eikenella corrodens*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 156 | <i>Eubacterium angustum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 157 | <i>Eubacterium barkeri</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 158 | <i>Eubacterium bifforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 159 | <i>Eubacterium brachy</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 160 | <i>Eubacterium budayi</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 161 | <i>Eubacterium combesii</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 162 | <i>Eubacterium contortum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 163 | <i>Eubacterium cylindroides</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 164 | <i>Eubacterium desmolans</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 165 | <i>Eubacterium dolichum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 166 | <i>Eubacterium fissicatena</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 167 | <i>Eubacterium hallii</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 168 | <i>Eubacterium limosum BGA</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 169 | <i>Eubacterium limosum BGB</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 170 | <i>Eubacterium moniliforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 171 | <i>Eubacterium multifforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 172 | <i>Eubacterium nodatum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 173 | <i>Eubacterium plautii</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 174 | <i>Eubacterium saburreum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 175 | <i>Eubacterium tortuosum</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 176 | <i>Eubacterium uniforme</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 177 | <i>Eubacterium yurii ss margaretae</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 178 | <i>Eubacterium yurii ss schtitka</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 179 | <i>Eubacterium yurii ss yurii</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 180 | <i>Facklamia hominis</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 181 | <i>Facklamia sourekii</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 182 | <i>Falcivibrio grandis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 183 | <i>Falcivibrio vaginalis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 184 | <i>Filifactor alocis*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 185 | <i>Filifactor villosus</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 186 | <i>Finegoldia magna</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 187 | <i>Fusobacterium gonidiaformans*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 188 | <i>Fusobacterium naviforme*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 189 | <i>Fusobacterium necrogenes*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 190 | <i>Fusobacterium necrophorum*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 191 | <i>Fusobacterium nucleatum ss animalis*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 192 | <i>Fusobacterium nucleatum ss fusiforme*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 193 | <i>Fusobacterium nucleatum ss nucleatum*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 194 | <i>Fusobacterium nucleatum ss polymorphum*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 195 | <i>Fusobacterium nucleatum ss vincentii*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 196 | <i>Fusobacterium russii*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 197 | <i>Fusobacterium simiae*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 198 | <i>Fusobacterium ulcerans*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 199 | <i>Fusobacterium varium*</i> | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 200 | <i>Gardnerella vaginalis</i> | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 201 | <i>Gemella bergeri</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 202 | <i>Gemella haemolysans</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 203 | <i>Gemella morbillorum</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 204 | <i>Granulicatella adiacens</i> | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |

*Indicates species that require a MicroStation Reader for identification.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|--------------|------|------|--------|-----|-------|
| 205 | Hallella seregens | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 206 | Lactobacillus acetotolerans | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 207 | Lactobacillus acidophilus BGA | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 208 | Lactobacillus acidophilus BGB | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 209 | Lactobacillus alimentarius | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 210 | Lactobacillus amylovorus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 211 | Lactobacillus aviarius ss araffinosus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 212 | Lactobacillus aviarius ss aviarius | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 213 | Lactobacillus bifermentans | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 214 | Lactobacillus brevis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 215 | Lactobacillus buchneri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 216 | Lactobacillus casei | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 217 | Lactobacillus cateniformis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 218 | Lactobacillus collinoides | AN GP-ROD | N/A | N | BUA | ANA | 26 |
| 219 | Lactobacillus coryniformis ss coryniformis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 220 | Lactobacillus coryniformis ss torquens | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 221 | Lactobacillus crispatus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 222 | Lactobacillus curvatus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 223 | Lactobacillus delbrueckii ss bulgaricus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 224 | Lactobacillus delbrueckii ss delbrueckii | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 225 | Lactobacillus delbrueckii ss lactis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 226 | Lactobacillus farciminis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 227 | Lactobacillus fermentum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 228 | Lactobacillus fructivorans | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 229 | Lactobacillus fructosus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 230 | Lactobacillus gasseri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 231 | Lactobacillus hamsteri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 232 | Lactobacillus helveticus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 233 | Lactobacillus hilgardii | AN GP-ROD | N/A | N | BUA | ANA | 26 |
| 234 | Lactobacillus jensenii | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 235 | Lactobacillus kefir | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 236 | Lactobacillus malefermentans | AN GP-ROD | N/A | N | BUA | ANA | 30 |
| 237 | Lactobacillus mali/sakei ss sakei | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 238 | Lactobacillus murinus/paracasei ss tolerans | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 239 | Lactobacillus oris/parabuchneri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 240 | Lactobacillus paracasei ss paracasei | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 241 | Lactobacillus pentosus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 242 | Lactobacillus plantarum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 243 | Lactobacillus reuteri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 244 | Lactobacillus rhamnosus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 245 | Lactobacillus salivarius ss salicinius | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 246 | Lactobacillus salivarius ss salivarius | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 247 | Lactobacillus sanfranciscensis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 248 | Lactobacillus suebicus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 249 | Lactobacillus vaginalis | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 250 | Lactococcus garvieae | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 251 | Lactococcus lactis ss cremoris | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 252 | Lactococcus lactis ss hordniae | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 253 | Lactococcus lactis ss lactis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 254 | Lactococcus plantarum | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 255 | Lactococcus raffinolactis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 256 | Leuconostoc carnosum | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 257 | Leuconostoc citreum | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 258 | Leuconostoc fallax | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 259 | Leuconostoc gelidum | AN GP-COCCUS | N/A | N | BUA | ANA | 26 |
| 260 | Leuconostoc lactis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 261 | Leuconostoc mesenteroides ss cremoris | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 262 | Leuconostoc mesenteroides ss dextranicum | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 263 | Leuconostoc mesenteroides ss mesenteroides | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 264 | Megamonas hypermegale | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 265 | Megasphaera cerevisiae | AN GN-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 266 | Megasphaera elsdenii | AN GN-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 267 | Micromonas micros | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 268 | Mobiluncus curtisii | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 269 | Mobiluncus mulieris | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 270 | Mongibacterium timidum | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 271 | Pectinatus cerevisiiphilus | AN GN-ROD | N/A | N | BUA | ANA | 30 |
| 272 | Pectinatus frisingensis | AN GN-ROD | N/A | N | BUA | ANA | 30 |
| 273 | Pediococcus acidilactici | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |

*Indicates species that require a MicroStation Reader for identification.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|--------------|------|------|--------|-----|-------|
| 274 | Pediococcus dextrinicus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 275 | Pediococcus parvulus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 276 | Pediococcus pentosaceus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 277 | Pediococcus urinaeaequi | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 278 | Peptostreptococcus anaerobius | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 279 | Peptostreptococcus asaccharolyticus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 280 | Peptostreptococcus barnesae | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 281 | Peptostreptococcus hydrogenalis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 282 | Peptostreptococcus indolicus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 283 | Peptostreptococcus lacrimalis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 284 | Peptostreptococcus lactolyticus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 285 | Peptostreptococcus prevotii | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 286 | Peptostreptococcus tetradius | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 287 | Peptostreptococcus vaginalis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 288 | Porphyromonas asaccharolytica | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 289 | Porphyromonas circumdentaria/endodontalis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 290 | Porphyromonas gingivalis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 291 | Porphyromonas levii | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 292 | Porphyromonas macacae | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 293 | Prevotella bivia | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 294 | Prevotella buccae | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 295 | Prevotella buccalis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 296 | Prevotella corporis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 297 | Prevotella denticola | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 298 | Prevotella disiens (GNR Kan-R setup) | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 299 | Prevotella disiens (Standard AN setup) | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 300 | Prevotella heparinolytica | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 301 | Prevotella intermedia | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 302 | Prevotella loescheii | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 303 | Prevotella nigrescens | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 304 | Prevotella oralis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 305 | Prevotella oris (GNR Kan-R setup) | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 306 | Prevotella oulora | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 307 | Prevotella species | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 308 | Propionibacterium acidipropionici* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 309 | Propionibacterium acnes* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 310 | Propionibacterium avidum* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 311 | Propionibacterium freudenreichii ss freudenreichii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 312 | Propionibacterium granulosum* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 313 | Propionibacterium jensenii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 314 | Propionibacterium lymphophilum* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 315 | Propionibacterium propionicus* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 316 | Propionibacterium thoenii* | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 317 | Pseudoramibacter alactolyticus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 318 | Rikenella microfus | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 319 | Ruminococcus gnavus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 320 | Ruminococcus hansenii | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 321 | Ruminococcus lactaris | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 322 | Ruminococcus productus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 323 | Ruminococcus torques | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 324 | Sarcina maxima | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 325 | Sarcina ventriculi | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 326 | Seibaldella termitidis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 327 | Selenomonas flueggei | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 328 | Selenomonas infelix | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 329 | Selenomonas noxia | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 330 | Selenomonas sputigena | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 331 | Slackia heliotrinireducens | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 332 | Sporolactobacillus inulinus | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 333 | Sporomusa acidovorans | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 334 | Sporomusa ovata | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 335 | Sporomusa sphaeroides | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 336 | Staphylococcus aureus ss anaerobius | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 337 | Staphylococcus saccharolyticus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 338 | Streptococcus alactolyticus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 339 | Streptococcus anginosus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 340 | Streptococcus bovis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 341 | Streptococcus macacae | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 342 | Streptococcus mutans | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |

*Indicates species that require a MicroStation Reader for identification.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|--------------|------|------|--------|-----|-------|
| 343 | Streptococcus pleomorphus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 344 | Streptococcus suis | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 345 | Streptococcus thermophilus | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 346 | Succinivibrio dextrinosolvens | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 347 | Sutterella wadsworthensis | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 348 | Tetragenococcus halophilus | AN GP-COCCUS | N/A | N | BUA | ANA | 30 |
| 349 | Tissierella praeacuta | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 350 | Veillonella atypica/caviae/dispar/parvula | AN GN-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 351 | Veillonella criceti | AN GN-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 352 | Veillonella ratti | AN GN-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 353 | Weissella confusa | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 354 | Weissella halotolerans/hellenica | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 355 | Weissella kandleri | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 356 | Weissella minor | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 357 | Weissella paramesenteroides | AN GP-COCCUS | N/A | N | BUA | ANA | 35-37 |
| 358 | Weissella viridescens | AN GP-ROD | N/A | N | BUA | ANA | 35-37 |
| 359 | Wolinella succinogenes | AN GN-ROD | N/A | N | BUA | ANA | 35-37 |
| 360 | Zymomonas mobilis ss mobilis* | AN GN-ROD | N/A | N | BUA | ANA | 30 |
| 361 | Zymophilus paucivorans/raffinosisivorans | AN GN-ROD | N/A | N | BUA | ANA | 30 |

**Indicates species that require a MicroStation Reader for identification.*

Yeasts

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|----|--------------------------------------|------|------|------|--------|-----|------|
| 1 | Arthroascus javanensis | N/A | N/A | N | BUY | Air | 26 |
| 2 | Bulleromyces albus | N/A | N/A | N | BUY | Air | 26 |
| 3 | Candida aaseri A | N/A | N/A | N | BUY | Air | 26 |
| 4 | Candida aaseri B | N/A | N/A | N | BUY | Air | 26 |
| 5 | Candida albicans | N/A | N/A | N | BUY | Air | 26 |
| 6 | Candida apicola | N/A | N/A | N | BUY | Air | 26 |
| 7 | Candida azyma | N/A | N/A | N | BUY | Air | 26 |
| 8 | Candida blankii | N/A | N/A | N | BUY | Air | 26 |
| 9 | Candida boidinii | N/A | N/A | N | BUY | Air | 26 |
| 10 | Candida bombi | N/A | N/A | N | BUY | Air | 26 |
| 11 | Candida cantarelli | N/A | N/A | N | BUY | Air | 26 |
| 12 | Candida cariosilignicola | N/A | N/A | N | BUY | Air | 26 |
| 13 | Candida castellii | N/A | N/A | N | BUY | Air | 26 |
| 14 | Candida catenulata | N/A | N/A | N | BUY | Air | 26 |
| 15 | Candida diddensiae | N/A | N/A | N | BUY | Air | 26 |
| 16 | Candida diversa | N/A | N/A | N | BUY | Air | 26 |
| 17 | Candida edax | N/A | N/A | N | BUY | Air | 26 |
| 18 | Candida entomophila | N/A | N/A | N | BUY | Air | 26 |
| 19 | Candida ergastensis | N/A | N/A | N | BUY | Air | 26 |
| 20 | Candida etchellsii | N/A | N/A | N | BUY | Air | 26 |
| 21 | Candida famata | N/A | N/A | N | BUY | Air | 26 |
| 22 | Candida fluvialis | N/A | N/A | N | BUY | Air | 26 |
| 23 | Candida freyschussii | N/A | N/A | N | BUY | Air | 26 |
| 24 | Candida friedrichii | N/A | N/A | N | BUY | Air | 26 |
| 25 | Candida fructus | N/A | N/A | N | BUY | Air | 26 |
| 26 | Candida fusiformata | N/A | N/A | N | BUY | Air | 26 |
| 27 | Candida galacta | N/A | N/A | N | BUY | Air | 26 |
| 28 | Candida geochares | N/A | N/A | N | BUY | Air | 26 |
| 29 | Candida glabrata | N/A | N/A | N | BUY | Air | 26 |
| 30 | Candida glabiosa | N/A | N/A | N | BUY | Air | 26 |
| 31 | Candida gropengiesseri | N/A | N/A | N | BUY | Air | 26 |
| 32 | Candida haemulonii | N/A | N/A | N | BUY | Air | 26 |
| 33 | Candida humilis | N/A | N/A | N | BUY | Air | 26 |
| 34 | Candida incommunis | N/A | N/A | N | BUY | Air | 26 |
| 35 | Candida insectalens | N/A | N/A | N | BUY | Air | 26 |
| 36 | Candida insectamans | N/A | N/A | N | BUY | Air | 26 |
| 37 | Candida insectorum | N/A | N/A | N | BUY | Air | 26 |
| 38 | Candida intermedia | N/A | N/A | N | BUY | Air | 26 |
| 39 | Candida ishiwadae | N/A | N/A | N | BUY | Air | 26 |
| 40 | Candida magnoliae | N/A | N/A | N | BUY | Air | 26 |
| 41 | Candida maltosa | N/A | N/A | N | BUY | Air | 26 |
| 42 | Candida maris | N/A | N/A | N | BUY | Air | 26 |
| 43 | Candida maritima | N/A | N/A | N | BUY | Air | 26 |
| 44 | Candida melibiosica | N/A | N/A | N | BUY | Air | 26 |
| 45 | Candida mogii | N/A | N/A | N | BUY | Air | 26 |
| 46 | Candida montana | N/A | N/A | N | BUY | Air | 26 |
| 47 | Candida multigemmis | N/A | N/A | N | BUY | Air | 26 |
| 48 | Candida musae | N/A | N/A | N | BUY | Air | 26 |
| 49 | Candida nemodendra | N/A | N/A | N | BUY | Air | 26 |
| 50 | Candida nitratophila | N/A | N/A | N | BUY | Air | 26 |
| 51 | Candida norvegica | N/A | N/A | N | BUY | Air | 26 |
| 52 | Candida oleophila | N/A | N/A | N | BUY | Air | 26 |
| 53 | Candida parapsilosis A | N/A | N/A | N | BUY | Air | 26 |
| 54 | Candida parapsilosis B | N/A | N/A | N | BUY | Air | 26 |
| 55 | Candida pararugosa | N/A | N/A | N | BUY | Air | 26 |
| 56 | Candida peltata | N/A | N/A | N | BUY | Air | 26 |
| 57 | Candida pinus | N/A | N/A | N | BUY | Air | 26 |
| 58 | Candida rugosa A | N/A | N/A | N | BUY | Air | 26 |
| 59 | Candida rugosa B | N/A | N/A | N | BUY | Air | 26 |
| 60 | Candida sake | N/A | N/A | N | BUY | Air | 26 |
| 61 | Candida salmanticensis | N/A | N/A | N | BUY | Air | 26 |
| 62 | Candida santamariae | N/A | N/A | N | BUY | Air | 26 |
| 63 | Candida savonica | N/A | N/A | N | BUY | Air | 26 |
| 64 | Candida shehatae | N/A | N/A | N | BUY | Air | 26 |
| 65 | Candida shehatae var shehatae | N/A | N/A | N | BUY | Air | 26 |

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|------|------|------|--------|-----|------|
| 66 | <i>Candida silvae</i> | N/A | N/A | N | BUY | Air | 26 |
| 67 | <i>Candida silvatica</i> | N/A | N/A | N | BUY | Air | 26 |
| 68 | <i>Candida solani</i> | N/A | N/A | N | BUY | Air | 26 |
| 69 | <i>Candida sonorensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 70 | <i>Candida sorbophila</i> | N/A | N/A | N | BUY | Air | 26 |
| 71 | <i>Candida sorboxylosa</i> | N/A | N/A | N | BUY | Air | 26 |
| 72 | <i>Candida spandovens</i> | N/A | N/A | N | BUY | Air | 26 |
| 73 | <i>Candida succiphila</i> | N/A | N/A | N | BUY | Air | 26 |
| 74 | <i>Candida tropicalis A</i> | N/A | N/A | N | BUY | Air | 26 |
| 75 | <i>Candida tropicalis B</i> | N/A | N/A | N | BUY | Air | 26 |
| 76 | <i>Candida vanderwaltii</i> | N/A | N/A | N | BUY | Air | 26 |
| 77 | <i>Candida vartiovaarai</i> | N/A | N/A | N | BUY | Air | 26 |
| 78 | <i>Candida versatilis</i> | N/A | N/A | N | BUY | Air | 26 |
| 79 | <i>Candida viswanathii</i> | N/A | N/A | N | BUY | Air | 26 |
| 80 | <i>Candida wickerhamii</i> | N/A | N/A | N | BUY | Air | 26 |
| 81 | <i>Candida zeylanoides</i> | N/A | N/A | N | BUY | Air | 26 |
| 82 | <i>Citeromyces matritensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 83 | <i>Clavispora lusitaniae</i> | N/A | N/A | N | BUY | Air | 26 |
| 84 | <i>Cryptococcus albidus</i> | N/A | N/A | N | BUY | Air | 26 |
| 85 | <i>Cryptococcus albidus var aerius</i> | N/A | N/A | N | BUY | Air | 26 |
| 86 | <i>Cryptococcus albidus var albidus</i> | N/A | N/A | N | BUY | Air | 26 |
| 87 | <i>Cryptococcus albidus var diffluens</i> | N/A | N/A | N | BUY | Air | 26 |
| 88 | <i>Cryptococcus amyloletus</i> | N/A | N/A | N | BUY | Air | 26 |
| 89 | <i>Cryptococcus curvatus A</i> | N/A | N/A | N | BUY | Air | 26 |
| 90 | <i>Cryptococcus curvatus B</i> | N/A | N/A | N | BUY | Air | 26 |
| 91 | <i>Cryptococcus dimennae</i> | N/A | N/A | N | BUY | Air | 26 |
| 92 | <i>Cryptococcus gastricus</i> | N/A | N/A | N | BUY | Air | 26 |
| 93 | <i>Cryptococcus kuetzingii</i> | N/A | N/A | N | BUY | Air | 26 |
| 94 | <i>Cryptococcus laurentii</i> | N/A | N/A | N | BUY | Air | 26 |
| 95 | <i>Cryptococcus luteolus</i> | N/A | N/A | N | BUY | Air | 26 |
| 96 | <i>Cryptococcus macerans</i> | N/A | N/A | N | BUY | Air | 26 |
| 97 | <i>Cryptococcus magnus</i> | N/A | N/A | N | BUY | Air | 26 |
| 98 | <i>Cryptococcus magnus var aerius</i> | N/A | N/A | N | BUY | Air | 26 |
| 99 | <i>Cryptococcus marinus</i> | N/A | N/A | N | BUY | Air | 26 |
| 100 | <i>Cryptococcus skinneri</i> | N/A | N/A | N | BUY | Air | 26 |
| 101 | <i>Cryptococcus terreus A</i> | N/A | N/A | N | BUY | Air | 26 |
| 102 | <i>Cryptococcus terreus B</i> | N/A | N/A | N | BUY | Air | 26 |
| 103 | <i>Cryptococcus tsukubaensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 104 | <i>Debaryomyces castellii</i> | N/A | N/A | N | BUY | Air | 26 |
| 105 | <i>Debaryomyces hansenii A</i> | N/A | N/A | N | BUY | Air | 26 |
| 106 | <i>Debaryomyces hansenii B</i> | N/A | N/A | N | BUY | Air | 26 |
| 107 | <i>Debaryomyces hansenii C</i> | N/A | N/A | N | BUY | Air | 26 |
| 108 | <i>Debaryomyces hansenii var fabryi</i> | N/A | N/A | N | BUY | Air | 26 |
| 109 | <i>Debaryomyces maramus</i> | N/A | N/A | N | BUY | Air | 26 |
| 110 | <i>Debaryomyces polymorphus</i> | N/A | N/A | N | BUY | Air | 26 |
| 111 | <i>Debaryomyces vanrijiae</i> | N/A | N/A | N | BUY | Air | 26 |
| 112 | <i>Dekkera anomala</i> | N/A | N/A | N | BUY | Air | 26 |
| 113 | <i>Dekkera bruxellensis A</i> | N/A | N/A | N | BUY | Air | 26 |
| 114 | <i>Dekkera bruxellensis B</i> | N/A | N/A | N | BUY | Air | 26 |
| 115 | <i>Dekkera custersiana</i> | N/A | N/A | N | BUY | Air | 26 |
| 116 | <i>Dekkera naardenensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 117 | <i>Dipodascus capitatus</i> | N/A | N/A | N | BUY | Air | 26 |
| 118 | <i>Dipodascus ovetensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 119 | <i>Endomyces fibuliger</i> | N/A | N/A | N | BUY | Air | 26 |
| 120 | <i>Endomycopsella vivi</i> | N/A | N/A | N | BUY | Air | 26 |
| 121 | <i>Eremothecium ashbyi</i> | N/A | N/A | N | BUY | Air | 26 |
| 122 | <i>Fellomyces fuzhouensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 123 | <i>Filobasidiella neoformans bacillisporus</i> | N/A | N/A | N | BUY | Air | 26 |
| 124 | <i>Filobasidiella neoformans neoformans A</i> | N/A | N/A | N | BUY | Air | 26 |
| 125 | <i>Filobasidiella neoformans neoformans B</i> | N/A | N/A | N | BUY | Air | 26 |
| 126 | <i>Filobasidium uniguttulatum</i> | N/A | N/A | N | BUY | Air | 26 |
| 127 | <i>Galactomyces geotrichum</i> | N/A | N/A | N | BUY | Air | 26 |
| 128 | <i>Geotrichum terrestre</i> | N/A | N/A | N | BUY | Air | 26 |
| 129 | <i>Guilliermondella selenospora</i> | N/A | N/A | N | BUY | Air | 26 |
| 130 | <i>Hanseniaspora guilliermondii/uvarum/valb</i> | N/A | N/A | N | BUY | Air | 26 |
| 131 | <i>Hanseniaspora occidentalis</i> | N/A | N/A | N | BUY | Air | 26 |
| 132 | <i>Hanseniaspora osmophila/vineae</i> | N/A | N/A | N | BUY | Air | 26 |
| 133 | <i>Hyphopichia burtonii A</i> | N/A | N/A | N | BUY | Air | 26 |
| 134 | <i>Hyphopichia burtonii B</i> | N/A | N/A | N | BUY | Air | 26 |

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|--|------|------|------|--------|-----|------|
| 135 | <i>Hyphopichia burtonii</i> C | N/A | N/A | N | BUY | Air | 26 |
| 136 | <i>Issatchenkia orientalis</i> | N/A | N/A | N | BUY | Air | 26 |
| 137 | <i>Issatchenkia scutulata</i> | N/A | N/A | N | BUY | Air | 26 |
| 138 | <i>Issatchenkia scutulata</i> var <i>exigua</i> | N/A | N/A | N | BUY | Air | 26 |
| 139 | <i>Issatchenkia scutulata</i> var <i>scutula</i> | N/A | N/A | N | BUY | Air | 26 |
| 140 | <i>Kluyveromyces delphensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 141 | <i>Kluyveromyces lactis</i> | N/A | N/A | N | BUY | Air | 26 |
| 142 | <i>Kluyveromyces lodderae</i> | N/A | N/A | N | BUY | Air | 26 |
| 143 | <i>Kluyveromyces marxianus</i> | N/A | N/A | N | BUY | Air | 26 |
| 144 | <i>Kluyveromyces thermotolerans</i> | N/A | N/A | N | BUY | Air | 26 |
| 145 | <i>Kluyveromyces wickerhamii</i> | N/A | N/A | N | BUY | Air | 26 |
| 146 | <i>Kurtzmanomyces nectairei</i> | N/A | N/A | N | BUY | Air | 26 |
| 147 | <i>Lodderomyces elongisporus</i> | N/A | N/A | N | BUY | Air | 26 |
| 148 | <i>Metschnikowia pulcherrima</i> | N/A | N/A | N | BUY | Air | 26 |
| 149 | <i>Metschnikowia reukauffii</i> | N/A | N/A | N | BUY | Air | 26 |
| 150 | <i>Metschnikowia zobellii</i> | N/A | N/A | N | BUY | Air | 26 |
| 151 | <i>Nadsonia fulvescens</i> | N/A | N/A | N | BUY | Air | 26 |
| 152 | <i>Pachysolen tannophilus</i> | N/A | N/A | N | BUY | Air | 26 |
| 153 | <i>Phaffia rhodozyma</i> | N/A | N/A | N | BUY | Air | 26 |
| 154 | <i>Pichia alni</i> | N/A | N/A | N | BUY | Air | 26 |
| 155 | <i>Pichia amenthionina</i> var <i>amethonina</i> | N/A | N/A | N | BUY | Air | 26 |
| 156 | <i>Pichia amenthionina</i> var <i>pachy</i> | N/A | N/A | N | BUY | Air | 26 |
| 157 | <i>Pichia amylophila/mississippiensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 158 | <i>Pichia angusta</i> | N/A | N/A | N | BUY | Air | 26 |
| 159 | <i>Pichia anomala</i> | N/A | N/A | N | BUY | Air | 26 |
| 160 | <i>Pichia bisporea</i> | N/A | N/A | N | BUY | Air | 26 |
| 161 | <i>Pichia canadensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 162 | <i>Pichia carsonii</i> | N/A | N/A | N | BUY | Air | 26 |
| 163 | <i>Pichia etchellsii</i> | N/A | N/A | N | BUY | Air | 26 |
| 164 | <i>Pichia fabianii</i> | N/A | N/A | N | BUY | Air | 26 |
| 165 | <i>Pichia farinosa/muscolicola</i> | N/A | N/A | N | BUY | Air | 26 |
| 166 | <i>Pichia fermentans</i> | N/A | N/A | N | BUY | Air | 26 |
| 167 | <i>Pichia fluxuum</i> | N/A | N/A | N | BUY | Air | 26 |
| 168 | <i>Pichia glucozyma/methanolica</i> | N/A | N/A | N | BUY | Air | 26 |
| 169 | <i>Pichia guilliermondii</i> A | N/A | N/A | N | BUY | Air | 26 |
| 170 | <i>Pichia guilliermondii</i> B | N/A | N/A | N | BUY | Air | 26 |
| 171 | <i>Pichia haplophila</i> | N/A | N/A | N | BUY | Air | 26 |
| 172 | <i>Pichia holstii</i> | N/A | N/A | N | BUY | Air | 26 |
| 173 | <i>Pichia jadinii</i> | N/A | N/A | N | BUY | Air | 26 |
| 174 | <i>Pichia kluyveri</i> | N/A | N/A | N | BUY | Air | 26 |
| 175 | <i>Pichia media</i> | N/A | N/A | N | BUY | Air | 26 |
| 176 | <i>Pichia membranaefaciens</i> | N/A | N/A | N | BUY | Air | 26 |
| 177 | <i>Pichia mexicana</i> | N/A | N/A | N | BUY | Air | 26 |
| 178 | <i>Pichia minuta</i> | N/A | N/A | N | BUY | Air | 26 |
| 179 | <i>Pichia muscolicola</i> | N/A | N/A | N | BUY | Air | 26 |
| 180 | <i>Pichia norvegensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 181 | <i>Pichia ohmeri</i> A | N/A | N/A | N | BUY | Air | 26 |
| 182 | <i>Pichia ohmeri</i> B | N/A | N/A | N | BUY | Air | 26 |
| 183 | <i>Pichia onychis</i> | N/A | N/A | N | BUY | Air | 26 |
| 184 | <i>Pichia opuntiae</i> | N/A | N/A | N | BUY | Air | 26 |
| 185 | <i>Pichia pastoris</i> | N/A | N/A | N | BUY | Air | 26 |
| 186 | <i>Pichia petersonii</i> | N/A | N/A | N | BUY | Air | 26 |
| 187 | <i>Pichia pijperi</i> | N/A | N/A | N | BUY | Air | 26 |
| 188 | <i>Pichia pini</i> | N/A | N/A | N | BUY | Air | 26 |
| 189 | <i>Pichia rabaulensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 190 | <i>Pichia rhodanensis</i> | N/A | N/A | N | BUY | Air | 26 |
| 191 | <i>Pichia silvicola</i> | N/A | N/A | N | BUY | Air | 26 |
| 192 | <i>Pichia spartinae</i> | N/A | N/A | N | BUY | Air | 26 |
| 193 | <i>Pichia stipitis</i> | N/A | N/A | N | BUY | Air | 26 |
| 194 | <i>Pichia subpelliculosa</i> | N/A | N/A | N | BUY | Air | 26 |
| 195 | <i>Pichia sydowiorum</i> | N/A | N/A | N | BUY | Air | 26 |
| 196 | <i>Pichia thermotolerans</i> | N/A | N/A | N | BUY | Air | 26 |
| 197 | <i>Pichia toletana</i> | N/A | N/A | N | BUY | Air | 26 |
| 198 | <i>Pichia trehalophila</i> | N/A | N/A | N | BUY | Air | 26 |
| 199 | <i>Pichia triangularis</i> | N/A | N/A | N | BUY | Air | 26 |
| 200 | <i>Rhodospiridium diobovatum</i> | N/A | N/A | N | BUY | Air | 26 |
| 201 | <i>Rhodospiridium sphaerocarpum</i> | N/A | N/A | N | BUY | Air | 26 |
| 202 | <i>Rhodospiridium toruloides</i> | N/A | N/A | N | BUY | Air | 26 |

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Test | Thio | Medium | Atm | Temp |
|-----|---|------|------|------|--------|-----|------|
| 203 | Rhodotorula acheniorum | N/A | N/A | N | BUY | Air | 26 |
| 204 | Rhodotorula acuta | N/A | N/A | N | BUY | Air | 26 |
| 205 | Rhodotorula araucariae | N/A | N/A | N | BUY | Air | 26 |
| 206 | Rhodotorula aurantiaca A | N/A | N/A | N | BUY | Air | 26 |
| 207 | Rhodotorula aurantiaca B | N/A | N/A | N | BUY | Air | 26 |
| 208 | Rhodotorula bacarum | N/A | N/A | N | BUY | Air | 26 |
| 209 | Rhodotorula glutinis | N/A | N/A | N | BUY | Air | 26 |
| 210 | Rhodotorula glutinis var glutinis | N/A | N/A | N | BUY | Air | 26 |
| 211 | Rhodotorula graminis | N/A | N/A | N | BUY | Air | 26 |
| 212 | Rhodotorula hylophila | N/A | N/A | N | BUY | Air | 26 |
| 213 | Rhodotorula minuta | N/A | N/A | N | BUY | Air | 26 |
| 214 | Rhodotorula minuta var minuta | N/A | N/A | N | BUY | Air | 26 |
| 215 | Rhodotorula mucilaginosa | N/A | N/A | N | BUY | Air | 26 |
| 216 | Rhodotorula muscorum | N/A | N/A | N | BUY | Air | 26 |
| 217 | Rhodotorula philyla | N/A | N/A | N | BUY | Air | 26 |
| 218 | Rhodotorula pustula | N/A | N/A | N | BUY | Air | 26 |
| 219 | Saccharomyces bayanus | N/A | N/A | N | BUY | Air | 26 |
| 220 | Saccharomyces boulardii | N/A | N/A | N | BUY | Air | 26 |
| 221 | Saccharomyces cerevisiae A/Tor.pretorien | N/A | N/A | N | BUY | Air | 26 |
| 222 | Saccharomyces cerevisiae B | N/A | N/A | N | BUY | Air | 26 |
| 223 | Saccharomyces dairensis | N/A | N/A | N | BUY | Air | 26 |
| 224 | Saccharomycodes ludwigii | N/A | N/A | N | BUY | Air | 26 |
| 225 | Saccharomycopsis capsularis | N/A | N/A | N | BUY | Air | 26 |
| 226 | Saturnospora dispora | N/A | N/A | N | BUY | Air | 26 |
| 227 | Schizoblastosporon starkeyi-henricii | N/A | N/A | N | BUY | Air | 26 |
| 228 | Schizosaccharomyces japonicus | N/A | N/A | N | BUY | Air | 26 |
| 229 | Schizosaccharomyces japonicus var japoni | N/A | N/A | N | BUY | Air | 26 |
| 230 | Schizosaccharomyces octosporus | N/A | N/A | N | BUY | Air | 26 |
| 231 | Schizosaccharomyces pombe | N/A | N/A | N | BUY | Air | 26 |
| 232 | Schizosaccharomyces pombe var malidevora | N/A | N/A | N | BUY | Air | 26 |
| 233 | Schwanniomyces occidentalis | N/A | N/A | N | BUY | Air | 26 |
| 234 | Sporidiobolus johnsonii A | N/A | N/A | N | BUY | Air | 26 |
| 235 | Sporidiobolus johnsonii B | N/A | N/A | N | BUY | Air | 26 |
| 236 | Sporidiobolus johnsonii C | N/A | N/A | N | BUY | Air | 26 |
| 237 | Sporidiobolus pararoseus A | N/A | N/A | N | BUY | Air | 26 |
| 238 | Sporidiobolus pararoseus B | N/A | N/A | N | BUY | Air | 26 |
| 239 | Sporobolomyces albo-rubescens | N/A | N/A | N | BUY | Air | 26 |
| 240 | Sporopachydermia cereana | N/A | N/A | N | BUY | Air | 26 |
| 241 | Sporopachydermia lactativora | N/A | N/A | N | BUY | Air | 26 |
| 242 | Stephanoascus ciferrii | N/A | N/A | N | BUY | Air | 26 |
| 243 | Sterigmatomyces elviae | N/A | N/A | N | BUY | Air | 26 |
| 244 | Sterigmatomyces halophilus | N/A | N/A | N | BUY | Air | 26 |
| 245 | Torulaspora delbrueckii | N/A | N/A | N | BUY | Air | 26 |
| 246 | Torulaspora globosa | N/A | N/A | N | BUY | Air | 26 |
| 247 | Trichosporon beigelii A | N/A | N/A | N | BUY | Air | 26 |
| 248 | Trichosporon beigelii B | N/A | N/A | N | BUY | Air | 26 |
| 249 | Trichosporon brassicae | N/A | N/A | N | BUY | Air | 26 |
| 250 | Trichosporon inkin | N/A | N/A | N | BUY | Air | 26 |
| 251 | Trigonopsis triangularis/variabilis | N/A | N/A | N | BUY | Air | 26 |
| 252 | Ustilago maydis | N/A | N/A | N | BUY | Air | 26 |
| 253 | Wickerhamiella domercqiae | N/A | N/A | N | BUY | Air | 26 |
| 254 | Williopsis californica | N/A | N/A | N | BUY | Air | 26 |
| 255 | Williopsis saturnus | N/A | N/A | N | BUY | Air | 26 |
| 256 | Williopsis saturnus var mrakii | N/A | N/A | N | BUY | Air | 26 |
| 257 | Williopsis saturnus var saturnus | N/A | N/A | N | BUY | Air | 26 |
| 258 | Wingea robertsiae | N/A | N/A | N | BUY | Air | 26 |
| 259 | Yarrowia lipolytica | N/A | N/A | N | BUY | Air | 26 |
| 260 | Zygoascus hellenicus | N/A | N/A | N | BUY | Air | 26 |
| 261 | Zygosaccharomyces bailii | N/A | N/A | N | BUY | Air | 26 |
| 262 | Zygosaccharomyces bisporus | N/A | N/A | N | BUY | Air | 26 |
| 263 | Zygosaccharomyces cidri | N/A | N/A | N | BUY | Air | 26 |
| 264 | Zygosaccharomyces fermentati | N/A | N/A | N | BUY | Air | 26 |
| 265 | Zygosaccharomyces florentinus | N/A | N/A | N | BUY | Air | 26 |
| 266 | Zygosaccharomyces microellipsoides/mraki | N/A | N/A | N | BUY | Air | 26 |
| 267 | Zygosaccharomyces rouxii | N/A | N/A | N | BUY | Air | 26 |

Filamentous Fungi (and select Yeasts)

| | Species Name | Type | Medium | Atm | Temp | Photo |
|----|--|-------------------|--------|-----|------|-------|
| 1 | Absidia californica J.Ellis & Hesselstine | FF-AIR | 2% ME | Air | 26 | Mm |
| 2 | Absidia corymbifera (Cohn) Sacc. & Trotter | FF-AIR, FOOD | 2% ME | Air | 26 | |
| 3 | Absidia clindrospora var clindrospora Hagen | FF-AIR | 2% ME | Air | 26 | m |
| 4 | Absidia glauca Hagen | FF-AIR | 2% ME | Air | 26 | Mm |
| 5 | Absidia spinosa var spinosa Lendner | FF-AIR | 2% ME | Air | 26 | |
| 6 | Acremonium atrogrisum (Panasenko) W. Gams | FF-AIR | 2% ME | Air | 26 | M |
| 7 | Acremonium bacillisporum (Onions & Barron) W. Gams | FF-AIR | 2% ME | Air | 26 | M |
| 8 | Acremonium berkeleyanum (P. Karsten) W. Gams | FF-FOOD | 2% ME | Air | 26 | M |
| 9 | Acremonium breve (Sukapure & Thirumalachar) W. Gams BGA | FF-AIR | 2% ME | Air | 26 | M |
| 10 | Acremonium breve (Sukapure & Thirumalachar) W. Gams BGB | FF-AIR | 2% ME | Air | 26 | M |
| 11 | Acremonium charticola (Lindau) W. Gams | FF-AIR, FOOD | 2% ME | Air | 26 | M |
| 12 | Acremonium furcatum (F. & V. Moreau) ex W. Gams | FF-AIR | 2% ME | Air | 26 | M |
| 13 | Acremonium fusidioides (Nicot) W. Gams | FF-AIR | 2% ME | Air | 26 | M |
| 14 | Acremonium hyalinulum | FF-AIR | 2% ME | Air | 26 | |
| 15 | Acremonium kiliense Gruetz | FF-AIR | 2% ME | Air | 26 | M |
| 16 | Acremonium murorum var felinum (Marchal) S. Hughes | FF-AIR | 2% ME | Air | 26 | M |
| 17 | Acremonium murorum var murorum (Corda) W. Gams | FF-AIR | 2% ME | Air | 26 | |
| 18 | Acremonium polychromum (v. Beyma) W. Gams | FF-AIR | 2% ME | Air | 26 | |
| 19 | Acremonium potronii Vuill. | FF-AIR | 2% ME | Air | 26 | M |
| 20 | Acremonium strictum W. Gams BGA | FF-AIR, FOOD | 2% ME | Air | 26 | M |
| 21 | Acremonium strictum W. Gams BGB | FF-AIR, FOOD | 2% ME | Air | 26 | M |
| 22 | Acrodontium crateriforme (van Beyma) de Hoog | FF-AIR | 2% ME | Air | 26 | M |
| 23 | Acrodontium griseum (Fassatiava') de Hoog | FF-AIR | 2% ME | Air | 26 | |
| 24 | Acrodontium salmoneum de Hoog BGA | FF-AIR | 2% ME | Air | 26 | M |
| 25 | Acrodontium salmoneum de Hoog BGB | FF-AIR | 2% ME | Air | 26 | M |
| 26 | Acrodontium simplex (Mangenot) de Hoog | FF-AIR | 2% ME | Air | 26 | M |
| 27 | Actinomucor elegans (Eidam) C.R. Benjamin & Hesselstine BGA | FF-AIR | 2% ME | Air | 26 | Mm |
| 28 | Actinomucor elegans (Eidam) C.R. Benjamin & Hesselstine BGB | FF-AIR | 2% ME | Air | 26 | Mm |
| 29 | Alternaria alternata (Fr.) Keissl | FF-AIR, FOOD | 2% ME | Air | 26 | Mm |
| 30 | Amorphotheca resinae D. Parbery | FF-AIR | 2% ME | Air | 26 | M |
| 31 | Aphanocladium album (Preuss) W. Gams | FF-AIR | 2% ME | Air | 26 | M |
| 32 | Apiospora montagnei Sacc. | FF-AIR | 2% ME | Air | 26 | M |
| 33 | Arthrinium phacospermum (Corda) M.B. Ellis | FF-AIR, FOOD | 2% ME | Air | 26 | M |
| 34 | Aspergillus aculeatus Izuka | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 35 | Aspergillus aeneus Sappa | FF-ASP | 2% ME | Air | 26 | M |
| 36 | Aspergillus asperescens Stolk | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 37 | Aspergillus aureolatus Muntanola-Cvctkovic & Bata BGA | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 38 | Aspergillus aureolatus Muntanola-Cvctkovic & Bata BGB | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 39 | Aspergillus auricomus (Guegen) Saito | FF-ASP | 2% ME | Air | 26 | Mm |
| 40 | Aspergillus avenaceus G. Smith | FF-ASP | 2% ME | Air | 26 | M |
| 41 | Aspergillus awamori Nakazawa | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 42 | Aspergillus brevipes G. Smith | FF-ASP | 2% ME | Air | 26 | M |
| 43 | Aspergillus caesiellus Saito | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 44 | Aspergillus caespitosus Raper & Thom | FF-ASP | 2% ME | Air | 26 | M |
| 45 | Aspergillus candidus Link BGA | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 46 | Aspergillus candidus Link BGB | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 47 | Aspergillus carbonarius (Bainier) Thom | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 48 | Aspergillus carneus (V. Tiegham) Blockwitz | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 49 | Aspergillus cervinus Massee | FF-ASP | 2% ME | Air | 26 | M |
| 50 | Aspergillus clavatus Desm. | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 51 | Aspergillus egyptiacus Moubasher & Moustafa | FF-ASP | 2% ME | Air | 26 | M |
| 52 | Aspergillus elegans Gasperini | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 53 | Aspergillus flavofurcatus Batista & Maia | FF-ASP | 2% ME | Air | 26 | Mm |
| 54 | Aspergillus flavus var flavus Link | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 55 | Aspergillus foetidus Thom & Raper | FF-ASP | 2% ME | Air | 26 | M |
| 56 | Aspergillus fresenii Subramanian | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 57 | Aspergillus fumigatus Fresen. | FF-AIR, FOOD, ASP | 2% ME | Air | 26 | Mm |
| 58 | Aspergillus furniculosus G. Smith | FF-ASP | 2% ME | Air | 26 | M |
| 59 | Aspergillus giganteus Wehmer | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 60 | Aspergillus japonicus var japonicus K. Saito | FF-AIR, ASP | 2% ME | Air | 26 | Mm |
| 61 | Aspergillus kanagawaensis Nehira BGA | FF-ASP | 2% ME | Air | 26 | M |
| 62 | Aspergillus kanagawaensis Nehira BGB | FF-ASP | 2% ME | Air | 26 | M |
| 63 | Aspergillus lanosus Kamal & Bhargava | FF-ASP | 2% ME | Air | 26 | M |
| 64 | Aspergillus malignus Lindt | FF-ASP | 2% ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 65 | <i>Aspergillus niger</i> v. Tiegham BGA | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 66 | <i>Aspergillus niger</i> v. Tiegham BGB | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 67 | <i>Aspergillus nutans</i> McLennan & Ducker | FF-ASP | 2%ME | Air | 26 | Mm |
| 68 | <i>Aspergillus ochraceus</i> K. Wilhelm BGA | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 69 | <i>Aspergillus ochraceus</i> K. Wilhelm BGB | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 70 | <i>Aspergillus oryzae</i> var <i>oryzae</i> (Ahlburg) Cohn | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 71 | <i>Aspergillus ostianus</i> Wehmer | FF-AIR, ASP | 2%ME | Air | 26 | Mm |
| 72 | <i>Aspergillus pallidus</i> Kamyschko | FF-ASP | 2%ME | Air | 26 | Mm |
| 73 | <i>Aspergillus parasiticus</i> Speare BGA | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 74 | <i>Aspergillus parasiticus</i> Speare BGB | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 75 | <i>Aspergillus parvulus</i> G. Sm. BGA | FF-ASP | 2%ME | Air | 26 | M |
| 76 | <i>Aspergillus parvulus</i> G. Sm. BGB | FF-ASP | 2%ME | Air | 26 | M |
| 77 | <i>Aspergillus penicilliformis</i> Kamyschko | FF-ASP | 2%ME | Air | 26 | M |
| 78 | <i>Aspergillus petrakii</i> Voros | FF-ASP | 2%ME | Air | 26 | M |
| 79 | <i>Aspergillus phoenicis</i> (Corda) Thom | FF-ASP | 2%ME | Air | 26 | M |
| 80 | <i>Aspergillus proliferans</i> G. Smith | FF-ASP | 2%ME | Air | 26 | M |
| 81 | <i>Aspergillus pulverulentus</i> (McAlp.) Thom | FF-ASP | 2%ME | Air | 26 | Mm |
| 82 | <i>Aspergillus puniceus</i> Kwon & Fennell | FF-AIR, ASP | 2%ME | Air | 26 | Mm |
| 83 | <i>Aspergillus raperi</i> Stolk & J.A. Meyer | FF-AIR, ASP | 2%ME | Air | 26 | M |
| 84 | <i>Aspergillus restrictus</i> G. Smith | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 85 | <i>Aspergillus sclerotiorum</i> Huber | FF-AIR, ASP | 2%ME | Air | 26 | Mm |
| 86 | <i>Aspergillus sepultus</i> Tuthill & Christensen | FF-ASP | 2%ME | Air | 26 | M |
| 87 | <i>Aspergillus sparsus</i> Raper & Thom | FF-ASP | 2%ME | Air | 26 | M |
| 88 | <i>Aspergillus speluneus</i> Raper & Fennell | FF-ASP | 2%ME | Air | 26 | Mm |
| 89 | <i>Aspergillus sydowii</i> (Brainier & Sartory) Thom & Church | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 90 | <i>Aspergillus tamaritii</i> Kita | FF-FOOD, ASP | 2%ME | Air | 26 | M |
| 91 | <i>Aspergillus terreus</i> var <i>africanus</i> Fennell & Rapier | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 92 | <i>Aspergillus terreus</i> var <i>terreus</i> Thom | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 93 | <i>Aspergillus terricola</i> var <i>americanus</i> Marchal | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | m |
| 94 | <i>Aspergillus terricola</i> var <i>terricola</i> Marchal | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | m |
| 95 | <i>Aspergillus ustus</i> (Bainier) Thom & Church | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 96 | <i>Aspergillus versicolor</i> (Vuill.) Tirab. | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 97 | <i>Aspergillus violaceofuscus</i> Gasperini | FF-ASP | 2%ME | Air | 26 | m |
| 98 | <i>Aspergillus wentii</i> Wermer | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 99 | <i>Aspergillus zonatus</i> Zonatus | FF-ASP | 2%ME | Air | 26 | M |
| 100 | <i>Aureobasidium pullulans</i> var <i>melanigenum</i> (de Barry) G. Arnaud | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 101 | <i>Aureobasidium pullulans</i> var <i>pullulans</i> (de Barry) G. Arnaud | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 102 | <i>Basipetospora halophila</i> (J.F.H. Beyma) Pitt & A.D. Hocking | FF-FOOD | 2%ME | Air | 26 | M |
| 103 | <i>Beauveria bassiana</i> | FF-AIR | 2%ME | Air | 26 | M |
| 104 | <i>Blastobotrys elegans</i> de Hoog et al. BGA | FF-AIR | 2%ME | Air | 26 | M |
| 105 | <i>Blastobotrys elegans</i> de Hoog et al. BGB | FF-AIR | 2%ME | Air | 26 | M |
| 106 | <i>Botryosporium longibrachiatum</i> (Oud.) Maire | FF-AIR | 2%ME | Air | 26 | M |
| 107 | <i>Botrytis aclada</i> Fresen. | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 108 | <i>Botrytis cinerea</i> Pers. BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 109 | <i>Botrytis cinerea</i> Pers. BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 110 | <i>Bulleromyces albus</i> | FF-AIR, YST | 2%ME | Air | 26 | Mm |
| 111 | <i>Byssochlamys fulva</i> Olliver & G. Smith BGA | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 112 | <i>Byssochlamys fulva</i> Olliver & G. Smith BGB | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 113 | <i>Byssochlamys nivea</i> Westling | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 114 | <i>Candida albicans</i> (Robin) Berkhout | FF-YST | 2%ME | Air | 26 | |
| 115 | <i>Candida guilliermondii</i> (Castellani) Langeron & Guerra Kodamaea ohmeri) (teleo. Kodamaea ohmeri) | FF-FOOD, YST | 2%ME | Air | 26 | |
| 116 | <i>Candida insectorum</i> D.B. Scott et al. | FF-YST | 2%ME | Air | 26 | |
| 117 | <i>Candida intermedia</i> (AS) (Ciferri & Ashford) Langeron & Guerra | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 118 | <i>Candida montana</i> S. Goto & Oguri | FF-YST | 2%ME | Air | 26 | |
| 119 | <i>Candida parapsilosis</i> (Ashford) Langeron & Talice | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 120 | <i>Candida rugosa</i> (H.W. Anderson) Diddens & Lodder | FF-YST | 2%ME | Air | 26 | |
| 121 | <i>Candida sake</i> (Saito & Oda) v. Uden & Buckley | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 122 | <i>Candida sorbophila</i> (Nakase) S.A. Meyer & Yarrow | FF-YST | 2%ME | Air | 26 | |
| 123 | <i>Candida zeylanoides</i> (Castellani) Langeron & Guerra | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 124 | <i>Cephalophora tropica</i> Thaxler | FF-AIR | 2%ME | Air | 26 | Mm |
| 125 | <i>Cerinosterus cyanesens</i> (de Hoog & de Vries) R.T. Moore | FF-AIR | 2%ME | Air | 26 | |
| 126 | <i>Chaetomium funicula</i> Cooke | FF-FOOD | 2%ME | Air | 26 | |
| 127 | <i>Chaetomium globosum</i> Kunze: Fries | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 128 | <i>Chaetosartorya stromatoides</i> Wiley & Simmons | FF-ASP | 2%ME | Air | 26 | M |
| 129 | <i>Chrysosporium indicum</i> (Randhawa & Sandhu) Garg (anam. <i>Aspergillus stromatoides</i>) | FF-AIR | 2%ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 130 | Chrysosporium inops J.W. Carmich. | FF-FOOD | 2%ME | Air | 26 | M |
| 131 | Chrysosporium keratinophilum D.Frey ex Charmichael | FF-FOOD | 2%ME | Air | 26 | |
| 132 | Chrysosporium merdarium (Link: Fries) Carmichael. | FF-AIR | 2%ME | Air | 26 | |
| 133 | Chrysosporium pannicola (Corda) van Oorschot & Stalpers | FF-AIR | 2%ME | Air | 26 | |
| 134 | Cladosporium cladosporioides (Fresen.) G.A. de Vries | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 135 | Cladosporium herbarum (Pers.) Link BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 136 | Cladosporium herbarum (Pers.) Link BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 137 | Cladosporium macrocarpum Preuss | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 138 | Cladosporium sphaerospermum Penz. BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 139 | Cladosporium sphaerospermum Penz. BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 140 | Cladosporium tenuissimum Cooke | FF-AIR | 2%ME | Air | 26 | M |
| 141 | Clonostachys rosea (Link: Fr.) Schroers, Samuels, Seifert & W. Gams | FF-AIR | 2%ME | Air | 26 | Mm |
| 142 | Colletotrichum acutatum Simmonds | FF-COL | 2%ME | Air | 26 | m |
| 143 | Colletotrichum capsici (H. Sydow) E. Butler & Bisby | FF-COL | 2%ME | Air | 26 | |
| 144 | Colletotrichum coccodes (Wallroth) Hughes | FF-COL | 2%ME | Air | 26 | |
| 145 | Colletotrichum coffeanum Noack sensu Hindorf | FF-COL | 2%ME | Air | 26 | |
| 146 | Colletotrichum crassipes (Spegazzini) von Arx | FF-COL | 2%ME | Air | 26 | |
| 147 | Colletotrichum dematium (Persoon: Fries) Grove | FF-COL | 2%ME | Air | 26 | |
| 148 | Colletotrichum destructivum O'Gara | FF-COL | 2%ME | Air | 26 | |
| 149 | Colletotrichum gloeosporioides (Penzig) Penzig & Saccardo | FF-FOOD, COL | 2%ME | Air | 26 | M |
| 150 | Colletotrichum lindenuhianum (Saccardo & Magnus) Briosi | FF-COL | 2%ME | Air | 26 | |
| 151 | Colletotrichum trichellum (Fries) Duke | FF-COL | 2%ME | Air | 26 | |
| 152 | Colletotrichum trifolii Bain & Essary | FF-COL | 2%ME | Air | 26 | |
| 153 | Colletotrichum truncatum (Schweinitz) Andrus & Moore BGA | FF-COL | 2%ME | Air | 26 | |
| 154 | Colletotrichum truncatum (Schweinitz) Andrus & Moore BGB | FF-COL | 2%ME | Air | 26 | |
| 155 | Cryptococcus albidus var albidus (Saito) Skinner | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 156 | Cryptococcus laurentii var laurentii (Kufferath) Skinner | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 157 | Cryptococcus terreus di Menna | FF-YST | 2%ME | Air | 26 | |
| 158 | Cunninghamella elegans Lendner | FF-AIR, FOOD | 2%ME | Air | 26 | |
| 159 | Curvularia lunata var lunata (Wakker) Boedijn BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 160 | Curvularia lunata var lunata (Wakker) Boedijn BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 161 | Debaryomyces castellii Capriotti | FF-YST | 2%ME | Air | 26 | |
| 162 | Debaryomyces hansenii var hansenii (Zopl) Lodder & Kreger-v.Rij BGA (anam. CANADA FAMATA) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 163 | Debaryomyces hansenii var hansenii (Zopl) Lodder & Kreger-v.Rij BGB (anam. CANADA FAMATA) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 164 | Debaryomyces hansenii var hansenii (Zopl) Lodder & Kreger-v.Rij BGC (anam. CANADA FAMATA) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 165 | Debaryomyces hansenii var hansenii (Zopl) Lodder & Kreger-v.Rij BGD (anam. CANADA FAMATA) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 166 | Debaryomyces occidentalis var occidentalis (Klocker) Kurtzman & Robnett | FF-YST | 2%ME | Air | 26 | |
| 167 | Debaryomyces polymorphus var polymorphus (Klocker) C.W. Price & Phaff | FF-YST | 2%ME | Air | 26 | |
| 168 | Dekkera bruxellensis van der Walt (anam. Brettanomyces bruxellensis) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 169 | Doratomyces microsporus (Sacc.) Morton & G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 170 | Doratomyces stemonitis (Pers.:Fr) Morton & G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 171 | Emericella fruticulosa (Raper & Fennel) Malloch & Cain (anam. Aspergillus triticans) | FF-ASP | 2%ME | Air | 26 | M |
| 172 | Emericella nidulans var nidulans (Eidam) Vuillemin (anam. Aspergillus nidulellus) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 173 | Emericella quadrilineata (Thom & Raper) C.R. Benjamin (anam. Aspergillus tetrazonus) | FF-ASP | 2%ME | Air | 26 | M |
| 174 | Emericella rugulosa (Thom & Raper) C.R. Benjamin (anam. Aspergillus rugulovalvus) | FF-ASP | 2%ME | Air | 26 | M |
| 175 | Emericella striata (J.N. Rai et al.) Malloch & Cain (anam. Aspergillus striatulus) | FF-ASP | 2%ME | Air | 26 | M |
| 176 | Emericella unguis Malloch & Cain BGA (anam. Aspergillus unguis) | FF-ASP | 2%ME | Air | 26 | M |
| 177 | Emericella unguis Malloch & Cain BGB (anam. Aspergillus unguis) | FF-ASP | 2%ME | Air | 26 | M |
| 178 | Emericella varicolor var varicolor Berk. & Br. In Berk. (anam. Aspergillus stellifer) | FF-ASP | 2%ME | Air | 26 | M |
| 179 | Emericella violacea (Fennell & Raper) Malloch & Cain (anam. Aspergillus violacobrunneus) | FF-ASP | 2%ME | Air | 26 | M |
| 180 | Endomyces fibuliger Lindner | FF-AIR, FOOD, YST | 2%ME | Air | 26 | Mm |
| 181 | Engyodontium album (Limber) de Hoog | FF-AIR | 2%ME | Air | 26 | M |
| 182 | Epicoccum nigrum (purapurascens) Link | FF-AIR | 2%ME | Air | 26 | M |
| 183 | Eremascus fertilis Stoppel | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 184 | Eupenicillium baarnense (V.Beyma)Stolk&Scott (anam. Penicillium turbatum) | FF-PEN | 2%ME | Air | 26 | M |
| 185 | Eupenicillium brefeldianum (Dodge)Stolk&Scott (anam. Penicillium dodgei) | FF-PEN | 2%ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|--|-------------------|--------|-----|------|-------|
| 186 | Eupenicillium cinnamopurpureum D.B. Scott & Stolk BGA (anam. <i>Penicillium cinnamopurpureum</i>) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 187 | Eupenicillium cinnamopurpureum D.B. Scott & Stolk BGB (anam. <i>Penicillium cinnamopurpureum</i>) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 188 | Eupenicillium crustaceum Ludwig (anam. <i>Penicillium nilense</i>) | FF-AIR, PEN | 2%ME | Air | 26 | M |
| 189 | Eupenicillium euglaucum (van Beyma) Stolk & Samson BGA (anam. <i>Penicillium citreonigrum</i>) | FF-PEN | 2%ME | Air | 26 | Mm |
| 190 | Eupenicillium euglaucum (van Beyma) Stolk & Samson BGB (anam. <i>Penicillium citreonigrum</i>) | FF-PEN | 2%ME | Air | 26 | Mm |
| 191 | Eupenicillium javanicum var javanicum (Van Beyma) Stolk & Scott (<i>P. simplicissimum</i> p.p.) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 192 | Eupenicillium lapidosum Scott & Stolk (anam. <i>Penicillium lapidosum</i>) | FF-PEN | 2%ME | Air | 26 | M |
| 193 | Eupenicillium meridianum Scott (anam. <i>Penicillium decumbens</i>) | FF-PEN | 2%ME | Air | 26 | M |
| 194 | Eupenicillium pinetorum Stolk (anam. <i>Penicillium fuscum</i>) | FF-PEN | 2%ME | Air | 26 | M |
| 195 | Eurotium amstelodami Mangin (anam. <i>Aspergillus hollandicus</i>) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 196 | Eurotium chevalieri L.Mangin BGA (anam. <i>Aspergillus equitis</i>) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 197 | Eurotium chevalieri L.Mangin BGB (anam. <i>Aspergillus equitis</i>) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 198 | Eurotium echinulatum Delacroix (anam. <i>Aspergillus echinulatus</i>) | FF-AIR, ASP | 2%ME | Air | 26 | M |
| 199 | Eurotium herbariorum (Wiggers:Fr.) Link (anam. <i>Aspergillus glaucus</i>) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 200 | Eurotium intermedium Blaser (anam. <i>Aspergillus intermedium</i>) | FF-ASP | 2%ME | Air | 26 | Mm |
| 201 | Eurotium repens | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 202 | Eurotium rubrum Konig et al. (anam. <i>Aspergillus rubrobrunneus</i>) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | Mm |
| 203 | Eurotium spiculosum Blaser BGA (anam. <i>Aspergillus spiculosum</i>) | FF-ASP | 2%ME | Air | 26 | |
| 204 | Eurotium spiculosum Blaser BGB (anam. <i>Aspergillus spiculosum</i>) | FF-ASP | 2%ME | Air | 26 | |
| 205 | Fennellia flavipes Wiley & Simmons (anam. <i>Aspergillus flavipes</i>) | FF-FOOD, ASP | 2%ME | Air | 26 | M |
| 206 | Fennellia nivea (Wiley & Simmons) Samson BGA (anam. <i>Aspergillus niveus</i>) | FF-FOOD, ASP | 2%ME | Air | 26 | M |
| 207 | Fennellia nivea (Wiley & Simmons) Samson BGB (anam. <i>Aspergillus niveus</i>) | FF-FOOD, ASP | 2%ME | Air | 26 | M |
| 208 | Fusarium acuminatum Ellis & Everh. BGA (teleo. <i>Gibberella acuminatum</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 209 | Fusarium acuminatum Ellis & Everh. BGB (teleo. <i>Gibberella acuminatum</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 210 | Fusarium annulatum Bugnicourt | FF-FUS | 2%ME | Air | 26 | M |
| 211 | Fusarium avenaceum (Corda:Fr.) Sacc. BGA (teleo. <i>Gibberella avenacea</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 212 | Fusarium avenaceum (Corda:Fr.) Sacc. BGB (teleo. <i>Gibberella avenacea</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 213 | Fusarium avenaceum (Corda:Fr.) Sacc. BGC (teleo. <i>Gibberella avenacea</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 214 | Fusarium avenaceum s sp aywerte Sangalang & L.W. Burgess | FF-FUS | 2%ME | Air | 26 | Mm |
| 215 | Fusarium avenaceum s sp nurragai Summerell & L.W. Burgess | FF-FUS | 2%ME | Air | 26 | Mm |
| 216 | Fusarium brevicatenulatum Nirenberg, O'Donnell, Kroschel & Andrianairo | FF-FUS | 2%ME | Air | 26 | M |
| 217 | Fusarium camptoceras Wollenw. & Reinking | FF-FUS | 2%ME | Air | 26 | M |
| 218 | Fusarium chlamydosporum var chlamydosporum Wollenw. & Reinking | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 219 | Fusarium ciliatum Link | FF-FUS | 2%ME | Air | 26 | M |
| 220 | Fusarium compactum (Wollenw.) Raillo | FF-FUS | 2%ME | Air | 26 | M |
| 221 | Fusarium crookwellense Burgess, Nelson & Toussoun BGA | FF-FUS | 2%ME | Air | 26 | M |
| 222 | Fusarium crookwellense Burgess, Nelson & Toussoun BGB | FF-FUS | 2%ME | Air | 26 | M |
| 223 | Fusarium crookwellense Burgess, Nelson & Toussoun BGC | FF-FUS | 2%ME | Air | 26 | M |
| 224 | Fusarium culmorum (W.G. Smith) Sacc. BGA | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 225 | Fusarium culmorum (W.G. Smith) Sacc. BGB | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 226 | Fusarium decemcellulare Brick (teleo. <i>Nectria rigidiuscula</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 227 | Fusarium equiseti (Corda) Sacc. (teleo. <i>Gibberella intricans</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 228 | Fusarium eumartii Carpenter (teleo. <i>Nectria haematococca</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 229 | Fusarium flocciferum Corda (teleo. <i>Gibberella heterochroma</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 230 | Fusarium fujikuroi Nirenberg (teleo. <i>Gibberella fujikuroi</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 231 | Fusarium globosum Rheeder, Marasas & Nelson | FF-FUS | 2%ME | Air | 26 | M |
| 232 | Fusarium graminearum Schwabe BGA (teleo. <i>Gibberella zeae</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 233 | Fusarium graminearum Schwabe BGB (teleo. <i>Gibberella zeae</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 234 | Fusarium graminearum Schwabe BGC (teleo. <i>Gibberella zeae</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 235 | Fusarium graminearum Schwabe BGD (teleo. <i>Gibberella zeae</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 236 | Fusarium graminearum Schwabe BGE (teleo. <i>Gibberella zeae</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 237 | Fusarium graminum Corda | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 238 | Fusarium heterosporum Nees BGA (teleo. <i>Gibberella gordonia</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 239 | Fusarium heterosporum Nees BGB (teleo. <i>Gibberella gordonia</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 240 | Fusarium heterosporum Nees BGC (teleo. <i>Gibberella gordonia</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 241 | Fusarium juruanum P. Henn. (teleo. <i>Nectria diplos</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 242 | Fusarium lactis Pirodda & Riboni | FF-FUS | 2%ME | Air | 26 | M |
| 243 | Fusarium lateritium var lateritium Nees (teleo. <i>Gibberella baccata</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 244 | Fusarium longipes Wollenw. & Reinking | FF-FOOD, FUS | 2%ME | Air | 26 | Mm |
| 245 | Fusarium melanochlorum (Caspary) Sacc. (teleo. <i>Nectria flavoviridis</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 246 | Fusarium merismoides VAR MERISMOIDES Corda | FF-AIR, FUS | 2%ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 247 | Fusarium oxysporum Schlect.:Fr. BGA | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 248 | Fusarium oxysporum Schlect.:Fr. BGB | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 249 | Fusarium pallidoroseum (Cooke) Sacc. | FF-FUS | 2%ME | Air | 26 | M |
| 250 | Fusarium poae (Peck) Wollenw. BGA | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 251 | Fusarium poae (Peck) Wollenw. BGB | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 252 | Fusarium proliferatum var proliferatum (Matsush.) Nirenberg BGA (teleo. <i>Gibberella intermedia</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 253 | Fusarium proliferatum var proliferatum (Matsush.) Nirenberg BGB (teleo. <i>Gibberella intermedia</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 254 | Fusarium pseudoanthophilum Nirenberg, O'Donnell & Mabatanhema | FF-FUS | 2%ME | Air | 26 | M |
| 255 | Fusarium pseudograminearum Aoki & O'Donnell (teleo. <i>Gibberella coronicola</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 256 | Fusarium ramigenum O'Donnell & Nirenberg | FF-FUS | 2%ME | Air | 26 | M |
| 257 | Fusarium redolens Wollenw. | FF-FUS | 2%ME | Air | 26 | M |
| 258 | Fusarium reticulatum Mont. (teleo. <i>Gibberella cyanea</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 259 | Fusarium robustum Gerlach | FF-FUS | 2%ME | Air | 26 | M |
| 260 | Fusarium sacchari (Butler) W. Gams BGA | FF-FOOD, FUS | 2%ME | Air | 26 | M |
| 261 | Fusarium sacchari (Butler) W. Gams BGB | FF-FOOD, FUS | 2%ME | Air | 26 | M |
| 262 | Fusarium sambucinum var sambucinum Fuckel (teleo. <i>Gibberella pulicaris</i>) | FF-FOOD, FUS | 2%ME | Air | 26 | Mm |
| 263 | Fusarium solani (Mart.) Sacc. BGA | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 264 | Fusarium solani (Mart.) Sacc. BGB | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | Mm |
| 265 | Fusarium sporotrichioides var minus Sherb. BGA | FF-AIR, FUS | 2%ME | Air | 26 | Mm |
| 266 | Fusarium sporotrichioides var minus Sherb. BGB | FF-AIR, FUS | 2%ME | Air | 26 | Mm |
| 267 | Fusarium sporotrichioides var sporotrichioides Sherb. BGA | FF-AIR, FUS | 2%ME | Air | 26 | Mm |
| 268 | Fusarium sporotrichioides var sporotrichioides Sherb. BGB | FF-AIR, FUS | 2%ME | Air | 26 | Mm |
| 269 | Fusarium stilboides Wollenw. BGA (teleo. <i>Gibberella stilboides</i>) | FF-FOOD, FUS | 2%ME | Air | 26 | Mm |
| 270 | Fusarium stilboides Wollenw. BGB (teleo. <i>Gibberella stilboides</i>) | FF-FOOD, FUS | 2%ME | Air | 26 | Mm |
| 271 | Fusarium stilboides Wollenw. BGC (teleo. <i>Gibberella stilboides</i>) | FF-FOOD, FUS | 2%ME | Air | 26 | Mm |
| 272 | Fusarium subglutinans (Wollenw. & Reinking) P.E. Nelson et al. (teleo. <i>Gibberella subglutinans</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 273 | Fusarium thapsinum Klittich, Leslie, Nelson & Marasas (teleo. <i>Gibberella thapsinum</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 274 | Fusarium torulosum (Berk. & Curt.) Nirenberg | FF-FUS | 2%ME | Air | 26 | M |
| 275 | Fusarium tricinctum (Corda) Sacc. BGA | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 276 | Fusarium tricinctum (Corda) Sacc. BGB | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 277 | Fusarium tricinctum (Corda) Sacc. BGC | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 278 | Fusarium tricinctum (Corda) Sacc. BGD | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 279 | Fusarium tumidum Scherbakoff (teleo. <i>Gibberella tumida</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 280 | Fusarium udum E. Butler | FF-FUS | 2%ME | Air | 26 | M |
| 281 | Fusarium venenatum Nirenberg | FF-FUS | 2%ME | Air | 26 | M |
| 282 | Fusarium ventricosum Appel & Wollenw. | FF-FUS | 2%ME | Air | 26 | M |
| 283 | Fusarium verticillioides (Sacc.) Nirenberg BGA (teleo. <i>Gibberella moniliformis</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 284 | Fusarium verticillioides (Sacc.) Nirenberg BGB (teleo. <i>Gibberella moniliformis</i>) | FF-AIR, FOOD, FUS | 2%ME | Air | 26 | M |
| 285 | Fusarium xylarioides Steyaert (teleo. <i>Gibberella xylarioides</i>) | FF-FUS | 2%ME | Air | 26 | M |
| 286 | Galactomyces geotrichum (SS)(Butler & Perterson) Redhead & Malloch BGA (anam. <i>Geotrichum candidum</i>) | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 287 | Galactomyces geotrichum (SS)(Butler & Perterson) Redhead & Malloch BGB (anam. <i>Geotrichum candidum</i>) | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 288 | Galactomyces geotrichum (SS)(Butler & Perterson) Redhead & Malloch BGC (anam. <i>Geotrichum candidum</i>) | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 289 | Galactomyces geotrichum (SS)(Butler & Perterson) Redhead & Malloch BGD (anam. <i>Geotrichum candidum</i>) | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 290 | Geomyces pannorum var pannorum (Link) Sigler & Carmichael BGA | FF-AIR | 2%ME | Air | 26 | M |
| 291 | Geomyces pannorum var pannorum (Link) Sigler & Carmichael BGB | FF-AIR | 2%ME | Air | 26 | M |
| 292 | Geomyces pulvereus Hocking & Pitt | FF-AIR | 2%ME | Air | 26 | |
| 293 | Geosmithia putterilli (Thom) Pitt | FF-FOOD, PEN | 2%ME | Air | 26 | M |
| 294 | Geotrichum klebahnii (Stautz) Morenz | FF-AIR | 2%ME | Air | 26 | |
| 295 | Geotrichum sericeum (Stautz) de Hoog et al. | FF-YST | 2%ME | Air | 26 | |
| 296 | Gibberella subglutinans (Edwards) Nelson, Toussoun & Marasas (anam. <i>Fusarium subglutinans</i>) | FF-FUS | 2%ME | Air | 26 | |
| 297 | Gliocladium catenulatum Gilman & Abbott | FF-AIR | 2%ME | Air | 26 | M |
| 298 | Gliocladium viride Matruchot | FF-AIR | 2%ME | Air | 26 | Mm |
| 299 | Gonatobotrys simplex Corda | FF-AIR | 2%ME | Air | 26 | M |
| 300 | Hanseniaspora uvarum (Niehaus) El-Taby Shehata et al. (anam. <i>Kloeckera apiculata</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|--|-------------------|--------|-----|------|-------|
| 301 | Hanseniaspora valbyensis Kloecker (anam. Kloeckera japonica) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 302 | Hemicarpenales paradoxus Sarbhoy & Elphick (anam. Aspergillus paradoxus) | FF-ASP | 2%ME | Air | 26 | M |
| 303 | Hortaea werneckii (Horta) Nishimura & Miyaji | FF-FOOD | 2%ME | Air | 26 | |
| 304 | Hyphopichia burtonii (Boidin et al.) Arx & Van der Walt (anam. Candida chodatii) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 305 | Hypocrea andinensis Samuels & O. Petrini | FF-TRI | 2%ME | Air | 26 | M |
| 306 | Hypocrea aurantia P. Hennings | FF-TRI | 2%ME | Air | 26 | M |
| 307 | Hypocrea aureoviridis Plowright & Cooke (anam. Trichoderma aureoviride) | FF-TRI | 2%ME | Air | 26 | |
| 308 | Hypocrea gelatinosa (Tode: Fr.) Fr. | FF-TRI | 2%ME | Air | 26 | M |
| 309 | Hypocrea jecorina Berk. & Broome BGA (anam. Trichoderma reesei) | FF-TRI | 2%ME | Air | 26 | M |
| 310 | Hypocrea jecorina Berk. & Broome BGB (anam. Trichoderma reesei) | FF-TRI | 2%ME | Air | 26 | M |
| 311 | Hypocrea nigricans (Imai) Doi | FF-TRI | 2%ME | Air | 26 | M |
| 312 | Hypocrea novaezelandiae Samuels & O. Petrini | FF-TRI | 2%ME | Air | 26 | M |
| 313 | Hypocrea orientalis Samuels & O. Petrini | FF-TRI | 2%ME | Air | 26 | M |
| 314 | Hypocrea pseudokoningii Samuels & O. Petrini | FF-TRI | 2%ME | Air | 26 | M |
| 315 | Hypocrea schweinitzii (Fr.) Saccardo (anam. Trichoderma citrinoviride) | FF-TRI | 2%ME | Air | 26 | M |
| 316 | Issatchenkia orientalis Kudryavtsev (anam. Candida krusei) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 317 | Issatchenkia scutulata (Phaff et al.) Kurtzman et al. | FF-YST | 2%ME | Air | 26 | |
| 318 | Kluyveromyces lactis var lactis (Dombrowski) Van der Walt (anam. Candida sphaerica) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 319 | | | | | | |
| 320 | Kluyveromyces marxianus (Hansen) Van der Walt BGA (anam. Candida kefyri) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 321 | Kluyveromyces marxianus (Hansen) Van der Walt BGB (anam. Candida kefyri) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 322 | Kluyveromyces marxianus (Hansen) Van der Walt BGC (anam. Candida kefyri) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 323 | Kluyveromyces marxianus (Hansen) Van der Walt BGD (anam. Candida kefyri) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 324 | Kluyveromyces marxianus (Hansen) Van der Walt BGE (anam. Candida kefyri) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 325 | Lecythophora hoffmannii (van.Beyma) W. Gams BGA | FF-FOOD | 2%ME | Air | 26 | M |
| 326 | Lecythophora hoffmannii (van.Beyma) W. Gams BGB | FF-FOOD | 2%ME | Air | 26 | M |
| 327 | Lecythophora mutabilis (van Beyma) W. Gams | FF-AIR | 2%ME | Air | 26 | M |
| 328 | Mariannaea elegans var elegans (Corda) G. Amand ex Samson | FF-AIR | 2%ME | Air | 26 | Mm |
| 329 | Mariannaea elegans var punicea Samson | FF-AIR | 2%ME | Air | 26 | M |
| 330 | Microdochium bolleyi (Sprague) de Hoog & Herm.-Nijhof | FF-AIR | 2%ME | Air | 26 | M |
| 331 | Microsphaeropsis olivacea (Bonorden) v. Hoehnel | FF-AIR | 2%ME | Air | 26 | M |
| 332 | Monascus ruber v. Tiegham | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 333 | Moniliella acetoabutens Stolk & Dakin | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 334 | Moniliella suaveolens var suaveolens (Lindner) v.Arxx | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 335 | Mucor circinelloides f sp circinelloides v. Tiegham | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 336 | Mucor hiemalis f sp hiemalis Wehmer | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 337 | Mucor mucedo Linn.: Fr. | FF-AIR | 2%ME | Air | 26 | M |
| 338 | Mucor plumbeus Bonord. | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 339 | Mucor racemosus f sp racemosus Fresen. | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 340 | Myceliophthora lutea Costantin | FF-AIR | 2%ME | Air | 26 | |
| 341 | Myrothecium roridum Tode:Fr. | FF-AIR | 2%ME | Air | 26 | M |
| 342 | Myrothecium verrucaria (Alb. & Schwein.:Fr.) Ditmar | FF-AIR | 2%ME | Air | 26 | M |
| 343 | Nectria aureofulva Cooke & Ellis (anam. Clonostachys sp.) | FF-FUS | 2%ME | Air | 26 | M |
| 344 | Nectria hematococca Berk. & Br. (anam. Fusarium eumartii) | FF-FUS | 2%ME | Air | 26 | |
| 345 | Nectria ochroleuca (Schweinitz) Berkeley | FF-FUS | 2%ME | Air | 26 | |
| 346 | Nectria sesquicilli Samuels (anam. Clonostachys sp.) | FF-FUS | 2%ME | Air | 26 | M |
| 347 | Neosartorya fischeri var fischeri (Wehmer) Malloch & Cain (anam. Aspergillus .fischerianus) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 348 | Neurospora crassa Schear & Dodge | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 349 | Neurospora sitophila Schear & Dodge | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 350 | Nigrospora oryzae (Berk. & Broome) Petch | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 351 | Nigrospora sphaerica (sacc.) Mason | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 352 | Ophiostoma piceae (Munch) H. & P. Sydow BGA | FF-AIR | 2%ME | Air | 26 | Mm |
| 353 | Ophiostoma piceae (Munch) H. & P. Sydow BGB | FF-AIR | 2%ME | Air | 26 | Mm |
| 354 | Ophiostoma ulmi (Buisman) Nannfeldt | FF-AIR | 2%ME | Air | 26 | M |
| 355 | Paecilomyces carneus (Duché & Heim) A. Brown & G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 356 | Paecilomyces farinosus (Holm:Fr.) A. Brown & G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 357 | Paecilomyces fumosoroseus (Wize) A. Brown & G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 358 | Paecilomyces lilacinus (Thom) Samson | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 359 | Paecilomyces marquandii (Masse) S. Hughes | FF-AIR | 2%ME | Air | 26 | M |
| 360 | Paecilomyces variotii (Bainier) BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|--|-------------------|--------|-----|------|-------|
| 361 | Paecilomyces variotii (Bainer) BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 362 | Penicillium adametzii Zaleski BGA | FF-PEN | 2%ME | Air | 26 | Mm |
| 363 | Penicillium adametzii Zaleski BGB | FF-PEN | 2%ME | Air | 26 | Mm |
| 364 | Penicillium aethiopicum Frisvad | FF-FOOD, PEN | 2%ME | Air | 26 | M |
| 365 | Penicillium allii Vincent & Pitt | FF-FOOD, PEN | 2%ME | Air | 26 | Mm |
| 366 | Penicillium argillaceum Stolk et al. | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 367 | Penicillium atramentosum Thom | FF-PEN | 2%ME | Air | 26 | M |
| 368 | Penicillium aurantiogriseum Dierckx BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 369 | Penicillium aurantiogriseum Dierckx BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 370 | Penicillium aurantiovirens Biourge | FF-PEN | 2%ME | Air | 26 | Mm |
| 371 | Penicillium bilaiae Chalabuda | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 372 | Penicillium brevicompactum Dierckx BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 373 | Penicillium brevicompactum Dierckx BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 374 | Penicillium brevicompactum Dierckx BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 375 | Penicillium camemberti Thom BGA | FF-FOOD, PEN | 2%ME | Air | 26 | Mm |
| 376 | Penicillium camemberti Thom BGB | FF-FOOD, PEN | 2%ME | Air | 26 | Mm |
| 377 | Penicillium canescens Sopp | FF-PEN | 2%ME | Air | 26 | Mm |
| 378 | Penicillium capsulatum Raper & Fennell | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 379 | Penicillium chermesinum Biourge | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 380 | Penicillium chrysogenum Thom | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 381 | Penicillium citreonigrum Dierckx BGA (teleo. Eupenicillium euglaucum) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 382 | Penicillium citreonigrum Dierckx BGB (teleo. Eupenicillium euglaucum) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 383 | Penicillium citreonigrum Dierckx BGC (teleo. Eupenicillium euglaucum) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 384 | Penicillium citrinum Thom | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 385 | Penicillium clavigerum Demelius | FF-PEN | 2%ME | Air | 26 | M |
| 386 | Penicillium commune Thom BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 387 | Penicillium commune Thom BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 388 | Penicillium corylophilum Dierckx | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 389 | Penicillium crateriforme Gilman & Abbott | FF-PEN | 2%ME | Air | 26 | |
| 390 | Penicillium crustosum Thom | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 391 | Penicillium cyclopium Westling | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 392 | Penicillium decumbens Thom (teleo. Eupenicillium meridianum) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 393 | Penicillium digitatum (Pers.:Fr.) Sacc. BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 394 | Penicillium digitatum (Pers.:Fr.) Sacc. BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 395 | Penicillium discolor Frisvad & Samson | FF-PEN | 2%ME | Air | 26 | M |
| 396 | Penicillium duclauxii Deluc. | FF-PEN | 2%ME | Air | 26 | M |
| 397 | Penicillium echinulatum Fassatiava | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 398 | Penicillium erythromellis Hocking | FF-PEN | 2%ME | Air | 26 | Mm |
| 399 | Penicillium expansum Link BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 400 | Penicillium expansum Link BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 401 | Penicillium expansum Link BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 402 | Penicillium expansum Link BGD | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 403 | Penicillium fellutanum Biourge | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 404 | Penicillium freii Frisvad & Samson BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 405 | Penicillium freii Frisvad & Samson BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 406 | Penicillium freii Frisvad & Samson BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 407 | Penicillium funiculosum Thom BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 408 | Penicillium funiculosum Thom BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 409 | Penicillium glabrum (Wehmer) Westling | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 410 | Penicillium glandicola var glandicola (Oud.) Seifert & Samson | FF-PEN | 2%ME | Air | 26 | M |
| 411 | Penicillium griseofulvum Dierckx | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 412 | Penicillium herquei Bainier & Sartory | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 413 | Penicillium hirsutum Dierckx | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 414 | Penicillium hordei Stolk | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 415 | Penicillium implicatum Biourge | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | m |
| 416 | Penicillium islandicum Sopp | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 417 | Penicillium italicum Wehmer | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 418 | Penicillium janczewskii K.M. Zaleski | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 419 | Penicillium janthienellum Biourge | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 420 | Penicillium janssenii Zaleski | FF-PEN | 2%ME | Air | 26 | M |
| 421 | Penicillium lagena (Delitsch) Stolk & Samson | FF-PEN | 2%ME | Air | 26 | Mm |
| 422 | Penicillium lanosum Westling | FF-PEN | 2%ME | Air | 26 | Mm |
| 423 | Penicillium lividum Westling | FF-PEN | 2%ME | Air | 26 | M |
| 424 | Penicillium madriti G.Smith | FF-PEN | 2%ME | Air | 26 | M |
| 425 | Penicillium melanoconidium (Frisvad) Frisvad & Samson BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 426 | Penicillium melanoconidium (Frisvad) Frisvad & Samson BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 427 | Penicillium melinii Thom | FF-PEN | 2%ME | Air | 26 | Mm |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 428 | Penicillium miczynski Zaleski | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 429 | Penicillium minioluteum Dierckx BGA | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 430 | Penicillium minioluteum Dierckx BGB | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 431 | Penicillium minioluteum Dierckx BGC | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 432 | Penicillium montanense M. Christensen & Backus | FF-PEN | 2%ME | Air | 26 | M |
| 433 | Penicillium nalgiovense Laxa | FF-FOOD, PEN | 2%ME | Air | 26 | M |
| 434 | Penicillium neoechinulatum (Frisvad et al.) Frisvad & Samson | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 435 | Penicillium oblatum Pitt & Hocking | FF-PEN | 2%ME | Air | 26 | Mm |
| 436 | Penicillium ochrochloron Biourge BGA | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 437 | Penicillium ochrochloron Biourge BGB | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 438 | Penicillium olsonii Bainier & Sartory BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 439 | Penicillium olsonii Bainier & Sartory BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 440 | Penicillium oxalicum Currie & Thom | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 441 | Penicillium paxilli Bainier | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 442 | Penicillium phoeniceum van.Beyma | FF-AIR, PEN | 2%ME | Air | 26 | M |
| 443 | Penicillium piceum Raper & Fennell | FF-AIR, PEN | 2%ME | Air | 26 | M |
| 444 | Penicillium pinophilum Hedge. BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 445 | Penicillium pinophilum Hedge. BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 446 | Penicillium pinophilum Hedge. BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 447 | Penicillium pinophilum Hedge. BGD | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 448 | Penicillium polonicum Zaleski BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 449 | Penicillium polonicum Zaleski BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 450 | Penicillium primulinum Pitt | FF-PEN | 2%ME | Air | 26 | Mm |
| 451 | Penicillium purpurescens (Sopp.) Biourge | FF-PEN | 2%ME | Air | 26 | M |
| 452 | Penicillium purpurogenum Stoll BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 453 | Penicillium purpurogenum Stoll BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 454 | Penicillium purpurogenum Stoll BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 455 | Penicillium purpurogenum var rubisclorotium Thom | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 456 | Penicillium raistrickii G.Smith | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 457 | Penicillium restrictum J.C. Gilman & E.V. Abbott BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 458 | Penicillium restrictum J.C. Gilman & E.V. Abbott BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 459 | Penicillium roqueforti Thom BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 460 | Penicillium roqueforti Thom BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 461 | Penicillium roqueforti Thom BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 462 | Penicillium roqueforti Thom BGD | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 463 | Penicillium roqueforti Thom BGE | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 464 | Penicillium rubrum Stoll BGA | FF-AIR, PEN | 2%ME | Air | 26 | M |
| 465 | Penicillium rubrum Stoll BGB | FF-AIR, PEN | 2%ME | Air | 26 | M |
| 466 | Penicillium rugulosum Thom BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 467 | Penicillium rugulosum Thom BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 468 | Penicillium rugulosum Thom BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 469 | Penicillium sclerotiorum Van.Beyma | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 470 | Penicillium siamense Manoch & Ramirez | FF- PEN | 2%ME | Air | 26 | Mm |
| 471 | Penicillium simplicissimum (Oudem.) Thom BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 472 | Penicillium simplicissimum (Oudem.) Thom BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 473 | Penicillium smithii Quintanilla | FF-PEN | 2%ME | Air | 26 | M |
| 474 | Penicillium solitum Westling BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 475 | Penicillium solitum Westling BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 476 | Penicillium solitum Westling BGC | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 477 | Penicillium steckii Zaleski BGA | FF-AIR, PEN | 2%ME | Air | 26 | m |
| 478 | Penicillium steckii Zaleski BGB | FF-AIR, PEN | 2%ME | Air | 26 | m |
| 479 | Penicillium thomii Maire BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 480 | Penicillium thomii Maire BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 481 | Penicillium tricolor Fresvad et al. | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 482 | Penicillium trubatum Westling BGA (teleo. Eupenicillium baarnese) | FF-PEN | 2%ME | Air | 26 | Mm |
| 483 | Penicillium trubatum Westling BGB (teleo. Eupenicillium baarnese) | FF-PEN | 2%ME | Air | 26 | Mm |
| 484 | Penicillium variabile Sopp BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 485 | Penicillium variabile Sopp BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 486 | Penicillium verrucosum var verrucosum Dierckx BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 487 | Penicillium verrucosum var verrucosum Dierckx BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 488 | Penicillium viridicatum Westling BGA | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 489 | Penicillium viridicatum Westling BGB | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | Mm |
| 490 | Penicillium vulpinum (Cooke & Massee) Seifert & Samson BGA | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 491 | Penicillium vulpinum (Cooke & Massee) Seifert & Samson BGB | FF-AIR, PEN | 2%ME | Air | 26 | Mm |
| 492 | Pestalotiopsis maculans (Corda) Nag Raj | FF-AIR | 2%ME | Air | 26 | |
| 493 | Petromyces alliaceus Malloch & Cain BGA (anam. Aspergillus alliaceous) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |
| 494 | Petromyces alliaceus Malloch & Cain BGB (anam. Aspergillus alliaceous) | FF-AIR, FOOD, ASP | 2%ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 495 | Phanerochaete chrysosporium Burdsall | FF-FOOD | 2%ME | Air | 26 | M |
| 496 | Phialophora fastigiata (Lagerb. & Melin) Conant | FF-AIR, FOOD | 2%ME | Air | 26 | |
| 497 | Phialophora malorum (Kidd & Beaumont) McColloch | FF-AIR | 2%ME | Air | 26 | M |
| 498 | Phialophora richardsiae (Nannf. & Melin) Conant | FF-AIR | 2%ME | Air | 26 | M |
| 499 | Phoma chrysanthemicola Hollos | FF-AIR | 2%ME | Air | 26 | |
| 500 | Phoma exigua var exigua Desm. | FF-AIR | 2%ME | Air | 26 | M |
| 501 | Phoma glomerata (Corda) Wollenw. & Hochapfel BGA | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 502 | Phoma glomerata (Corda) Wollenw. & Hochapfel BGB | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 503 | Phoma herbarum Westend. | FF-AIR | 2%ME | Air | 26 | M |
| 504 | Phoma macrostoma var macrostoma Mont. | FF-AIR | 2%ME | Air | 26 | M |
| 505 | Phoma multirostrata (Mathur et al.) Dorenbosch & Boerema | FF-AIR | 2%ME | Air | 26 | |
| 506 | Phoma septicalis Boerema | FF-AIR | 2%ME | Air | 26 | M |
| 507 | Phoma sorghina (Sacc.) Boerema | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 508 | Pichia anomala (Sydow & Sydow) Kurtzman BGA (anam. <i>Candida pelliculosa</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 509 | Pichia anomala (Sydow & Sydow) Kurtzman BGB (anam. <i>Candida pelliculosa</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 510 | Pichia fermentans Lodder BGA (anam. <i>Candida lambica</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 511 | Pichia fermentans Lodder BGB (anam. <i>Candida lambica</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 512 | Pichia guilliermondii Wickernam (anam. <i>Candida guilliermondii</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 513 | Pichia haplophila Shifrine & Phaff | FF-YST | 2%ME | Air | 26 | |
| 514 | Pichia membranifaciens Hansen (anam. <i>Candida valida</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 515 | Pichia pastoris (Guilliermond) Phaff | FF-YST | 2%ME | Air | 26 | |
| 516 | Pichia petersonii (Wickerham) Kurtzman | FF-YST | 2%ME | Air | 26 | |
| 517 | Pichia spartinae Ahearn et al. | FF-YST | 2%ME | Air | 26 | |
| 518 | Pithomyces chartarum (Berk. & Curtis) M.B. Ellis | FF-AIR | 2%ME | Air | 26 | |
| 519 | Pithomyces sacchari (Speg.) M.B. Ellis | FF-AIR | 2%ME | Air | 26 | M |
| 520 | Rhinochadiella atrovirens Nannf. BGA | FF-AIR | 2%ME | Air | 26 | M |
| 521 | Rhinochadiella atrovirens Nannf. BGB | FF-AIR | 2%ME | Air | 26 | M |
| 522 | Rhizomucor pusillus (Lindt) Schipper | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 523 | Rhizopus microsporus var microsporus v. Tiegh | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 524 | Rhizopus oligosporus Saito | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 525 | Rhizopus oryzae Went & Prins. Geerl. | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 526 | Rhizopus sexualis var sexualis (G. Sm.) Callen | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 527 | Rhizopus stolonifer var stolonifer (EHRENB.:FR.) LINDNER BGA | FF-AIR, FOOD | 2%ME | AIR | 26 | Mm |
| 528 | Rhizopus stolonifer var stolonifer (Ehrenb.:Fr.) Lindner BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 529 | Rhodotorula aurantica (Saito) Lodder | FF-AIR, YST | 2%ME | Air | 26 | Mm |
| 530 | Rhodotorula minuta (Saito) F.C. Harrison | FF-YST | 2%ME | Air | 26 | |
| 531 | Rhodotorula mucilaginosa (AS) (A.Jörg) F.C. Harrison BGA | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 532 | Rhodotorula mucilaginosa (AS) (A.Jörg) F.C. Harrison BGB | FF-AIR, FOOD, YST | 2%ME | Air | 26 | M |
| 533 | Rhodotorula muscorum (di Menna) von Arx & Weijman | FF-YST | 2%ME | Air | 26 | |
| 534 | Saccharomyces bayanus Sacc. | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 535 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGA (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 536 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGB (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 537 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGC (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 538 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGD (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 539 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGE (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 540 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGF (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 541 | Saccharomyces cerevisiae Meyen ex E.C. Hanson BGG (anam. <i>Candida robusta</i>) | FF-FOOD, YST | 2%ME | Air | 26 | M |
| 542 | Saccharomyces exiguus (SS) Reess ex Hanson (anam. <i>Candida holmii</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 543 | Schizosaccharomyces octosporus Beijerinck (<i>Octosporomyces octosporus</i>) | FF-FOOD, YST | 2%ME | Air | 26 | m |
| 544 | Schizosaccharomyces pombe Lindler | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 545 | Scopulariopsis asperula (Sacc.) Hughes | FF-AIR | 2%ME | Air | 26 | M |
| 546 | Scopulariopsis brevicaulis (Sacc.) Bainier BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 547 | Scopulariopsis brevicaulis (Sacc.) Bainier BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 548 | Scopulariopsis brumptii Salvanet-Duval | FF-AIR | 2%ME | Air | 26 | M |
| 549 | Scopulariopsis candida (Guéguen) Vuill. | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 550 | Scopulariopsis chartarum (G. Smith) Morton & G. Smith BGA | FF-AIR | 2%ME | Air | 26 | M |
| 551 | Scopulariopsis chartarum (G. Smith) Morton & G. Smith BGB | FF-AIR | 2%ME | Air | 26 | M |
| 552 | Scopulariopsis fusca Zach BGA | FF-FOOD | 2%ME | Air | 26 | Mm |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|-------------------|--------|-----|------|-------|
| 553 | Scopulariopsis fusca Zach BGB | FF-FOOD | 2%ME | Air | 26 | Mm |
| 554 | Sistotrema brinkmannii (Bresadola) John Eriksson | FF-AIR | 2%ME | Air | 26 | |
| 551 | Sistotrema raduloides (Karsten) Donk | FF-AIR | 2%ME | Air | 26 | |
| 552 | Sporobolomyces roseus Kluyver & Van Niel | FF-AIR, YST | 2%ME | Air | 26 | M |
| 553 | Sporopachydermia ceriana Rodrigues de Miranda | FF-YST | 2%ME | Air | 26 | |
| 554 | Sporothrix schenckii var schenckii Kektoen & Perkins | FF-AIR | 2%ME | Air | 26 | M |
| 555 | Stachybotrys bisbyi (Srinavesan) Barron | FF-AIR | 2%ME | Air | 26 | |
| 556 | Stachybotrys chartarum (Ehrenb.) Hughes | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 557 | Stachybotrys cylindrospora C.N. Jensen | FF-AIR | 2%ME | Air | 26 | Mm |
| 558 | Stachybotrys echinulata (Rivolta) G. Smith | FF-AIR | 2%ME | Air | 26 | M |
| 559 | Sterigmatomyces elviae Sonck & Yarrow | FF-AIR, FOOD, YST | 2%ME | Air | 26 | |
| 560 | Sterigmatomyces halophilus Fell | FF-AIR, YST | 2%ME | Air | 26 | |
| 561 | Syncephalastrum racemosum J.Schrot. | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 562 | Talaromyces bacillisporus (Swift) C.R. Benjamin (anam. <i>Penicillium bacillisporum</i>) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 563 | Talaromyces flavus var flavus (Klocker) Stolk & Samson (anam. <i>Penicillium vermiculatum</i>) | FF-AIR, FOOD, PEN | 2%ME | Air | 26 | M |
| 564 | Talaromyces macrosporus (Stolk & Samson) Frisvad et al. (anam. <i>Penicillium dangesdii</i>) | FF-FOOD, PEN | 2%ME | Air | 26 | Mm |
| 565 | Thamnidium elegans Link | FF-AIR, FOOD | 2%ME | Air | 26 | M |
| 566 | Tilletiopsis albescens Gokhale | FF-AIR | 2%ME | Air | 26 | M |
| 567 | Torula herbarum Link: Fr. | FF-AIR | 2%ME | Air | 26 | M |
| 568 | Torulaspora delbrueckii (Lindner) Lindner (anam. <i>Candida colliculosa</i>) | FF-FOOD, YST | 2%ME | Air | 26 | Mm |
| 569 | Trichoderma asperellum | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | M |
| 570 | Trichoderma atroviride Karsten | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 571 | Trichoderma aureoviride Rifai (teleo. <i>Hypocrea aureoviridis</i>) | FF-TRI | 2%ME | Air | 26 | Mm |
| 572 | Trichoderma citrinoviride Bissett BGA (teleo. <i>Hypocrea schweinitzii</i>) | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 573 | Trichoderma citrinoviride Bissett BGB (teleo. <i>Hypocrea schweinitzii</i>) | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 574 | Trichoderma citrinoviride Bissett BGC (teleo. <i>Hypocrea schweinitzii</i>) | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 575 | Trichoderma crassum Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 576 | Trichoderma fasciculatum Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 577 | Trichoderma fertile Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 578 | Trichoderma ghanense Y. Doi, Y. Abe & J. Sugiyama | FF-TRI | 2%ME | Air | 26 | M |
| 579 | Trichoderma hamatum (Bon.) Bainier | FF-AIR, TRI | 2%ME | Air | 26 | M |
| 580 | Trichoderma harzianum Rifai BGA | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 581 | Trichoderma harzianum Rifai BGB | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 582 | Trichoderma harzianum Rifai BGC | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 583 | Trichoderma koningii Oud. (teleo. <i>Hypocrea koningii</i>) | FF-AIR, TRI | 2%ME | Air | 26 | Mm |
| 584 | Trichoderma longibrachiatum Rifai BGA | FF-AIR, TRI | 2%ME | Air | 26 | Mm |
| 585 | Trichoderma longibrachiatum Rifai BGB | FF-AIR, TRI | 2%ME | Air | 26 | Mm |
| 586 | Trichoderma longibrachiatum Rifai BGC | FF-AIR, TRI | 2%ME | Air | 26 | Mm |
| 587 | Trichoderma longibrachiatum Rifai BGD | FF-AIR, TRI | 2%ME | Air | 26 | Mm |
| 588 | Trichoderma minutisporum Bissett | FF-AIR, TRI | 2%ME | Air | 26 | M |
| 589 | Trichoderma oblongisporum Bissett | FF-TRI | 2%ME | Air | 26 | |
| 590 | Trichoderma piluliferum Webster & Rifai (teleo. <i>Hypocrea pilulifera</i>) | FF-TRI | 2%ME | Air | 26 | M |
| 591 | Trichoderma polysporum (Link:Fr.) Rifai | FF-AIR, TRI | 2%ME | Air | 26 | M |
| 592 | Trichoderma pseudokoningii Rifai (teleo. <i>Hypocrea pseudokoningii</i>) | FF-TRI | 2%ME | Air | 26 | Mm |
| 593 | Trichoderma reesei Simmons BGA (teleo. <i>Hypocrea jecorina</i>) | FF-TRI | 2%ME | Air | 26 | Mm |
| 594 | Trichoderma reesei Simmons BGB (teleo. <i>Hypocrea jecorina</i>) | FF-TRI | 2%ME | Air | 26 | Mm |
| 595 | Trichoderma saturnisporum Hammill | FF-TRI | 2%ME | Air | 26 | Mm |
| 596 | Trichoderma spirale Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 597 | Trichoderma strictipile Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 598 | Trichoderma strigosum Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 599 | Trichoderma tomentosum Bissett | FF-TRI | 2%ME | Air | 26 | M |
| 600 | Trichoderma virens (Miller et al.)v.Ar. BGA | FF-TRI | 2%ME | Air | 26 | Mm |
| 601 | Trichoderma virens (Miller et al.)v.Ar. BGB | FF-TRI | 2%ME | Air | 26 | Mm |
| 602 | Trichoderma viride Pers.:Fr. | FF-AIR, FOOD, TRI | 2%ME | Air | 26 | Mm |
| 603 | Trichothecium roseum (Pers.) Link | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 604 | Trichurus spiralis Hasselbring BGA | FF-AIR | 2%ME | Air | 26 | M |
| 605 | Trichurus spiralis Hasselbring BGB | FF-AIR | 2%ME | Air | 26 | M |
| 606 | Tritirachium oryzae (Vincens) de Hoog | FF-AIR | 2%ME | Air | 26 | M |
| 607 | Ulocladium atrum Preuss | FF-AIR | 2%ME | Air | 26 | M |
| 608 | Ulocladium botrytis Preuss BGA | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 609 | Ulocladium botrytis Preuss BGB | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 610 | Ulocladium chartarum (Preuss) Simmons | FF-AIR, FOOD | 2%ME | Air | 26 | Mm |
| 611 | Ustilago maydis (de Candolle) Corda | FF-AIR, YST | 2%ME | Air | 26 | Mm |
| 612 | Verticillium lecanii (Zimmermann) Viegas | FF-AIR | 2%ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 4: Database Species Lists and Their Characteristics

| | Species Name | Type | Medium | Atm | Temp | Photo |
|-----|---|--------------|--------|-----|------|-------|
| 613 | Wallemia sebi (Fr.) v. Arx | FF-AIR, FOOD | 2% ME | Air | 26 | M |
| 614 | Williopsis saturnus var saturnus (Klocker) Zender | FF-YST | 2% ME | Air | 26 | |
| 615 | Yarrowia lipolytica (SS) (Wickerham et al.) v.d. Walt & v. Arx (anam. <i>Candida lipolytica</i>) | FF-FOOD, YST | 2% ME | Air | 26 | Mm |
| 616 | Zygorhynchus moelleri Vuillemin | FF-AIR | 2% ME | Air | 26 | |
| 617 | Zygosaccharomyces bailii (Lindner) Guillerm | FF-FOOD, YST | 2% ME | Air | 26 | Mm |
| 618 | Zygosaccharomyces rouxii (Boutroux) Yarrow | FF-FOOD, YST | 2% ME | Air | 26 | Mm |
| 619 | Zygosporium mycophilum (Vuillemin) Saccardo | FF-AIR | 2% ME | Air | 26 | M |

M Indicates that a Macroscopic photo is present in Biolog's photo library.

m Indicates that a Microscopic photo is present in Biolog's photo library.

Appendix 5: Program Printouts

Manual printout (sample)

| | | | | | | | | | | | | | | | |
|-----|---|---|-----------------|-----|-----|-----|-----|-----|-----|-----|-----|------|------|-------|--------|
| 1→ | Program | : Biolog MicroLog3 4.20 | | | | | | | | | | | | | |
| 2→ | Read From File | : C:\Biolog420\DataFiles\GNDemo2a.D4C | | | | | | | | | | | | | |
| 3→ | Edit Status | : OK | | | | | | | | | | | | | |
| 4→ | Read Time | : Jul 03 2001 12:07 | Creator : PG | | | | | | | | | | | | |
| 5→ | Modified Date | : Jul 09 2001 07:44 | Modifier: lisa | | | | | | | | | | | | |
| 6→ | Unrestricted Access? | : No | | | | | | | | | | | | | |
| 7→ | Plate Number | : 5 | | | | | | | | | | | | | |
| 8→ | Incubation Time | : 16-24 | | | | | | | | | | | | | |
| | Sample Number | : 5 | Plate Type: GN2 | | | | | | | | | | | | |
| | Strain Type | : GN-ENT | | | | | | | | | | | | | |
| | Strain Number | : 555 | | | | | | | | | | | | | |
| | Strain Name | : Cedecia | | | | | | | | | | | | | |
| | Other | : | | | | | | | | | | | | | |
| 9→ | Data Input Mode | : File | | | | | | | | | | | | | |
| | Parent File | : C:\Biolog420\DataFiles\GNDemo2.D4C | | | | | | | | | | | | | |
| 10→ | Number +/- Reactions | : 43/ 0/ 53 | | | | | | | | | | | | | |
| 11→ | Database To Search | : MicroLog | | | | | | | | | | | | | |
| | Data Base(s) Searched | : C:\BIOLOG420\Databases\GN601.KID | | | | | | | | | | | | | |
| | Key | : <X>: positive; <X-: mismatched positive; X: negative; X+: mismatched negative {X}: borderline; -X: less than A1 well | | | | | | | | | | | | | |
| | Color | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | | |
| | A | - | - | <+> | - | - | - | <+> | <+> | - | - | <+> | <+> | | |
| | B | - | <+> | - | <+> | <+> | <+> | <+> | - | - | - | <+> | <+> | | |
| | C | - | <+> | <+> | - | - | <+> | <+> | <+> | - | <+> | <+> | <+> | | |
| 12→ | D | - | - | <+> | <+> | <+> | - | <+> | - | - | - | - | - | | |
| | E | - | - | - | - | - | <+> | - | - | - | - | - | <+> | | |
| | F | - | <+> | - | - | - | <+> | <+> | <+> | <+> | <+> | <+> | - | | |
| | G | - | <+> | - | - | - | - | - | <+> | <+> | - | - | - | | |
| | H | - | <+> | <+> | <+> | - | - | - | - | <+> | <+> | <+> | <+> | | |
| 13→ | => Species ID: CEDECEA NETERI <= | | | | | | | | | | | | | | |
| | Species | | | | | | | | | | | PROB | SIM | DIST | TYPE |
| | =>1) CEDECEA NETERI | | | | | | | | | | | 100 | 0.61 | 5.96 | GN-ENT |
| | 2) CEDECEA DAVISAE | | | | | | | | | | | 0 | 0.00 | 10.71 | GN-ENT |
| | 3) KLEBSIELLA PLANTICOLA/ORNITHINOLYTICA | | | | | | | | | | | 0 | 0.00 | 10.97 | GN-ENT |
| | 4) KLEBSIELLA PLANTICOLA | | | | | | | | | | | 0 | 0.00 | 12.94 | GN-ENT |
| 14→ | 5) ENTEROBACTER NIMIPRESSURALIS | | | | | | | | | | | 0 | 0.00 | 13.35 | GN-ENT |
| | 6) ENTEROBACTER ASBURIAE | | | | | | | | | | | 0 | 0.00 | 13.43 | GN-ENT |
| | 7) ENTEROBACTER AEROGENES (KLB. MOBILIS) | | | | | | | | | | | 0 | 0.00 | 13.70 | GN-ENT |
| | 8) ENTEROBACTER CLOACAE | | | | | | | | | | | 0 | 0.00 | 14.04 | GN-ENT |
| | 9) CITROBACTER KOSERI | | | | | | | | | | | 0 | 0.00 | 14.38 | GN-ENT |
| | 10) EWINGELLA AMERICANA | | | | | | | | | | | 0 | 0.00 | 14.44 | GN-ENT |
| | Other) | | | | | | | | | | | | | | |
| | Print Time = | Jul 09 2001 07:50 | | | | | | | | | | | | | |

Manual Printout Key

1. Software version
2. Name of file
3. Data file record Edit status: OK, Out, Hold, Atypical
4. Time MicroPlate was read
5. Time file was modified
6. Does user have Unrestricted Access to software?
7. Number of MicroPlate
8. Sample identifiers
9. Data generation mode (file)
10. Number of positive, borderline, and negative reactions
11. Name of database searched
12. Positive (<+>), borderline ({/}), and negative reactions (-). A negative (-) well with a plus (+) on the right side is a mismatch. The well read as negative, but at least 80% of the strains of that organism tested in our database were positive for that well. A plus (+) with a < on the left and a minus (–) on the right is a mismatch, but the well read as positive, while the strains in our database were mostly negative for that well.
13. Species identification/ #1 identification
14. #2-9 identifications

Reader printout (sample)

| | | | |
|-----|--|---|------------------------|
| 1→ | Program | : Biolog MicroLog3 4.20 | |
| 2→ | Read From File | : C:\Biolog420\DataFiles\GNDemo1a.D4C | |
| | Edit Status | : OK | |
| 3→ | Threshold Mode | : Automatic: Color: 35/86 | |
| | 590/750 Filters Used | : 6 / 5 | |
| 4→ | Read Time | : Feb 12 1998 08:17 | Creator : AFB |
| 5→ | Modified Date | : Jul 03 2001 14:51 | Modifier: lisa |
| 6→ | Unrestricted Access? | : No | |
| 7→ | Plate Number | : 1 | |
| | Incubation Time | : 16-24 | |
| | Sample Number | : 1 | Plate Type: GN2 |
| | Strain Type | : GN-NENT OXI+ | |
| 8→ | Strain Number | : 3901 | |
| | Strain Name | : AER.SAL SS SAL | |
| | Other | : | |
| 9→ | Data Input Mode | : File | |
| | Parent File | : C:\Biolog420\DataFiles\GNDemo1.D4C | |
| 10→ | Number +/- Reactions | : 25 / 3 / 68 | |
| 11→ | Database To Search | : MicroLog | |
| | Data Base(s) Searched | : C:\BIOLOG420\Databases\GN601.KID | |
| | Key | : <X>: positive; <X-: mismatched positive; X: negative; X+: mismatched negative {X}: borderline; -X: less than A1 well | |
| | Color | 1 2 3 4 5 6 7 8 9 10 11 12 | |
| | A | 0 -9 < 286> < 284> < 374> < 209> -14 < 313> -5 -6 -3 10 | |
| | B | -11 < 347> 1 < 165> 7 < 353> -2 -9 -3 < 325> < 94> 5 | |
| 12→ | C | 11 -1 { 60} 3 4 1 -3 5 20 -7 < 228> 3 | |
| | D | 0 0 4 -9 -1 3 1 -2 -1 -4 -2 -4 | |
| | E | -5 2 0 19 -2 0 1 -2 2 -2 -6 < 278> | |
| | F | 17 < 108> 0 { 73} 4 < 151> < 202> < 300> < 228> < 253> 21 < 101> | |
| | G | 6 1 -3 11 -5 -2 10 < 261> < 242> -3 -5 2 | |
| | H | { 54} < 399> < 329> < 214> -3 2 9 9 < 306> 11 1 2 | |
| 13→ | => Species ID: AEROMONAS SALMONICIDA SS SALMONICIDA <= | | |
| | Species | PROB | SIM DIST TYPE |
| | =>1) AEROMONAS SALMONICIDA SS SALMONICIDA | 100 | 0.99 0.09 GN-NENT OXI+ |
| | 2) AEROMONAS ENCHELEIA | 0 | 0.00 6.69 GN-NENT OXI+ |
| | 3) AEROMONAS HYDROPHILA-LIKE DNA GROUP 2 | 0 | 0.00 8.26 GN-NENT OXI+ |
| | 4) AEROMONAS EUCRENOPHILA DNA GROUP 6 | 0 | 0.00 8.49 GN-NENT OXI+ |
| | 5) VIBRIO FURNISSII | 0 | 0.00 8.75 GN-NENT OXI+ |
| 14→ | 6) AEROMONAS ALLOSACCHAROPHILA | 0 | 0.00 9.07 GN-NENT OXI+ |
| | 7) VIBRIO HARVEYI | 0 | 0.00 9.08 GN-NENT OXI+ |
| | 8) VIBRIO FLUVIALIS | 0 | 0.00 9.80 GN-NENT OXI+ |
| | 9) AEROMONAS VERONII/SOBRIA DNA GROUP 8 | 0 | 0.00 9.83 GN-NENT OXI- |
| | 10) VIBRIO ALGINOLYTICUS | 0 | 0.00 9.89 GN-NENT OXI- |
| | Other) | | |

Reader and Saved File Modes Printout Key

1. Software version
2. If the printout is from a saved file, this entry lists file name, position of the entry in the file, and edit status.
3. Threshold mode. This listing may say “Manual” or “Automatic.” In Automatic mode, the software automatically calculates a threshold cutoff for negative and positive reactions, indicated by numerical OD values. Any values between these two numbers are considered borderline. The threshold can change from MicroPlate to MicroPlate, depending on the intensity of the positive reactions and the clarity of the negative reactions. In Manual Threshold mode, the user has changed the threshold cutoffs to match the visual interpretation of the plate.
4. Time MicroPlate was read
5. Time file was modified
6. Does user have Unrestricted Access to software?
7. Number of MicroPlate
8. Sample identifiers
9. Operating mode (file)
10. Number of positive, borderline, and negative reactions
11. Name of database searched
12. Numbers seen in each well is the percent change data of the OD compared to well A1 (A1 is always 0). $DWD (dual\ wavelength\ data) = (590nm-750nm)x - (590nm-750nm)A1 \times 1000$, where x = any well. Positive (<>), borderline ({ }), and negative reactions (no brackets). Numbers with no brackets and a plus (+) on the right side are mismatches. The well read as negative, but at least 80% of the strains of that organism tested in our database were positive for that well. Numbers with a < on the left and a minus (–) on the right are also mismatches, but the well read as positive, while the strains in our database were mostly negative for that well. Numbers with a minus (-) to the left of the well value mean that the well’s OD was less than the OD of the control well.
13. Species Identification/#1 identification
14. #2-9 identifications

Fungal printout (sample)

| | | |
|-----|---|---|
| 1→ | Program | : Biolog MicroLog3 4.20 |
| 2→ | Read From File | : C:\Biolog420\DataFiles\FID2.D4C |
| | Save To File | : |
| 3→ | Threshold Mode | : Automatic: Color: 47/145; Turbidity: 53/147 |
| 4→ | 490/750 Filters Used | : 3 / 5 |
| 5→ | Read Time | : Jun 01 3898 13:03 Creator : AFB |
| 6→ | Modified Date | : --- Modifier: |
| 7→ | Unrestricted Access? | : No |
| 8→ | Plate Number | : 1 |
| | Incubation Time | : 72 HR |
| | Sample Number | : 1 Plate Type: FF |
| | Strain Type | : Aspergillus |
| | Strain Number | : FBF 137 |
| 9→ | Strain Name | : ASP.WEN |
| | Other | : |
| 10→ | Data Input Mode | : File |
| | Parent File | : C:\Biolog420\DataFiles\Ffid.D4C |
| 11→ | Number +/- Reactions | : 93 / 26 / 73; Color: 47 / 10 / 39; Turbidity: 46 / 16 / 34 |
| | Database To Search | : MicroLog |
| 12→ | Data Base(s) Searched | : deprecated database format |
| | Key | : <X>: positive; <X-: mismatched positive; X: negative; X+: mismatched negative {X}: borderline; -X: less than A1 well |
| | Color | 1..... 2..... 3..... 4..... 5..... 6..... 7..... 8..... 9..... 10..... 11..... 12..... |
| | A | 0 < 270> -13 <1292> -7 < 187> < 169> -6 < 460> < 164> { 143} < 310> |
| | B | -15 -25 < 353> < 372> < 870> -38 33 <1333> < 465> <1119> -12 < 749> |
| | C | { 64} -40 2 < 991> { 59} < 487> <1125> -7 -13 7 < 543> 17 |
| | D | < 684> < 796> < 299> 46 3 8 0 < 764> < 209> { 57} < 232> < 253> |
| | E | 32 < 199> -8 < 753> < 181> { 107} < 625> -16 < 179> < 315> 42 < 557> |
| | F | < 705> -25 < 849> -37 -7 < 733> { 52} -52 -12 21 <1487> <1232> |
| | G | <1571> { 68} < 277> < 166> { 112} -4 9 <1610> < 573> <1753> < 267> <1090> |
| | H | { 83} 19 31 <1100> <1501> { 77} 40 2 32 -13 -2 -14 |
| 13→ | Turbidity | |
| | A | 0 < 421> 0 < 322> 6 < 259> { 136} 34 < 478> < 192> < 254> < 413> |
| | B | < 189> < 197> < 458> < 210> <1040> -11 { 107} < 282> < 506> < 674> 19 < 997> |
| | C | { 61} 26 24 < 307> { 75} < 195> < 218> 34 29 38 < 550> 32 |
| | D | < 346> < 669> < 981> { 78} 29 2 48 < 611> < 230> { 118} < 429> < 325> |
| | E | 27 { 82} 11 <1102> < 160> < 162> <2475> 13 < 192> < 218> { 75} < 653> |
| | F | < 278> { 62} < 174> 13 34 < 163> 35 14 32 34 < 189> < 725> |
| | G | < 229> < 224> { 134} { 86} < 370> 19 3 < 206> { 118} < 240> { 75} { 123} |
| | H | 42 -2 { 94} < 989> < 181> 8 13 5 29 { 59} -10 8 |
| 14→ | => Species ID: Aspergillus wentii Wehmer <= | |
| | Species | PROB SIM DIST TYPE |
| | => 1) Aspergillus wentii Wehmer | 100 0.70 4.50 ASP |
| | 2) Aspergillus sepultus Tuthill & Christensen | 0 0.00 8.69 ASP |
| | 3) Emericella quadrilineata (Thom & Raper) C.R. Benjamin | 0 0.00 9.60 ASP |
| | 4) Aspergillus carbonarius (Bainier) Thom | 0 0.00 10.29 ASP |
| 15→ | 5) Aspergillus carneus (v. Tiegham) Blockwitz | 0 0.00 11.31 ASP |
| | 6) Aspergillus niger v. Tiegham BGB | 0 0.00 11.88 ASP |
| | 7) Emericella violacea (Fennell & Raper) Malloch & Cain | 0 0.00 11.96 ASP |
| | 8) Aspergillus phoenicis (Corda) Thom | 0 0.00 12.02 ASP |
| | 9) Aspergillus ostianus Wehmer | 0 0.00 13.19 ASP |
| | 10) Aspergillus sparsus Raper & Thom | 0 0.00 13.23 ASP |
| | Other) | |
| | Print Time = Jul 09 2001 08:00 | |

FF Printout Key

1. Software version
2. If the printout is from a saved file, this entry lists file name, position of the entry in the file, and edit status.
3. Threshold mode. This listing may say “Manual” or “Automatic.” In Automatic mode, the software automatically calculates a threshold cutoff for negative and positive reactions, indicated by numerical OD values. Any values between these two numbers are considered borderline. The threshold can change from MicroPlate to MicroPlate, depending on the intensity of the positive reactions and the clarity of the negative reactions. Two thresholds are given for the FF MicroPlate: color (for 490 nm readings) and turbidity (for 750 nm readings).
4. Filters and their positions
5. Time MicroPlate was read
6. Time file was modified
7. Does user have Unrestricted Access to software?
8. Number of MicroPlate
9. Sample identifiers
10. Operating mode (file)
11. Number of positive, borderline, and negative reactions and listed for color and turbidity
12. Name of database searched
13. OD data for each well (in 8x12 format)
 - For color data: The numbers seen in each well is the percent change data of the optical density of the difference between 490nm and 750 nm compared to the A1 control well (A1 will always be 0). The Dual Wavelength Data formula = $DWD = (490nm - 750nm)_x - (490nm - 750nm)_{A1} \times 1000$, where x = any well. The turbidity is subtracted out to obtain the OD contributed by the color
 - For the Turbidity Data: The numbers seen in each well is the percent change data of the optical density of the 750 nm readings compared to the A1 control well (A1 will always be 0). The % Change formula = $\% \text{ Change } 750nm = (750nm_x - 750nm_{A1}) \times 1600$, where x = any well*
 - Positive (<>), borderline ({}), and negative reactions (no brackets). Numbers with no brackets and a plus (+) on the right side are mismatches. The well read as negative, but at least 80% of the strains of that organism tested in our database were positive for that well. Numbers with a < on the left and a minus (–) on the right are also mismatches, but the well read as positive, while the strains in our database were mostly negative for that well. Numbers with a minus (-) to the left of the well value mean that the well’s OD was less than the OD of the control well.
14. Species Identification/#1 identification
15. #2-9 identifications

*We use a larger multiplication factor for the turbidity to have the OD readings on the same scale as 490 readings.

Print ranked data headers printout (sample)

Print Ranked Data Headers

DATA HEADERS

Current Date = 7/9/01 10:13:27 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNPden.D4C
 Index File = C:\Biolog420\DataFiles\GNPden.ISC
 Sort File = C:\Biolog420\DataFiles\GNPden.SSC
 Organism Type = 5 : GP-COCCUS
 Day Index = 0

2/ 2) Enterococcus

| (F #) | OK MK Created | Plate | Sample | Hr | Plate | Strain Type | Strain Name | Strain # | Other | ID Result |
|-------|---------------|-------------|--------|----|-----------|-------------|--------------|----------|-------|------------|
| (4) | OK | Jan 19 2001 | ? | 7 | 24 HR GP2 | GP-COCCUS | Enterococcus | 777 | | Species ID |
| (5) | OK | Jan 19 2001 | ? | 8 | 24 HR GP2 | GP-COCCUS | Enterococcus | 888 | | Species ID |
| (6) | OK | Jan 19 2001 | ? | 9 | 24 HR GP2 | GP-COCCUS | Enterococcus | 999 | | Species ID |

DATA HEADERS

Current Date = 7/9/01 10:13:27 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNPden.D4C
 Index File = C:\Biolog420\DataFiles\GNPden.ISC
 Sort File = C:\Biolog420\DataFiles\GNPden.SSC
 Organism Type = 1 : GN-ENT
 Day Index = 0

1/ 2) Cedecia

| (F #) | OK MK Created | Plate | Sample | Hr | Plate | Strain Type | Strain Name | Strain # | Other | ID Result |
|-------|---------------|-------------|--------|----|-----------|-------------|-------------|----------|-------|------------|
| (1) | OK | Jan 19 2001 | ? | 4 | 24 HR GN2 | GN-ENT | Cedecia | 444 | | Species ID |
| (2) | OK | Jan 19 2001 | ? | 5 | 24 HR GN2 | GN-ENT | Cedecia | 555 | | Species ID |
| (3) | OK | Jan 19 2001 | ? | 6 | 24 HR GN2 | GN-ENT | Cedecia | 666 | | Species ID |

Print species counts (sample)

Print Species Counts

SPECIES COUNT

Current Date = 7/9/01 10:13:27 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNGPden.D4C
 Index File = C:\Biolog420\DataFiles\GNGPden.ISC
 Sort File = C:\Biolog420\DataFiles\GNGPden.SSC
 Organism Type = 1 : GN-ENT

1/ 2) Cedecia

#Plates = 3

#STR = 3

total plates = 3
 total strains = 3
 total species = 1

SPECIES COUNT

Current Date = 7/9/01 10:13:27 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNGPden.D4C
 Index File = C:\Biolog420\DataFiles\GNGPden.ISC
 Sort File = C:\Biolog420\DataFiles\GNGPden.SSC
 Organism Type = 5 : GP-COCCUS

2/ 2) Enterococcus

#Plates = 3

#STR = 3

total plates = 3
 total strains = 3
 total species = 1

Print species/plate dendrograms, page 1 (sample)

Print Species / Plate Dendrograms

SPECIES / PLATE DENDROGRAM

Current Date = 7/9/01 10:09:49 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\BIOLOG420\DataFiles\GNPden.D4C
 Out File = C:\BIOLOG420\tree2.drm
 Ref Dist = 10
 Comment = #5 = GP-COCCUS: Species / Plate = Enterococcus

```

.....|.....|.....|
/---( 0) A 1 ( 3) Enterococcus SA:7 ST:777
/-----+---( 1) A 2 ( 4) Enterococcus SA:8 ST:888
|
|-----+---( 2) B 3 ( 5) Enterococcus SA:9 ST:999
|.....|.....|.....|

```

TOTAL DATA = 1 # Entries = 13

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---------------------------|-----|-----|-----|-----|-----|-----|-----|-----|----|-----|-----|-----|
| A | 0 | 0 | 17 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 67 | 0 | 67 | 0 | 67 | 0 | 100 | 0 |
| C | 67 | 0 | 67 | 67 | 67 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 67 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 100 | 0 | 0 | 0 | 0 | 67 | 67 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 33 | 0 | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 33 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 67 |
| H | 100 | 67 | 33 | 67 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| avg % = 93.575 # pos = 13 | | | | | | | | | | | | |

TREE GROUP A: # Entries = 2

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 0 | 25 | 100 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 50 | 0 | 100 | 0 | 100 | 0 | 100 | 0 |
| C | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 100 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 50 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 50 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 100 |
| H | 100 | 100 | 50 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| avg % = 97.656 # pos = 27 | | | | | | | | | | | | |

TREE GROUP B: # Entries = 1

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|----------------------------|-----|-----|-----|-----|-----|-----|-----|-----|---|-----|-----|-----|
| A | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 100 | 0 |
| C | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 0 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| H | 100 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| avg % = 100.000 # pos = 14 | | | | | | | | | | | | |

Print species/plate dendrograms, page 2 (sample)

Print Species / Plate Dendrograms

SPECIES / PLATE DENDROGRAM:

Current Date = 7/10/01 9:41:08 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNGPden.D4C
 Out File = C:\BIOLOG420\tree2.prn
 Ref Dist = 10
 Comment = #1 = GN-ENT: Species / Plate = Cedecia

```
|.....|.....|.....|.....|
                        /----- ( 0)  A 1 ( 0) Cedecia SA:4 ST:444
-----+-----+---+ ( 1)  A 2 ( 1) Cedecia SA:5 ST:555
                        \-- ( 2)  A 3 ( 2) Cedecia SA:6 ST:666
|.....|.....|.....|.....|
```

TOTAL DATA = : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 33 | 100 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 100 | 100 |
| B | 0 | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 100 | 100 |
| C | 0 | 100 | 100 | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 100 | 100 |
| D | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 33 | 0 | 0 | 0 | 0 |
| E | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 100 |
| F | 100 | 33 | 33 | 0 | 100 | 100 | 100 | 100 | 100 | 100 | 0 | 0 |
| G | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 67 | 0 | 0 |
| H | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 100 | 67 |

avg % = 97.569 # pos = 41

TREE GROUP A : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 33 | 100 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 100 | 100 |
| B | 0 | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 100 | 100 |
| C | 0 | 100 | 100 | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 100 | 100 |
| D | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 33 | 0 | 0 | 0 | 0 |
| E | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 100 |
| F | 100 | 33 | 33 | 0 | 100 | 100 | 100 | 100 | 100 | 100 | 0 | 0 |
| G | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 67 | 0 | 0 |
| H | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 100 | 67 |

avg % = 97.569 # pos = 41

Print species/plate dendrograms, page 3 (sample)

SPECIES / PLATE DENDROGRAM

Current Date = 7/10/01 9:49:44 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNPden.D4C
 Out File = C:\BIOLOG420\tree2.prn
 Ref Dist = 10
 Comment = #1 = GN-ENT: Species / Plate = Cedecia

```

|.....|.....|.....|.....|
                        /-----( 0)  A 1 ( 0) Cedecia SA:4 ST:444
-----+-----+-----( 1)  A 2 ( 1) Cedecia SA:5 ST:555
                        \-( 2)  A 3 ( 2) Cedecia SA:6 ST:666
|.....|.....|.....|.....|

```

TOTAL DATA = : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 33 | 100 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 100 | 100 |
| B | 0 | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 100 | 100 |
| C | 0 | 100 | 100 | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 100 | 100 |
| D | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 33 | 0 | 0 | 0 | 0 |
| E | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 100 |
| F | 100 | 33 | 33 | 0 | 100 | 100 | 100 | 100 | 100 | 100 | 0 | 0 |
| G | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 67 | 0 | 0 |
| H | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 100 | 67 |

avg % = 97.569 # pos = 41

TREE GROUP A : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 33 | 100 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 100 | 100 |
| B | 0 | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 100 | 100 |
| C | 0 | 100 | 100 | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 100 | 100 |
| D | 0 | 0 | 100 | 100 | 100 | 0 | 100 | 33 | 0 | 0 | 0 | 0 |
| E | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 100 |
| F | 100 | 33 | 33 | 0 | 100 | 100 | 100 | 100 | 100 | 100 | 0 | 0 |
| G | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 67 | 0 | 0 |
| H | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 100 | 67 |

avg % = 97.569 # pos = 41

Print species/plate dendrograms, page 4 (sample)

SPECIES / PLATE DENDROGRAM

Current Date = 7/10/01 9:41:09 AM
 Program = C:\BIOLOG420\ML3_420.EXE
 Data File = C:\Biolog420\DataFiles\GNPden.D4C
 Out File = C:\BIOLOG420\tree2.prn
 Ref Dist = 10
 Comment = #5 = GP-COCCUS: Species / Plate = Enterococcus

```

|.....|.....|.....|.....|
|------( 0)  A 1 ( 3) Enterococcus SA:7 ST:777
|-----+---( 1)  A 2 ( 4) Enterococcus SA:8 ST:888
|-----+
|-----+------( 2)  B 3 ( 5) Enterococcus SA:9 ST:999
|.....|.....|.....|.....|

```

TOTAL DATA = : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|----|-----|-----|-----|
| A | 0 | 0 | 17 | 67 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 67 | 0 | 67 | 0 | 67 | 0 | 100 | 0 |
| C | 67 | 0 | 67 | 67 | 67 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 67 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 100 | 0 | 0 | 0 | 0 | 67 | 67 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 33 | 0 | 0 | 67 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 33 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 67 |
| H | 100 | 67 | 33 | 67 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

avg % = 93.576 # pos = 13

TREE GROUP A : # Entries = 2

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| A | 0 | 0 | 25 | 100 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 50 | 0 | 100 | 0 | 100 | 0 | 100 | 0 |
| C | 100 | 0 | 100 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 100 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 100 | 0 | 0 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 50 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 50 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 100 |
| H | 100 | 100 | 50 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

avg % = 97.656 # pos = 27

TREE GROUP B : # Entries = 3

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|---|-----|-----|-----|-----|-----|-----|---|-----|-----|-----|
| A | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 100 | 100 | 0 |
| B | 0 | 0 | 100 | 100 | 100 | 0 | 0 | 0 | 0 | 0 | 100 | 0 |
| C | 0 | 0 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 |
| D | 0 | 0 | 0 | 100 | 0 | 0 | 100 | 100 | 0 | 0 | 0 | 100 |
| E | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| F | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| G | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |
| H | 100 | 0 | 0 | 0 | 100 | 0 | 0 | 0 | 0 | 0 | 0 | 0 |

avg % = 66.667 # pos = 14

Appendix 6: Dangerous Pathogen (DP) Database

Biolog makes a supplemental database for species falling into a category called Dangerous Pathogens (e.g., the organisms that generally require BL2-BL3 handling protocols). If your laboratory is testing for these pathogens, you will know in advance if any samples are suspected of containing these organisms. Follow the instructions in the body of this user guide. As appropriate, also follow the special procedures listed here.

Safety Considerations

To comply with government safety regulations¹ use Biosafety Level 2 (BL2) methods, containment equipment, and facilities. Immunization of laboratory personnel may be required. It is recommended that Biosafety Level 3 (BL3) practices, containment equipment, and facilities are used for work involving production volumes or concentration of cultures, and for activities that have a high potential for aerosol production. In these facilities, immunization is recommended for all persons working with the agent, all persons working in the same laboratory room where the cultures are handled, and all persons working with infected animals.

Special Procedures for Dangerous Pathogens

For the most part, you will prepare, read, and interpret these samples exactly as described throughout this manual. There are, however, a few exceptions. Table 13-1 lists the special procedures required for dangerous pathogens.

¹ Biosafety in Microbiological and Biomedical Laboratories, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, and National Institutes of Health, Fourth Edition, May 1999.

| Organism | Special Procedures |
|---|--|
| <i>Francisella tularensis</i> | Set up as a Gram-negative Fastidious organism Do not add thioglycolate to inoculating fluid |
| <i>Yersinia pestis</i> | Set up as an Gram-negative Enteric organism Do not add thioglycolate to inoculating fluid Grow only on BUG + Blood (not TSA + Blood) |
| <i>Burkholderia mallei</i> and <i>Burkholderia pseudomallei</i> | Set up as a Gram-negative Non-Enteric organism Grow only on BUG + Blood (not TSA + Blood) Incubate at 30° C |
| <i>Brucella mellitensis</i> (Use this procedure for <i>B. abortus</i> , <i>B. canis</i> , <i>B. suis</i> and other species that have been reclassified as <i>B. mellitensis</i> .) | Set up as a Gram-negative Fastidious organism Do not add thioglycolate to inoculating fluid Multiple chocolate plates are needed to obtain proper inoculum |

TABLE 12-1: PROCEDURE EXCEPTIONS FOR DANGEROUS PATHOGENS