

P.aeruginosa biofilm media

	Final concentration	10L	8L	6L	5L	4L	3L	2L	1L	0.5L
FB										
Milli-Q H2O		9 L	7.2 L	5.4 L	4.5 L	3.6 L	2.7 L	1.8 L	0.9 L	0.45 L
1 M MgCl2	1mM	10 mL	8 mL	6 mL	5 mL	4 mL	3 mL	2 mL	1 mL	0.5 mL
1 M CaCl2	0.1mM	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
Trace metals	-	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
autoclave										
A10		1 L	0.8 L	0.6 L	0.5 L	0.4 L	0.3 L	0.2 L	0.1 L	0.05 L
Carbon source										
20% glucose	0.3mM	2.7ml	2.16 mL	1.62 mL	1.35 mL	1.08 mL	810 uL	540 uL	270 uL	135 uL

S.epi biofilm media

		10L	8L	6L	5L	4L	3L	2L	1L	0.5L
<i>P.aeruginosa</i> biofilm media										
TSB	3%(V/V)	300 mL	240 mL	180 mL	150 mL	120 mL	90 mL	60 mL	30 mL	15 mL

E.coli biofilm media

		10L	8L	6L	5L	4L	3L	2L	1L	0.5L
FB										
Milli-Q H2O		9 L	7.2 L	5.4 L	4.5 L	3.6 L	2.7 L	1.8 L	0.9 L	0.45 L
1 M MgCl2	1mM	10 mL	8 mL	6 mL	5 mL	4 mL	3 mL	2 mL	1 mL	0.5 mL
1 M CaCl2	0.1mM	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
10.000x Fe-EDTA	3.7mg/L	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
autoclave										
A10		1 L	0.8 L	0.6 L	0.5 L	0.4 L	0.3 L	0.2 L	0.1 L	0.05 L
Carbon source										
glucose	1mM	9mL	8.1mL	5.4mL	4.5mL	3.6mL	2.7mL	1.8mL	0.9mL	0.45mL
1.0 mg/ml thiamine	1µg/ml	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
10 mg/ml proline	10µg/ml	1 mL	0.8 mL	0.6 mL	0.5 mL	0.4 mL	0.3 mL	0.2 mL	0.1 mL	0.05 mL
10 mg/ml uridine*	25µg/ml	2.5mL	2mL	1.5mL	1.25mL	1mL	0.75mL	0.5mL	0.25mL	0.125mL

(* , uridine makes E. coli K12 wildtypes as MG1655 or W3110 grow better because they are starving from pyrimidine)