


Confocal Microscopy and Atomic Force Microscopy (AFM) of biofilms


A very brief primer...




Fundamentals of Confocal Microscopy

Based on a conventional fluorescence microscope


Fluorescent Microscope




Confocal Microscope



3D reconstruction



Pseudomonas putida
cells mixed with
Acinetobacter cells in a
microbial biofilm





Benefits of Confocal Microscopy


- Reduced blurring of the image from light scattering
- Increased effective resolution
- Improved signal to noise ratio
- Clear examination of thick specimens
- Z-axis scanning (3D-reconstruction possible)
- Magnification can be adjusted electronically
- X-Y resolution: $\sim 200\text{-}250\text{ nm}$ (Ernst Abbe)

Disadvantages of Confocal Microscopy

- Requires fluorescent samples
- Uses laser illumination (expensive, few wavelengths)
- Instrument expensive to acquire and run
- Z-resolution typically $> 500\text{ nm}$








New developments in confocal microscopy


- **MP (multi-photon) or two-photon confocal microscopy:** Two or more photons from a long wavelength illumination at a time excites a fluorophore. High resolution in Z-axis (practically down to ~ 200 nm)
- **White Laser (Leica):** A continuous wave "white" laser (tunable from 470-670 nm, with up to 8 lines simultaneously)
- **STED (Stimulated Emission Depletion)(Leica):** A new method where fluorescence is depleted around the area of interest (see next slide)
- **Superresolution structured microscopy (Zeiss):** A method where images are rotated and combined to create a moiré pattern which is deconvolved to create a high res image.
- **Photo-Activated localization microscopy (Zeiss):** Sequential illumination and localization of fluorophores combined with computational reconstruction of high res images.
- **Raman-confocal microscopy (Leica, under development):** Confocal microscope combined with raman spectroscopy. Enables localized determination of [changes in] concentrations of metabolites etc.

STED (Stimulated Emission Depletion)


- In a Leica TCS STED microscope the sample is illuminated by two pulsed laser beams, tightly synchronized.
- The 635 nm wavelength excites the fluorophores of the sample the same way a conventional confocal system does. The excitation laser pulses are directly followed by a ring shaped illumination of a Ti:Sapphire Infrared laser (730-780 nm).
- This pulse inhibits/depletes the fluorescence in the outer regions of the illuminated spot.
- The result: A smaller fluorescence spot that allows much more accurate scanning than with other methods using focused light. X-Y res: <90 nm



Sample STED image




Confocal image




STED image




Atomic Force Microscopy (AFM)



- A method to record the topographic property of a surface
- Physical interaction with the surface is necessary



The principle







Benefits of AFM

- Extremely high resolution (down to atomic level few Å, but typically 10-20 nm)
- No staining required
- Possible to measure attractive/repulsive forces

Disadvantages of AFM

- Extremely high resolution (very low "focal depth")
- Small image area
- Sample must be "flat" and tightly fixed
- Very difficult to work in humid or wet environments
- Need direct contact to sample
- Imaging depends on tip shape
- Slow

