

GEN III Database Species and Their Characteristics

Gram-Negative Aerobic Bacteria

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
1	<i>Achromobacter cholinophagum</i>	GN-NENT	A	BUG+B	Air	33
2	<i>Achromobacter denitrificans</i>	GN-NENT	A	BUG+B	Air	33
3	<i>Achromobacter insolitus</i>	GN-NENT	A	BUG+B	Air	33
4	<i>Achromobacter piechaudii</i>	GN-NENT	A	BUG+B	Air	33
5	<i>Achromobacter spanius</i>	GN-NENT	A	BUG+B	Air	33
6	<i>Achromobacter xylosoxidans</i> ss <i>xylosoxidans</i>	GN-NENT	A	BUG+B	Air	33
7	<i>Acidovorax avenae</i> ss <i>avenae</i>	GN-NENT	A	BUG	Air	30
8	<i>Acidovorax avenae</i> ss <i>cattleyae</i>	GN-NENT	A	BUG	Air	30
9	<i>Acidovorax avenae</i> ss <i>citrulli</i>	GN-NENT	A	BUG	Air	30
10	<i>Acidovorax delafieldii</i>	GN-NENT	A	BUG	Air	30
11	<i>Acidovorax facilis</i>	GN-NENT	A	BUG	Air	30
12	<i>Acidovorax konjaci</i>	GN-NENT	A	BUG	Air	30
13	<i>Acidovorax temperans</i>	GN-NENT	A	BUG	Air	30
14	<i>Acinetobacter baumannii</i>	GN-NENT	A	BUG+B	Air	33
15	<i>Acinetobacter baylyi</i>	GN-NENT	A	BUG+B	Air	33
16	<i>Acinetobacter bouvetii</i>	GN-NENT	A	BUG+B	Air	33
17	<i>Acinetobacter calcoaceticus</i>	GN-NENT	A	BUG+B	Air	33
18	<i>Acinetobacter genospecies 10</i>	GN-NENT	A	BUG+B	Air	33
19	<i>Acinetobacter genospecies 11</i>	GN-NENT	A	BUG+B	Air	33
20	<i>Acinetobacter genospecies 13TU</i>	GN-NENT	A	BUG+B	Air	33
21	<i>Acinetobacter genospecies 14TU</i>	GN-NENT	A	BUG+B	Air	33
22	<i>Acinetobacter genospecies 15TU</i>	GN-NENT	A	BUG+B	Air	33
23	<i>Acinetobacter genospecies 3</i>	GN-NENT	A	BUG+B	Air	33
24	<i>Acinetobacter genospecies 6</i>	GN-NENT	A	BUG+B	Air	33
25	<i>Acinetobacter gerneri</i>	GN-NENT	A	BUG+B	Air	33
26	<i>Acinetobacter grimontii</i>	GN-NENT	A	BUG+B	Air	33
27	<i>Acinetobacter haemolyticus</i>	GN-NENT	A	BUG+B	Air	33
28	<i>Acinetobacter johnsonii</i>	GN-NENT	A	BUG+B	Air	33
29	<i>Acinetobacter junii</i>	GN-NENT	A	BUG+B	Air	33
30	<i>Acinetobacter lwoffii</i>	GN-NENT	A	BUG+B	Air	33
31	<i>Acinetobacter parvus</i>	GN-NENT	A	BUG+B	Air	33
32	<i>Acinetobacter radioresistens</i>	GN-NENT	A	BUG+B	Air	33
33	<i>Acinetobacter schindleri</i>	GN-NENT	A	BUG+B	Air	33
34	<i>Acinetobacter tandoii</i>	GN-NENT	A	BUG+B	Air	33
35	<i>Acinetobacter tjernbergiae</i>	GN-NENT	A	BUG+B	Air	33
36	<i>Acinetobacter townneri</i>	GN-NENT	A	BUG+B	Air	33
37	<i>Acinetobacter ursingii</i>	GN-NENT	A	BUG+B	Air	33
38	<i>Acinetobacter venetianus</i>	GN-NENT	A	BUG+B	Air	33
39	<i>Actinobacillus capsulatus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
40	<i>Actinobacillus equuli</i> ss <i>equuli</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
41	<i>Actinobacillus hominis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
42	<i>Actinobacillus indolicus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
43	<i>Actinobacillus lignieresii</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
44	<i>Actinobacillus minor</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
45	<i>Actinobacillus muris</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
46	<i>Actinobacillus pleuropneumoniae</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
47	<i>Actinobacillus porcinus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
48	<i>Actinobacillus porcitisillarum</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
49	<i>Actinobacillus rossii</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
50	<i>Actinobacillus seminis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
51	<i>Actinobacillus succinogenes</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
52	<i>Actinobacillus suis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
53	<i>Actinobacillus ureae</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
54	<i>Aeromonas allosaccharophila</i>	GN-NENT	A	BUG+B	Air	33
55	<i>Aeromonas bestiarum</i>	GN-NENT	A	BUG+B	Air	33
56	<i>Aeromonas caviae</i> DNA group 4	GN-NENT	A	BUG+B	Air	33
57	<i>Aeromonas</i> DNA group 11	GN-NENT	A	BUG+B	Air	33
58	<i>Aeromonas encheleia</i>	GN-NENT	A	BUG+B	Air	33
59	<i>Aeromonas enteropelogenes</i> DNA 13	GN-NENT	A	BUG+B	Air	33
60	<i>Aeromonas eucrenophila</i> DNA group 6	GN-NENT	A	BUG+B	Air	33
61	<i>Aeromonas hydrophila</i> DNA group 1	GN-NENT	A	BUG+B	Air	33
62	<i>Aeromonas hydrophila</i> -like DNA group 2	GN-NENT	A	BUG+B	Air	33
63	<i>Aeromonas hydrophila</i> -like DNA group 3	GN-NENT	A	BUG+B	Air	33
64	<i>Aeromonas ichthiosmia</i>	GN-NENT	A	BUG+B	Air	33
65	<i>Aeromonas jandaei</i> DNA group 9	GN-NENT	A	BUG+B	Air	33
66	<i>Aeromonas media</i> DNA group 5B	GN-NENT	A	BUG+B	Air	33
67	<i>Aeromonas media</i> -like DNA group 5A	GN-NENT	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
68	<i>Aeromonas popoffii</i>	GN-NENT	A	BUG+B	Air	33
69	<i>Aeromonas salmonicida</i> ss achromogenes	GN-NENT	A	BUG+B	Air	33
70	<i>Aeromonas salmonicida</i> ss masoucida	GN-NENT	A	BUG+B	Air	33
71	<i>Aeromonas salmonicida</i> ss pectinolytica	GN-NENT	A	BUG+B	Air	33
72	<i>Aeromonas salmonicida</i> ss salmonicida	GN-NENT	A	BUG+B	Air	33
73	<i>Aeromonas schubertii</i> DNA group 12	GN-NENT	A	BUG+B	Air	33
74	<i>Aeromonas sobria</i> DNA group 7	GN-NENT	A	BUG+B	Air	33
75	<i>Aeromonas veronii</i> DNA group 10	GN-NENT	A	BUG+B	Air	33
76	<i>Aeromonas veronii/sobria</i> DNA group 8	GN-NENT	A	BUG+B	Air	33
77	<i>Aggregatibacter actinomycetemcomitans</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
78	<i>Aggregatibacter aphrophilus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
79	<i>Aggregatibacter segnis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
80	<i>Alcaligenes defragrans</i>	GN-NENT	A	BUG+B	Air	33
81	<i>Alcaligenes faecalis</i> ss faecalis	GN-NENT	A	BUG+B	Air	33
82	<i>Alysiella crassa</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
83	<i>Alysiella filiformis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
84	<i>Aminobacter aminovorans</i>	GN-NENT	A	BUG+B	Air	30-33
85	<i>Aquaspirillum peregrinum</i> ss integrum	GN-NENT	A	BUG+B	Air	30-33
86	<i>Aquaspirillum peregrinum</i> ss peregrinum	GN-NENT	A	BUG+B	Air	30-33
87	<i>Aquaspirillum putridiconchylum</i>	GN-NENT	A	BUG+B	Air	30-33
88	<i>Avibacterium avium</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
89	<i>Avibacterium gallinarum</i> (formerly Pasteurella)	GN-FAS	A	BUG+B	6.5% CO2	33-37
90	<i>Avibacterium paragallinarum</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
91	<i>Avibacterium volantium</i> (formerly Pasteurella)	GN-FAS	A	BUG+B	6.5% CO2	33-37
92	<i>Azospirillum brasilense</i>	GN-NENT	A	BUG+B	Air	30-33
93	<i>Bergeriella denitrificans</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
94	<i>Bergeyella zoohelcum</i>	GN-NENT	A	BUG+B	Air	30-33
95	<i>Bibersteinia trehalosi</i> (formerly Pasteurella)	GN-FAS	C2	BUG+B	6.5% CO2	33-37
96	<i>Bordetella avium</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
97	<i>Bordetella bronchiseptica</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
98	<i>Bordetella hinzii</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
99	<i>Bordetella holmesii</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
100	<i>Bordetella parapertussis</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
101	<i>Bordetella petrii</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
102	<i>Bordetella trematum</i>	GN-FAS	C2	BUG+B	6.5% CO2	33-37
103	<i>Brachymonas denitrificans</i>	GN-NENT	A	BUG+B	Air	33
104	<i>Bradyrhizobium japonicum</i>	GN-NENT	A	BUG+B	Air	33
105	<i>Brenneria salicis</i>	GN-ENT	A	BUG	Air	30
106	<i>Brevundimonas diminuta</i>	GN-NENT	A	BUG+B	Air	33
107	<i>Brevundimonas vesicularis</i>	GN-NENT	A	BUG+B	Air	33
108	<i>Budvicia aquatica</i>	GN-ENT	A	BUG+B	Air	33
109	<i>Burkholderia ambifaria</i>	GN-NENT	A	BUG+B	Air	33
110	<i>Burkholderia andropogonis</i>	GN-NENT	A	BUG+B	Air	33
111	<i>Burkholderia anthina/caribensis</i>	GN-NENT	A	BUG+B	Air	33
112	<i>Burkholderia caledonica</i>	GN-NENT	A	BUG+B	Air	33
113	<i>Burkholderia caryophylli</i>	GN-NENT	A	BUG+B	Air	33
114	<i>Burkholderia cepacia</i>	GN-NENT	A	BUG+B	Air	33
115	<i>Burkholderia fungorum</i>	GN-NENT	A	BUG+B	Air	33
116	<i>Burkholderia gladioli</i>	GN-NENT	A	BUG+B	Air	33
117	<i>Burkholderia glathei</i>	GN-NENT	A	BUG+B	Air	33
118	<i>Burkholderia glumae</i>	GN-NENT	A	BUG+B	Air	33
119	<i>Burkholderia graminis</i>	GN-NENT	A	BUG+B	Air	33
120	<i>Burkholderia kururiensis</i>	GN-NENT	A	BUG+B	Air	33
121	<i>Burkholderia multivorans</i>	GN-NENT	A	BUG+B	Air	33
122	<i>Burkholderia phenazinium</i>	GN-NENT	A	BUG+B	Air	33
123	<i>Burkholderia plantarii</i>	GN-NENT	A	BUG+B	Air	33
124	<i>Burkholderia pyrrocinia</i>	GN-NENT	A	BUG+B	Air	33
125	<i>Burkholderia stabilis</i>	GN-NENT	A	BUG+B	Air	33
126	<i>Burkholderia terricola</i>	GN-NENT	A	BUG+B	Air	33
127	<i>Burkholderia vietnamiensis</i>	GN-NENT	A	BUG+B	Air	33
128	<i>Buttiauxella agrestis</i>	GN-ENT	A	BUG+B	Air	33
129	<i>Buttiauxella brennerae</i>	GN-ENT	A	BUG+B	Air	33
130	<i>Buttiauxella ferrugutiae</i>	GN-ENT	A	BUG+B	Air	33
131	<i>Buttiauxella gaviniae</i>	GN-ENT	A	BUG+B	Air	33
132	<i>Buttiauxella izardii</i>	GN-ENT	A	BUG+B	Air	33
133	<i>Buttiauxella noackiae</i>	GN-ENT	A	BUG+B	Air	33
134	<i>Buttiauxella warmboldiae</i>	GN-ENT	A	BUG+B	Air	33
135	<i>Capnocytophaga canimorsus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
136	<i>Capnocytophaga cynodegmi</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
137	<i>Capnocytophaga gingivalis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
138	<i>Capnocytophaga granulosa</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
139	<i>Capnocytophaga haemolytica</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium	Medium	
		Type	Protocol		Atm	Temp
140	Capnocytophaga ochracea/sputigena	GN-FAS	C2	CHOC	6.5% CO2	33-37
141	Cardiobacterium hominis	GN-FAS	C2	CHOC	6.5% CO2	33-37
142	Cardiobacterium valvarum	GN-FAS	C2	CHOC	6.5% CO2	33-37
143	CDC group DF-3 (Capnocytophaga)	GN-FAS	C2	CHOC	6.5% CO2	33-37
144	CDC group EF-4 (Neisseria)	GN-FAS	C2	CHOC	6.5% CO2	33-37
145	CDC group II-E subgroup A	GN-NENT	A	BUG+B	Air	33
146	CDC group II-E subgroup B	GN-NENT	A	BUG+B	Air	33
147	CDC group II-H	GN-NENT	A	BUG+B	Air	33
148	Cedecea davisae	GN-ENT	A	BUG+B	Air	33
149	Cedecea lapagei	GN-ENT	A	BUG+B	Air	33
150	Cedecea neteri	GN-ENT	A	BUG+B	Air	33
151	Chromobacterium violaceum	GN-NENT	A	BUG+B	Air	33
152	Chryseobacterium balustinum	GN-NENT	A	BUG+B	Air	33
153	Chryseobacterium gleum /indologenes	GN-NENT	A	BUG+B	Air	33
154	Chryseobacterium indoltheticum	GN-NENT	A	BUG+B	Air	33
155	Chryseobacterium scophthalmum	GN-NENT	A	BUG+B	Air	33
156	Citrobacter amalonaticus	GN-ENT	A	BUG+B	Air	33
157	Citrobacter braakii	GN-ENT	A	BUG+B	Air	33
158	Citrobacter farmeri	GN-ENT	A	BUG+B	Air	33
159	Citrobacter freundii	GN-ENT	A	BUG+B	Air	33
160	Citrobacter gillenii	GN-ENT	A	BUG+B	Air	33
161	Citrobacter koseri/youngae	GN-ENT	A	BUG+B	Air	33
162	Citrobacter murlinae	GN-ENT	A	BUG+B	Air	33
163	Citrobacter rodentium	GN-ENT	A	BUG+B	Air	33
164	Citrobacter sedlakii	GN-ENT	A	BUG+B	Air	33
165	Citrobacter werkmanii	GN-ENT	A	BUG+B	Air	33
166	Cloacibacterium normanense	GN-NENT	A	BUG+B	Air	33
167	Collimonas fungivorans	GN-NENT	A	BUG+B	Air	33
168	Comamonas aquatica	GN-NENT	A	BUG+B	Air	33
169	Comamonas denitrificans	GN-NENT	A	BUG+B	Air	33
170	Comamonas kerstersii	GN-NENT	A	BUG+B	Air	33
171	Comamonas terrigena	GN-NENT	A	BUG+B	Air	33
172	Comamonas testosteroni	GN-NENT	A	BUG+B	Air	33
173	Conchiformibius steedae	GN-FAS	C2	CHOC	6.5% CO2	33-37
174	Cronobacter dublinensis ss dublinensis	GN-ENT	A	BUG+B	Air	33
175	Cronobacter dublinensis ss lactaridi	GN-ENT	A	BUG+B	Air	33
176	Cronobacter dublinensis ss lausaneensis	GN-ENT	A	BUG+B	Air	33
177	Cronobacter genomospecies 1	GN-ENT	A	BUG+B	Air	33
178	Cronobacter malonaticus	GN-ENT	A	BUG+B	Air	33
179	Cronobacter muytjensii	GN-ENT	A	BUG+B	Air	33
180	Cronobacter sakazakii	GN-ENT	A	BUG+B	Air	33
181	Cronobacter turicensis	GN-ENT	A	BUG+B	Air	33
182	Cupriavidus campinensis	GN-NENT	A	BUG+B	Air	33
183	Cupriavidus gilardii	GN-NENT	A	BUG+B	Air	33
184	Cupriavidus necator	GN-NENT	A	BUG+B	Air	33
185	Cupriavidus pauculus	GN-NENT	A	BUG+B	Air	33
186	Cupriavidus taiwanensis	GN-NENT	A	BUG+B	Air	33
187	Delftia acidovorans	GN-NENT	A	BUG+B	Air	33
188	Delftia tsuruhatensis	GN-NENT	A	BUG+B	Air	33
189	Dickeya chrysanthemi	GN-ENT	A	BUG+B	Air	33
190	Dysgonomonas capnocytophagoides	GN-FAS	C2	CHOC	6.5% CO2	33-37
191	Edwardsiella hoshinae	GN-ENT	A	BUG+B	Air	33
192	Edwardsiella ictaluri	GN-ENT	A	BUG	Air	30
193	Edwardsiella tarda	GN-ENT	A	BUG+B	Air	33
194	Eikenella corrodens	GN-FAS	C2	CHOC	6.5% CO2	33-37
195	Elizabethkingia meningoseptica	GN-NENT	A	BUG+B	Air	33
196	Empedobacter brevis	GN-NENT	A	BUG+B	Air	33
197	Ensifer meliloti	GN-NENT	A	BUG	Air	30
198	Enterobacter aerogenes (Klebsiella mobilis)	GN-ENT	A	BUG+B	Air	33
199	Enterobacter amnigenus	GN-ENT	A	BUG+B	Air	33
200	Enterobacter asburiae	GN-ENT	A	BUG+B	Air	33
201	Enterobacter cancerogenus	GN-ENT	A	BUG+B	Air	33
202	Enterobacter cloacae ss cloacae	GN-ENT	A	BUG+B	Air	33
203	Enterobacter cloacae ss dissolvens	GN-ENT	A	BUG+B	Air	33
204	Enterobacter cowanii	GN-ENT	A	BUG+B	Air	33
205	Enterobacter gergoviae	GN-ENT	A	BUG+B	Air	33
206	Enterobacter helveticus	GN-ENT	A	BUG+B	Air	33
207	Enterobacter hormaechei	GN-ENT	A	BUG+B	Air	33
208	Enterobacter kobei	GN-ENT	A	BUG+B	Air	33
209	Enterobacter nimipressuralis	GN-ENT	A	BUG+B	Air	33
210	Enterobacter pulveris	GN-ENT	A	BUG+B	Air	33
211	Enterobacter turicensis	GN-ENT	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
212	<i>Erwinia amylovora</i>	GN-ENT	A	BUG	Air	30
213	<i>Escherichia albertii</i>	GN-ENT	A	BUG+B	Air	33
214	<i>Escherichia blattae</i>	GN-ENT	A	BUG+B	Air	33
215	<i>Escherichia coli</i>	GN-ENT	A	BUG+B	Air	33
216	<i>Escherichia coli</i> inactive	GN-ENT	A	BUG+B	Air	33
217	<i>Escherichia coli</i> O157:H7	GN-ENT	A	BUG+B	Air	33
218	<i>Escherichia fergusonii</i>	GN-ENT	A	BUG+B	Air	33
219	<i>Escherichia hermannii</i>	GN-ENT	A	BUG+B	Air	33
220	<i>Escherichia vulneris</i>	GN-ENT	A	BUG+B	Air	33
221	<i>Ewingella americana</i>	GN-ENT	A	BUG+B	Air	33
222	<i>Flavimonas oryzae</i>	GN-NENT	A	BUG+B	Air	33
223	<i>Flavobacterium aquatile</i>	GN-NENT	A	BUG+B	Air	33
224	<i>Flavobacterium columnare</i>	GN-NENT	A	BUG+B	Air	33
225	<i>Flavobacterium flevense</i>	GN-NENT	A	BUG+B	Air	33
226	<i>Flavobacterium hydatis</i>	GN-NENT	B	BUG+B	Air	26
227	<i>Flavobacterium johnsoniae</i>	GN-NENT	A	BUG+B	Air	33
228	<i>Flavobacterium mizutaii</i>	GN-NENT	A	BUG+B	Air	33
229	<i>Flavobacterium mizutaii</i> -like (CDC group II-I)	GN-NENT	A	BUG+B	Air	33
230	<i>Flavobacterium resinovorum</i>	GN-NENT	A	BUG+B	Air	33
231	<i>Flavobacterium tirrenicum</i>	GN-NENT	A	BUG+B	Air	33
232	<i>Francisella philomiragia</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
233	<i>Gallibacterium anatis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
234	<i>Gallibacterium genospecies 1</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
235	<i>Grimontia hollisae</i>	GN-NENT	A	BUG+B	Air	26
236	<i>Haemophilus aegyptius</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
237	<i>Haemophilus aphrophilus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
238	<i>Haemophilus ducreyi</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
239	<i>Haemophilus felis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
240	<i>Haemophilus haemoglobinophilus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
241	<i>Haemophilus haemolyticus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
242	<i>Haemophilus influenzae</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
243	<i>Haemophilus intermedius</i> ss <i>intermedius</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
244	<i>Haemophilus paracuniculus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
245	<i>Haemophilus parahaemolyticus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
246	<i>Haemophilus parainfluenzae</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
247	<i>Haemophilus paraphrohaemolyticus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
248	<i>Haemophilus parasuis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
249	<i>Hafnia alvei</i>	GN-ENT	A	BUG+B	Air	33
250	<i>Herbaspirillum frisingense</i>	GN-NENT	A	BUG+B	Air	33
251	<i>Herbaspirillum huttiense</i>	GN-NENT	A	BUG+B	Air	33
252	<i>Herbaspirillum rubrisubalbicans</i>	GN-NENT	A	BUG+B	Air	33
253	<i>Herbaspirillum seropedicae</i>	GN-NENT	A	BUG+B	Air	33
254	<i>Histophilus somnus</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
255	<i>Hydrogenophaga palleronii</i>	GN-NENT	A	BUG+B	Air	26
256	<i>Hydrogenophaga taeniospiralis</i>	GN-NENT	B	BUG+B	Air	26
257	<i>Inquilinus limosus</i> A	GN-NENT	A	BUG+B	Air	33
258	<i>Inquilinus limosus</i> B	GN-NENT	A	BUG+B	Air	33
259	<i>Iodobacter fluviatilis</i>	GN-NENT	A	BUG+B	Air	33
260	<i>Janthinobacterium lividum</i>	GN-NENT	A	BUG+B	Air	26
261	<i>Kerstersia gyiorum</i>	GN-NENT	A	BUG+B	Air	33
262	<i>Kingella denitrificans</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
263	<i>Kingella kingae</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
264	<i>Kingella oralis</i>	GN-FAS	C2	CHOC	6.5% CO2	33-37
265	<i>Klebsiella oxytoca</i>	GN-ENT	A	BUG+B	Air	33
266	<i>Klebsiella pneumoniae</i> ss <i>ozaenae</i>	GN-ENT	A	BUG+B	Air	33
267	<i>Klebsiella pneumoniae</i> ss <i>pneumoniae</i>	GN-ENT	A	BUG+B	Air	33
268	<i>Klebsiella pneumoniae</i> ss <i>rhinoscleromatis</i>	GN-ENT	A	BUG+B	Air	33
269	<i>Kluyvera ascorbata</i>	GN-ENT	A	BUG+B	Air	33
270	<i>Kluyvera cryocrescens</i>	GN-ENT	A	BUG+B	Air	33
271	<i>Kluyvera georgiana</i>	GN-ENT	A	BUG+B	Air	33
272	<i>Kluyvera intermedia</i>	GN-ENT	A	BUG+B	Air	33
273	<i>Lampropedia hyalina</i>	GN-NENT	A	BUG+B	Air	33
274	<i>Lautropia mirabilis</i>	GN-NENT	A	BUG+B	Air	33
275	<i>Leclercia adecarboxylata</i>	GN-ENT	A	BUG+B	Air	33
276	<i>Leminorella grimontii</i>	GN-ENT	A	BUG+B	Air	33
277	<i>Leminorella richardii</i>	GN-ENT	A	BUG+B	Air	33
278	<i>Listonella anguillarum</i>	GN-NENT	A	BUG+B	Air	33
279	<i>Listonella pelagia</i>	GN-NENT	A	BUG+B	Air	33
280	<i>Malikia spinosa</i>	GN-NENT	A	BUG+B	Air	33
281	<i>Mannheimia glucosida</i>	GN-NENT	A	BUG+B	Air	33
282	<i>Mannheimia granulomatis</i>	GN-NENT	A	BUG+B	Air	33
283	<i>Mannheimia haemolytica</i>	GN-NENT	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
284	Mannheimia ruminalis	GN-NENT	A	BUG+B	Air	33
285	Mannheimia varigena	GN-NENT	A	BUG+B	Air	33
286	Massilia timonae	GN-NENT	A	BUG+B	Air	33
287	Methylobacterium extorquens	GN-NENT	C2	R2A	Air	30
288	Methylobacterium fujisawaense	GN-NENT	C2	R2A	Air	30
289	Methylobacterium radiotolerans	GN-NENT	C2	R2A	Air	30
290	Methylobacterium rhodinum	GN-NENT	C2	R2A	Air	30
291	Microvirgula aerodenitrificans	GN-NENT	A	BUG+B	Air	33
292	Moellerella wisconsensis	GN-ENT	A	BUG+B	Air	33
293	Moraxella bovis	GN-FAS	C2	CHOC	6.5% CO2	33-37
294	Moraxella canis	GN-FAS	C2	CHOC	6.5% CO2	33-37
295	Moraxella caprae	GN-FAS	C2	CHOC	6.5% CO2	33-37
296	Moraxella catarrhalis	GN-FAS	C2	CHOC	6.5% CO2	33-37
297	Moraxella caviae	GN-FAS	C2	CHOC	6.5% CO2	33-37
298	Moraxella cuniculi	GN-FAS	C2	CHOC	6.5% CO2	33-37
299	Moraxella equi	GN-FAS	C2	CHOC	6.5% CO2	33-37
300	Moraxella lacunata	GN-FAS	C2	CHOC	6.5% CO2	33-37
301	Moraxella lincolnii	GN-FAS	C2	CHOC	6.5% CO2	33-37
302	Moraxella nonliquefaciens	GN-FAS	C2	CHOC	6.5% CO2	33-37
303	Moraxella osloensis	GN-FAS	C2	CHOC	6.5% CO2	33-37
304	Moraxella ovis	GN-FAS	C2	CHOC	6.5% CO2	33-37
305	Morganella morganii ss morganii	GN-ENT	A	BUG+B	Air	33
306	Morganella morganii ss sibonii	GN-ENT	A	BUG+B	Air	33
307	Myroides odoratimimus	GN-NENT	A	BUG+B	Air	33
308	Myroides odoratus	GN-NENT	A	BUG+B	Air	33
309	Neisseria animalis	GN-FAS	C2	CHOC	6.5% CO2	33-37
310	Neisseria animaloris	GN-FAS	C2	CHOC	6.5% CO2	33-37
311	Neisseria canis	GN-FAS	C2	CHOC	6.5% CO2	33-37
312	Neisseria cinerea	GN-FAS	C2	CHOC	6.5% CO2	33-37
313	Neisseria elongata ss elongata	GN-FAS	C2	CHOC	6.5% CO2	33-37
314	Neisseria flava	GN-FAS	C2	CHOC	6.5% CO2	33-37
315	Neisseria flavescens	GN-FAS	C2	CHOC	6.5% CO2	33-37
316	Neisseria gonorrhoeae	GN-FAS	C2	CHOC	6.5% CO2	33-37
317	Neisseria lactamica	GN-FAS	C2	CHOC	6.5% CO2	33-37
318	Neisseria meningitidis	GN-FAS	C2	CHOC	6.5% CO2	33-37
319	Neisseria mucosa	GN-FAS	C2	CHOC	6.5% CO2	33-37
320	Neisseria perflava	GN-FAS	C2	CHOC	6.5% CO2	33-37
321	Neisseria sicca	GN-FAS	C2	CHOC	6.5% CO2	33-37
322	Neisseria subflava	GN-FAS	C2	CHOC	6.5% CO2	33-37
323	Neisseria weaveri	GN-FAS	C2	CHOC	6.5% CO2	33-37
324	Obesumbacterium proteus	GN-ENT	A	BUG+B	Air	33
325	Obesumbacterium proteus biogroup 2	GN-ENT	A	BUG+B	Air	33
326	Oceanimonas doudoroffii	GN-NENT	A	BUG+B	Air	33
327	Ochrobactrum anthropi	GN-NENT	A	BUG+B	Air	33
328	Ochrobactrum grignonense	GN-NENT	A	BUG+B	Air	33
329	Ochrobactrum intermedium	GN-NENT	A	BUG+B	Air	33
330	Ochrobactrum tritici	GN-NENT	A	BUG+B	Air	33
331	Oligella ureolytica	GN-FAS	C2	CHOC	6.5% CO2	33-37
332	Oligella urethralis	GN-FAS	C2	CHOC	6.5% CO2	33-37
333	Ornithobacterium rhinotracheale	GN-NENT	C2	BUG+B	Air	33
334	Pandoraea apista	GN-NENT	A	BUG+B	Air	33
335	Pandoraea norimbergensis	GN-NENT	A	BUG+B	Air	33
336	Pandoraea pnomenusa	GN-NENT	A	BUG+B	Air	33
337	Pandoraea pulmonicola	GN-NENT	A	BUG+B	Air	33
338	Pandoraea sputorum	GN-NENT	A	BUG+B	Air	33
339	Pantoea agglomerans	GN-ENT	A	BUG+B	Air	33
340	Pantoea agglomerans bgp 2	GN-ENT	A	BUG+B	Air	33
341	Pantoea agglomerans bgp 3	GN-ENT	A	BUG+B	Air	33
342	Pantoea agglomerans bgp 4	GN-ENT	A	BUG+B	Air	33
343	Pantoea agglomerans bgp 5	GN-ENT	A	BUG+B	Air	33
344	Pantoea agglomerans bgp 6	GN-ENT	A	BUG	Air	30
345	Pantoea agglomerans bgp 7	GN-ENT	A	BUG+B	Air	33
346	Pantoea citrea	GN-ENT	A	BUG	Air	30
347	Pantoea dispersa	GN-ENT	A	BUG	Air	30
348	Pantoea punctata	GN-ENT	A	BUG	Air	30
349	Pantoea terrea	GN-ENT	A	BUG	Air	30
350	Paracoccus yeei	GN-NENT	A	BUG+B	Air	33
351	Pasteurella aerogenes	GN-NENT	A	BUG+B	Air	33-37
352	Pasteurella bettyae	GN-NENT	A	BUG+B	Air	33-37
353	Pasteurella caballi	GN-NENT	A	BUG+B	Air	33-37
354	Pasteurella canis/stomatitis	GN-NENT	A	BUG+B	Air	33-37

Organism

Medium

GEN III Database Species and Their Characteristics

	Species Name	Type	Protocol	Medium	Atm	Temp
355	<i>Pasteurella dagmatis</i>	GN-NENT	A	BUG+B	Air	33-37
356	<i>Pasteurella langaensis</i>	GN-NENT	A	BUG+B	Air	33-37
357	<i>Pasteurella lymphangitidis</i>	GN-NENT	C2	BUG+B	Air	33-37
358	<i>Pasteurella mairii</i>	GN-NENT	C2	BUG+B	Air	33-37
359	<i>Pasteurella multocida</i> ss <i>multocida</i>	GN-NENT	A	BUG+B	Air	33-37
360	<i>Pasteurella pneumotropica</i>	GN-NENT	A	BUG+B	Air	33-37
361	<i>Pasteurella testudinis</i>	GN-NENT	A	BUG+B	Air	33-37
362	<i>Pectobacterium atrosepticum</i>	GN-ENT	A	BUG	Air	30
363	<i>Pectobacterium betavascularum</i>	GN-ENT	A	BUG	Air	30
364	<i>Pectobacterium carotovorum</i> ss <i>carotovorum</i>	GN-ENT	A	BUG	Air	30
365	<i>Pectobacterium cypripedii</i>	GN-ENT	A	BUG	Air	30
366	<i>Pedobacter heparinus</i>	GN-NENT	A	BUG	Air	30
367	<i>Photobacterium angustum</i>	GN-NENT	A	BUG+B	Air	26
368	<i>Photobacterium damsela</i> ss <i>damsela</i>	GN-NENT	A	BUG+B	Air	30
369	<i>Photobacterium leiognathi</i>	GN-NENT	A	BUG+B	Air	33
370	<i>Photobacterium luminescens</i> ss <i>luminescens</i>	GN-ENT	A	BUG	Air	30
371	<i>Phyllobacterium myrsinacearum</i>	GN-NENT	A	BUG+B	Air	33
372	<i>Plesiomonas shigelloides</i>	GN-NENT	A	BUG+B	Air	33
373	<i>Pragia fontium</i>	GN-ENT	A	BUG+B	Air	33
374	<i>Proteus hauseri</i>	GN-ENT	A	BUG+B	Air	33
375	<i>Proteus mirabilis</i>	GN-ENT	A	BUG+B	Air	33
376	<i>Proteus myxofaciens</i>	GN-ENT	A	BUG+B	Air	33
377	<i>Proteus penneri/vulgaris</i>	GN-ENT	A	BUG+B	Air	33
378	<i>Providencia alcalifaciens</i>	GN-ENT	A	BUG+B	Air	33
379	<i>Providencia heimbachae</i>	GN-ENT	A	BUG+B	Air	33
380	<i>Providencia rettgeri</i>	GN-ENT	A	BUG+B	Air	33
381	<i>Providencia rustigianii</i>	GN-ENT	A	BUG+B	Air	33
382	<i>Providencia stuartii</i>	GN-ENT	A	BUG+B	Air	33
383	<i>Pseudomonas aeruginosa</i>	GN-NENT	A	BUG+B	Air	33
384	<i>Pseudomonas agarici</i>	GN-NENT	A	BUG+B	Air	33
385	<i>Pseudomonas alcaligenes</i>	GN-NENT	A	BUG+B	Air	33
386	<i>Pseudomonas asplenii</i>	GN-NENT	A	BUG+B	Air	33
387	<i>Pseudomonas bathycetes</i>	GN-NENT	A	BUG+B	Air	33
388	<i>Pseudomonas caricapapayae</i>	GN-NENT	A	BUG	Air	30
389	<i>Pseudomonas chlororaphis</i> ss <i>aurantiaca</i>	GN-NENT	A	BUG+B	Air	33
390	<i>Pseudomonas cichorii</i>	GN-NENT	A	BUG	Air	30
391	<i>Pseudomonas citronellolis</i>	GN-NENT	A	BUG+B	Air	33
392	<i>Pseudomonas corrugata</i>	GN-NENT	A	BUG+B	Air	33
393	<i>Pseudomonas fluorescens</i>	GN-NENT	A	BUG+B	Air	33
394	<i>Pseudomonas fluorescens</i> biotype A	GN-NENT	A	BUG+B	Air	33
395	<i>Pseudomonas fluorescens</i> biotype C	GN-NENT	A	BUG+B	Air	33
396	<i>Pseudomonas fluorescens</i> biotype F	GN-NENT	A	BUG+B	Air	33
397	<i>Pseudomonas fluorescens</i> biotype G	GN-NENT	A	BUG+B	Air	33
398	<i>Pseudomonas fragi</i>	GN-NENT	A	BUG+B	Air	33
399	<i>Pseudomonas fulva</i>	GN-NENT	A	BUG+B	Air	33
400	<i>Pseudomonas fuscovaginae</i>	GN-NENT	A	BUG+B	Air	33
401	<i>Pseudomonas lundensis</i>	GN-NENT	A	BUG+B	Air	33
402	<i>Pseudomonas luteola</i>	GN-NENT	A	BUG+B	Air	33
403	<i>Pseudomonas maculicola</i>	GN-NENT	A	BUG+B	Air	33
404	<i>Pseudomonas marginalis</i>	GN-NENT	A	BUG+B	Air	33
405	<i>Pseudomonas mendocina</i>	GN-NENT	A	BUG+B	Air	33
406	<i>Pseudomonas mephitica</i>	GN-NENT	A	BUG+B	Air	33
407	<i>Pseudomonas mucidolens</i>	GN-NENT	A	BUG+B	Air	33
408	<i>Pseudomonas nitroreducens/azelaica</i>	GN-NENT	A	BUG+B	Air	33
409	<i>Pseudomonas oleovorans</i>	GN-NENT	A	BUG+B	Air	33
410	<i>Pseudomonas pertucinogena</i>	GN-NENT	A	BUG+B	Air	33
411	<i>Pseudomonas pseudoalcaligenes</i> ss <i>pseudoalcaligenes</i>	GN-NENT	A	BUG+B	Air	33
412	<i>Pseudomonas putida</i>	GN-NENT	A	BUG+B	Air	33
413	<i>Pseudomonas putida</i> biotype A	GN-NENT	A	BUG+B	Air	33
414	<i>Pseudomonas putida</i> biotype B	GN-NENT	A	BUG+B	Air	33
415	<i>Pseudomonas resinovorans</i>	GN-NENT	A	BUG+B	Air	33
416	<i>Pseudomonas savastanoi</i> pv <i>fraxini</i>	GN-NENT	A	BUG	Air	30
417	<i>Pseudomonas savastanoi</i> pv <i>nerii</i>	GN-NENT	A	BUG	Air	30
418	<i>Pseudomonas straminea</i>	GN-NENT	A	BUG+B	Air	33
419	<i>Pseudomonas stutzeri</i>	GN-NENT	A	BUG+B	Air	33
420	<i>Pseudomonas synxantha</i>	GN-NENT	A	BUG+B	Air	33
421	<i>Pseudomonas syringae</i> pv <i>aceris</i>	GN-NENT	A	BUG	Air	30
422	<i>Pseudomonas syringae</i> pv <i>antirrhini</i>	GN-NENT	A	BUG	Air	30
423	<i>Pseudomonas syringae</i> pv <i>aptata</i>	GN-NENT	A	BUG	Air	30
424	<i>Pseudomonas syringae</i> pv <i>atrofaciens</i>	GN-NENT	A	BUG	Air	30
425	<i>Pseudomonas syringae</i> pv <i>delphinii</i>	GN-NENT	A	BUG	Air	30
426	<i>Pseudomonas syringae</i> pv <i>glycinea</i>	GN-NENT	A	BUG	Air	30
427	<i>Pseudomonas syringae</i> pv <i>helianthi</i>	GN-NENT	A	BUG	Air	30

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
428	Pseudomonas syringae pv lachrymans	GN-NENT	A	BUG	Air	30
429	Pseudomonas syringae pv phaseolicola	GN-NENT	A	BUG	Air	30
430	Pseudomonas syringae pv pisi	GN-NENT	A	BUG	Air	30
431	Pseudomonas syringae pv porri	GN-NENT	A	BUG	Air	30
432	Pseudomonas syringae pv primulae	GN-NENT	A	BUG	Air	30
433	Pseudomonas syringae pv syringae	GN-NENT	A	BUG	Air	30
434	Pseudomonas syringae pv tabaci A	GN-NENT	A	BUG	Air	30
435	Pseudomonas syringae pv tabaci B	GN-NENT	A	BUG	Air	30
436	Pseudomonas syringae pv tomato	GN-NENT	A	BUG	Air	30
437	Pseudomonas taetrolens	GN-NENT	A	BUG+B	Air	33
438	Pseudomonas tolaasii	GN-NENT	A	BUG+B	Air	33
439	Pseudomonas viridiflava	GN-NENT	A	BUG	Air	30
440	Pseudomonas viridilivida	GN-NENT	A	BUG	Air	30
441	Psychrobacter immobilis	GN-NENT	A	BUG+B	Air	33
442	Psychrobacter phenylpyruvicus	GN-NENT	A	BUG+B	Air	33
443	Rahnella aquatilis	GN-ENT	A	BUG+B	Air	33
444	Ralstonia insidiosa	GN-NENT	A	BUG+B	Air	33
445	Ralstonia mannitolilytica	GN-NENT	A	BUG+B	Air	33
446	Ralstonia pickettii	GN-NENT	A	BUG+B	Air	33
447	Ralstonia solanacearum	GN-NENT	A	BUG+B	Air	33
448	Raoultella planticola/ornithinolytica	GN-ENT	A	BUG+B	Air	33
449	Raoultella terrigena	GN-ENT	A	BUG+B	Air	33
450	Rhizobium radiobacter	GN-NENT	A	BUG	Air	30
451	Rhizobium rhizogenes	GN-NENT	A	BUG	Air	30
452	Rhizobium vitis (bv 3)	GN-NENT	A	BUG	Air	30
453	Riemerella anatipestifer	GN-NENT	A	BUG+B	Air	33
454	Riemerella columbina	GN-NENT	A	BUG+B	Air	33
455	Roseomonas cervicalis	GN-NENT	A	BUG+B	Air	33
456	Roseomonas genomospecies 4	GN-NENT	A	BUG+B	Air	33
457	Roseomonas genomospecies 5	GN-NENT	A	BUG+B	Air	33
458	Roseomonas gilardii ss gilardii	GN-NENT	A	BUG+B	Air	33
459	Salmonella bongori (gp 5)	GN-ENT	A	BUG+B	Air	33
460	Salmonella enterica (gp 1)	GN-ENT	A	BUG+B	Air	33
461	Salmonella enterica (gp 1) choleraesuis	GN-ENT	A	BUG+B	Air	33
462	Salmonella enterica (gp 1) gallinarum	GN-ENT	A	BUG+B	Air	33
463	Salmonella enterica (gp 1) paratyphi A	GN-ENT	A	BUG+B	Air	33
464	Salmonella enterica (gp 1) paratyphi B	GN-ENT	A	BUG+B	Air	33
465	Salmonella enterica (gp 1) paratyphi C	GN-ENT	A	BUG+B	Air	33
466	Salmonella enterica (gp 1) pullorum	GN-ENT	A	BUG+B	Air	33
467	Salmonella enterica (gp 1) typhi	GN-ENT	A	BUG+B	Air	33
468	Salmonella enterica (gp 1) typhimurium	GN-ENT	A	BUG+B	Air	33
469	Salmonella enterica (gp 2) ss salamae	GN-ENT	A	BUG+B	Air	33
470	Salmonella enterica (gp 3A) ss arizonae	GN-ENT	A	BUG+B	Air	35
471	Salmonella enterica (gp 3B) ss diarizonae	GN-ENT	A	BUG+B	Air	33
472	Salmonella enterica (gp 4) ss houtenae	GN-ENT	A	BUG+B	Air	33
473	Salmonella enterica (gp 6) ss indica	GN-ENT	A	BUG+B	Air	33
474	Serpens flexibilis	GN-NENT	A	BUG+B	Air	33
475	Serratia entomophila	GN-ENT	A	BUG+B	Air	33
476	Serratia ficaria	GN-ENT	A	BUG+B	Air	33
477	Serratia fonticola	GN-ENT	A	BUG+B	Air	33
478	Serratia liquefaciens/grimesii	GN-ENT	A	BUG+B	Air	33
479	Serratia marcescens ss marcescens	GN-ENT	A	BUG+B	Air	33
480	Serratia odorifera	GN-ENT	A	BUG+B	Air	33
481	Serratia plymuthica	GN-ENT	A	BUG+B	Air	33
482	Serratia proteamaculans	GN-ENT	A	BUG+B	Air	33
483	Serratia rubidaea	GN-ENT	A	BUG+B	Air	33
484	Shewanella algae	GN-NENT	A	BUG+B	Air	33
485	Shewanella putrefaciens A	GN-NENT	A	BUG+B	Air	33
486	Shewanella putrefaciens B	GN-NENT	A	BUG+B	Air	33
487	Shigella boydii	GN-ENT	A	BUG+B	Air	33
488	Shigella dysenteriae	GN-ENT	A	BUG+B	Air	33
489	Shigella flexneri	GN-ENT	A	BUG+B	Air	33
490	Shigella sonnei	GN-ENT	A	BUG+B	Air	33
491	Simonsiella muelleri	GN-FAS	C2	CHOC	6.5% CO2	33-37
492	Simplicispira metamorpha	GN-NENT	A	BUG+B	Air	33
493	Sphingobacterium multivorum	GN-NENT	A	BUG+B	Air	33
494	Sphingobacterium spiritovorum	GN-NENT	A	BUG+B	Air	33
495	Sphingobacterium thalpophilum	GN-NENT	A	BUG+B	Air	33
496	Sphingomonas adhaesiva	GN-NENT	A	BUG+B	Air	33
497	Sphingomonas capsulata	GN-NENT	A	BUG+B	Air	33
498	Sphingomonas macrogoltabidus	GN-NENT	A	BUG+B	Air	33
499	Sphingomonas parapaucimobilis	GN-NENT	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium	Medium	
		Type	Protocol		Atm	Temp
500	Sphingomonas paucimobilis A	GN-NENT	A	BUG+B	Air	33
501	Sphingomonas paucimobilis B	GN-NENT	A	BUG+B	Air	33
502	Sphingomonas sanguinis	GN-NENT	A	BUG+B	Air	33
503	Sphingomonas terrae	GN-NENT	A	BUG+B	Air	33
504	Sphingomonas yanoikuyae	GN-NENT	A	BUG+B	Air	33
505	Stenotrophomonas acidaminiphila	GN-NENT	A	BUG+B	Air	33
506	Stenotrophomonas maltophilia	GN-NENT	A	BUG+B	Air	33
507	Stenotrophomonas rhizophila	GN-NENT	A	BUG+B	Air	33
508	Suttonella indologenes	GN-FAS	C2	CHOC	6.5% CO2	33-37
509	Tatumella ptyseos	GN-ENT	A	BUG+B	Air	33
510	Taylorella equigenitalis	GN-FAS	C2	CHOC	6.5% CO2	33-37
511	Terrimonas ferruginea	GN-NENT	A	BUG+B	Air	33
512	Trabulsiella guamensis	GN-ENT	A	BUG+B	Air	33
513	Variovorax paradoxus	GN-NENT	A	BUG+B	Air	33
514	Vibrio aerogenes	GN-NENT	A	BUG+B	Air	33
515	Vibrio aestuarianus	GN-NENT	A	BUG+B	Air	33
516	Vibrio alginolyticus	GN-NENT	A	BUG+B	Air	33
517	Vibrio campbellii	GN-NENT	A	BUG+B	Air	33
518	Vibrio cholerae	GN-NENT	B	BUG+B	Air	33
519	Vibrio cholerae non O1	GN-NENT	B	BUG+B	Air	33
520	Vibrio cholerae O1/ATCC 25870	GN-NENT	B	BUG+B	Air	33
521	Vibrio cholerae O1/O139	GN-NENT	B	BUG+B	Air	33
522	Vibrio cincinnatiensis	GN-NENT	A	BUG+B	Air	33
523	Vibrio diazotrophicus	GN-NENT	A	BUG+B	Air	33
524	Vibrio fluvialis	GN-NENT	A	BUG+B	Air	33
525	Vibrio furnissii	GN-NENT	A	BUG+B	Air	33
526	Vibrio harveyi	GN-NENT	A	BUG+B	Air	33
527	Vibrio mediterranei	GN-NENT	A	BUG+B	Air	33
528	Vibrio metschnikovii	GN-NENT	A	BUG+B	Air	33
529	Vibrio mimicus	GN-NENT	A	BUG+B	Air	33
530	Vibrio natriegens	GN-NENT	A	BUG+B	Air	33
531	Vibrio ordalii	GN-NENT	A	BUG+B	Air	33
532	Vibrio parahaemolyticus	GN-NENT	A	BUG+B	Air	33
533	Vibrio proteolyticus	GN-NENT	A	BUG+B	Air	33
534	Vibrio splendidus	GN-NENT	A	BUG+B	Air	26
535	Vibrio tubiashii	GN-NENT	A	BUG+B	Air	33
536	Vibrio vulnificus	GN-NENT	A	BUG+B	Air	33
537	Vogesella indigofera	GN-NENT	A	BUG+B	Air	33
538	Weeksella virosa	GN-NENT	A	BUG+B	Air	33
539	Xanthomonas campestris pv begoniae A	GN-NENT	A	BUG	Air	30
540	Xanthomonas campestris pv begoniae B	GN-NENT	A	BUG	Air	30
541	Xanthomonas campestris pv campestris	GN-NENT	A	BUG	Air	30
542	Xanthomonas campestris pv carotae	GN-NENT	A	BUG	Air	30
543	Xanthomonas campestris pv dieffenbachiae	GN-NENT	A	BUG	Air	30
544	Xanthomonas campestris pv juglandis	GN-NENT	A	BUG	Air	30
545	Xanthomonas campestris pv malvacearum	GN-NENT	A	BUG	Air	30
546	Xanthomonas campestris pv pelargonii	GN-NENT	A	BUG	Air	30
547	Xanthomonas campestris pv phaseoli	GN-NENT	A	BUG	Air	30
548	Xanthomonas campestris pv poinsetticola	GN-NENT	A	BUG	Air	30
549	Xanthomonas campestris pv raphani	GN-NENT	A	BUG	Air	30
550	Xanthomonas campestris pv translucens	GN-NENT	A	BUG	Air	30
551	Xanthomonas campestris pv vesicatoria	GN-NENT	A	BUG	Air	30
552	Xenorhabdus bovienii	GN-ENT	A	BUG	Air	30
553	Xenorhabdus nematophila	GN-ENT	A	BUG	Air	30
554	Yersinia aldovae	GN-ENT	A	BUG+B	Air	33
555	Yersinia bercovieri	GN-ENT	A	BUG+B	Air	33
556	Yersinia enterocolitica ss enterocolitica	GN-ENT	A	BUG+B	Air	33
557	Yersinia fredericksonii	GN-ENT	A	BUG+B	Air	33
558	Yersinia intermedia	GN-ENT	A	BUG+B	Air	33
559	Yersinia kristensenii	GN-ENT	A	BUG+B	Air	33
560	Yersinia mollaretii	GN-ENT	A	BUG+B	Air	33
561	Yersinia pseudotuberculosis	GN-ENT	A	BUG+B	Air	33
562	Yersinia rohdei	GN-ENT	A	BUG+B	Air	33
563	Yersinia ruckeri	GN-ENT	A	BUG+B	Air	33
564	Yokenella regensburgei	GN-ENT	A	BUG+B	Air	33

Gram-Positive Aerobic Bacteria

	Species Name	Organism		Medium	Medium	
		Type	Protocol		Atm	Temp
1	<i>Actinomyces bovis</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
2	<i>Actinomyces canis</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
3	<i>Actinomyces hordeovulneris</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
4	<i>Actinomyces hyovaginalis</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
5	<i>Actinomyces naeslundii</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
6	<i>Actinomyces neuii ss anitratus</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
7	<i>Actinomyces neuii ss neuii</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
8	<i>Actinomyces odontolyticus</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
9	<i>Aerococcus christensenii</i>	GP-COCCUS	C2	BUG+B	6.5% CO2	33-37
10	<i>Aerococcus sanguinicola</i>	GP-COCCUS	C2	BUG+B	6.5% CO2	33-37
11	<i>Aerococcus urinae</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
12	<i>Aerococcus urinaeequi</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
13	<i>Aerococcus urinaehominis</i>	GP-COCCUS	C2	BUG+B	6.5% CO2	33-37
14	<i>Aerococcus viridans</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
15	<i>Alloiooccus otitis</i>	GP-COCCUS	C2	BUG+B	6.5% CO2	33-37
16	<i>Aneurinibacillus aneurinilyticus</i>	GP-ROD SB	A	BUG+B	Air	33
17	<i>Arcanobacterium bernardiae</i> (CDC.2)	GP-ROD	C1	BUG+B	6.5% CO2	33-37
18	<i>Arcanobacterium haemolyticum</i>	GP-ROD	C1	BUG+B	6.5% CO2	33-37
19	<i>Arcanobacterium hippocoleae</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
20	<i>Arcanobacterium phocae</i>	GP-ROD	C1	BUG+B	6.5% CO2	33-37
21	<i>Arcanobacterium pluranimalium</i>	GP-ROD	C2	BUG+B	6.5% CO2	33-37
22	<i>Arcanobacterium pyogenes</i>	GP-ROD	C1	BUG+B	6.5% CO2	33-37
23	<i>Arthrobacter cummingsii</i>	GP-ROD	A	BUG+B	Air	33
24	<i>Arthrobacter globiformis</i>	GP-ROD	A	BUG+B	Air	33
25	<i>Arthrobacter histidinovorans</i>	GP-ROD	A	BUG+B	Air	33
26	<i>Arthrobacter ilicis</i>	GP-ROD	A	BUG+B	Air	33
27	<i>Arthrobacter nicotianae</i>	GP-ROD	A	BUG+B	Air	33
28	<i>Arthrobacter woluwensis</i>	GP-ROD	A	BUG+B	Air	33
29	<i>Arthrobacter xylinum</i>	GP-ROD	A	BUG+B	Air	33
30	<i>Bacillus amyloliquefaciens</i>	GP-ROD SB	A	BUG+B	Air	33
31	<i>Bacillus atrophaeus</i>	GP-ROD SB	A	BUG+B	Air	33
32	<i>Bacillus badius</i>	GP-ROD SB	A	BUG+B	Air	33
33	<i>Bacillus cereus/thuringiensis</i>	GP-ROD SB	B	BUG+B	Air	33
34	<i>Bacillus circulans</i>	GP-ROD SB	A	BUG+B	Air	33
35	<i>Bacillus clausii</i>	GP-ROD SB	A	BUG+B	Air	33
36	<i>Bacillus coagulans</i>	GP-ROD SB	A	BUG+B	Air	33
37	<i>Bacillus cohnii</i>	GP-ROD SB	B	BUG+B	Air	33
38	<i>Bacillus firmus</i>	GP-ROD SB	A	BUG+B	Air	33
39	<i>Bacillus flexus</i>	GP-ROD SB	B	BUG+B	Air	33
40	<i>Bacillus gibsonii</i>	GP-ROD SB	A	BUG+B	Air	33
41	<i>Bacillus halodurans</i>	GP-ROD SB	A	BUG+B	Air	33
42	<i>Bacillus horikoshii</i>	GP-ROD SB	A	BUG+B	Air	33
43	<i>Bacillus humi</i>	GP-ROD SB	A	BUG+B	Air	33
44	<i>Bacillus lentus</i>	GP-ROD SB	A	BUG+B	Air	26
45	<i>Bacillus licheniformis</i>	GP-ROD SB	A	BUG+B	Air	33
46	<i>Bacillus maroccanus</i>	GP-ROD SB	B	BUG+B	Air	33
47	<i>Bacillus megaterium</i>	GP-ROD SB	B	BUG+B	Air	33
48	<i>Bacillus mycoides</i>	GP-ROD SB	B	BUG+B	Air	33
49	<i>Bacillus niacini</i>	GP-ROD SB	B	BUG+B	Air	33
50	<i>Bacillus oleronius</i>	GP-ROD SB	A	BUG+B	Air	33
51	<i>Bacillus pseudofirmus</i>	GP-ROD SB	A	BUG+B	Air	33
52	<i>Bacillus pseudomyoides</i>	GP-ROD SB	A	BUG+B	Air	33
53	<i>Bacillus psychrodurans</i>	GP-ROD SB	A	BUG+B	Air	33
54	<i>Bacillus psychrosaccharolyticus</i>	GP-ROD SB	A	BUG+B	Air	26
55	<i>Bacillus psychrotolerans</i>	GP-ROD SB	A	BUG+B	Air	33
56	<i>Bacillus pumilus</i>	GP-ROD SB	A	BUG+B	Air	33
57	<i>Bacillus silvestris</i>	GP-ROD SB	A	BUG+B	Air	33
58	<i>Bacillus simplex</i>	GP-ROD SB	B	BUG+B	Air	33
59	<i>Bacillus subtilis ss spizizenii</i>	GP-ROD SB	A	BUG+B	Air	33
60	<i>Bacillus subtilis ss subtilis</i>	GP-ROD SB	A	BUG+B	Air	33
61	<i>Bacillus weihenstephanensis</i>	GP-ROD SB	A	BUG+B	Air	33
62	<i>Brevibacillus agri</i>	GP-ROD SB	A	BUG+B	Air	33
63	<i>Brevibacillus borstelensis</i>	GP-ROD SB	A	BUG+B	Air	33
64	<i>Brevibacillus brevis</i>	GP-ROD SB	A	BUG+B	Air	33
65	<i>Brevibacillus choshinensis</i>	GP-ROD SB	B	BUG+B	Air	33
66	<i>Brevibacillus invocatus</i>	GP-ROD SB	A	BUG+B	Air	33
67	<i>Brevibacillus laterosporus</i>	GP-ROD SB	A	BUG+B	Air	33
68	<i>Brevibacillus parabravis</i>	GP-ROD SB	A	BUG+B	Air	33
69	<i>Brevibacterium casei</i>	GP-ROD	A	BUG+B	Air	33
70	<i>Brevibacterium epidermidis</i>	GP-ROD	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
71	Brevibacterium frigoritolerans	GP-ROD	B	BUG+B	Air	33
72	Brevibacterium halotolerans	GP-ROD	A	BUG+B	Air	33
73	Brevibacterium linens	GP-ROD	A	BUG+B	Air	33
74	Brevibacterium mcbrellneri	GP-ROD	A	BUG+B	Air	33
75	Brevibacterium otitidis	GP-ROD	A	BUG+B	Air	33
76	Brevibacterium paucivorans	GP-ROD	A	BUG+B	Air	33
77	Brochothrix campestris	GP-ROD	A	BUG+B	Air	30
78	Brochothrix thermosphacta	GP-ROD	A	BUG+B	Air	30
79	Carnobacterium alterfunditum	GP-ROD	C2	BUG+B	Air	30
80	Carnobacterium divergens	GP-ROD	A	BUG+B	Air	30
81	Carnobacterium gallinarum	GP-ROD	A	BUG+B	Air	30
82	Carnobacterium inhibens	GP-ROD	C2	BUG+B	Air	30
83	Carnobacterium maltaromaticum	GP-ROD	A	BUG+B	Air	30
84	Carnobacterium mobile	GP-ROD	A	BUG+B	Air	30
85	Cellulomonas biazotea	GP-ROD	A	BUG+B	Air	30-33
86	Cellulomonas cellasea	GP-ROD	C1	BUG+B	Air	30-33
87	Cellulomonas fimi	GP-ROD	A	BUG+B	Air	30-33
88	Cellulomonas flavigena	GP-ROD	A	BUG+B	Air	30-33
89	Cellulomonas gelida	GP-ROD	A	BUG+B	Air	30-33
90	Cellulomonas hominis (CDC.A-3)	GP-ROD	A	BUG+B	Air	30-33
91	Cellulomonas uda	GP-ROD	A	BUG+B	Air	30-33
92	Cellulosimicrobium cellulans	GP-ROD	A	BUG+B	Air	30-33
93	Clavibacter michiganensis ss insidiosus	GP-ROD	A	BUG	Air	30
94	Clavibacter michiganensis ss michiganensis	GP-ROD	A	BUG	Air	30
95	Clavibacter michiganensis ss nebraskensis	GP-ROD	A	BUG	Air	30
96	Clavibacter michiganensis ss sepedonicus	GP-ROD	A	BUG	Air	30
97	Clavibacter michiganensis ss tessellarius	GP-ROD	A	BUG	Air	30
98	Corynebacterium accolens	GP-ROD	A	BUG+B	Air	33
99	Corynebacterium afermentans ss afermentans(CDC.ANF-1)	GP-ROD	A	BUG+B	Air	33
100	Corynebacterium afermentans ss lipophilum	GP-ROD	A	BUG+B	Air	33
101	Corynebacterium ammoniagenes	GP-ROD	A	BUG+B	Air	33
102	Corynebacterium amycolatum (CDC.F-2)	GP-ROD	A	BUG+B	Air	33
103	Corynebacterium appendicis	GP-ROD	A	BUG+B	Air	33
104	Corynebacterium argensoratense	GP-ROD	A	BUG+B	Air	33
105	Corynebacterium aurimucosum	GP-ROD	A	BUG+B	Air	33
106	Corynebacterium auris	GP-ROD	A	BUG+B	Air	33
107	Corynebacterium auriscanis	GP-ROD	A	BUG+B	Air	33
108	Corynebacterium bovis	GP-ROD	A	BUG+B	Air	33
109	Corynebacterium callunae	GP-ROD	A	BUG+B	Air	33
110	Corynebacterium camporealensis	GP-ROD	A	BUG+B	Air	33
111	Corynebacterium capitovis	GP-ROD	A	BUG+B	Air	33
112	Corynebacterium confusum	GP-ROD	A	BUG+B	Air	33
113	Corynebacterium coyleae	GP-ROD	A	BUG+B	Air	33
114	Corynebacterium cystitidis	GP-ROD	C2	BUG+B	Air	33
115	Corynebacterium diphtheriae	GP-ROD	A	BUG+B	Air	33
116	Corynebacterium durum	GP-ROD	A	BUG+B	Air	33
117	Corynebacterium efficiens	GP-ROD	A	BUG+B	Air	33
118	Corynebacterium falsenii	GP-ROD	A	BUG+B	Air	33
119	Corynebacterium felinum	GP-ROD	A	BUG+B	Air	33
120	Corynebacterium flavescens	GP-ROD	A	BUG+B	Air	33
121	Corynebacterium freneyi	GP-ROD	A	BUG+B	Air	33
122	Corynebacterium glucuronolyticum	GP-ROD	A	BUG+B	Air	33
123	Corynebacterium glutamicum	GP-ROD	A	BUG+B	Air	33
124	Corynebacterium imitans	GP-ROD	A	BUG+B	Air	33
125	Corynebacterium jeikeium	GP-ROD	A	BUG+B	Air	33
126	Corynebacterium kroppenstedtii	GP-ROD	C2	BUG+B	Air	33
127	Corynebacterium kutscheri	GP-ROD	A	BUG+B	Air	33
128	Corynebacterium lipophiloflavum	GP-ROD	C1	BUG+B	Air	33
129	Corynebacterium macginleyi	GP-ROD	A	BUG+B	Air	33
130	Corynebacterium mastitidis	GP-ROD	A	BUG+B	Air	33
131	Corynebacterium matruchotii	GP-ROD	A	BUG+B	Air	33
132	Corynebacterium minutissimum	GP-ROD	A	BUG+B	Air	33
133	Corynebacterium mucifaciens	GP-ROD	A	BUG+B	Air	33
134	Corynebacterium mycetoides	GP-ROD	A	BUG+B	Air	33
135	Corynebacterium phocae	GP-ROD	A	BUG+B	Air	33
136	Corynebacterium pilosum	GP-ROD	C2	BUG+B	Air	33
137	Corynebacterium pseudodiphtheriticum/propinquum (CDC.ANF-3)	GP-ROD	A	BUG+B	Air	33
138	Corynebacterium pseudotuberculosis	GP-ROD	A	BUG+B	Air	33
139	Corynebacterium renale	GP-ROD	A	BUG+B	Air	33
140	Corynebacterium riegliei	GP-ROD	A	BUG+B	Air	33
141	Corynebacterium seminale	GP-ROD	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium	Medium	
		Type	Test		Atm	Temp
142	<i>Corynebacterium simulans</i>	GP-ROD	A	BUG+B	Air	33
143	<i>Corynebacterium singulare</i>	GP-ROD	A	BUG+B	Air	33
144	<i>Corynebacterium spheioscorum</i>	GP-ROD	A	BUG+B	Air	33
145	<i>Corynebacterium</i> spp. (CDC.G)	GP-ROD	A	BUG+B	Air	33
146	<i>Corynebacterium striatum</i> (CDC.I-1)	GP-ROD	A	BUG+B	Air	33
147	<i>Corynebacterium terpenotabidum</i>	GP-ROD	A	BUG+B	Air	33
148	<i>Corynebacterium testudinoris</i>	GP-ROD	A	BUG+B	Air	33
149	<i>Corynebacterium thomssenii</i>	GP-ROD	C1	BUG+B	Air	33
150	<i>Corynebacterium ulcerans</i>	GP-ROD	A	BUG+B	Air	33
151	<i>Corynebacterium urealyticum</i>	GP-ROD	A	BUG+B	Air	33
152	<i>Corynebacterium variabile</i>	GP-ROD	A	BUG+B	Air	33
153	<i>Corynebacterium vitaeruminis</i>	GP-ROD	A	BUG+B	Air	33
154	<i>Corynebacterium xerosis</i> (GPC)	GP-COCCUS	A	BUG+B	Air	33
155	<i>Curtobacterium albidum</i>	GP-ROD	A	BUG	Air	30
156	<i>Curtobacterium citreum</i>	GP-ROD	A	BUG	Air	30
157	<i>Curtobacterium flaccumfaciens</i>	GP-ROD	A	BUG	Air	30
158	<i>Curtobacterium luteum</i>	GP-ROD	A	BUG	Air	30
159	<i>Curtobacterium pusillum</i>	GP-ROD	A	BUG	Air	30
160	<i>Deinococcus proteolyticus</i>	GP-COCCUS	B	BUG+B	Air	30
161	<i>Deinococcus radiodurans</i>	GP-COCCUS	A	BUG+B	Air	30
162	<i>Deinococcus radiopugnans</i>	GP-COCCUS	A	BUG+B	Air	30
163	<i>Dermabacter hominis</i>	GP-ROD	A	BUG+B	Air	33
164	<i>Dermaococcus nishinomiyaensis</i>	GP-COCCUS	A	BUG+B	Air	33
165	<i>Dietzia maris</i>	GP-ROD	A	BUG+B	Air	33
166	<i>Dolosicoccus paucivorans</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
167	<i>Dolosigranulum pigrum</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
168	<i>Enterococcus asini</i>	GP-COCCUS	A	BUG+B	Air	33
169	<i>Enterococcus avium</i>	GP-COCCUS	A	BUG+B	Air	33
170	<i>Enterococcus casseliflavus</i>	GP-COCCUS	A	BUG+B	Air	33
171	<i>Enterococcus cecorum</i>	GP-COCCUS	A	BUG+B	Air	33
172	<i>Enterococcus columbae</i>	GP-COCCUS	A	BUG+B	Air	33
173	<i>Enterococcus dispar</i>	GP-COCCUS	A	BUG+B	Air	33
174	<i>Enterococcus durans</i>	GP-COCCUS	A	BUG+B	Air	33
175	<i>Enterococcus faecalis</i>	GP-COCCUS	A	BUG+B	Air	33
176	<i>Enterococcus faecium</i>	GP-COCCUS	A	BUG+B	Air	33
177	<i>Enterococcus gallinarum</i>	GP-COCCUS	A	BUG+B	Air	33
178	<i>Enterococcus gilvus</i>	GP-COCCUS	A	BUG+B	Air	33
179	<i>Enterococcus haemoperoxidus</i>	GP-COCCUS	A	BUG+B	Air	33
180	<i>Enterococcus hirae</i>	GP-COCCUS	A	BUG+B	Air	33
181	<i>Enterococcus malodoratus</i>	GP-COCCUS	A	BUG+B	Air	33
182	<i>Enterococcus mundtii</i>	GP-COCCUS	A	BUG+B	Air	33
183	<i>Enterococcus pallens</i>	GP-COCCUS	A	BUG+B	Air	33
184	<i>Enterococcus pseudoavium</i>	GP-COCCUS	C1	BUG+B	Air	33
185	<i>Enterococcus raffinosus</i>	GP-COCCUS	A	BUG+B	Air	33
186	<i>Enterococcus ratti</i>	GP-COCCUS	C2	BUG+B	Air	33
187	<i>Enterococcus saccharolyticus</i>	GP-COCCUS	A	BUG+B	Air	33
188	<i>Enterococcus sulfureus</i>	GP-COCCUS	A	BUG+B	Air	33
189	<i>Enterococcus villorum</i>	GP-COCCUS	A	BUG+B	Air	33
190	<i>Eremococcus coleocola</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
191	<i>Erysipelothrix rhusiopathiae/tonsillarum</i>	GP-ROD	C2	BUG+B	6.5%CO2	33-37
192	<i>Exiguobacterium acetylicum</i>	GP-ROD	A	BUG+B	Air	33
193	<i>Exiguobacterium aurantiacum</i>	GP-ROD	A	BUG+B	Air	33
194	<i>Exiguobacterium undae</i>	GP-ROD	A	BUG+B	Air	33
195	<i>Gardnerella vaginalis</i>	GP-ROD	C2	CHOC	6.5%CO2	33-37
196	<i>Gemella bergeri</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
197	<i>Gemella cuniculi</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
198	<i>Gemella haemolysans/morbillorum</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
199	<i>Gemella palaticanis</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
200	<i>Gemella sanguinis</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
201	<i>Geobacillus stearothermophilus</i>	GP-ROD SB	A	BUG+B	Air	55
202	<i>Globicatella sanguinis</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
203	<i>Gordonia aichiensis</i>	GP-ROD	A	BUG+B	Air	30
204	<i>Gordonia bronchialis</i>	GP-ROD	A	BUG+B	Air	30
205	<i>Gordonia rubripertincta</i>	GP-ROD	A	BUG+B	Air	30
206	<i>Gordonia sputi</i>	GP-ROD	A	BUG+B	Air	30
207	<i>Gordonia terrae</i>	GP-ROD	A	BUG+B	Air	30
208	<i>Granulicatell adiacens</i>	GP-COCCUS	C2	BUG+B	6.5%CO2	33-37
209	<i>Granulicatell elegans</i>	GP-COCCUS	C2	BUG+B	6.5%CO2	33-37
210	<i>Helcococcus kunzii</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
211	<i>Helcococcus ovis</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
212	<i>Ignavigranum rouffiae</i>	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
213	<i>Jonesia denitrificans</i>	GP-ROD	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium	Medium	
		Type	Protocol		Atm	Temp
214	Kocuria kristinae	GP-COCCUS	A	BUG+B	Air	33
215	Kocuria rhizophila	GP-COCCUS	A	BUG+B	Air	33
216	Kocuria rosea	GP-COCCUS	A	BUG+B	Air	33
217	Kocuria varians	GP-COCCUS	A	BUG+B	Air	33
218	Kurthia gibsonii	GP-ROD	A	BUG+B	Air	33
219	Kurthia sibirica	GP-ROD	A	BUG+B	Air	33
220	Kurthia zopfii	GP-ROD	A	BUG+B	Air	33
221	Kytococcus schroeteri	GP-COCCUS	A	BUG+B	Air	33
222	Kytococcus sedentarius	GP-COCCUS	A	BUG+B	Air	33
223	Lactobacillus acidophilus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
224	Lactobacillus alimentarius	GP-ROD	C2	BUG+B	6.5%CO2	33-37
225	Lactobacillus amylophilus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
226	Lactobacillus animalis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
227	Lactobacillus aviarius ss araffinosus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
228	Lactobacillus bifermentans	GP-ROD	C1	BUG+B	6.5%CO2	33-37
229	Lactobacillus brevis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
230	Lactobacillus buchneri	GP-ROD	C2	BUG+B	6.5%CO2	33-37
231	Lactobacillus casei	GP-ROD	C1	BUG+B	6.5%CO2	33-37
232	Lactobacillus coryniformis ss coryniformis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
233	Lactobacillus coryniformis ss torquens	GP-ROD	C2	BUG+B	6.5%CO2	33-37
234	Lactobacillus crispatus	GP-ROD	C2	BUG+B	6.5%CO2	33-37
235	Lactobacillus curvatus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
236	Lactobacillus delbrueckii ss delbrueckii	GP-ROD	C2	BUG+B	6.5%CO2	33-37
237	Lactobacillus delbrueckii ss lactis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
238	Lactobacillus farciminis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
239	Lactobacillus fermentum	GP-ROD	C2	BUG+B	6.5%CO2	33-37
240	Lactobacillus fructivorans	GP-ROD	C1	BUG+B	6.5%CO2	33-37
241	Lactobacillus gasseri	GP-ROD	C1	BUG+B	6.5%CO2	33-37
242	Lactobacillus hamsteri	GP-ROD	C2	BUG+B	6.5%CO2	33-37
243	Lactobacillus helveticus	GP-ROD	C2	BUG+B	6.5%CO2	33-37
244	Lactobacillus jensenii	GP-ROD	C2	BUG+B	6.5%CO2	33-37
245	Lactobacillus kefir	GP-ROD	C2	BUG+B	6.5%CO2	33-37
246	Lactobacillus mali	GP-ROD	C1	BUG+B	6.5%CO2	33-37
247	Lactobacillus murinus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
248	Lactobacillus oris	GP-ROD	C1	BUG+B	6.5%CO2	33-37
249	Lactobacillus parabuchneri	GP-ROD	C2	BUG+B	6.5%CO2	33-37
250	Lactobacillus paracasei ss tolerans	GP-ROD	C2	BUG+B	6.5%CO2	33-37
251	Lactobacillus pentosus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
252	Lactobacillus plantarum	GP-ROD	C1	BUG+B	6.5%CO2	33-37
253	Lactobacillus reuteri	GP-ROD	C1	BUG+B	6.5%CO2	33-37
254	Lactobacillus rhamnosus	GP-ROD	C1	BUG+B	6.5%CO2	33-37
255	Lactobacillus sake	GP-ROD	C2	BUG+B	6.5%CO2	33-37
256	Lactobacillus salivarius ss salicinius	GP-ROD	C1	BUG+B	6.5%CO2	33-37
257	Lactobacillus salivarius ss salivarius	GP-ROD	C1	BUG+B	6.5%CO2	33-37
258	Lactobacillus vaginalis	GP-ROD	C1	BUG+B	6.5%CO2	33-37
259	Lactococcus garvieae	GP-COCCUS	A	BUG+B	6.5%CO2	33-37
260	Lactococcus lactis ss cremoris	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
261	Lactococcus lactis ss hordniae	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
262	Lactococcus lactis ss lactis	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
263	Lactococcus plantarum	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
264	Lactococcus raffinolactis	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
265	Leifsonia aquatica	GP-ROD	A	BUG+B	Air	33
266	Leifsonia poae	GP-ROD	A	BUG+B	Air	33
267	Leuconostoc carnosum	GP-COCCUS	C2	BUG+B	Air	26
268	Leuconostoc citreum	GP-COCCUS	C1	BUG+B	Air	33
269	Leuconostoc fallax	GP-COCCUS	A	BUG+B	Air	33
270	Leuconostoc gelidum	GP-COCCUS	A	BUG+B	Air	33
271	Leuconostoc lactis	GP-COCCUS	C1	BUG+B	Air	33
272	Leuconostoc mesenteroides ss cremoris	GP-COCCUS	A	BUG+B	Air	33
273	Leuconostoc mesenteroides ss dextranicum	GP-COCCUS	A	BUG+B	Air	33
274	Leuconostoc mesenteroides ss mesenteroides	GP-COCCUS	A	BUG+B	Air	33
275	Listeria grayi	GP-ROD	A	BUG+B	Air	33
276	Listeria innocua	GP-ROD	A	BUG+B	Air	33
277	Listeria ivanovii ss ivanovii	GP-ROD	A	BUG+B	Air	33
278	Listeria ivanovii ss londoniensis	GP-ROD	A	BUG+B	Air	33
279	Listeria monocytogenes	GP-ROD	A	BUG+B	Air	33
280	Listeria seeligeri	GP-ROD	A	BUG+B	Air	33
281	Listeria welshimeri	GP-ROD	A	BUG+B	Air	33
282	Lysinibacillus fusiformis	GP-ROD SB	A	BUG+B	Air	33
283	Lysinibacillus sphaericus	GP-ROD SB	A	BUG+B	Air	33
284	Macroccoccus bovicus	GP-COCCUS	A	BUG+B	Air	33
285	Macroccoccus carouelicus	GP-COCCUS	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
286	Macrococcus caseolyticus	GP-COCCUS	A	BUG+B	Air	33
287	Macrococcus equiperdus	GP-COCCUS	A	BUG+B	Air	33
288	Microbacterium arborescens	GP-ROD	A	BUG+B	Air	33
289	Microbacterium barkeri	GP-ROD	A	BUG+B	Air	33
290	Microbacterium dextranolyticum	GP-ROD	A	BUG+B	Air	33
291	Microbacterium esteraromaticum	GP-ROD	A	BUG+B	Air	33
292	Microbacterium flavescens	GP-ROD	C2	BUG+B	Air	33
293	Microbacterium imperiale	GP-ROD	A	BUG+B	Air	33
294	Microbacterium lacticum	GP-ROD	A	BUG+B	Air	33
295	Microbacterium laevaniformans	GP-ROD	A	BUG+B	Air	33
296	Microbacterium liquefaciens	GP-ROD	A	BUG+B	Air	33
297	Microbacterium maritimum	GP-ROD	A	BUG+B	Air	33
298	Microbacterium resistens	GP-ROD	A	BUG+B	Air	33
299	Microbacterium saperdae	GP-ROD	A	BUG+B	Air	33
300	Microbacterium spp. (CDC.A-4)	GP-ROD	A	BUG+B	Air	33
301	Microbacterium spp. (CDC.A-5)	GP-ROD	A	BUG+B	Air	33
302	Microbacterium terregens	GP-ROD	A	BUG+B	Air	33
303	Microbacterium testaceum	GP-ROD	A	BUG+B	Air	33
304	Micrococcus diversus	GP-COCCUS	A	BUG+B	Air	33
305	Micrococcus luteus A	GP-COCCUS	A	BUG+B	Air	33
306	Micrococcus luteus B	GP-COCCUS	A	BUG+B	Air	33
307	Micrococcus lylae	GP-COCCUS	A	BUG+B	Air	33
308	Nocardia abscessus	GP-ROD	A	BUG+B	Air	33
309	Nocardia asteroides	GP-ROD	A	BUG+B	Air	33
310	Nocardia brasiliensis	GP-ROD	A	BUG+B	Air	33
311	Nocardia brevicatena	GP-ROD	A	BUG+B	Air	33
312	Nocardia farcinica	GP-ROD	A	BUG+B	Air	33
313	Nocardia nova	GP-ROD	A	BUG+B	Air	33
314	Nocardia otitidiscaviarum	GP-ROD	B	BUG+B	Air	33
315	Nocardia trasvalensis	GP-ROD	A	BUG+B	Air	33
316	Oerskovia turbata	GP-ROD	A	BUG+B	Air	33
317	Paenibacillus alvei	GP-ROD SB	A	BUG+B	Air	33
318	Paenibacillus amylolyticus	GP-ROD SB	A	BUG+B	Air	33
319	Paenibacillus azotofixans	GP-ROD SB	A	BUG+B	Air	33
320	Paenibacillus glucanolyticus	GP-ROD SB	A	BUG+B	Air	33
321	Paenibacillus larvae	GP-ROD SB	A	BUG+B	Air	33
322	Paenibacillus macerans	GP-ROD SB	A	BUG+B	Air	33
323	Paenibacillus pabuli	GP-ROD SB	A	BUG+B	Air	33
324	Paenibacillus peoriae	GP-ROD SB	A	BUG+B	Air	33
325	Paenibacillus polymyxa	GP-ROD SB	A	BUG+B	Air	33
326	Paenibacillus popilliae	GP-ROD SB	A	BUG+B	Air	33
327	Paenibacillus thiaminolyticus	GP-ROD SB	A	BUG+B	Air	33
328	Paenibacillus validus	GP-ROD SB	A	BUG+B	Air	33
329	Pediococcus acidilactici	GP-COCCUS	C2	BUG+B	6.5%CO2	33-37
330	Pediococcus parvulus	GP-COCCUS	C2	BUG+B	6.5%CO2	33-37
331	Pediococcus pentosaceus	GP-COCCUS	C1	BUG+B	6.5%CO2	33-37
332	Rathayibacter rathayi	GP-ROD	A	BUG	Air	30
333	Rathayibacter tritici	GP-ROD	A	BUG	Air	30
334	Renibacterium salmoninarum	GP-ROD	A	BUG+B	Air	33
335	Rhodococcus australis	GP-ROD	A	BUG+B	Air	30
336	Rhodococcus coprophilus	GP-ROD	A	BUG+B	Air	30
337	Rhodococcus corynebacterioides	GP-ROD	A	BUG+B	Air	30
338	Rhodococcus equi	GP-ROD	A	BUG+B	Air	30
339	Rhodococcus erythropolis	GP-ROD	A	BUG+B	Air	30
340	Rhodococcus fascians	GP-ROD	A	BUG+B	Air	30
341	Rhodococcus globululus	GP-ROD	A	BUG+B	Air	30
342	Rhodococcus pyridinivorans	GP-ROD	A	BUG+B	Air	30
343	Rhodococcus rhodnii	GP-ROD	A	BUG+B	Air	30
344	Rhodococcus rhodochrous	GP-ROD	A	BUG+B	Air	30
345	Rhodococcus ruber	GP-ROD	A	BUG+B	Air	30
346	Rhodococcus wratislaviensis	GP-ROD	A	BUG+B	Air	30
347	Rothia amarae	GP-ROD	A	BUG+B	Air	33
348	Rothia dentocariosa	GP-ROD	A	BUG+B	Air	33
349	Rothia mucilaginosa	GP-ROD	A	BUG+B	Air	33
350	Rothia nasimurium	GP-ROD	A	BUG+B	Air	33
351	Sanguibacter inulinus	GP-ROD	A	BUG+B	Air	33
352	Sanguibacter keddiei	GP-ROD	A	BUG+B	Air	33
353	Sanguibacter suarezi	GP-ROD	A	BUG+B	Air	33
354	Sporolactobacillus laevolacticus	GP-ROD SB	A	BUG+B	Air	33
355	Sporosarcina globispora	GP-ROD SB	A	BUG+B	Air	26
356	Sporosarcina pasteurii	GP-ROD SB	A	BUG+B	Air	33
357	Sporosarcina ureae	GP-ROD SB	A	BUG+B	Air	33

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium		
		Type	Protocol	Medium	Atm	Temp
358	<i>Staphylococcus arlettae</i>	GP-COCCUS	A	BUG+B	Air	33
359	<i>Staphylococcus aureus</i> ss <i>aureus</i>	GP-COCCUS	A	BUG+B	Air	33
360	<i>Staphylococcus auricularis</i>	GP-COCCUS	A	BUG+B	Air	33
361	<i>Staphylococcus capitis</i> ss <i>capitis</i>	GP-COCCUS	A	BUG+B	Air	33
362	<i>Staphylococcus capitis</i> ss <i>ureolyticus</i>	GP-COCCUS	A	BUG+B	Air	33
363	<i>Staphylococcus caprae</i>	GP-COCCUS	A	BUG+B	Air	33
364	<i>Staphylococcus carnosus</i> ss <i>carnosus</i>	GP-COCCUS	A	BUG+B	Air	33
365	<i>Staphylococcus carnosus</i> ss <i>utilis</i>	GP-COCCUS	A	BUG+B	Air	33
366	<i>Staphylococcus chromogenes</i>	GP-COCCUS	A	BUG+B	Air	33
367	<i>Staphylococcus cohnii</i> ss <i>cohnii</i>	GP-COCCUS	A	BUG+B	Air	33
368	<i>Staphylococcus cohnii</i> ss <i>urealyticum</i>	GP-COCCUS	A	BUG+B	Air	33
369	<i>Staphylococcus condimentii</i>	GP-COCCUS	A	BUG+B	Air	33
370	<i>Staphylococcus delphini</i>	GP-COCCUS	A	BUG+B	Air	33
371	<i>Staphylococcus epidermidis</i>	GP-COCCUS	A	BUG+B	Air	33
372	<i>Staphylococcus equorum</i> ss <i>equorum</i>	GP-COCCUS	A	BUG+B	Air	33
373	<i>Staphylococcus felis</i>	GP-COCCUS	A	BUG+B	Air	33
374	<i>Staphylococcus fleurettii</i>	GP-COCCUS	A	BUG+B	Air	33
375	<i>Staphylococcus gallinarum</i>	GP-COCCUS	A	BUG+B	Air	33
376	<i>Staphylococcus haemolyticus</i>	GP-COCCUS	A	BUG+B	Air	33
377	<i>Staphylococcus hominis</i> ss <i>hominis</i>	GP-COCCUS	A	BUG+B	Air	33
378	<i>Staphylococcus hominis</i> ss <i>novobiosepticus</i>	GP-COCCUS	A	BUG+B	Air	33
379	<i>Staphylococcus hyicus</i>	GP-COCCUS	A	BUG+B	Air	33
380	<i>Staphylococcus intermedius</i>	GP-COCCUS	A	BUG+B	Air	33
381	<i>Staphylococcus kloosii</i>	GP-COCCUS	A	BUG+B	Air	33
382	<i>Staphylococcus lentus</i>	GP-COCCUS	A	BUG+B	Air	33
383	<i>Staphylococcus lugdunensis</i>	GP-COCCUS	A	BUG+B	Air	33
384	<i>Staphylococcus lutrae</i>	GP-COCCUS	A	BUG+B	Air	35
385	<i>Staphylococcus muscae</i>	GP-COCCUS	A	BUG+B	Air	33
386	<i>Staphylococcus pasteurii</i>	GP-COCCUS	A	BUG+B	Air	33
387	<i>Staphylococcus piscifermentans</i>	GP-COCCUS	A	BUG+B	Air	33
388	<i>Staphylococcus saprophyticus</i> ss <i>bovis</i>	GP-COCCUS	A	BUG+B	Air	33
389	<i>Staphylococcus saprophyticus</i> ss <i>saprophyticus</i>	GP-COCCUS	A	BUG+B	Air	33
390	<i>Staphylococcus schleiferi</i> ss <i>coagulans</i>	GP-COCCUS	A	BUG+B	Air	33
391	<i>Staphylococcus schleiferi</i> ss <i>schleiferi</i>	GP-COCCUS	A	BUG+B	Air	33
392	<i>Staphylococcus sciuri</i> ss <i>carnaticus</i>	GP-COCCUS	A	BUG+B	Air	33
393	<i>Staphylococcus sciuri</i> ss <i>rodentium</i>	GP-COCCUS	A	BUG+B	Air	33
394	<i>Staphylococcus sciuri</i> ss <i>sciuri</i>	GP-COCCUS	A	BUG+B	Air	33
395	<i>Staphylococcus simulans</i>	GP-COCCUS	A	BUG+B	Air	33
396	<i>Staphylococcus succinus</i> ss <i>succinus</i>	GP-COCCUS	A	BUG+B	Air	33
397	<i>Staphylococcus vitulinus</i>	GP-COCCUS	A	BUG+B	Air	33
398	<i>Staphylococcus warneri</i>	GP-COCCUS	A	BUG+B	Air	33
399	<i>Staphylococcus xylosus</i>	GP-COCCUS	A	BUG+B	Air	33
400	<i>Streptococcus acidominimus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
401	<i>Streptococcus agalactiae</i> (GP B)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
402	<i>Streptococcus alactolyticus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
403	<i>Streptococcus anginosus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
404	<i>Streptococcus australis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
405	<i>Streptococcus bovis</i> (GP D)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
406	<i>Streptococcus canis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
407	<i>Streptococcus constellatus</i> ss <i>constellatus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
408	<i>Streptococcus constellatus</i> ss <i>pharyngis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
409	<i>Streptococcus criceti</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
410	<i>Streptococcus cristatus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
411	<i>Streptococcus didelphis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
412	<i>Streptococcus downei</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
413	<i>Streptococcus dysgalactiae</i> ss <i>dysgalactiae</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
414	<i>Streptococcus dysgalactiae</i> ss <i>equisimilis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
415	<i>Streptococcus entericus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
416	<i>Streptococcus equi</i> ss <i>equi</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
417	<i>Streptococcus equi</i> ss <i>zoepidemicus</i> (GP C)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
418	<i>Streptococcus equinus</i> (GP D)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
419	<i>Streptococcus ferus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
420	<i>Streptococcus gallinaceus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
421	<i>Streptococcus gallolyticus</i> ss <i>gallolyticus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
422	<i>Streptococcus gallolyticus</i> ss <i>macedonicus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
423	<i>Streptococcus gallolyticus</i> ss <i>pasteurianus</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
424	<i>Streptococcus gordonii</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
425	<i>Streptococcus halichoeri</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
426	<i>Streptococcus hyointestinalis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
427	<i>Streptococcus hyovaginalis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
428	<i>Streptococcus infantis</i>	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37

GEN III Database Species and Their Characteristics

	Species Name	Organism		Medium	Medium	
		Type	Protocol		Atm	Temp
429	Streptococcus iniae	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
430	Streptococcus intermedius	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
431	Streptococcus lutetiensis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
432	Streptococcus macacae	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
433	Streptococcus marimammalium	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
434	Streptococcus mitis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
435	Streptococcus mutans/ratti	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
436	Streptococcus oralis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
437	Streptococcus orisratti	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
438	Streptococcus ovis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
439	Streptococcus parasanguinis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
440	Streptococcus parauberis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
441	Streptococcus peroris	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
442	Streptococcus phocae	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
443	Streptococcus pluranimalium	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
444	Streptococcus pneumoniae	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
445	Streptococcus porcinus	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
446	Streptococcus pseudopneumoniae	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
447	Streptococcus pyogenes (GP A)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
448	Streptococcus salivarius	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
449	Streptococcus sanguinis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
450	Streptococcus sinensis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
451	Streptococcus sobrinus	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
452	Streptococcus suis (GP R)	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
453	Streptococcus suis serogroup 1/2	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
454	Streptococcus suis serogroup 4/6	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
455	Streptococcus suis serogroup 5	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
456	Streptococcus suis serogroup 7	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
457	Streptococcus thermophilus	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
458	Streptococcus thoraltensis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
459	Streptococcus uberis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
460	Streptococcus vestibularis	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
461	Tetragenococcus halophilus	GP-COCCUS	C1	BUG+B	6.5% CO2	33-37
462	Tetragenococcus solitarius	GP-COCCUS	C2	BUG+B	6.5% CO2	33-37
463	Tsukamurella inchoensis	GP-ROD	A	BUG+B	Air	30
464	Tsukamurella paurometabola	GP-ROD	B	BUG+B	Air	30
465	Tsukamurella pulmonis	GP-ROD	A	BUG+B	Air	30
466	Tsukamurella tyrosinosolvens	GP-ROD	A	BUG+B	Air	30
467	Turicella otitidis	GP-ROD	A	BUG+B	Air	33
468	Vagococcus fessus	GP-COCCUS	A	BUG+B	6.5% CO2	33-37
469	Vagococcus fluvialis	GP-COCCUS	A	BUG+B	6.5% CO2	33-37
470	Vagococcus lutrae	GP-COCCUS	A	BUG+B	6.5% CO2	33-37
471	Vagococcus salmoninarum	GP-COCCUS	A	BUG+B	Air	26
472	Virgibacillus pantothenicus	GP-ROD SB	A	BUG+B	Air	33
473	Virgibacillus proomii	GP-ROD SB	A	BUG+B	Air	33
474	Weissella confusa	GP-ROD	C1	BUG+B	6.5% CO2	33-37
475	Weissella halotolerans	GP-ROD	C1	BUG+B	6.5% CO2	33-37
476	Weissella minor	GP-ROD	C1	BUG+B	6.5% CO2	33-37
477	Weissella viridescens	GP-ROD	C1	BUG+B	6.5% CO2	33-37
478	Williamsia deligens	GP-ROD	A	BUG+B	Air	33
479	Williamsia maris	GP-ROD	A	BUG+B	Air	26
480	Williamsia muralis	GP-ROD	A	BUG+B	Air	33