
Microbial Diagnostics Exercise



Rasmus Lykke Marvig

Infection Microbiology Group
Department of Systems Biology
Technical University of Denmark

Microbial diagnostics: Identification and characterization of bacteria in a medical scenario


Learning objectives:

A student who has met the objectives of the course will be able to:

- Perform biofilm flow cell setup and operation: Manipulation and monitoring of growing biofilms, antibiotic treatment of biofilms, viability assessment, on-line measurements of oxygen in biofilms, harvesting of biofilms.
- Perform biofilm micro titer tray assays.
- Use advanced optical microscopy techniques: Confocal microscopy and fluorescence microscopy, time-lapse (confocal) microscopy.
- Use atomic force microscopy (AFM) for nano scale imaging.
- Use flow cytometry for single cell analysis.
- Understand chip technology as a tool for gene analysis.
- Perform mathematical image processing using COMSTAT analysis.
- Perform digital image preparation and presentation (Imaris).
- Present biofilm experiment results for an audience.

Genetic typing methods...

Outline of the ArrayTube Chip from Clondiaq Chip Technologies GmbH




P. aeruginosa array of 14 x 14 features

M	C-45	C-46_1	C-47-1_2	SG17M-1	SG17M-8	C-specific-1	M
	C-islid.spec.1	C-islid.spec.4	C-islid.spec.5	C-islid.spec.6			
M	pKL-1	pKL-3	TB-C47-3	TB-C47-3_2	PAPI-1_pili-chap	PAPI-1_lum-bin_pro	
M	pKLC102-unknown	pKLC102-adhesin	pKLC102-metab				
	PyovRec-I	PyovRec-IIa	PyovRec-IIb	PyovRec-III	PyovRec-Fpv_B	LES	
	PA0636	PA0722			PAGI-1-1	PA0980	
	PA0728	PA2185	fla-islid-1	fla-islid-2_orfA	tRNA-Pro-islid	PAPI-2_rec.tr	
	PA2221	PA2835	fla-islid-2_orfI	fla-islid-2_orfJ	tRNA-Pro-islid	PAPI-2-XF1753	
	ampC_7 wt	ampC_7 mut	flhC a	flhC b	exoS2-1	exoU/PAPI-2	
	ampC_4 wt	ampC_4 mut	ampC_5 wt	ampC_5 mut	ampC_6 wt	ampC_6 mut	
	oprI b wt	oprI b mut	ampC_1 wt	ampC_1 mut	ampC_3 wt	ampC_3 mut	
	citS-1 wt	citS-1mut	citS-2 wt	citS-2 mut	oprI a wt	oprI a mut	
	flhC 1 wt	flhC 1 mut	flhC 2 wt	flhC 2 mut	alkB2 wt	alkB2 mut	
M	oprI wt	oprI mut	oprL a wt	oprL a mut	oprL b wt	oprL b mut	M



Gene islands and variable genes
(e.g. pathogenicity factors such as *exoS* and *exoU*)

SNPs within conserved genes

Flow chart of *P. aeruginosa* genotyping with the low-resolution microarray




The AT chip is able to discriminate different clone types of *P. aeruginosa*




Detection of specific SNPs and presence/absence of genes measures the genotype which can be easily compared and used to define clone types

Chip-images



Genotype information




Clone types

 **DK1**

 **DK2**


Motivating Case story – Molecular evolution of
Pseudomonas aeruginosa in cystic fibrosis patients

Epidemiology of cystic fibrosis patients can be monitored over time by genotyping bacterial isolates (each color indicate specific genotypes)





The *P. aeruginosa* DK2 clone type is a successful colonizer of the airways of cystic fibrosis patients

Genetic adaptation to host airways



The DK2 clone has disseminated through a cohort of more than 40 cystic fibrosis patients over a period of 38 years ($\approx 200,000$ bacterial generations)

The genomes of 55 isolates of the DK2 clone type was sequenced to study the evolution of this bacterial pathogen




How predict the genetic relationship of clones?

Building a maximum-parsimonious phylogenetic tree

Genotype of clones

Genetic relationship




?

WT MUT

How predict the genetic relationship of clones?


Building a maximum-parsimonious phylogenetic tree

Genotype of clones




WT MUT


Genetic relationship



The DK2 genome is reduced during long-term infections (up to 8%)
No horizontal acquisition of novel DNA was observed



Case story - Parallel adaptive evolution: Multiple independent mutational events within the domain V of 23s rRNA



Conclusions


Genotyping of specific SNPs and variable genomics regions can be used to identify the clone type of infecting bacteria (low resolution)

Whole-genome sequencing can detect genomic differences between strains of the same clone type (high resolution)

Genotype-phenotype relationships can be established by co-comparison of genotypes and phenotypes

Diagnostics Exercise Mission

Patient suffering from chronic infections in sinuses, jaw and lungs
We need to make a microbial diagnosis of sampled *P. aeruginosa* clones




Exercise questions:


What is the genotype?

How is the antibiogram?

Overview of Microbial Diagnostics Exercise



Questions?



Acknowledgements

Infection Microbiology Group

Søren Molin

Lars Jelsbak

Morten Sommer

Martin Holm Rau

Mette Søndergaard

Trine Markussen

Søren Damkiær

Nicholas Jochumsen

Collaboration/Funding

Copenhagen University Hospital

/ Niels Høiby

/ Helle Johansen

The Lundbeck Foundation

Danish Council for Independent Research