

Practical information

- Experiments and talks
 - Short 10 min intro every day before the first lab session
- Work in the lab
 - Lab coats
 - Waste
 - P2 rules (for tagged *P. aeruginosa*). If anyone is pregnant let us know
- Talks every day
- Presentations at the last day
 - Treat data/results on the way so that you “just” have to put everything together in the end
- The manual
 - Schedule for the different experiments
 - Flowchart of the different experiments
 - List of speakers
 - Laboratory manuals and protocols
 - The flow systems
 - Additional protocols
 - Accompanying papers
- Valuables
- Teams
- The Coffee machine

Working in P2 classified lab

- Only max. 3 groups inside P2 at a time due to limited work space.
- Introduction to the P2 lab is carried out by the supervisor.
- In general:
- Leave your own lab coat, course material, protocols, pencils, etc. outside the P2 lab.
- Wear the lab coats and shoe covers which are provided in the P2 lab.
- Wear protective gloves when working.
- All waste, syringes etc. go into special waste-containers (metal buckets, yellow-cap boxes).
- Wash your hands before leaving the P2 lab.


The exercises

- **5 different biofilm systems will be set up:**
- System 1: Development and architecture of *E. coli* biofilms
- System 2: FISH in biofilms(*P. aeruginosa* and *Acinetobacter* variants)
- System 3: Differentiation, structure development and colistin tolerance in *P. aeruginosa* biofilms
- System 4: *P. aeruginosa* mutants structure development and resistance
- System 5: Free Exercise
- Characterization of *Pseudomonas aeruginosa* from chronic ear infections in dogs
- Microtiter assays of different strains in connection to system 1 and 4
- Additionally:
 - The biofilm protocol
 - The FISH protocol
 - Adhesion assay
 - Conjugation/plasmid transfer
 - MIC determinations
 - Cell sorting
 - Quorum Sensing assays

The Flow Chamber System

**The standard set-up
at DTU for hydro-
dynamic biofilm
development.**

**The flow-cells are
mountable directly
on the Confocal
Microscopes for in
situ investigations
(Zeiss LSM510
or Leica DMRXA)**



Biofilm system 1:

***E. coli* biofilm development, the affect of F pilus**

**Isogenic F plasmid
derivatives:**

***traA*: No pilus synthesis**

***traD*: No DNA transfer**

F⁺

FtraA

FtraD

Images for COMSTAT at 2 time points during the next days

Biofilm system 3:

- Impact of type IV pili-driven motility on structural biofilm formation

P. aeruginosa wild-type Yfp + *P. aeruginosa pilA* Cfp

- Tolerance-development towards Colistin in biofilms

P. aeruginosa wild-type Gfp + Colistin + Propidium Iodide