

**EuroBiofilms Workshop 3 - Molecular diagnostics of
biofilm bacteria**

Microfluidic Systems and Confocal Microscopy.

FLOW CELL SYSTEMS

Paul Stoodley

Lecturer Microbial Tribology

Affiliate Assoc. Prof. Microbiology and Immunology,

Drexel University

Medical College, Pittsburgh Campus, USA

School of Engineering Sciences

University of Southampton



Workshop Objectives

- Introduction to importance of shear and biofilm accumulation
- Identify experimental objectives
- Hydrodynamic parameters relevant to biofilms
- Selection of biofilm flow cells
- Growth systems



Biofilms – Wide Range of Shear

Marine Biofouling

Hulls
Structures

Host tissue


Infective endocarditis
CF lung

Nosocomial

Catheters
Dental water lines

Public Health

Drinking water



Biofilms in Industry

Problems

Product contamination



Public health / spoilage

Membrane plugging

Corrosion

Energy losses in pipelines

SRB oil souring




Experimental Systems for Quantifying Biofilm Attachment, Accumulation and Detachment in Flow



Experimental Design


- Identify Relevant Parameters

Laboratory World



Highly controlled
Investigate factors in isolation
Relevance?


Real World



Highly relevant
Little control
Interpretation?







Imaging Biofilms



1st Choice: Emphasis on Imaging or Sampling

Imaging: 2 basic systems “once through” or
“recirculating”

Sampling: Robbins device, annular reactor
(rototorque)



PS8

Flow Cell / Reactor Design

Hydrodynamics

Magnification range – scales of observation

Reactor wall thickness may limit high definition high power microscopy

Surface of interest – transparent or opaque?

Sampling required?



Experimental System

**Digital time-lapse
image capture**

Camera


**Computer +
framestore board**

Flow cell

**Mixing
chamber**

Heat tape


Microscope




PS10

Flow Cells for Biofilm Research


Microscope and image analysis *in situ*



“Once-through” System




Recirculating System




Flow cell operating requirements

High resolution microscopy – requires flat surface and a transparent window



Flow cell operating requirements

High resolution microscopy – requires flat surface and a transparent window



Channel wall width limited by objective working distance, limits rigidity



Where to get flow cells

BioSurface Technologies (<http://cu.imt.net/~mitbst/flowcell.html>)

Fluxion

Stovall (<http://www.slscience.com/flowcell.html>)

Make your own



national Centre for Advanced Tribology at Southampton (nCATS) University of Southampton
<http://www.soton.ac.uk/nCATS>

PS16

Microscopy

Bright field (phase / interference contrast)

Epi-fluorescent

Reflected (interference contrast)

Confocal (3D fluorescent / transmitted)

SEM/TEM



Stains

DNA - AO, PI, EB, SYBR Green

Activity - CTC

Viability - Live dead

Species identification - FISH (16S RNA)

Gene expression - GFP transcriptional or constitutive reporters

EPS - Lectins

Proteins

Check out Molecular Probes <http://www.probes.com/>



Free Image Analysis Packages

NIH Image (Mac) Image (Mac or PC)

rsb.info.nih.gov/nih-image/

(also lists commercial packages)

Scion Image (PC version of NIH image)

www.scioncorp.com/

ImageTool (PC)

macorb.uthscsa.edu/dig/itdesc.html



What can we measure?

Microscopically

Surface area coverage

Thickness

Length dimensions of biofilm structures

Count objects

Brightness

Rates of change

Molecular tools - GFP expression

Material properties

External monitoring


Pressure drop

Effluent cell count

Biofilm cell count (end point)



Thickness



$$\text{Real thickness} = 1.6 \times \text{apparent thickness}$$



Recirculating System

The residence time (θ) can be controlled by the nutrient flow rate (Q_n) according to $\theta = V/Q_n$, where V = the volume of the mixing chamber *plus* the recirculation loop.

The flow rate in the flow cells can be adjusted independently of the nutrient flow rate so that much higher flow rates can be achieved in the flow cells without using impractical volumes of media.

Q_n can be adjusted so no planktonically dividing population by making the dilution rate (D) greater than the growth rate (μ)




$$D = Q_n/V$$

$$\mu = \ln(2)/Dt$$

PS22

Hydrodynamic characterization



Hydrodynamic Parameters for Closed Channel

Flow rate (Q) ml/min or m³/s


Flow velocity (u_{ave}) = Q/CSA (m/s)

Channel cross sectional area, CSA = hw (m²)

Wetted perimeter, WP = 2w + 2h (m)

Reynolds number = $ul_{ch}\rho/\eta$

$L_{ch} = 4D_h = 4CSA/WP$



η absolute viscosity

ρ density

ν kinematic viscosity = η / ρ





Hydrodynamic Parameters for Closed Channel II

Reynolds number

<500 laminar (low shear)

500 – 1500 transition

> 1500 turbulent (high shear)



η absolute viscosity


ρ density

ν kinematic viscosity = η/ρ


PS25

Pressure drop as fluid moves down a tube due to friction

$$f = \frac{\Delta P \times D_h}{2l_f \rho_w u_{(ave)}^2}$$




Biofilm Flow Cell Hydrodynamics



Friction factor and shear stress

$$f = \frac{16}{\text{Re}} \quad \text{Laminar flow}$$

$$f = 0.0791/\text{Re}^{0.25} \quad \text{Blasius equation for turbulent flow}$$



Shear stress

$$\tau_w = \frac{f \rho_w u_{(ave)}^2}{2} \quad \text{Turbulent flow (empirical)}$$

$$\tau_w = \frac{4\eta u_{(max)}}{D_h} \quad \text{Laminar flow (analytical solution)}$$


Square or circular tube $u_{max} = 2 u_{ave}$

Flat plate $u_{max} = 3/2 u_{ave}$


D_h = channel depth



ATTACHMENT ASSAYS



316 Stainless Steel





Fluoropolymer

P. aeruginosa PAO1


PS30

Growth Assays




PS31

Liquid Flow in Biofilm Channels




$$\tau_w = \frac{\eta du}{dz}$$




du/dz can be found using particle velocimetry with confocal microscopy



Adhesion Assays






Shear induced detachment




Cell Signaling Experiments

- Expression of GFP was only observed in the biofilms after flow was turned off.



P.aeruginosa pMH509
lasB::GFP

5-day old biofilm,
Q=1ml/min



Thank You to BioSurface Technologies for Supplying Flow Cell Materials

Questions?

